



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6044 Cornerstone Court W, Ste E
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet
EZH2 (Y641N) TR-FRET Assay Kit
Catalog #52078
Size: 384 reactions

DESCRIPTION: The *EZH2 (Y641N) TR-FRET Assay Kit* is designed to measure activity of the mutant EZH2 complex [EZH2 (Y641N)/EED/SUZ12/RbAp48/AEBP] in a homogeneous 384 reaction format. This FRET-based assay requires no time-consuming washing steps, making it especially suitable for high throughput screening applications. The *EZH2 (Y641N) TR-FRET Assay Kit* comes in a convenient format, with histone H3 peptide substrate, a Tb-labeled antibody against methylated K27 residue of Histone H3, S-adenosylmethionine, methyltransferase assay buffer, TR-FRET detection buffer, dye-labeled acceptor, and purified EZH2 (Y641N) complex for 384 enzyme reactions. The key to the EZH2 (Y641N) Assay Kit is a highly specific antibody that recognizes methylated Histone H3K27. With this kit, only three simple steps on a microtiter plate are required for methyltransferase detection. First, S-adenosylmethionine is incubated with a sample containing assay buffer and methyltransferase enzyme for 4 hours. Next, antibody to methylated H3K27 is added. Finally, dye-labeled acceptor is added followed by fluorescence detection.

COMPONENTS:

Catalog #	Component	Amount	Storage	<i>Avoid freeze/ thaw cycles!</i>
51028	EZH2 (Y641N)/EED/SUZ12/RbAp48/AEBP2	100 µg	-80°C	
52120	20 µM S-adenosylmethionine	500 µl	-80°C	
52089	Tb-labeled antibody	5 µl	-80°C	
	Biotinylated histone H3 peptide substrate*	500 µl	-80°C	
52170-A	4x HMT Assay Buffer 2A	3 ml	-20°C	
	Dye-labeled acceptor	20 µl	-80°C	
	TR-FRET Detection Buffer	4 ml	-20°C	
Fisher 07-200-330	White, Nonbinding Corning, low volume, microtiter plate	1	Room temp.	

**Resuspend in 500 µl of distilled water.*

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

Fluorescent microplate reader capable of measuring Time Resolved Fluorescence Resonance Energy Transfer (TR-FRET)
Adjustable micropipettor and sterile tips

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6044 Cornerstone Court W, Ste E
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

APPLICATIONS: Great for screening small molecular inhibitors for drug discovery and HTS applications.

STABILITY: At least 6 months from date of receipt when stored as directed.

REFERENCES:

1. Dillon, S.C., *et al.* *Genome Biology* 2005; **6**:227.
2. Morin, R.D., *et al.* *Nat Genet.* 2010, **42**(2):181.

ASSAY PROTOCOL:

All samples and controls should be tested in duplicate.

Step 1:

- 1) Thaw **S-adenosylmethionine** on ice. Upon first thaw, briefly spin tube containing **S-adenosylmethionine** to recover full contents of the tube. Aliquot **S-adenosylmethionine** into single use aliquots. Store remaining **S-adenosylmethionine** in aliquots at -80°C immediately. *Note: S-adenosylmethionine is very sensitive to freeze/thaw cycles. Avoid multiple freeze-thaw cycles.*
- 2) Prepare the master mixture: N wells × (2.5 µl **4x HMT Assay Buffer 2A** + 1 µl **Histone Substrate** + 1 µl **20 µM S-adenosylmethionine**)

	Blank	Substrate Control	Positive Control	Test Inhibitor
4x HMT assay buffer 2A	2.5 µl	2.5 µl	2.5 µl	2.5 µl
20 µM S-adenosylmethionine	1 µl	–	1 µl	1 µl
Histone substrate	1 µl	1 µl	1 µl	1 µl
1x HMT assay buffer 2A	2.5 µl	–	–	–
H ₂ O	–	1 µl	–	–
Test Inhibitor/Activator	–	–	–	3 µl
Inhibitor buffer (no inhibitor)	3 µl	3 µl	3 µl	–
EZH2(Y641N) (100 ng/µl)	–	2.5 µl	2.5 µl	2.5 µl
Total	10 µl	10 µl	10 µl	10 µl

- 3) Add 3.5 µl of master mixture to each well designated for the “Positive Control”, “Test Inhibitor”, and “Blank”. For the “Substrate Control”, add 2.5 µl **4x HMT Assay Buffer 2A** + 1 µl **Histone Substrate** + 1 µl distilled water.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6044 Cornerstone Court W, Ste E
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

- 4) Add 3 μ l of inhibitor solution to each well designated "Test Inhibitor". For the "Positive Control", "Substrate Control" and "Blank", add 3 μ l of the same solution without inhibitor (inhibitor buffer).
- 5) Dilute one part **4x HMT Assay Buffer 2A** with 3 parts distilled water (4-fold dilution) to make **1x HMT assay buffer 2A**. Make only a sufficient quantity needed for the assay; store remaining stock solution in aliquots at -20°C. Add 2.5 μ l of **1x HMT assay buffer 2A** to the well designated "Blank".
- 6) Thaw **EZH2 (Y641N) enzyme** on ice. Upon first thaw, briefly spin tube containing enzyme to recover full contents of the tube. Aliquot **EZH2 (Y641N) enzyme** into single use aliquots. Store remaining undiluted enzyme in aliquots at -80°C immediately. *Note: **EZH2 (Y641N) enzyme** is very sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted enzyme.*
- 7) Dilute **EZH2(Y641N) enzyme** in **1x HMT assay buffer 2A** at 100 ng/ μ l (250 ng/2.5 μ l). Keep diluted enzyme on ice until use. Discard any unused diluted enzyme after use.
- 8) Initiate reaction by adding 2.5 μ l of diluted **EZH2(Y641N)** prepared as described above to wells labeled "Positive Control", "Test Inhibitor", and "Substrate Control". Incubate at room temperature for 4 hours.
- 9) Cover the plate with a plate sealer if necessary.

Step 2:

- 1) Thaw **TR-FRET Detection Buffer** on ice.
- 2) Dilute **Tb-labeled antibody** 400-fold with **TR-FRET Detection Buffer**.
- 3) Add 5 μ l to every well. Incubate 30 minutes at room temperature with slow shaking.

Step 3:

- 1) Dilute **Dye-labeled acceptor** 100-fold with **TR-FRET Detection Buffer**.
- 2) Add 5 μ l to every well. Incubate for 30 min. at room temperature with slow shaking.
(Alternatively, dilute Tb-labeled antibody (1:800) and Dye-labeled acceptor (1:200) with TR-FRET Detection Buffer in one step. Add 10 μ l of Antibody/Acceptor mixture per well and incubate 1 hour.)

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6044 Cornerstone Court W, Ste E
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

- 3) Read the fluorescent intensity using a microtiter-plate reader capable of measuring TR-FRET.

Instrument Settings

Reading Mode	Time Resolved
Excitation Wavelength	340±20 nm
Emission Wavelength	620±10 nm
Lag Time	60 µs
Integration Time	500 µs
Excitation Wavelength	340±20 nm
Emission Wavelength	665±10 nm
Lag Time	60 µs
Integration Time	500 µs

CALCULATING RESULTS:

Two sequential measurements should be conducted. Tb-donor emission should be measured at 620 nm followed by dye-acceptor emission at 665 nm. Data analysis is performed using the TR-FRET ratio (665 nm emission/620 nm emission).

When percentage activity is calculated, the FRET value from the negative control (Blank or Substrate Control) can be set as zero percent activity and the FRET value from the positive control can be set as one hundred percent activity.

$$\% \text{ Activity} = \frac{\text{FRET}_s - \text{FRET}_{\text{neg}}}{\text{FRET}_p - \text{FRET}_{\text{neg}}} \times 100\%$$

Where FRET_s = Sample FRET, FRET_{neg} = negative control FRET, and FRET_p = Positive control FRET.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**


Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6044 Cornerstone Court W, Ste E
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Example of Assay Results:



EZH2(Y641N) enzyme activity, measured using the *EZH2 (Y641N) TR-FRET Assay Kit*, BPS Bioscience #52078. Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at info@bpsbioscience.com.

RELATED PRODUCTS

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
EZH2 (Y641F)/EED/SUZ12/RbAp48/AEBP2	51017	20 µg
EZH2 (Y641C)/EED/SUZ12/RbAp48/AEBP2	51029	20 µg
EZH2 (Y641N)/EED/SUZ12/RbAp48/AEBP2	51028	20 µg
EZH2 (Y641S)/EED/SUZ12/RbAp48/AEBP2	51013	20 µg
EZH2 (Y641H)/EED/SUZ12/RbAp48/AEBP2	51011	20 µg
EZH2/EED/SUZ12/RbAp48/AEBP2	51004	50 µg
EZH2 Chemiluminescent Assay Kit	52009L	96 rxns
EZH2 (Y641F) Chemiluminescent Assay Kit	52075	96 rxns
EZH2 (Y641N) Chemiluminescent Assay Kit	52076	96 rxns
EZH2 (A677G) Chemiluminescent Assay Kit	52077	96 rxns
EZH2 Homogeneous Assay Kit	52059	384 rxns

Note: Dye-labeled acceptor is a product of Cisbio Bioassays.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com