



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet

MAPK/ERK Signaling Pathway SRE Reporter – HEK293 Cell Line Catalog #: 60406

Description

The MAPK/ERK signaling pathway is a major participant in the regulation of cell growth and differentiation. It can be activated by various extracellular stimuli including mitogens, growth factors, and cytokines. Upon stimulation, MEK1/2 phosphorylates and activates ERK1/2. The activated ERK translocates to the nucleus where it phosphorylates and activates transcription factors. The TCFs (Ternary Complex Factors), including Elk1, are among the best-characterized transcription factor substrates of ERK. When phosphorylated by ERK, Elk1 forms a complex with Serum Response Factor (SRF) and binds to Serum Response Element (SRE), resulting in the expression of numerous mitogen-inducible genes.

The SRE Reporter – HEK293 cell line contains a firefly luciferase gene under the control of SRE responsive elements stably integrated into HEK293 cells, resulting in an ERK pathway-responsive reporter cell line. This cell line is validated for the response to the stimulation of EGF or serum and to the treatment of inhibitors of ERK signaling pathway.

Application

- Monitor MAPK/ERK signaling pathway activity and SRF-mediated activity.
- Screen for activators or inhibitors of the MAPK/ERK signaling pathway.

Format

Each vial contains $\sim 1.5 \times 10^6$ cells in 1 ml of 10% DMSO.

Mycoplasma Testing

The cell line has been screened using the PCR-based VenorGeM Mycoplasma Detection kit (Sigma-Aldrich) to confirm the absence of Mycoplasma species.

Storage

Immediately upon receipt, store in liquid nitrogen.

General culture conditions

Thaw Medium 1 (BPS Cat. #60187): MEM medium (Hyclone #SH30024.01) supplemented with 10% FBS (Invitrogen #26140-079), 1% non-essential amino acids (Hyclone #SH30238.01), 1 mM Na pyruvate (Hyclone #SH30239.01), 1% Penicillin/Streptomycin (Hyclone SV30010.01)

Growth Medium 1B (BPS Cat. #79531): Thaw Medium 1 (BPS Cat. #60187) plus 400 μ g/ml of Geneticin (Invitrogen #11811031)

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Cells should be grown at 37°C with 5% CO₂ using Growth Medium 1B.

It may be necessary to adjust the percentage of CO₂ in the incubator depending on the NaHCO₃ level in the basal medium.

To thaw the cells, it is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath, transfer to a tube containing 10 ml of Thaw Medium 1 (**no Geneticin**), spin down cells, and resuspend cells in pre-warmed Thaw Medium 1 (**no Geneticin**). Transfer resuspended cells to a T25 flask and culture at 37°C in a CO₂ incubator. At first passage, switch to Growth Medium 1B (**contains Geneticin**). Cells should be split before they reach complete confluence.

To passage the cells, rinse cells with phosphate buffered saline (PBS), detach cells from culture vessel with 0.05% Trypsin/EDTA, add Growth Medium 1B and transfer to a tube. Spin down cells, resuspend cells, and seed appropriate aliquots of cell suspension into new culture vessels. Subcultivation ratio: 1:10 to 1:20 weekly.

To freeze down the cells, rinse cells with phosphate buffered saline (PBS), and detach cells from culture vessel with Trypsin/EDTA. Add Growth Medium 1B and transfer to a tube, spin down cells, and resuspend in freezing medium (10% DMSO + 90% FBS). Place at -80°C overnight and place in liquid nitrogen the next day. Alternatively, vials may be placed directly in liquid nitrogen.

Functional Validation and Assay Performance

The following assays are designed for a 96-well plate format. To perform assay in different tissue culture formats, the cell number and reagent volumes should be scaled appropriately.

Materials Required but Not Supplied

- Thaw Medium 1 (BPS Cat. #60187)
- Growth Medium 1B (BPS Cat. #79531)
- Recombinant human EGF (BPS bioscience # 90201-1)
- U0126 (BPS Bioscience #27012): inhibitor of ERK pathway (MEK inhibitor). Prepare stock solution of U0126 in DMSO.
- Assay Medium 1B (BPS Bioscience #79617): MEM medium (Hyclone #SH30024.01) supplemented with 0.5% FBS (Thermo Fisher, #26140079), 1% non-essential amino acids (Hyclone #SH30238.01), 1 mM Na pyruvate (Hyclone #SH30239.01), 1% Penicillin/Streptomycin (Hyclone SV30010.01).
- 96-well tissue culture plate or 96-well tissue culture-treated white clear-bottom assay plate
- ONE-Step™ Luciferase Assay System (BPS, Cat. #60690)
- Luminometer

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

A. Response of SRE Reporter – HEK293 cells to EGF or serum

1. Harvest SRE Reporter – HEK293 cells from culture in Growth Medium 1B and seed cells into the white clear-bottom 96-well microplate at a density of ~ 30,000 cells per well in 100 µl of Thaw Medium 1.
2. Incubate cells at 37°C in a CO₂ incubator overnight.
3. The next day, carefully remove the medium from wells. Add 45 µl Assay Medium 1B to wells.

Assay Medium 1B (BPS Bioscience #79617): MEM medium (Hyclone #SH30024.01) supplemented with 0.5% FBS (Thermo Fisher, #26140079), 1% non-essential amino acids (Hyclone #SH30238.01), 1 mM Na pyruvate (Hyclone #SH30239.01), 1% Penicillin/Streptomycin (Hyclone SV30010.01).

4. Incubate the plate at 37°C in a CO₂ incubator for 18 to 24 hours.
5. The next day, add 5 µl of FBS or threefold serial dilution of human EGF in Assay Medium 1B to stimulated wells.
Add 5 µl of Assay Medium 1B to the unstimulated control wells.
Add 50 µl of Assay Medium 1B to cell-free control wells (for determining background luminescence).
Set up each treatment in at least triplicate.
6. Incubate the plate at 37°C in a CO₂ incubator for ~ 6 hours.
7. Perform luciferase assay using ONE-Step™ Luciferase Assay System according to the protocol provided: Add 100 µl of ONE-Step™ Luciferase reagent per well and rock at room temperature for ~15 minutes. Measure luminescence using a luminometer.
If using luciferase reagents from other vendors, follow the manufacturer's assay protocol.
8. Data Analysis: Subtract the average background luminescence (cell-free control wells) from the luminescence reading of all wells.
The fold induction of SRE luciferase reporter expression = background-subtracted luminescence of stimulated well / average background-subtracted luminescence of unstimulated control wells

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.
To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com

Figure 1. EGF or serum induced the expression of SRE reporter in SRE Reporter – HEK293.





Figure 2. Dose response of SRE Reporter – HEK293 cells to EGF. The results were shown as fold induction of SRE luciferase reporter expression.



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.
To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
 San Diego, CA 92121
 Tel: 1.858.829.3082
 Fax: 1.858.481.8694
 Email: info@bpsbioscience.com

B. Inhibition of EGF-induced reporter activity by inhibitor of ERK signaling pathway in SRE Reporter – HEK293 cells

1. Harvest SRE Reporter – HEK293 cells from culture in Growth Medium 1B and seed cells into the white clear-bottom 96-well microplate at a density of 30,000 cells per well in 100 µl of Thaw Medium 1.
2. Incubate cells at 37°C in a CO₂ incubator for overnight.
3. The next day, dilute the inhibitor (U0126) stock in Assay Medium 1B. Carefully remove the medium from wells and add 45 µl of diluted inhibitor in Assay Medium 1B to the wells. The final concentration of DMSO in Assay Medium 1B can be up to 0.5%.
 Add 45 µl of Assay Medium 1B with same concentration of DMSO without inhibitor to inhibitor control wells.
 Add 45 µl of Assay Medium 1B with DMSO to cell-free control wells (for determining background luminescence).
4. Incubate the plate at 37°C in a CO₂ incubator for 18-24 hours.
5. The next day, add 5 µl of diluted human EGF in Assay Medium 1B to stimulated wells (with and without inhibitor) (final [EGF] = 10 ng/ml).
 Add 5 µl of Assay Medium 1B to the unstimulated control wells (cells without inhibitor and EGF treatment for determining the basal activity).
 Add 5 µl of Assay Medium 1B to cell-free control wells.
 Set up each treatment in at least triplicate.

Treatment Reference Guide

	Stimulated Wells		Unstimulated Control Wells	Cell-free Control Wells
	With inhibitor	Without inhibitor (control well)		
Step 3	45 µl diluted inhibitor in Assay Medium 1B	45 µl Assay Medium 1B with DMSO only	45 µl Assay Medium 1B with DMSO only	45 µl Assay Medium 1B with DMSO only
Step 5	5 µl EGF in Assay Medium 1B (final [EGF] = 10 ng/ml)	5 µl EGF in Assay Medium 1B (final [EGF] = 10 ng/ml)	5 µl Assay Medium 1B	5 µl Assay Medium 1B

6. Incubate the plate at 37°C in a CO₂ incubator for ~6 hours.
7. Perform luciferase assay using ONE-Step™ Luciferase Assay System according to the protocol provided: Add 100 µl of ONE-Step™ Luciferase reagent per well and rock at room temperature for ~15 minutes. Measure luminescence using a luminometer.


OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.
 To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
 Or you can Email us at: info@bpsbioscience.com
 Please visit our website at: www.bpsbioscience.com

If using other luciferase reagents from other vendors follow the manufacturer's assay protocol.


8. Data Analysis: Obtain background-subtracted luminescence by subtracting average background luminescence (cell-free control wells) from luminescence reading of all wells.

Figure 3. Inhibition of EGF-induced reporter activity by ERK pathway inhibitor in SRE Reporter – HEK293 cells

3a. U0126 blocked EGF-induced SRE reporter activity.



3b. U0126 inhibition dose response curve. The results were shown as percentage of luminescence. The background-subtracted luminescence of cells stimulated with EGF in the absence of U0126 was set at 100%.



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.
 To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
 Or you can Email us at: info@bpsbioscience.com
 Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

References

1. Wong, K.K. (2009) Recent developments in anti-cancer agents targeting the Ras/Raf/MEK/ERK pathway. *Recent Pat. Anticancer Drug Discov.* **4(1)**:28-35.
2. Treisman, R. (1992) The serum response element. *Trends Biochem. Sci.* **17(10)**: 423-426.

Related Products

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
SRE Reporter Kit (MAPK/ERK Signaling Pathway)	60511	500 reactions
FOXO Reporter Kit (PI3K/AKT Pathway)	60643	500 reactions
U0126	27012	5 mg
EGF, human	90201-1	100 µg
EGF, human	90201-2	500 µg
EGF, mouse	90200-1	100 µg
EGF, mouse	90200-2	500 µg
ERK1	40055	10 µg
ERK2	40299	10 µg
ONE-Step™ Luciferase Assay System	60690-1	10 ml
ONE-Step™ Luciferase Assay System	60690-2	100 ml

License Disclosure

Purchase of this cell line grants you with a 10-year license to use this cell line in your immediate laboratory, for research use only. This license does not permit you to share, distribute, sell, sublicense, or otherwise make the cell line available for use to other laboratories, departments, research institutions, hospitals, universities, or biotech companies. The license does not permit use of this cell line in humans or for therapeutic or drug use. The license does not permit modification of the cell line in any way. Inappropriate use or distribution of this cell line will result in revocation of the license and result in an immediate cease of sales and distribution of BPS products to your laboratory. BPS does not warrant the suitability of the cell line for any particular use, and does not accept any liability in connection with the handling or use of the cell line. Modifications of this cell line, transfer to another facility, or commercial use of the cells may require a separate license and additional fees; contact sales@bpsbioscience.com for details. Publications using this cell line should reference BPS Bioscience, Inc., San Diego.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com