



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Description

Adeno-Associated Virus serotype 9 (AAV9) is one of the most promising serotypes for gene therapy applications. AAV9 transduces a wide range of tissue types, including cardiac and skeletal muscle, liver, pancreas, and eye tissue. AAV8 and AAV9 have recently been used to correct disease-causing mutations and improve muscle function in mouse models of Duchenne muscular dystrophy. AAV9 has significantly lower seroprevalence in the human population than other AAV serotypes, making it a desirable candidate for therapeutic applications.

These AAV particles constitutively express the firefly (*Photinus pyralis*) luciferase and enhanced Green Fluorescent Protein (eGFP) genes connected via a T2A linker, under the control of a CMV promoter. The T2A self-cleaving peptide (derived from *Thosea asigna* virus 2A) leads to the efficient cleavage of the transcript and expression of luciferase and eGFP as two separate proteins.

Application(s)

- Use as a positive control for transduction
- Optimize transduction assays and track protein expression over time

Serotype

Wild-type AAV Serotype 9

Formulation

AAV9 was produced in HEK293-AAV cells and is supplied in PBS-MK (PBS Magnesium-Potassium) buffer containing 0.01% Pluronic F68.

Purification

The purity of the AAV particles was confirmed to be greater than 90% by staining with One-Step Lumitein™ UV Protein Gel Stain (Biotium #21005-1L). The purity will vary with each lot; the exact value will be provided with each shipment.




Figure 1. Purified AAV9 Luciferase-eGFP particles.

Staining of a 4-20% SDS-PAGE gel. The protein ladder is in lane 1, and 2×10^9 GC (genome copy number) of AAV9 is shown in lane 2. AAV viral proteins VP1, VP2, and VP3 are labeled.

Titer

Two vials (50 μ l x 2) of AAV at a titer $\geq 1 \times 10^{12}$ TU/ml. The titer is determined by qPCR and will vary with each lot; the exact value will be provided with each shipment.

Storage

AAV is shipped with dry ice. For long-term storage, it is recommended to store AAV at -80°C . Avoid repeated freeze-thaw cycles. Titers can drop significantly with each freeze-thaw cycle.

Biosafety

Recombinant AAV is inherently replication-deficient and not known to cause any human diseases. Additionally, following transduction, AAV vectors exist episomally and do not integrate into or disrupt the host cell's genome. AAV requires the use of a Biosafety Level 1 facility. BPS Bioscience recommends following all local, federal, state, and institutional regulations and using all the appropriate safety precautions.


Validation Data

Figure 2. Transduction of HEK293 cells using AAV9 Luciferase-eGFP particles.

1×10^5 cells/well were transduced in a 6-well plate with AAV9 Luciferase-eGFP at an MOI of 2×10^4 . After 72 hours of transduction, eGFP expression in the target cells was observed under a fluorescence microscope. eGFP expression was stable over time and still observed 30 days after transduction.




Figure 3. Luciferase activity of HEK293 cells transduced by AAV9 Luciferase-eGFP- particles. 1 x 10⁵ cells/well were transduced in a 6-well plate with AAV9 Luciferase-eGFP at an MOI of 2 x 10⁴. After 72 hours of transduction, transduced cells or parental HEK293 cells were seeded in a 96-well plate at a density of 2 x 10⁴ cells/well, and luciferase activity was measured using the ONE-Step™ luciferase assay system (BPS Bioscience #60690).

Troubleshooting Guide

Visit bpsbioscience.com/lentivirus-faq for detailed troubleshooting instructions. For all further questions, please email support@bpsbioscience.com.

Related Products

<i>Products</i>	<i>Catalog #</i>	<i>Size</i>
AAV1 ZsGreen	78443	50 µl x 2
AAV2 ZsGreen	78444	50 µl x 2
AAV5 ZsGreen	78447	50 µl x 2
AAV8 ZsGreen	78449	50 µl x 2
AAV9 ZsGreen	78450	50 µl x 2
AAV-DJ ZsGreen	78442	50 µl x 2
AAV5-Luciferase-eGFP	78645	50 µl x 2
AAV6-Luciferase-eGFP	78646	50 µl x 2
AAV8-Luciferase-eGFP	78647	50 µl x 2
AAV-DJ Luciferase-eGFP	78460	50 µl x 2