



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Description

The TNKS1 (also known as PARP5a) Homogeneous Assay Kit is designed to measure TNKS1 activity for screening and profiling applications. TNKS1 is known to catalyze the NAD-dependent addition of poly(ADP-ribose) to its substrate. The TNKS1 Homogeneous Assay Kit comes in a convenient AlphaLISA® format, with biotinylated histone substrate, recombinant antibody (ADP-ribose binding reagent 1), PP-01 assay buffer, and purified TNKS1 for 384 enzyme reactions. The key to the TNKS1 Homogeneous Assay Kit is a highly specific antibody that recognizes PARylated substrate. With this kit, only three simple steps are required for TNKS1 reactions. First, a sample containing TNKS1 enzyme is incubated with biotinylated substrate and NAD for one hour. Next, acceptor beads and primary antibody are added, then donor beads, followed by reading the Alpha-counts.




Figure 1: TNKS1 (PARP5a) Homogenous Assay Kit schematic

## Application(s)

- Screen for molecules that inhibit PARP1 activity.

## Supplied Materials

Catalog #	Name	Amount	Storage	
80504	TNKS1 (PARP5a)*	2 µg	-80°C	<b>Avoid multiple freeze/thaw cycles</b>
	5X PP-01 assay buffer	2 x 1 ml	-80°C	
	Biotinylated histone substrate**	500 rxns	-80°C	
78311	ADP-ribose binding reagent 1	10 µl	-80°C	
	NAD+ (750 µM)	400 µl	-80°C	
52301	4X Detection buffer 1	2 ml	-80°C	

\*The initial concentration of enzyme is lot-specific and will be indicated on the tube containing the protein.

\*\*Reconstitute in 500 µl H<sub>2</sub>O.

**Materials Required but Not Supplied**

Name	Catalog #
AlphaLISA anti-rabbit IgG acceptor beads	Perkin Elmer #AL104C
AlphaScreen Streptavidin-conjugated donor beads	Perkin Elmer #6760002S
Optiplate -384	Perkin Elmer #6007290
AlphaScreen microplate reader	
0.5 M DTT	

**Storage Conditions**

This assay kit will perform optimally for up to 6 months from date of receipt when the materials are stored as directed. **Avoid multiple freeze/thaw cycles!**

**Safety**

This product is for research purposes only and not for human or therapeutic use. This product should be considered hazardous and is harmful by inhalation, in contact with skin, eyes, clothing, and if swallowed. If contact occurs, wash thoroughly.

**Contraindications**

The TNKS1 Homogenous Assay Kit is compatible with up to 1% final DMSO concentration.

We recommend preparing the inhibitor in no higher than 5% DMSO solution in buffer and using 3  $\mu$ l per well.

Avoid green and blue dyes that absorb light in the AlphaScreen signal emission range (520-620 nm), such as Trypan Blue. Avoid the use of the potent singlet oxygen quenchers such as sodium azide ( $\text{NaN}_3$ ) or metal ions ( $\text{Fe}^{2+}$ ,  $\text{Fe}^{3+}$ ,  $\text{Cu}^{2+}$ ,  $\text{Zn}^{2+}$  and  $\text{Ni}^{2+}$ ). The presence of >1% RPMI 1640 culture medium leads to a signal reduction due to the presence of excess biotin and iron in this medium. MEM, which lacks these components, does not affect AlphaScreen assays.

**Assay Protocol**

- All samples and controls should be performed in triplicates
- The assay should include a “Blank”, a “Positive control”, and a “Negative control”

*Reagent Preparation*

1. Reconstitute the biotinylated histone substrate in 500  $\mu$ l of distilled water. Aliquot biotinylated histone substrate into single use aliquots. Store aliquots at  $-80^\circ\text{C}$ .
2. Thaw 5X PP-01 assay buffer. Add 0.5 M DTT to 5X PP-01 assay buffer for a final concentration of 10 mM DTT.

*Note: DTT should be added just before use. Prepare only enough DTT-containing buffer as required for the assay. Store the remaining assay buffer at  $-20^\circ\text{C}$ .*

3. Prepare 1X PP-01 assay buffer by adding 1 part of 5X PP-01 assay buffer (**with DTT**) to 4 parts distilled water.
4. Prepare the compound solution.
  - a. If the compound is dissolved in DMSO, make a 100-fold higher concentration of the compound in DMSO than the highest concentration you want to test in the assay. Then dilute 20-fold in water.

*Note: To run an  $IC_{50}$  or test lower concentrations of the compound, prepare serial dilutions using water containing 5% DMSO, so the final concentration of DMSO will be 1% in all samples.*

- b. If the compound is soluble in water, prepare a solution of the compound in water that is 5-fold higher than the final assay concentration.

#### Reaction Setup

1. Make Master Mix: N wells  $\times$  (1  $\mu$ l of **Biotinylated histone substrate** + 1  $\mu$ l of **NAD<sup>+</sup> (750  $\mu$ M)** + 2  $\mu$ l of **5x stock assay buffer and** + 3  $\mu$ l of **water**).
2. Add 7  $\mu$ l of Master Mix to each well.
3. To the wells designated as “Blank”, add 5  $\mu$ l of **1X PP-01 assay buffer** and 3  $\mu$ l of **diluent solution without inhibitor** (for example DMSO 5%).

Component	Blank
Master Mix	7 $\mu$ l
Diluent solution* (no inhibitor)	3 $\mu$ l
1X PP-01 assay buffer	5 $\mu$ l
<b>Total</b>	<b>15 <math>\mu</math>l</b>

*\*The diluent solution contains the water with the same concentration of solvent (e.g., DMSO) as the test compound solution.*

4. Thaw TNKS1 on ice. Briefly spin the tube containing enzyme to recover the full content of the tube.
5. Calculate the amount of enzyme required for your assay and dilute TNKS1 in 1X PP-01 assay buffer (**with DTT**) to 1 ng/ $\mu$ l (the final amount of TNKS1 in the assay will be 5 ng/rxn). Keep the diluted enzyme on ice until use. Discard any unused diluted enzyme after use.

*Note: Aliquot the remaining undiluted TNKS1 enzyme into single use aliquots and store at  $-80^{\circ}\text{C}$ .*



*TNKS1 is sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted enzyme.*

6. Add 3  $\mu$ l of inhibitor solution to each well designated “Test Inhibitor”. To “Positive Control” add 3  $\mu$ l of the diluent solution without inhibitor.

7. Initiate the reaction by adding 5  $\mu\text{l}$  of diluted **TNKS1** to the wells labeled “Positive Control” and “Test Inhibitor”.

Component	Test Sample	Positive Control
Master Mix	7 $\mu\text{l}$	7 $\mu\text{l}$
Test compound	3 $\mu\text{l}$	–
Diluent solution* (no inhibitor)	–	3 $\mu\text{l}$
TNKS1 (1 ng/ $\mu\text{l}$ )	5 $\mu\text{l}$	5 $\mu\text{l}$
<b>Total</b>	<b>15 <math>\mu\text{l}</math></b>	<b>15 <math>\mu\text{l}</math></b>

\*The diluent contains the same concentration of solvent (e.g., DMSO) as the test compound solution.

8. Once the reaction mixture (15  $\mu\text{l}$ ) is added to 384-well plate, incubate at room temperature for 1 hour with slow shaking.

#### Reaction detection and reading results



**Protect your samples from direct exposure to light. Photobleaching will occur.**

1. Prepare the 1X Detection buffer by adding 1-part 4X Detection buffer to 3 parts distilled water.
2. Prepare a mixture in 1X Detection buffer containing:
  - a. anti-Rabbit Acceptor beads (Perkin Elmer #AL104C) diluted 250-fold
  - b. ADP-ribose binding reagent 1 diluted 400-fold

*Example: For 100 wells, prepare 1 ml of 1X Detection buffer and add 4  $\mu\text{l}$  of anti-Rabbit Acceptor beads as well as 2.5  $\mu\text{l}$  of ADP-ribose binding reagent 1.*

3. Add 10  $\mu\text{l}$  of mixture from step 2 to each well then briefly shake plate.
4. Dilute the Streptavidin-conjugated donor beads (PE #6760002S) 250-fold with 1x Detection buffer and add 10  $\mu\text{l}$  per well.

*Example: For 100 wells, prepare 1 ml of 1X Detection buffer and add 4  $\mu\text{l}$  Streptavidin-conjugated donor beads.*

5. Incubate for at least 5-15 min at room temperature and read Alpha-counts.

## Example Results




Figure 2: Inhibition of TNKS1 (PARP5a) activity by XAV939 (Cayman Chemical), measured using the TNKS1 Homogenous Assay Kit (BPS Bioscience #78489). Data shown is representative. For lot-specific information, please contact BPS Bioscience, Inc. at [support@bpsbioscience.com](mailto:support@bpsbioscience.com).

### Troubleshooting Guide

Visit [bpsbioscience.com/assay-kits-faq](https://bpsbioscience.com/assay-kits-faq) for detailed troubleshooting instructions. For all further questions, please email [support@bpsbioscience.com](mailto:support@bpsbioscience.com)

**Related Products**

Products	Catalog #	Size
PARP1 Homogeneous Assay Kit	78438	384 rxns
PARP1 Chemiluminescent Assay Kit	80551	96 rxns
PARP1 Chemiluminescent Assay Kit (384-well)	80569	384 rxns
PARP1 Colorimetric Assay Kit	80580	96 rxns
PARP1, FLAG-Avi-tag Recombinant	80521	20 µg
PARP2 Homogeneous Assay Kit	80702	384 rxns
PARP2 Chemiluminescent Assay Kit	80552	96 rxns
PARP2 Colorimetric Assay Kit	80581	96 rxns
PARP3 Homogeneous Assay Kit	78491	384 rxns
PARP3 Chemiluminescent Assay Kit	80553	96 rxns
PARP6 Chemiluminescent Assay Kit	80556	96 rxns
PARP10 Chemiluminescent Assay Kit	80560	96 rxns
PARP10, FLAG-Strep-Tag Recombinant	80522	10 µg
PARP11 Homogeneous Assay Kit	78492	384 rxns
PARP11 Chemiluminescent Assay Kit	80561	96 rxns
PARP14 Chemiluminescent Assay Kit	80568	96 rxns
PARP15 Chemiluminescent Assay Kit	80567	96 rxns
PARPtrap™ Assay Kit for PARP1	80584	96 rxns
PARPtrap™ Assay Kit for PARP2	78296	96 rxns
PARPtrap™ Combo Assay Kit for PARP1 and PARP2	78317	384 rxns
TNKS1 (PARP5A) Chemiluminescent Assay Kit	78405	96 rxns
TNKS1 Histone Ribosylation Assay Kit (Biotin-labeled NAD+)	80579	384 rxns
TNKS1 Histone Ribosylation Colorimetric Assay Kit	80582	96 rxns
TNKS2 Homogeneous Assay Kit	78490	384 rxns
TNKS2 (PARP5B) Chemiluminescent Assay Kit	78406	96 rxns
TNKS2 Histone Ribosylation Assay Kit (Biotin-labeled NAD+)	80572	384 rxns
TNKS2 Histone Ribosylation Colorimetric Assay Kit	80583	96 rxns
Set of PARP Inhibitors (8 x 50 µl)	78318	8 x 50 µl
Streptavidin-HRP (For PARP & Cytokine Assay Kits)	80611	100 µl