



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Description

The TR β -GAL4 Luciferase Reporter HEK293 Cell Line is a HEK293 cell line expressing firefly luciferase under the control of the GAL4 upstream activation sequence (UAS) with constitutive expression of human thyroid receptor β ligand binding domain (TR β LBD, amino acids 173-461) fused to the DNA binding domain (DBD) of GAL4 (GAL4 DBD). This system allows specific detection of thyroid hormone-induced activation of the thyroid receptor β with low cross-reactivity from other nuclear receptors. This cell line has been validated by stimulation with triiodothyronine (T-3).

Background

Thyroid hormones play an important role in growth, development, and metabolism. These are mediated by two different thyroid hormone receptor (TR) isoforms, TR α and TR β . Thyroid receptors are transcriptional factors that control various genes by interacting with specific co-activators, co-repressors and DNA sequences. In humans, TR α is the predominant form of the thyroid receptor in the heart, brain, and bone while TR β is mainly expressed in the liver, kidney and brain. Interestingly, it has been shown that the patients with non-alcoholic steatohepatitis (NASH) display lower TR β expression in the liver, and more importantly, TR β agonist treatment decreased liver steatosis and circulating lipids as well as showed metabolic benefits in preclinical models of diabetes and obesity.

Application(s)

Screen for TR β agonists.

Materials Provided

| Components | Format |
|-------------------------|---|
| 2 vials of frozen cells | Each vial contains $\geq 1 \times 10^6$ cells in 1 ml of Cell Freezing Medium (BPS Bioscience #79796) |

Parental Cell Line

HEK293, Human Embryonic Kidney, epithelial-like cells, adherent.

Mycoplasma Testing

The cell line has been screened to confirm the absence of Mycoplasma species.

Materials Required but Not Supplied



These materials are not supplied with this cell line but are necessary for cell culture and cellular assays. BPS Bioscience reagents systems are validated and optimized for use with this cell line and are highly recommended for best results. Media components are provided in the Media Formulations section.

Materials Required for Cell Culture

| Name | Ordering Information |
|------------------|---------------------------------------|
| Thaw Medium 1 | BPS Bioscience #60187 |
| Growth Medium 1M | BPS Bioscience #79723 |

Materials Required for Cellular Assay

| Name | Ordering Information |
|--|---------------------------------------|
| 3, 3', 5-Triiodo-L-thyronine (T-3) | Cayman #16028 |
| Assay Medium 6B | BPS Bioscience #82202 |
| 96-well tissue culture treated white clear-bottom assay plate | Corning #3610 |
| 384-well tissue culture treated white clear-bottom assay plate | PerkinElmer (#6007680) |
| ONE-Step™ Luciferase Assay System | BPS Bioscience #60690 |
| Luminometer | |

Storage Conditions

Cells will arrive upon dry ice and should immediately be thawed or stored in liquid nitrogen upon receipt. Do not use a -80°C freezer for long term storage. Contact technical support at support@bpsbioscience.com if the cells are not frozen in dry ice upon arrival.

Media Formulations

For best results, the use of validated and optimized media from BPS Bioscience is *highly recommended*. Note that using similar but not BPS Bioscience validated reagents can result in suboptimal performance.



Note: Thaw Media does *not* contain selective antibiotics. However, Growth Media *does* contain selective antibiotics, which are used to maintain selective pressure on the cell population expressing the gene of interest.

Cells should be grown at 37°C with 5% CO₂. BPS Bioscience's cell lines are stable for at least 10 passages when grown under proper conditions.

Media Required for Cell Culture

Thaw Medium 1 (BPS Bioscience #60187):

MEM medium supplemented with 10% FBS, 1% non-essential amino acids, 1 mM Na pyruvate, and 1% Penicillin/Streptomycin

Growth Medium 1M (BPS Bioscience #79723):

MEM medium supplemented with 10% FBS, 1% non-essential amino acids, 1 mM Na pyruvate, 1% Penicillin/Streptomycin, 400 µg/ml of G418, and 0.5 µg/ml of Puromycin.

Media Required for Functional Cellular Assay

Assay Medium:

DMEM (without phenol red) medium supplemented with 2% charcoal/dextran treated FBS, 1% Gluta-Max™, and 1% Penicillin/Streptomycin.

Cell Culture Protocol

Note: HEK293 cells are derived from human material and thus the use of adequate safety precautions is recommended.

Cell Thawing

1. Retrieve a cell vial from liquid nitrogen storage. Keep on dry ice until ready to thaw.

2. When ready to thaw, swirl the vial of frozen cells for approximately 60 seconds in a 37°C water bath. Once cells are thawed (it may be slightly faster or slower than 60 seconds), quickly transfer the entire content of the vial to an empty 50 ml conical tube.

Note: Leaving the cells in the water bath at 37°C for too long will result in rapid loss of viability.

3. Using a 10 ml serological pipette, slowly add 10 ml of pre-warmed Thaw Medium 1 to the conical tube containing the cells. Thaw Medium 1 should be added dropwise while gently rocking the conical tube to permit gentle mixing and avoid osmotic shock.
4. Immediately spin down the cells at 300 x *g* for 5 minutes, remove the medium and resuspend the cells in 5 ml of pre-warmed Thaw Medium 1.
5. Transfer the resuspended cells to a T25 flask or T75 flask and incubate at 37°C in a 5% CO₂ incubator.
6. After 24 hours of culture, check for cell attachment and viability. Change medium to fresh Thaw Medium 1 and continue growing in a 5% CO₂ incubator at 37°C until the cells are ready to passage.
7. Cells should be passaged before they are fully confluent. At first passage and subsequent passages, use Growth Medium 1M.

Cell Passage

1. Aspirate the medium, wash the cells with phosphate buffered saline (PBS) without Ca²⁺/Mg²⁺, and detach the cells from the culture vessel with 0.05% Trypsin/EDTA.
2. Once the cells have detached, add Growth Medium 1M and transfer to a tube.
3. Spin down cells at 300 x *g* for 5 minutes, remove the medium and resuspend the cells in Growth Medium 1M.
4. Seed into new culture vessels at the recommended sub-cultivation ratio of 1:10 once a week.

Cell Freezing

1. Aspirate the medium, wash the cells with phosphate buffered saline (PBS) without Ca²⁺/Mg²⁺ and detach the cells from the culture vessel with 0.05% Trypsin/EDTA.
2. Once the cells have detached, add Growth Medium 1M and count the cells.
3. Spin down the cells at 300 x *g* for 5 minutes, remove the medium and resuspend the cells in 4°C Cell Freezing Medium (BPS Bioscience #79796) at 1~2 x 10⁶ cells/ml.
4. Dispense 1 ml of cell suspension into each cryogenic vial. Place the vials in an insulated container for slow cooling and store at -80°C overnight.
5. Transfer the vials to liquid nitrogen the next day for long term storage.



Note: It is recommended to expand the cells and freeze down at least 10 vials of cells at an early passage for future use.

Validation Data

- The following assay was designed for a 96-well/384-well format. To perform the assay in different tissue culture formats, the cell number and reagent volume should be scaled appropriately.
- All conditions should be performed in triplicate.
- The assay should include “Stimulated Cells”, “Background Control” and “Unstimulated Control” conditions.

A. Dose Response of TRβ-GAL4 Luciferase Reporter HEK293 Cell Line to T-3 (96-well)

1. Harvest TRβ-GAL4 Luciferase Reporter HEK293 cells from culture in Growth Medium 1M and seed cells at a density of 20,000 ~ 40,000 cells per well in 90 μl of Assay Medium 6B into a white clear-bottom 96-well microplate. Leave empty wells as cell-free control wells (“Background Control”).
2. Prepare a serial dilution of T-3 in Assay Medium 6B at 10x the final testing concentrations (10 μl/well).
3. Add 10 μl of diluted T-3 to the “Stimulated Cells” wells.
4. Add 10 μl of Assay Medium 6B to the “Unstimulated Control” wells.
5. Add 100 μl of Assay Medium 6B to “Background Control” wells (cell-free wells).
6. Incubate at 37°C with 5% CO₂ for 16 ~ 24 hours.
7. Add 100 μl of ONE-Step™ Luciferase reagent per well.
8. Incubate at Room Temperature (RT) for ~10 minutes.
9. Measure luminescence using a luminometer.
10. The “Background Control” luminescence value should be subtracted from all readings.
11. Data Analysis: Subtract the average background luminescence from the luminescence reading of all other wells. The fold induction of luciferase reporter expression is the average background-subtracted luminescence of stimulated wells divided by the average background-subtracted luminescence of unstimulated control wells.

$$\text{Fold induction} = \frac{\text{Luminescence of Stimulated Wells} - \text{avg. background}}{\text{Avg. Luminescence of Unstimulated Wells} - \text{avg. background}}$$




Figure 1. Dose response curve of TR β -GAL4 Luciferase Reporter HEK293 Cell Line to T-3 (96-well format).

TR β -GAL4 Luciferase Reporter HEK293 cells were treated with increasing concentrations of T-3 in a 96-well plate format. Luciferase activity was measured using ONE-Step™ Luciferase Assay System. The results are shown as fold induction of luciferase reporter expression in relation to the activity of cells without treatment (unstimulated control) and represent four independently performed experiments.

B. Dose Response of TR β -GAL4 Luciferase Reporter HEK293 Cell Line to T-3 (384-well)

1. Harvest TR β -GAL4 Luciferase Reporter HEK293 cells from culture in Growth Medium 1M and seed cells at a density of ~8,000 cells per well in 25 μ l of Assay Medium 6B into a white clear-bottom 96-well microplate. Leave empty wells as cell-free control wells ("Background Control").
2. Prepare a serial dilution of T-3 in Assay Medium 6B at 2x the final testing concentrations (25 μ l/well).
3. Add 25 μ l of diluted T-3 to the "Stimulated Cells" wells.
4. Add 25 μ l of Assay Medium 6B to the "Unstimulated Control" wells.
5. Add 50 μ l of Assay Medium 6B to "Background Control" wells (cell-free wells).
6. Incubate at 37°C with 5% CO₂ for 16 ~ 24 hours.
7. Add 50 μ l of ONE-Step™ Luciferase reagent per well.
8. Incubate at Room Temperature (RT) for ~10 minutes.
9. Measure luminescence using a luminometer.
10. The "Background Control" luminescence value should be subtracted from all readings.
11. Data Analysis: Subtract the average background luminescence from the luminescence reading of all other wells. The fold induction of luciferase reporter expression is the average background-subtracted luminescence of stimulated wells divided by the average background-subtracted luminescence of unstimulated control wells.

$$\text{Fold induction} = \frac{\text{Luminescence of Stimulated Wells} - \text{avg. background}}{\text{Avg. Luminescence of Unstimulated Wells} - \text{avg. background}}$$




Figure 2. Dose response curve of TR β -GAL4 Luciferase Reporter HEK293 Cell Line to T-3 (384-well format).

TR β -GAL4 Luciferase Reporter HEK293 cells were treated with increasing concentrations of T-3 in a 384-well plate. Luciferase activity was measured using the ONE-Step™ Luciferase Assay System. The results are shown as fold induction of luciferase reporter expression in relation to the activity of cells without treatment (unstimulated control) and represent a comparison of two different plate formats.

License Disclosure

Visit bpsbioscience.com/license for the label license and other key information about this product.

Troubleshooting Guide

Visit bpsbioscience.com/cell-line-faq for detailed troubleshooting instructions. For all further questions, please email support@bpsbioscience.com.

References:

Paguio A, *et al.*, 2010 *Curr Chem Genomics*. 4: 43-49.

Related Products

| Products | Catalog # | Size |
|---|-----------|-----------------|
| Transfection Collection™: GAL4 Transient Pack Glucocorticoid Receptor Pathway | 79265 | 100 reactions |
| GAL4 Luciferase Reporter HEK293 Cell Line | 60656 | 2 vials |
| GAL4 Reporter Kit (Glucocorticoid Receptor Pathway) | 60522 | 2 vials |
| GAL4 DBR-GR Lentivirus | 78632 | 500 μ l x 2 |
| GR-GAL Luciferase Reporter Jurkat Cell Line | 78632 | 2 vials |

Version 11172023