



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.829.3082  
**Fax:** 1.858.481.8694  
**Email:** [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

## **Data Sheet**

### **SLC5A5 – HEK293 Recombinant Cell line**

**Catalog #:90333**

#### **Description**

Stable recombinant HEK293 cell line expressing human SLC5A5 (UniProt # Q92911) with a C-terminal FLAG-tag. SLC5A5 is a sodium-dependent iodide transporter and is a member of the sodium glucose co-transporter family. The expression of SLC5A5 is confirmed by Western blotting. This cell line can be used to evaluate agents that may affect iodide transport.

#### **Background**

Iodine transport is essential for the production of thyroid hormone (T4). Sodium Iodide Symporter SLC5A5 (NIS) is a sodium-dependent transporter of iodide into thyroid follicular cells and other peripheral tissues to support T4 synthesis [1]. Lack of iodide transport can result in hypothyroidism and can manifest itself in several disease states like goiter formation and weight gain. SLC5A5 also plays an important role in hyperthyroidism (an excessive production of T4). Expression and activity of SCL5A5 is altered in several forms of cancer. Inhibitors of SLC5A5 have potential therapeutic benefit for the regulation of iodine overload induced by amiodarone or other thyromimetics that can result in hyperthyroidism [2].

#### **Applications**

- Drug compound screening
- Functional assays

#### **Format**

Each vial contains  $1.5 \times 10^6$  cells in 1 ml of 10% DMSO

#### **Host cell**

HEK293 cells

#### **Recommended Storage**

Immediately upon receipt, store in liquid nitrogen.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.829.3082  
**Fax:** 1.858.481.8694  
**Email:** [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

## Functional Validation and Assay Performance

The following assays are designed for a 96-well format. To perform the assay in different tissue culture formats, the cell number and reagent volumes should be scaled appropriately.

### Culture Conditions

**Thaw Medium 1 (BPS Cat. #60187):** MEM medium (Hyclone #SH30024.01) supplemented with 10% FBS (Invitrogen #26140-079), 1% non-essential amino acids (Hyclone #SH30238.01), 1 mM Na pyruvate (Hyclone #SH30239.01), 1% Penicillin/Streptomycin (Hyclone SV30010.01)

**Growth Medium 1B (BPS Cat. #79531):** Thaw Medium 1 (BPS Cat. #60187) plus 400 µg/ml of Geneticin (Life Technologies #11811031).

Cells should be grown at 37°C with 5% CO<sub>2</sub> using Growth Medium 1B.

It is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath, transfer to a tube containing 10 ml of Thaw Medium 1 (**no Geneticin**), spin down cells, and resuspend cells in pre-warmed Thaw Medium 1 (**no Geneticin**). Transfer resuspended cells to a T25 flask and culture in a CO<sub>2</sub> incubator at 37°C. At first passage, switch to Growth Medium 1B. Cells should be split before they reach complete confluence.

To passage the cells, rinse cells with phosphate buffered saline (PBS), detach cells from culture vessel using 0.05% Trypsin/EDTA, and add Growth Medium 1B. Transfer to a tube, spin down cells, resuspend cells and seed appropriate aliquots of cell suspension into new culture vessels. Subcultivation ration: 1:10 to 1:20 weekly.

### Materials Required but Not Supplied

- Assay Medium: Opti-MEM I (Life technologies #31985-062), 0.5% FBS, 1% non-essential amino acids, 1 mM Na pyruvate, and 1% Pen/Strep
- 96-well tissue culture plate coated with poly-D-lysine  
[We used 0.1 mg/ml poly-D-lysine (Sigma #P6407) in H<sub>2</sub>O for 1 hour at RT.]
- Hanks Buffered Salt Solution (HBSS, Hyclone #SH30588.01)
- Reagents for the Ce/As reaction to determine iodide concentration as reported previously [3]
  - 4x Ammonium cerium (IV) sulfate solution (4x Ce, 42 mM)
  - 4x Sodium arsenite solution (4x As, 96 mM)
- Iodide standards (0- 1 µM NaI, Sigma #383112)
- Plate reader for measuring absorbance at 420 nm

### Transport activity of – HEK293 stably expressing SLC5A5

1. Harvest SLC5A5 – HEK293 cells from culture in Growth Medium 1B and seed cells at a density of ~ 25,000 cells per well into a white clear-bottom 96-well microplate in

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.829.3082  
**Fax:** 1.858.481.8694  
**Email:** [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

100  $\mu$ l of Assay Medium. Make sure to leave two columns of the plate (16 wells total) empty for the iodide standards.

2. Incubate cells at 37°C with 5% CO<sub>2</sub> overnight.
3. The next day, carefully remove the medium from the wells. Wash cells once with Hanks Buffered Salt Solution (HBSS, Hyclone #SH30588.01) and then add 80  $\mu$ l HBSS to each well.
4. Pretreat cells with vehicle, ITB5 (Chembridge #5461570) or any test inhibitors by adding 10  $\mu$ l of 10x stock solutions to each well and incubating for 30 minutes at RT. Ideally, it is best to keep the final DMSO concentrations at 0.1% or less.
5. Begin the uptake assay by adding 10  $\mu$ l of 100  $\mu$ M NaI (Sigma #383112) diluted in HBSS (final NaI concentration of 10  $\mu$ M) and incubate for 1h at RT.
6. Uptake should then be terminated by carefully washing twice with 200  $\mu$ l cold HBSS, making sure not to disrupt the cell layer.
7. Accumulated iodide can be measured using the Ce/As method as previously described [3].

#### **Protocol for Ce/As Method**

- Dilute 4x Ce and 4x As solutions to 1x with H<sub>2</sub>O immediately before reaction
  - After the final wash with HBSS (step 6), add 100  $\mu$ l H<sub>2</sub>O to wells.
    - Prepare and add 100  $\mu$ l iodide standards (S0- S7) to empty columns in duplicate (0, 100, 250, 400, 550, 700, 850, 1000 nM respectively).
  - Add 100  $\mu$ l 1x Ce solution to each well
  - Add 100  $\mu$ l 1x As solution to each well
  - Allow reaction to continue for 20 min at RT
  - Measure absorbance at 420 nm (A<sub>420</sub>)
8. Data Analysis: log [A<sub>420</sub>] for samples was normalized to the iodide standards to determine pmol iodide accumulated per well. Data was analyzed with GraphPad Prism.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.


To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)


Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)

### Sample Data


**Figure 1. Concentration dependent uptake of Iodide by SLC5A5 - HEK293 cells.** WT or SLC5A5 – HEK293 cells were incubated with the indicated concentrations of NaI (0-100  $\mu$ M) for 1h at RT. Iodide uptake was measured using the Ce/As protocol described above and normalized to freshly made iodide standards.




**Figure 2. Time dependent uptake of Iodide by SLC5A5 - HEK293 cells.** WT or SLC5A5 – HEK293 cells were incubated with 5  $\mu$ M NaI for the indicated times (0-240 min) at RT. Additions of NaI were staggered to allow termination of the uptake reaction at the 240 min. timepoint. Iodide uptake was measured using the Ce/As protocol described above and normalized to freshly made iodide standards.



**Figure 3. ITB5 dependent inhibition of Iodide uptake by SLC5A5 - HEK293 cells.** WT or SLC5A5 – HEK293 cells were pre with the indicated concentrations of the known SLC5A5 inhibitor ITB5 (0-30  $\mu$ M, Chembridge # 5461570) for 30 minutes at RT. Subsequently, 5  $\mu$ M NaI (final concentration) was added to the reaction and incubated for 1h at RT. Iodide uptake was measured using the Ce/As protocol described above and normalized to freshly made iodide standards.



**Figure 4. Western Blot of the SLC5A5 expressing cells.** Western Blot of WT HEK293 cells (lane 1) and SLC5A5 - HEK293 cells (lane 2) stained with rabbit anti-FLAG (Sigma # F7425) M: molecular weight marker.





6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.829.3082  
**Fax:** 1.858.481.8694  
**Email:** [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

## Vector and Sequence

MEAVETGERPTFGAWDYGVFALMLLVSTGIGLWVGLARGGQRSAEDFFTGGRRRLAALPVGLS  
LSASFMSAVQVLGVPSEAYRYGLKFLWMCLGQLLNSVLTALLFMPVYRGLTSTYEYLEMRF  
SRAVRLCGTLQYIVATMLYTGIVYAPALILNQVTGLDIWASLLSTGIICTFYTAVGGMKAVVWTD  
VFQVVVMLSGFWVVLARGVMLVGGPRQVLTLAQNHSRINLMDFNPDPRSRYTFWTFVVGGL  
VWLSMYGVNQAQVQRYVACRTEKQAKLALLINQVGLFLIVSSAACCGIVMFVYTDGCDPLLLGR  
ISAPDQYMPLLVLDIFEDLPGVPLFLACAYSGTLSTASTSINAMA AVTVEDLIKPRRLSLAPRKL  
VIISKGLSLIYGSACTVAALSSLLGGGVLQGSFTVMGVISGPLLGAFILGMFLPACNTPGVLGGL  
GAGLALS L WVALGATLYPPSEQTMRVLPSSAARCVALS VNASGLLDPALLPANDSSRAPSSGM  
DASRPALADSFYAISYLYYGALGTLTTVLCGALISCLTGPTKRSTLAPGLLWWDLARQTASVAP  
KEEVAILDNDLVKGPEELPTGNKKPPGFLPTNEDRLFFLGQKELEGAGSWTPCVGHGGRDQ  
QETNLAADYKDDDDK

## Stability

The cell line has been demonstrated to be stable for at least seven continuous passages. For optimal results, it is recommended to use the cells prior to the 7<sup>th</sup> passage. Upon receipt, amplify the cells in culture and make several frozen aliquots for future use.

## References

1. Dohán et al., (2003) The sodium/iodide Symporter (NIS): characterization, regulation, and medical significance. *Endocrine Reviews* **24(1)**:48-77.
2. Lindenthal et al., (2009) Characterization of NIS inhibitors. *Journal of Endocrinology* **200**: 357–365.
3. Waltz et al., (2010) A nonradioactive iodide uptake assay for sodium iodide symporter function. *Analytical Biochemistry* **396(1)**: 91-5.

## License Disclosure

Purchase of this cell line grants you with a 10-year license to use this cell line in your immediate laboratory, for research use only. This license does not permit you to share, distribute, sell, sublicense, or otherwise make the cell line available for use to other laboratories, departments, research institutions, hospitals, universities, or biotech companies. The license does not permit use of this cell line in humans or for therapeutic or drug use. The license does not permit modification of the cell line in any way. Inappropriate use or distribution of this cell line will result in revocation of the license and result in an immediate cease of sales and distribution of BPS products to your laboratory. BPS does not warrant the suitability of the cell line for any particular use, and does not accept any liability in connection with the handling or use of the cell line. Modifications of this cell line, transfer to another facility, or commercial use of the cells may require a separate license and additional fees; contact [sales@bpsbioscience.com](mailto:sales@bpsbioscience.com) for details. Publications using this cell line should reference BPS Bioscience, Inc., San Diego.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)