



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC Handels GmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



# Product Information



## 8-Isoprostane Affinity Column

Item No. 401111

### IMPORTANT

Water used to prepare all EIA reagents and buffers must be deionized and free of trace organic contaminants ('UltraPure'). Use activated carbon filter cartridges or other organic scavengers. Glass distilled water (even if double distilled), HPLC-grade water, and sterile water (for injections) are not adequate for EIA. UltraPure Water is available for purchase from Cayman (Item No. 400000).

### Product Description

The column contains 0.5 ml affinity sorbent (mouse anti-8-isoprostane antibody covalently bound to Sepharose 4B) and has a binding capacity of 5 ng 8-isoprostane (10 ng 8-isoprostane/ml sorbent). This column will be stable for at least two years when stored at 4°C. Storing the columns frozen is not recommended. Be certain that the column is stored in an upright position.

There are two general uses of this affinity column. The first, and most common, is for rapid purification of 8-isoprostane from biological samples for subsequent analysis for 8-isoprostane by EIA. The second use is for removal of 8-isoprostane from biological samples, particularly when 8-isoprostane interferes with assaying a different target molecule.

*NOTE: This affinity columns may be used to isolate free 8-isoprostane from samples. If membrane-bound or total 8-isoprostane measurement is desired, please see the 8-Isoprostane EIA Kit booklet (Item No. 516351).*

The columns are intended for single use only. We do not recommend reuse.

### Sample Preparation

All samples must be free of particulates and precipitates to avoid plugging the column. This can be achieved either by filtration or by centrifugation. All samples must be at approximately neutral pH (6.5-7.5). Urine samples should be centrifuged briefly to remove sediment and may be applied directly to the column. Plasma samples should be centrifuged briefly and diluted 1:5 with Column Buffer (see below).

### Reagent Preparation

Prepare the following reagents for use with the affinity sorbent.

#### 1. Eicosanoid Affinity Column Buffer

Prepare a 0.1 M phosphate buffer solution by combining 13.3 g potassium phosphate (dibasic) 3.22 g potassium phosphate (monobasic) 0.5 g sodium azide, and 29.2 g sodium chloride. Dilute to a total volume of 1.0 liter with UltraPure Water. The pH of this buffer will be 7.4. This buffer may be purchased as a 5X concentrated buffer (Item No. 400220).

#### 2. Eicosanoid Affinity Column Elution Solution

Prepare a solution containing 95% absolute ethanol and 5% UltraPure Water. This solution may be purchased (Item No. 400230).

#### 3. EIA Buffer

Combine 13.3 g potassium phosphate (dibasic) 3.2 g potassium phosphate (monobasic) 23.4 g sodium chloride, 370 mg tetrasodium EDTA, 100 mg sodium azide, and 1 g bovine serum albumin (Sigma A7030 or equivalent) in a final volume of 1,000 ml. Stir at room temperature until completely dissolved. The pH will be 7.0-7.4. This buffer may be purchased as a 10X concentrated buffer (Item No. 400060).

**WARNING: THIS PRODUCT IS FOR LABORATORY RESEARCH ONLY: NOT FOR ADMINISTRATION TO HUMANS. NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

#### SAFETY DATA

This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. This information contains some, but not all, of the information required for the safe and proper use of this material. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

Cayman Chemical Company makes **no warranty or guarantee** of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman warrants only to the original customer that the material will meet our specifications at the time of delivery.

Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have any obligation or liability, whether in tort (including negligence) or in contract, for any direct, indirect, incidental or consequential damages, even if Cayman is informed about their possible existence.

This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees.

Buyer's exclusive remedy and Cayman's sole liability hereunder shall be limited to a refund of the purchase price, or at Cayman's option, the replacement, at no cost to Buyer, of all material that does not meet our specifications.

Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

For further details, please refer to our Warranty and Limitation of Remedy located on our website and in our catalog.

Copyright Cayman Chemical Company, 06/12/2014

### Cayman Chemical

#### Mailing address

1180 E. Ellsworth Road  
Ann Arbor, MI  
48108 USA

#### Phone

(800) 364-9897  
(734) 971-3335

#### Fax

(734) 971-3640

#### E-Mail

custserv@caymanchem.com

#### Web

www.caymanchem.com

# Product Information



## Protocol for Purification or Removal of 8-Isoprostane


1. Remove the top cap from the column and then gently remove the bottom plug. Be certain to remove the top cap first to avoid air bubbles being drawn into the settled sorbent.
2. Allow the storage solution to pass through the sorbent and then wash the column with 2 ml of Column Buffer. Repeat wash.
3. Apply the sample to the column and allow the entire sample to pass through the sorbent.
4. Wash the column with 2 ml Column Buffer, followed by 2 ml UltraPure Water. Allow all the water to pass through the sorbent. Discard both of these washes.
5. Place a test tube capable of holding at least 5 ml under the column. Elute the 8-isoprostane from the column by adding 2 ml Elution Solution and allowing it to pass through the sorbent. Repeat with an additional 2 ml Elution Solution. If the analysis cannot be performed immediately, store the sample in the Elution Solution at -80°C; it will be stable for at least one year.
6. Evaporate the Elution Solution to dryness either by vacuum centrifugation or under a stream of dry nitrogen. It is important that all the organic solvent be removed as even trace quantities may adversely affect an immunoassay.
7. Immediately dissolve the 8-isoprostane in buffer or solvent appropriate for your application. If you are assaying the sample with one of our 8-Isoprostane EIA Kits (Item No. 516351 or Item No. 516360), dissolve the sample in EIA Buffer. The amount of EIA Buffer depends on the original sample volume and the expected concentration of 8-isoprostane in the sample. A dilute sample may be concentrated by dissolving the residue in a smaller volume of EIA Buffer than the original sample volume.

## Recovery

Urine and plasma recoveries average >90%.

## Sample Data

The data shown in the figure below was generated by purifying 8-isoprostane from urine or plasma using the affinity columns. Urine and plasma were spiked with the indicated amounts of 8-isoprostane and purifications were performed for each concentration. Samples were analyzed by EIA.



## Related Product

For a list of related products please visit: [www.caymanchem.com/catalog/401111](http://www.caymanchem.com/catalog/401111)

## Cayman Chemical

**Mailing address**  
1180 E. Ellsworth Road  
Ann Arbor, MI  
48108 USA

**Phone**  
(800) 364-9897  
(734) 971-3335

**Fax**  
(734) 971-3640

**E-Mail**  
[custserv@caymanchem.com](mailto:custserv@caymanchem.com)

**Web**  
[www.caymanchem.com](http://www.caymanchem.com)