Duplica^{RealTime} Next FACTOR V H1299R Kit

REV. EER041032_IFU_REV.02D_ENITA

REF : EER041032-32 tests

Instructions For Use

INTENDED USE

Duplica^{RealTime} Next FACTOR V H1299R kit is an *in vitro* nucleic acid amplification test for the detection of the nucleotide substitution A4070G of the blood coagulation Factor V gene in human genomic DNA, extracted from peripheral whole blood samples collected in EDTA.

INTRODUCTION

Thrombophilia is a multigenic disorder caused by inherited and acquired defects. Individuals with multiple defects have a highly increased risk of development of Deep Venous Thrombosis (DVT). The A4070G (FV1299 His-Arg) polymorphism of the Factor V gene (FV; this allele is a part of R2 haplotype) in exon 13 was recently shown to influence circulating FV levels and to contribute to the Activated Protein C (APC) resistance phenotype. Heterozygous of both FV G1691A and FV A4070G substitutions is associated with a significant increasing risk of developing DVT (3-4 times if compared with the single FV R506Q substitution).

PRINCIPLE OF THE TEST

Duplica^{RealTime} Next FACTOR V H1299R kit is designed to identify the nucleotide substitution A4070G in the human Factor V gene. Two reaction mixes are provided for the amplification:

- **AMPLIFICATION MIX**, containing Hot Start Taq DNA polymerase, nucleotides, MgCl₂ and buffer.
- **OLIGO MIX**, containing primers and fluorogenic probes.

Duplica^{RealTime} Next FACTOR V H1299R kit is based on specific recognition and amplification of target sequences by PCR, and the simultaneous detection of the accumulation of PCR amplification products by fluorescent DNA probes. In particular, the probe designed to detect the Wild Type allele carries the fluorophore FAM (6-carboxy-fluorescein) at the 5' end, while the probe detecting the Mutated allele, is labelled with the fluorophore HEX (hexa-chloro-fluorescein). Both the probes have a non fluorescent black quencher at the 3' end. If excited, the whole probe does not emit fluorescence, since the proximity of the quencher to the reporter prevents the emission of the fluorescence from the reporter (quenching effect).

REAGENTS PROVIDED

Each kit contains enough reagents to perform 32 tests when used in 4 analytical sessions with 5 **samples**, 1 **Wild Type control (Control 1, C1)**, 1 **Mutated control (Control 2, C2)** and one **Reaction Blank (BM)** each.

Kit Components

<i>Reagent</i>	<i>Color Code</i>	<i>Storage (range, °C)</i>	<i>Volume (µl)</i>	<i>Quantity (tubes)</i>
Oligo Mix (OM)*	Green Cap	-22÷-18	400	1
Amplification Mix (AM)	Blue Cap	-22÷-18	400	1
Control 1 (C1, Wild Type)	Red Cap	-22÷-18	>50	1
Control 2 (C2, Mutated)	Yellow Cap	-22÷-18	>50	1
Reaction Blank (BM)	White Cap	-22÷-18	>50	1

***protect the tube from direct light**

STORAGE AND HANDLING

All reagents must be stored at **-22÷-18°C** and can be used until the expiry date printed on the labels. Do not freeze and thaw the products more than six times.

MATERIAL REQUIRED BUT NOT PROVIDED

- Extraction kit for DNA purification (refer to the specific handbook's section)
- Optical tubes or microplate for Real Time PCR
- Disposable powder-free gloves and laboratory coat
- Variable volumes pipettes (5-20 µl, 20-200 µl and 100-1000 µl)
- Disposable RNase/DNase-free tips with aerosol barriers
- Tube racks
- Desktop centrifuge
- PCR box
- Refrigerator
- Deep-freezer
- Thermal cycler for Real Time PCR

The kit has been optimized to be used on DX®/CFX (Bio-Rad) and Applied Biosystems® 7500 (Life Technologies™) thermal cyclers. Other makes and models should be fully tested and evaluated for optimal performance by the user before reporting results.

The equipments should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure an optimal performance.

PRECAUTIONS AND WARNINGS

- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood
- If required, Clonit offers the necessary technical support for the correct use of the kit.
- Carefully read this Instruction for Use before using the kit
- Do not use the reagents after the expiry date
- Thaw and carefully mix the reagents of the kit before use
- Do not mix the reagents from different lots of the product
- Use calibrated and regularly checked pipettes and instrumentation only
- Use dedicated laboratory equipments. Change gloves frequently
- Periodically wipe the working area with 0,5% hypochlorite
- Use powder-free gloves. Do not leave fingerprints on optical caps. Do not write on caps as this may cause an interference with fluorescent detection
- Materials containing or potentially containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood
- In case of damaged package, contact the technical support before using the kit
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet (MSDS, available on request)
- The kit reagents, individual protective equipments, used materials, biological samples and test residuals must be disposed in accordance with local regulations
- Patient Drug treatment may interfere with the final result of the molecular biology analysis

OPERATING PROCEDURE

a) DNA purification

For genomic DNA purification Clonit recommends to use:

-Duplica Blood DNA kit (ref. EDI002250) and Duplica NA Body Fluid kit (ref. EDI004200) with Duplica® PREP Automatic Extractor (ref. EDI001) for automated purification.

For manual extraction/purification:

-Fasst DNA Releaser (ref. EMR057050) for **fresh** peripheral blood (i.e. stored for up to 24 hours at 2÷8°C).
-Spin DNA Purification Kit (ref. EMR061050) for **frozen** blood or stored at 2÷8°C for more than 24 hours.

Other extraction reagents and methods should be fully tested and evaluated for optimal performance by the user before reporting results.

Attention! Use blood samples in EDTA anti-coagulant solution.

b) Thermal cycler Setup

Refer to the specific handbook of the equipment used to set the thermal profile indicated in the **Thermal Profile Table**. We recommend to switch on the instrument and to set the thermal profile before preparing the reaction mix.

Attention! The Thermal Profile is platform-specific: be sure to set the correct one.

Bio-Rad DX®/CFX platform

- Select *Create a new Experiment*
- In the *Protocol* section select *Create New* to set the thermal profile
- In *Plate Editor*: set FAM (Wild Type genotype, Allele 1) and HEX (Mutated genotype, Allele 2) as *Fluorophores* and define the sample name

Thermal Profile Table (specific for DX®/CFX)

Time	Temperature	Cycles
10 min	95°C	1
20 sec	95°C	35
60 sec	60°C	Fluorescence Acquisition

Applied Biosystems® 7500 platform

- Select *Create a new Experiment*
- In the *Experiment Properties* select *Quantitation-Standard curve* as experiment type and *Taqman® Reagents* as reagent type
- Select the Run Mode (*7500* or *7500 Fast*) and the Ramp Speed (*Standard*)
- In the *Setup* menu select *Run Method* to set the thermal profile
- In the *Setup* menu select *Plate Setup-Define Targets and Samples* to select Samples and assign VIC target to MUT (Mutated genotype) and FAM target to WT (Wild-Type genotype). Leave the default value (*NFQ-MGB*) as quencher type
- In *Define Samples* insert the sample name
- Select the used wells in *Plate Setup-Assign Targets and Samples*
- Select *None* as passive reference

Thermal Profile Table (specific for Applied Biosystems® 7500)

Time	Temperature	Cycles
10 min	95°C	1
20 sec	95°C	35
60 sec	62°C	Fluorescence Acquisition

c) Preparation of PCR mix

The total reaction volume is **25 µl**.

For each experiment prepare a PCR mix for the **2 controls (C1 and C2)**, 1 **Reaction Blank (BM)** and **n+1** samples. The reagents of the PCR mix have to be mixed as indicated in the table below:

REAGENT	VOLUME (µl)
Amplification mix	10
Oligo mix	10
Extracted DNA	5

The PCR mix has to be freshly prepared every time.

After its preparation, aliquot **20 µl of Master Mix** in the tubes or in the microplates well for PCR then add in each tube/well **5 µl** (100-250 ng/reaction) from the **extracted DNA** or **control DNA**, place in order the tubes/microplate in the instrument and start the program of amplification. At the end of the program remove the tubes/microplate from the thermal cycler.

d) ANALYSIS and INTERPRETATION of the RESULTS

d.1) Pre-analysis settings

DX/CFX platform

-In *Plate Editor*: set FAM (wild-type genotype, Allele 1) and HEX (mutated genotype, Allele 2) as *Fluorophores* and define the sample name.

ABI 7500 platform

-In the *Setup* menu select *Plate Setup-Define Targets and Samples* to select Samples and to assign VIC target to MUT (Mutated genotype) and FAM target to WT (Wild Type genotype)

-In *Define Samples* insert the sample name

-Select the used wells in *Plate Setup-Assign Targets and Samples*

-Select *None* as passive reference

d.2) Data analysis and Interpretation of the Results

DX/CFX platform

-Select *C(t) Determination Mode* and *Regression* in the *Settings* menu to display the relative scatter plot

-In the *View/Edit Plate* menu set the reaction blank as *NTC*

-In the *Allelic Discrimination* menu select *RFU* and flag the *Normalize Data* in the *Display Mode*


-For both Fluorophores set manually the Threshold to 0,15

The scatter plot displays four quadrants in which each spot represents the two fluorescent signals for a single sample:

Condition in which no signal is detected indicates PCR inhibition.

Control 1 results positive in FAM and negative in HEX/VIC, whereas **Control 2** results positive in HEX/VIC and negative in FAM.

If all these conditions have been met, the run is valid and it is possible to analyse the data; otherwise the run is not valid. It is responsibility of the user to validate the run by checking that these conditions have been met.


The corresponding genotype of the samples is reported in a table displayed close to the Scatter Plot:

Bottom Left Quadrant	Reaction Blank/No Amplification
Bottom Right Quadrant	WT Samples (<i>Allele 1</i>)
Upper Left Quadrant	MUT Samples (<i>Allele 2</i>)
Upper Right Quadrant	Heterozygote Samples (<i>Both Alleles</i>)

ABI 7500 platform

-Select *Analysis Settings* in the *Analysis* menu to manually set the threshold to 35000 and flag the *Automatic Baseline*

-Select *Analyze* in the *Analysis* menu to visualize the amplification graph of the targets for each sample:

-Select *View Well Table* to display the *C_T* values for all samples.

Refer to the following table for a correct interpretation of the results:

<i>C_T</i> value range ^	Sample suitability	Result	Genotype
<i>C_T</i> FAM and <i>C_T</i> VIC >30	Unsuitable	Invalid*	-
<i>C_T</i> FAM <30 and <i>C_T</i> VIC - <i>C_T</i> FAM >5	Valid	Valid	Homozygote Wild Type (WT)
<i>C_T</i> VIC <30 and <i>C_T</i> FAM - <i>C_T</i> VIC >5	Valid	Valid	Homozygote Mutated (MUT)
<i>C_T</i> FAM or <i>C_T</i> VIC <30 and <i>C_T</i> VIC - <i>C_T</i> FAM ≤2.5	Valid	Valid	Heterozygote
<i>C_T</i> FAM or <i>C_T</i> VIC <30 and 2.5< <i>C_T</i> VIC - <i>C_T</i> FAM ≤5	Unsuitable	Equivocal*	Undetermined

^ If the *C_T* value is “undetermined”, assume ***C_T*=35**.

° | ***C_T* VIC - *C_T* FAM**| indicates the absolute value of *C_T* difference.

* **Repeat Extraction and Analysis. Invalid:** Amplification not present or inefficient; problem with the extraction (i.e. DNA loss, presence of inhibitors in the sample or starting material too poor). **Equivocal:** Inefficient amplification or extraction (i.e. inhibitors in the sample or starting material too high) do not allow the correct genotyping of the sample.

-Refer to the following table for Control 1 and Control 2 suitability:

Control	FAM	HEX/VIC	<i>C_T</i> value
1	positive	negative	22-25
2	negative	positive	22-25

If all these conditions have been met, the run is valid and it is possible to analyse the data; otherwise the run is not valid. It is responsibility of the user to validate the run by checking that these conditions have been met.

TROUBLESHOOTING

Problem 1: Weak or no signal of unknown samples.

1. The PCR was inhibited:
 - Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.
2. Pipetting error due to omitted reagents or samples:
 - Repeat the analysis starting from the PCR.
3. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.
4. Very low starting amount and/or low purity of genomic DNA. Improper DNA extraction:
 - Repeat the DNA purification.
5. Wrong channel/filter was chosen. The PCR conditions did not comply with the instructions:
 - Check the PCR conditions and select the fluorescence channels reported in the protocol for the Unknown Sample detection.

Problem 2: Weak or no signal of the Wild Type Control C1 and/or Mutated Control C2.

1. The PCR conditions did not comply with the instructions:
 - Check the amplification protocol and select the fluorescence channel reported in the manual.
2. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.

Problem 3: Any signal with Reaction Blank BM.

1. Contamination during PCR preparation procedure. All samples results are INVALID:
 - Decontaminate the working area and all instruments.
 - Pipette the controls C1 and C2 at last.
 - Repeat the PCR preparation using a new set of reagents.

Problem 4: HEX signal in Wild Type Control C1 and/or FAM signal in Mutated Control C2.

1. Contamination during PCR preparation procedure. All samples results are INVALID:
 - Decontaminate the working area and all instruments.
 - Use tips with aerosol barriers only. Change tips between tubes/wells.
 - Repeat the PCR preparation using a new set of reagents.

Problem 5: Fluorescence intensity varies.

1. The PCR Master Mix is not well prepared:
 - Carefully repeat the PCR preparation procedure.
2. Air bubbles trapped in the PCR tubes:
 - Check the presence of air bubbles before starting a new run.

Problem 6: Absence of any fluorescent signal.

1. Verify the performance of the thermal cycler:
 - Calibrate the equipment.
2. Deterioration of dyes and/or primers. The storage conditions did not comply with the instructions:
 - Check storage conditions.
 - Check the expiry date of the kit.

Problem 7: The thermal cycler gives an error message.

1. Refer to the Real Time PCR instrument user manual or contact the local technical support of the Real Time PCR instrument company.

Problem 8: The kit reagents left out of the storage range temperature.

1. These reagents must be stored as indicated for a proper execution of the test. The performance of the product is not guaranteed if the reagents have not been properly stored.

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CLONIT S.r.l.