Duplica^{RealTime} Beta-Fibrinogen 455G>A Genotyping Kit

REV. EER049050_IFU_REV.01D_ENITA

REF : EER049050-50 tests

Instructions For Use

INTENDED USE

Duplica^{RealTime} Beta-Fibrinogen 455G>A Genotyping Kit is an *in vitro* nucleic acid amplification test for the detection of the allelic variant 455G>A in the promoter region of the Human Beta-Fibrinogen gene, in genomic DNA purified from peripheral blood collected in EDTA.

INTRODUCTION

Increased plasma fibrinogen levels have been identified as a risk indicator for myocardial infarction, stroke, and thrombosis. Both environmental and genetic factors make an important contribution to plasma fibrinogen levels in humans. Fibrinogen (Factor I) is a glycoprotein, synthesized in the liver, that participates in the coagulation cascade. Fibrinogen converts to fibrin in the presence of thrombin. Elevated fibrinogen levels have been associated with increased risk for atherosclerosis and deep venous thrombosis as well as cardiovascular diseases. Fibrinogen is composed of three pairs of polypeptide chains (named α , β and γ) which are linked by disulfide bonds. These three chains are encoded by three different genes (α , β and γ , respectively), that are located on chromosome 4.

Among genetic factors, a common G $\rightarrow$ A polymorphism at position - 455 in the promoter region of the Beta-Fibrinogen gene has been associated with elevated fibrinogen levels (for the AA homozygous carriers).

PRINCIPLE OF THE TEST

Duplica^{RealTime} Beta-Fibrinogen 455G>A Genotyping has been designed to detect the 455G>A polymorphism in the promoter region of the Beta-Fibrinogen gene.

The reagents for the amplification reaction are ready to use and provided in two reaction mixes:

- **AMPLIFICATION MIX**: containing Hot Start Taq DNA polymerase, nucleotides, MgCl₂ and buffer.
- **OLIGO MIX**: containing primers and fluorogenic probes.

Duplica^{RealTime} Beta-Fibrinogen 455G>A Genotyping Kit is based on specific recognition and amplification of target sequences by PCR, and the simultaneous detection of the accumulation of PCR amplification products by fluorescent DNA probes. In particular, the probe designed to detect the Wild Type allele carries the fluorophore FAM (6-carboxy-fluorescein) at the 5' end, while the probe detecting the Mutated allele, is labelled with the fluorophore HEX (hexa-chloro-fluorescein). Both the probes have a non fluorescent black quencher at the 3' end. If excited, the whole probe does not emit fluorescence, since the proximity of the quencher to the reporter prevents the emission of the fluorescence from the reporter (quenching effect).

REAGENTS PROVIDED

The kit has been designed to perform 50 reactions, this format allows to perform 2 runs, each of which includes: **25 samples, 1 Wild Type Control (Control 1, C1), 1 Mutated Control (Control 2, C2) and 1 Reaction Blank (BM)** each.

Kit Components

<i>Reagent</i>	<i>Color Code</i>	<i>Storage (range, °C)</i>	<i>Volume (µl)</i>	<i>Quantity (tubes)</i>
Oligo Mix (OM)*	Green Cap	-22 ÷ -18	600	1
Amplification Mix (AM)	Blue Cap	-22 ÷ -18	600	1
Control 1 (C1, Wild Type)	Red Cap	-22 ÷ -18	>50	1
Control 2 (C2, Mutated)	Yellow Cap	-22 ÷ -18	>50	1
Reaction Blank (BM)	White Cap	-22 ÷ -18	>50	1

*** protect the tube from direct light**

STORAGE AND HANDLING

All reagents must be stored at **-22÷-18°C** and can be used until the expiry date printed on the labels. Do not freeze and thaw the products more than six times.

MATERIAL REQUIRED BUT NOT PROVIDED

- Extraction kit for DNA purification (refer to specific handbook's section)
- Disposable powder-free gloves and laboratory coat
- Variable volumes pipettes (5-20 µl, 20-200 µl, 100-1000 µl)
- Disposable DNase/RNase-free tips with aerosol barriers
- Tube racks
- Desktop centrifuge
- PCR box
- Disposable optical polypropylene tubes for Real Time PCR
- Refrigerator
- Deep-freezer
- Thermal cycler for Real Time PCR **Rotor-Gene®Q** (Qiagen)

All the equipments should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure optimal performance.

PRECAUTIONS AND WARNINGS

- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood
- If required, Clonit offers the necessary technical support for the correct use of the kit
- Carefully read this Instruction for Use before using the kit
- Do not use the reagents after the expiry date
- Thaw and carefully mix the reagents of the kit before use
- Do not mix the reagents from different lots of the product
- Use calibrated and regularly checked pipettes and instrumentation only
- Use dedicated laboratory equipments. Change gloves frequently
- Periodically wipe the working area with 0,5% hypochlorite
- Use powder-free gloves. Do not leave fingerprints or write on optical part of the tubes/microplates
- Materials containing or potentially containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood
- In case of damaged package, contact the technical support before using the kit
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet of the product (MSDS, available on request)
- The kit reagents, individual protective equipments, used materials, biological samples and test residuals must be disposed of according to local regulations
- Patient Drug treatment may interfere with the final result of the molecular biology analysis

OPERATING PROCEDURE

a) DNA purification

For Genomic DNA purification from whole blood samples (200µl extraction volume), Clonit recommends to use:

- the automated platform **QIAasymphony®SP** (Qiagen). For the extraction procedure follow the protocol reported in the user manual of instrumentation (**QIAasymphony DNA mini kit, Qiagen**), starting from 200 µl of whole blood and eluting in 200 µl of Elution Buffer.

- the automated platform **Duplica® PREP Automatic Extractor** (ref. EDI001) to be used with Duplica Blood DNA kit (ref. EDI002250) and Duplica NA Body Fluid kit (ref. EDI004200).

For manual extraction/purification:

- Fasts DNA Releaser (ref. EMR057050) for **fresh** peripheral blood (i.e. stored for up to 24 hours at 2÷8°C).

- Spin DNA Purification Kit (ref. EMR061050) for **frozen** blood or stored at 2÷8°C for more than 24 hours.

Other extraction reagents and methods should be fully tested and evaluated for optimal performance by the user before reporting results.

Attention! Use blood samples in EDTA anti-coagulant solution.

Conservation and quality of purified DNA

The purified genomic DNA should be kept at 2÷8°C if the analysis is performed within 24-72 h, in case of a longer storage period is better to store it at -20°C (or -80°C).

b) Rotor-Gene®Q Thermalcycler setup

- Open the software and select *Advanced* on the box *New Run*

- Select a *new template* in *Empty Run* or a pre-existing one

- Select the Rotor Type of your instrument and then *Next*

- Identify the operator, the sample and the final volume of the amplification reaction (**25 µl**) and then *Next*

- Select *Edit Profile* in the *New Run Wizard* and set up the correct Thermal Profile as indicated in the table below

- Select *Gain Optimisation* in the *New Run Wizard*

- In *Channel Settings* scroll menu, select *Acquiring Channels* and click *Add*

- Set the following parameters for channels Green and Yellow: tube position "1", target sample range "5 FI to 10 FI", acceptable gain range " -10 to 10"

- Flag the option *Perform Optimisation before 1st acquisition*

- Close the window and select *Next* and *Start Run*.

N.B.: before starting the run is recommended to save the file as "Template". In this way it is possible to save the Thermal Profile and Gain Optimisation settings and recall them in subsequent runs.

Thermal Profile.

Time	Temperature	Cycles
10 min	95°C	1
20 sec	95°C	45
30 sec	60°C	Fluorescence Acquisition Green and Yellow channels

c) Preparation of PCR mix

For each experiment prepare a Master Mix considering **2 controls (C1 and C2)**, **1 Reaction Blank (BM)** and **n+1 samples**. The reagents for the Master Mix have to be mixed as indicated in the table below:

REAGENT	VOLUME (µl)
Amplification mix	10
Oligo mix	10
Extracted DNA	5

After its preparation, aliquot **20 µl of Master Mix** in the PCR tubes, then add in each tube **5 µl of extracted DNA or control DNA**, place in order the tubes in the instrument and start the program of amplification. At the end of the program remove the tubes from the thermalcycler.

d) ANALYSIS and INTERPRETATION of RESULTS

The fluorescence in each channel indicates the hybridisation of the allelic specific probes: **Channel Green** for **Wild Type Allele probe (fluorophore FAM)** and **Channel Yellow** for **Mutated Allele probe (fluorophore HEX)**. If a sample shows fluorescence in **Green (FAM)**, the sample has the **Wild Type Allele**. If a sample shows fluorescence in **Yellow (HEX)**, the sample has the **Mutated Allele**.

The results are interpreted through the analysis of the Threshold Cycle (C_T) for each sample using the **"Allelic Discrimination"** function.

The threshold cycle is defined as the cycle at which the amplification curve crosses the threshold (Threshold line). During the first few cycles of PCR amplification the background fluorescence is calculated. Subsequently, it will be used for the calculation of the fluorescence baseline.

For the interpretation of results on Rotor-Gene®Q platform refer to the following indications:

- Select *Analysis* from the toolbar menu then click *Allelic Discrimination*
- Select the 2 channels: Cycling A.Green e Cycling A.Yellow by keeping pressed the SHIFT key on the keyboard
- Click *Show*

In this way, for each sample, the fluorescence data combined are shown (**Fig.1**): the **smooth line** represents the fluorescence signal in **Green (Wild Type Allele)** and the **line with circles** the fluorescence signal in **Yellow (Mutated Allele)**. Before setting the threshold value for genotyping of the samples, it is appropriate to optimize the fluorescence data using the options of normalization of the signal and background reduction of fluorescence offered by the software. The use of these functions allows to avoid any non-specific signals.


Fig.1

The background fluorescence depends on several factors not related to the kit, and its intensity may differ between different sessions of analysis. For this reason the main functions of Rotor-Gene®Q software are shown below, leaving to the operator the possibility to vary, within the limits described, the specific values of these functions (**Fig.2**).

- **Dynamic Tube:** this feature allows the normalization of fluorescence data. The level of the background average is calculated for each sample from cycle 1 to the output cycle (C_T). This function is set by default and is recommended to not disable it.
- **Slope Correct:** the background fluorescence of each sample increases progressively during the reaction cycles. This is mainly due to the chemistry used (Taqman probes) and can result in a growth of the basal fluorescence lines which could be translated into a non-specific signal if it crosses the threshold line. If this increase (slope) is observed, it is recommended to correct it by activating the *Slope Correct*.
- **Ignore first:** the fluorescence signal of the first reaction cycles may not be representative of the rest of the run. For this reason, in some cases is possible to get better results if the first few cycles are ignored. By selecting this function up to first 10 cycles of reaction can be skipped. To enable this feature, click *Ignore First* and indicate the

number of cycles to be ignored by selecting a number between 1 and 10.

- **Outlier removal – NTC Threshold:** This function allows to discriminate the real fluorescence signals from background noise. Clicking on *Outlier Removal* it is possible to select, for the NTC-Threshold function, a value comprised between 0 and 30%. For example, by setting 15% all the signals with an increase of fluorescence less than 15% compared to the maximum increase of fluorescence detected, will be considered non-specific and neutralized by the software. It is recommended not to use the underlying function "Reaction Efficiency Threshold" to avoid that some samples will be excluded from genotyping.

Once removed the background fluorescence it is possible to set the threshold values on the basis of the batch related *Data Sheet* attached to this manual (Fig.2).


Fig.2

Once the threshold is set, check that **Control 1** is detected in **Green channel** and not detected in **Yellow channel**, whereas **Control 2** is detected in **Yellow channel** and not detected in **Green channel**. The **Reaction Blank** must not be detected in any channel.

If these three conditions have been met, the run is valid and it is possible to analyse the data; otherwise the run is not valid. It is responsibility of the user to validate the run.

Therefore, if only a **FAM/Green signal** is detected the sample is **Homozygous Wild Type**, whereas if only a **HEX/Yellow signal** is detected the sample is **Homozygous Mutated**. Finally, if both **FAM/Green** and **HEX/Yellow** are detected the sample is **Heterozygous**.

Condition in which no signal is detected for samples or controls indicates PCR inhibition. In this case, refer to troubleshooting.

Result Interpretation Table

Green channel	Yellow channel	Results
Detected	Not detected	Wild Type
Not detected	Detected	Mutated
Detected	Detected	Heterozygous
Not detected	Not detected	Inhibition

Selecting all the samples it is possible to generate a report of all the data (**Fig.3**)


Fig.3

To generate a Report:

- Select all the samples and controls
- Click *Reports*, select *Allelic Discrimination Analysis* and then *Show*

Additionally, **Scatter Plot Analysis** can be used to visualize the results. In this analysis each sample is identified by a spot representing the end-point fluorescence values for the two dyes (**Fig.4**).

- In the *Analysis* menu select *Other*
- Select *Scatter Graph Analysis*, highlight both *Green* and *Yellow* channels and press *Show* to display the scatter plot. The scatter plot is divided in four areas corresponding to samples classified as Wild Type, Mutant, Heterozygous.
- In the *Scatter Analysis graph* window that appears select with the mouse the area corresponding to a specific Genotype and name it as show in the figure below (*Green* → *Wild-type*; *Yellow* → *Mutated*; and *Green/Yellow* → *Heterozygous*)


Fig.4

TROUBLESHOOTING

Problem 1: Weak or no signal of unknown samples.

1. The PCR was inhibited:
 - Make sure to follow to the manufacturer's instructions for DNA purification
2. The reagents storage conditions did not comply with the instructions:
 - Check the storage conditions.
3. Improper DNA extraction:
 - Repeat analysis starting from the DNA extraction stage
4. Very low starting amount and/or low purity of genomic DNA. Improper DNA extraction or sample storage:
 - Check sample (blood and extracted DNA) storage conditions
 - Repeat DNA purification
5. Wrong channel/filter was chosen:
 - Check fluorescence channels reported in this manual.

Problem 2: Weak or no signal of the Wild Type Control (C1) and/or Mutated Control (C2).

1. The PCR conditions did not comply with the instructions:
 - Check the amplification protocol and select the fluorescence channel reported in the manual.
2. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.

Problem 3: Signal in Yellow channel in C1 and/or in Green channel in C2.

1. Contamination during DNA extraction procedure. All samples results are INVALID:
 - Decontaminate all surfaces and instruments with sodium hypochlorite (0,5%) and ethanol (70%) or special DNA decontamination reagents
 - In case of manual procedure: Use only filter tips during the extraction procedure. Change tips between tubes
 - Repeat the DNA extraction with the new set of reagents.

Problem 4: Any signal with Reaction Blank (BM).

1. Contamination during PCR preparation procedure. All samples results are INVALID:
 - Decontaminate all surfaces and instruments with sodium hypochlorite (0,5%) and ethanol (70%) or special DNA decontamination reagents
 - Pipette the Positive controls at last
 - Repeat the PCR preparation with the new set of reagents.

Problem 5: Fluorescence intensity varies.

1. The PCR Master Mix is not well prepared:
 - Carefully repeat the PCR preparation procedure
2. Air bubbles trapped in the PCR tubes:
 - Check the presence of air bubbles before starting a new run.

Problem 6: Absence of any fluorescent signal.

1. Verify the performance of the thermal cycler:
 - Calibrate the equipment.
2. Deterioration of dyes and/or primers. The storage conditions did not comply with the instructions:
 - Check the storage conditions
 - Check the expiry date of the kit.

Problem 7: The thermal cycler gives an error message.

1. Refer to the Real Time PCR instrument user manual or contact the local technical support.

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CLONIT S.r.l.