**SZABO
SCANDIC**

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic

Factor V Leiden (FV)

REF: RT-50

Determination of the G1691A mutation
on Factor V Leiden in Real Time PCR

INTRODUCTION AND PURPOSE OF USE

Heredity thrombophilia is a genetically determined tendency to develop thromboembolic disease at a young age with recurrent character. Environmental risk factors (smoking, alcohol, stress, diet, ...) and genetic mutations, such as the Leiden mutation localized on the Factor V gene, interact resulting in the expression of this disease.

The Factor V of Leiden (FV) Kit is a qualitative test that allows the allelic discrimination by means of *Real Time PCR*, of Leiden mutation (G1691A) associated to Hereditary Trombophilia.

The Procedure allows the amplification of Wild-type alleles and mutated alleles of Factor V Leiden. Allelic discrimination is performed making a scatter plot of mutated allele's fluorescence (Y axis) versus wild-type allele's fluorescence (X axis), discriminating in this way the three possible genotype: Homozigote Wild-Type, Homozygote Mutated and Heterozygote Mutated.

The analysis of the results is made by an instrument of Real Time PCR, composed by a thermal cycler with a system of fluorescence detection.

CONTENT

The kit contains reagents enough to perform 48 amplification tests:

Quantity	Description
R1 3 x 220 µl	Amplification mMix dNTPs, Tris-HCl, KCl, MgCl ₂ , Taq Polymerase, Nuclease-free water (Blue cap)
R2 3 x 130 µl	Factor V probes Mix Upstream primer, downstream primer, Target probes (FAM for Wild Type and VIC for Mutant), Nuclease-free water (Green cap)
R3 3 x 35 µl	Positive Control Homozygote Wild-Type cloned DNA corresponding to FV Wild-Type gene.
R4 3 x 35 µl	Positive Control Homozygote Mutated cloned DNA corresponding to FV Mutated gene.
R5 1 x 30 µl	Negative Control (Yellow cap)

Instruction for use: ST. RT50-ENG.7

MATERIALS AND STRUMENTATION REQUIRED BUT NOT SUPPLIED

Disposable latex powder-free gloves or similar material;
Bench microcentrifuge (12,000 - 14,000 rpm);
Micropipettes and Sterile tips with aerosol filter;
Vortex;
Plastic materials (microplate and optical adhesive cover);
Heat block (only for extraction);
Dry block shaker for 1.5ml conical tubes
Magnetic rack for 1.5ml conical tubes
EZ1 ADV XL DSP DNA Blood Card (ref. 9018702)

Reagents

The Factor V of Leiden kit was developed and validated to be used with the following extraction method:

Manual Extraction

Ref. 51304/51306

QIAamp DNA mini kit. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 50/250 samples. (QIAGEN)

Automatic extraction

Ref. 62124

EZ1 DSP DNA Blood kit. The kit allows the automatic DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 48 samples. (QIAGEN)

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10629800 - **VERSANT® Sample Preparation 1.2 Reagents kit box**
1. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 96 samples.

10629801 - **VERSANT® Sample Preparation 1.2 Reagents kit box**
2. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 96 samples.

Strumentation

Automatic Extraction

Ref. 9001492

EZ1 Advanced XL Robotic Workstation. (QIAGEN) for the automatic nucleic acid purification up to 14 samples simultaneously (QIAGEN).

Real Time PCR

The Factor V of Leiden kit was developed and validated to be used with the following real time PCR instruments:

- 7500 Fast from Lifetechnologies
- StepOne Plus from Lifetechnologies
- Rotor-Gene Q MDx from QIAGEN
- Versant kPCR AD from Siemens or Stratagene MX3005P/MX3000P
- CFX96 Real Time PCR System from BioRad

SAMPLES AND STORAGE

The Factor V Leiden (FV) system must be used with extracted DNA from the following biological samples: Whole Blood EDTA. Collected samples must be shipped and stored at +2 - +8°C and used within 3 days from the collected data. Store the sample at -20°C if it is used after 3 days.

PRECAUTIONS USE

This kit is for *in vitro* diagnostic (IVD), for professional use only and not for in vivo use.

After reconstitution, the amplification master mix must be used in one time (16 reactions). Do not re-frost already used material. Repeat thawing and freezing of reagents (more than twice) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.

At all times follow Good Laboratory Practice (GLP) guidelines.

Wear protective clothing such as laboratory coats and disposable gloves while assaying samples.

Avoid any contact between hands and eyes or nose during specimens collection and testing.

Handle and dispose all used materials into appropriate bio-hazard waste containers. It should be discarded according to local law. Keep separated the extraction and the reagents preparation.

Never pipette solutions by mouth.

Avoid the air bubbles during the master mix dispensing. Eliminate them before starting amplification.

Wash hands carefully after handling samples and reagents.

Do not mix reagents from different lots.

It is not infectious and hazardous for the health (see Material Safety data Sheet – MSDS).

Do not eat, drink or smoke in the area where specimens and kit reagents are handled.

Read carefully the instructions notice before using this test.

Do not use beyond the expiration date which appears on the package label.

Do not use a test from a damaged protective wrapper.

LIMIT OF THE METHOD

Use only extracted DNA from **Whole blood collected in EDTA**.

Do not use contaminated DNA with extracted mucoproteins or haemoglobin; the latter inhibit the amplification reaction of nucleic acids and may cause invalid results.

The extreme sensitivity of gene amplification may cause false positives due to cross-contamination between samples and controls. Therefore, you should:

- physically separate all the products and reagents used for amplification reactions from those used for other reactions, as well as from post-amplification products;
- use tips with filters to prevent cross-contamination between samples;
- use disposable gloves and change them frequently;
- carefully open test tubes to prevent aerosol formation;
- close every test tube before opening another one.

The proper functioning of the amplification mix depends on the correct collection, correct transportation, correct storage and correct preparation of a biological sample.

As with any diagnostic device, the results obtained with this product must be interpreted taking in consideration all the clinical data and other laboratory tests done on the patient.

As with any diagnostic device, with this product there is a residual risk of obtaining invalid, false positives or false negatives results.

STORAGE AND STABILITY

Store the product Factor V Leiden (FV) at **-20°C**.

The **Factor V of Leiden** kit is shipped on dry ice. The kit components should be frozen.

If one or more components are not frozen upon receipt or if the tubes have been compromised during transport, contact Clint srl for assistance.

An intact and well stored product has a stability of 18 months from the date of production. Do not use beyond the expiration date which appears on the package label.

Repeat thawing and freezing of reagents (more than twice) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.

ANALYTICAL PROCEDURE

Human DNA Extraction

Manual Extraction

Ref. 51304/51306 - **QIAamp DNA mini kit** (QIAGEN).

Follow the instructions inside the kit **QIAamp DNA Mini Kit**. Elute the sample in 50 µl of buffer AE.

Automatic extraction

Ref. 62124 - **EZ1 DSP DNA Blood kit** on **EZ1 Advanced XL instrument**.

Follow the instructions inside the kit **EZ1 DSP DNA Blood kit**. Start from 200 µl of sample and elute it in 50 µl of buffer AE.

Manual extraction (SIEMENS)

Ref. 10629800 - **VERSANT® Sample Preparation 1.2 Reagents kit box 1**.

Ref. 10629801 - **VERSANT® Sample Preparation 1.2 Reagents kit box 2**.

Follow the instructions supplied by Siemens and elute it in 70 µl of Elution buffer. Transfer 55 µl of eluted sample to an appropriately size tube.

Sample can be stored at -20°C.

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	60° C 1 min

Select the annealing / extension from the thermal profile and click on "Acquiring A to cycling".

In the next window, select yellow from the available channels and add it to acquiring channel along with the green channel and click "OK". In the next window click on "OK" and then click "Next".

click on "Edit Gain" button and set the following values for each channel:

Reporter	Gain
Green	8
Yellow	2

To begin the course, click on the button "Start Run". You can save the model before you begin your run by clicking on "Save Template".

After clicking on the button "Start Run" window appears "Save As". The stroke can be saved in the desired position by the user.

Once the run started, the window "Edit Samples" allows you to set the name of samples and controls in the positions in which they were loaded on the instrument.

Select the locations where they were positioned the Wild Type and Mutant controls designate them as **FV WT positive CTR and FV Mutant Positive CTR**. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analyzed. Select "Positive Controls".

Select the location where you placed the Negative Control and name it as **Negative Control**. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analyzed. Select "Negative Controls".

Select the location of each sample and enter the name or code of the patient. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analyzed. Select "Unknown".

At the end of the operation click "OK" in the "edit samples" and wait until the end of the race for the analysis (see "Interpretation of Results").

Analysis (see "Interpretation of Results").

Versant kPCR AD or Stratagene MX3005P/MX3000PMX3000

Turn on the instrument, the computer and start the control software. Turn on the light at least 20 minutes before starting a new experiment. You can click on the lamp image for turning on the light from toolbar or you can select "Lamp On" from menu "Instruments".

Verify the correct setting of fluorescent reporters gains: In the setting menu choose and then "Filter set gain setting".

Reporter	Gain
FAM	8
JOE/HEX	2
ROX	1

From the major screen of the software it will open a window "New option - Select experiment/Project type": select "Allelic Discrimination/SNPs Real Time". The software automatically open the window plate setup.

Choose a plate zone where you can put the positive control Factor V Wild Type and select in the bar on the right from the menu.

Well type:	Collect	Fluorescent Data:	Reference	Replicate
Pos. Control FAM		FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window "well information" here you can set the name of the calibrator: FV-WT

It's possible to set the name of the dye near the name of the analyzed target:

FAM	HEX
FV-WT	FV-MUT

Choose a plate zone where you can put the positive control Factor V mutated and select in the bar on the right from the menu.

Well type:	Collect	Fluorescent Data:	Reference	Replicate

<tbl_r cells="5" ix="2" maxcspan="1" maxrspan="1" usedcols="

QUALITATIVE ANALYSIS

Lifetechnologies 7500 Fast, StepOne Plus.

Scatter Plot Analysis

At the end of the amplification reaction, the software automatically shows the obtained results in the "Allelic Discrimination Plot".


The Lifetechnologies 7500Fast/StepOne Plus instrument automatically perform the genotyping of unknown samples by comparing controls versus homozygous wild-type and homozygous mutated contained in the kit.

A proper results' analysis needs a correct settings of the instrumentation. For this aim, set:

Baseline fluorescence level from cycle 6;

	FV-WT - FAM Threshold	FV-MUT - VIC Threshold
7500 Fast	0.2	0.35
StepOne Plus	0.2	0.35

It is recommended to verify correct placement of each individual sample on the scatter plot. For viewing the report containing all data obtained during the analysis, click the sheet "view table well".

The position of the negative control in other plot region could be a contamination of the reaction mix.
Be sure that in the negative control not be increasing of specific fluorescence by examining targets (FAM and VIC).

ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose it is necessary a correct setting of the software

Set the level of background fluorescence (Baseline) from cycle 6; Set the following threshold:

	FV-WT - FAM Threshold	FV-MUT - VIC Threshold
7500 Fast	0.25	0.25
StepOne Plus	0.25	0.25

Export data to Excel and set the formula for each sample and control:

Allele2 Ct (FV-MUT) – Allele1 Ct (FV-WT)

See paragraph "INTERPRETATION OF RESULTS"

ROTOR - GENE Q MDx

Scatter Plot Analysis

At the end of the amplification reaction, click on Analysis from RotorGene Software. In the Analysis window select Allelic Discrimination sheet, click on Cycling A green-Cycling A yellow using CTRL button and click "show". The amplification plot will appear.

Set in the "CT calculation – Threshold" the following value:

	FV-WT - FAM Threshold	FV-MUT - VIC Threshold
RotorGene-Q	0.3	0.3

Click genotypes button and set:

	Reacting Channel	Reacting Channel
Wild type	cycling A (green)	cycling A (green)
Heterozygous	cycling A (green)	cycling A (yellow)*
None		
Mutant		cycling A (yellow)*

In the Analysis window select "Scatter graph analisys" sheet, click on Cycling A green-Cycling A yellow and click "show". The scatter plot will appear with the 2 clouds: Wild Type (bottom right) and Mutant (high left). Is it possible to print the final report.

ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

At the end of the PCR run open the "Analysis" window. Select the "Quantification" sheet and click on "cycling A (green)".

Select from the menu "Dynamic Tube" and subsequently "Slope correct". Set the correct setting of the threshold in the space provided "CT calculation – Threshold".

	FV-WT – Green (FAM) Threshold
RotorGene-Q	0.05

Open the "Analysis" window. Select the "Quantification" sheet and click on "cycling A (yellow)".

Select from the menu "Dynamic Tube" and subsequently "Slope correct". Set the correct setting of the threshold in the space provided "CT calculation – Threshold".

	FV-MUT – Yellow (VIC) Threshold
RotorGene-Q	0.05

Export data to Excel, save the file as "Excel Analysis Sheet" and enter the following formula for each sample and control:

Yellow Ct (FV-MUT) – Green Ct (FV-WT)

See paragraph "INTERPRETATION OF RESULTS"

Versant kPCR AD or Stratagene MX3005P/MX3000PMX3000P

Scatter Plot Analysis
Click on button "Analysis" in the toolbar. The software immediatley open sheet "Analysis Term Setting". Activate the keys FAM and HEX in below screen and select the samples and disposable controls.

On the window "Analysis Term Setting" open the sheet "Results". Select in the right screen the voice "Amplification Plot" in the area "Area to Analysis". Set the values in the area "Threshold Fluorescence" the values:

	FV-WT - FAM Threshold	FV-MUT - VIC Threshold
Kpcr Versant	0.4	0.5

Select in the right screen the voice "Dual Color Scatter Plot" in the area "Area to Analysis". Select "Fluorescence" in the area "Display value for" and choose from menu "Rlast/Rfirst".

Select in the area "Allele Association":

A Allele: FV Wild Type

B Allele: FV Mutated

At the end it will be possible obtain the detailed account, clicking "Text Report" in the area "Area to Analysis".

Only for Versant kPCR AD from Siemens or Stratagene MX3005P/MX3000PMX3000P

If the software can't discriminate the wild type and mutated positive controls click "Show all genotypes". It will be now possible to modify the genotyping windows, select the samples placed near the positive wild type and identify them as wild type. Select the samples placed near the positive mutated and identify them as mutated. Select the samples placed in an intermediate position between wild type positive control and mutated positive control and identify them as Heterozygote for the Leiden Mutation.

ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

Click "Analysis" in the toolbar.

Click the sheet "Results", and choose the analysis "Amplification plot". Check the correct setting of the threshold in the window "Threshold fluorescence" and set the following values:

	FV-WT - FAM Threshold	FV-MUT - HEX Threshold
MX3000P – MX3005P	0.1	0.1

From the Text Report window you can export the results by clicking on the main menu: file, export

Export data to Excel and set the formula for each sample and control:

HEX Ct (FV-MUT) – FAM Ct (FV-WT)

See paragraph "INTERPRETATION OF RESULTS"

CFX96 real time PCR system

Scatter Plot Analysis
At the end of the PCR, select the "Allelic Discrimination" sheet. On the bottom of the screen, set "Selected fluorophores": X = FAM and Y = VIC. Choose RFU from "Display Mode" and "Normalize data". Set the following Threshold:

	FV-WT - FAM Vertical Threshold	FV-MUT - VIC Horizontal Threshold
CFX 96	0.4	0.3

You can export the report pushing the paper block figure on the top of the screen.

ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

At the end of the PCR, select the sheet "Quantitation". Analyzing a fluorophore at a time, open in the menu "Setting" the voice Baseline Threshold and set the following values of the Threshold in the "User Defined"

	FV-WT - FAM Threshold	FV-MUT - VIC Threshold
CFX 96	200	200

You can export the report by clicking on the figure notepad at the top of the screen.

Export data to Excel, save the file as "Excel Analysis Sheet" and enter the following formula for each sample and control:

VIC Ct (FV-MUT) – FAM Ct (FV-WT)

See paragraph "INTERPRETATION OF RESULTS"

INTERPRETATION OF RESULTS

The use of positive and negative control in each amplification session allows to verify the correct functioning of the amplification mix and the absence of any contamination.

The instrument software is able to analyze the fluorescences that are emitted by the specific probe for FV Wild-Type (FAM) and by the specific probe for the FV mutated (VIC/HEX).

Scatter Plot Analysis

A regular functioning of the amplification mix can be verified analyzing the correct position of positive controls and negative controls on the scatter plot.

- Positive control Homozygote Wild-Type: horizontal position on X axis (down on the right) Ct < 30
- positive Control homozgote Mutated: vertical position on Y axis (up on the left) Ct < 30
- negative Control: placed at the origin of cartesian plane (down on the left)

Genotyping tests (allelic discrimination) are endpoint experiments: fluorescence data are collected at the end of the reaction (Post PCR Read) and subtracted to initial read fluorescence (Pre PCR Read).

The software makes a scatter plot with obtained results: Y axis is the normalized fluorescence of Mutated Allele, while X axis shows the normalized fluorescence of Wild-Type Allele. The diagnosis obtained with the comparison between unknown samples and Homozygote Wild-Type and Homozygote Mutate, given by the kit. Selecting the controls, we obtain their disposition on the scatter plot, depending of their relation between fluorescence emitted by two probes FAM (Wild-Type) and VIC/HEX (Mutated).

ΔCt Analysis

By analyzing ΔCt you can identify the correct genotype of the sample being analyzed by performing subtraction of the assigned Mutated allele Ct to assigned wild-type allele Ct. The genotype is determined by following the table below:

Genotype	ΔCt (Ct Mutated - Ct Wild Type)
Wild Type	ΔCt > 2
Mutated	ΔCt < -2
Heterozygous	-2 > ΔCt < 2

In the amplification of each positive control (wild-type and Mutated), the Ct values of the allele-specific probe are used to validate the assay.

Make sure that the fluorescence emitted by the amplification allele has a correct Ct < 27.

If the result of the allele-specific amplification of each control has a Ct > 27 or undetermined the session can not be considered valid and must be repeated.

In the amplification of each sample, the Ct values of the allele-specific probe are used to validate the assay from the extraction process up to the stage of detection.

Make sure that the fluorescence emitted by the allele-specific amplification of the sample, identified after analysis, has not a Ct > 27.

Sample	Wild Type Allele (FAM)	Mutated Allele (VIC)	Assay and Genotype
Wild Type	Ct < 27	Not relevant	Valid
Mutated	Not relevant	Ct < 27	Valid
Heterozygous	Ct < 27	Ct < 27	Valid

If a sample has a Ct > 27 means that there are problems in the extraction phase or in the amplification and therefore could be assigned to a sample genotype wrong. Repeat sample.

Sample	Wild Type Allele (FAM)	Mutated Allele (

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