# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation

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[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## MTHFR C677T

### REF: RT-52 Determination of the C677T mutation on MTHFR in Real Time PCR

#### INTRODUCTION AND PURPOSE OF USE

Hereditary thrombophilia is a genetically determined tendency to develop thromboembolic disease at a young age with recurrent character. Environmental risk factors (smoking, alcohol, stress, diet, ...) and genetic mutations, such as the C677T mutation localized on the MTHFR gene, interact resulting in the expression of this disease.

The MTHFR C677T Kit is a qualitative test that allows the allelic discrimination, by means of *Real Time PCR*, of C677T mutation associated to Hereditary Trombophilia.

The Procedure allows the amplification of Wild-type alleles and mutated alleles of MTHFR C677T. Allelic discrimination is performed making a scatter plot of mutated allele's fluorescence (Y axis) versus wild-type allele's fluorescence (X axis), discriminating in this way the three possible genotype: Homozigote Wild-Type, Homozigote Mutated and Heterozigote Mutated.

The analysis of the results is made by an instrument of Real Time PCR, composed by a thermal cycler with a system of fluorescence detection.

#### CONTENT

The kit contains reagents enough to perform 48 amplification tests:

Quantity	Description
R1 3 x 220 µl	<b>Amplification mMix</b> dNTPs, Tris-HCl, KCl, MgCl <sub>2</sub> , Taq Polymerase, Nuclease-free water (Blu cap)
R2 3 x 130 µl	<b>MTHFR C677T probes Mix</b> Upstream primer, downstream primer, Target probes (FAM for Wild Type and VIC for Mutant), Nuclease-free water (Green cap)
R3 3 x 35 µl	<b>Positive Control Homozygote Wild-Type</b> cloned DNA corresponding to MTHFR C677T Wild-Type gene.
R4 3 x 35 µl	<b>Positive Control Homozygote Mutated</b> cloned DNA corresponding to MTHFR C677T Mutated gene.
R5 1 x 30 µl	<b>Negative Control</b> (Yellow cap)

Instruction for use: **ST. RT52-ENG.7**

#### MATERIALS AND STRUMENTATION REQUIRED BUT NOT SUPPLIED

Disposable latex powder-free gloves or similar material;  
Bench microcentrifuge (12,000 - 14,000 rpm);  
Micropipettes and Sterile tips with aerosol filter;  
Vortex;  
Plastic materials (microplate and optical adhesive cover);  
Heat block (only for extraction)  
Dry block shacker for 1.5ml conical tubes  
Magnetic rack for 1.5ml conical tubes  
*EZ1 ADV XL DSP DNA* Blood Card (ref. 9018702)

#### Reagents

The MTHFR C677T kit was developed and validated to be used with the following extraction method:

#### Manual Extraction

Ref. 51304/51306

*QIAmp DNA mini kit*. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 50/250 samples. (QIAGEN)

#### Automatic extraction

Ref. 62124

*EZ1 DSP DNA Blood kit*. The kit allows the automatic DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 48 samples. (QIAGEN)

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10629800 - *VERSANT® Sample Preparation 1.2 Reagents kit box 1*. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 96 samples.

10629801 - *VERSANT® Sample Preparation 1.2 Reagents kit box 2*. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 96 samples.

#### Strumentation

##### Automatic Extraction

Ref. 9001492.

*EZ1 Advanced XL*. Robotic Workstation. (QIAGEN) for the automatic nucleic acid purification up to 14 samples simultaneously (QIAGEN).

#### Real Time PCR

The MTHFR C677T kit was developed and validated to be used with the following real time PCR instruments:

- 7500 *Fast* from Lifetechnologies
- StepOne Plus* from Lifetechnologies
- Rotor-Gene Q MDx* from QIAGEN
- Versant kPCR AD* from Siemens or *Stratagene MX3005P/MX3000P*
- CFX96 Real Time PCR System* from BioRad

#### SAMPLES AND STORAGE

The MTHFR C677T system must be used with extracted DNA from the following biological samples: Whole Blood EDTA. Collected samples must be shipped and stored at +2 - +8°C and used within 3 days from the collected data.  
Store the sample at -20°C if it is used after 3 days.

#### PRECAUTIONS USE

This kit is for *in vitro* diagnostic (IVD), for professional use only and not for *in vivo* use.

**After reconstitution, the amplification master mix must be used in one time (16 reactions). Do not re-frost already used material. Repeat thawing and freezing of reagents (more than twice) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.**

At all times follow Good Laboratory Practice (GLP) guidelines.

Wear protective clothing such as laboratory coats and disposable gloves while assaying samples.

Avoid any contact between hands and eyes or nose during specimens collection and testing.

Handle and dispose all used materials into appropriate bio-hazard waste containers. It should be discarded according to local law.

Keep separated the extraction and the reagents preparation.

Never pipette solutions by mouth.

Avoid the air bubbles during the master mix dispensing. Eliminate them before starting amplification.

Wash hands carefully after handling samples and reagents.

Do not mix reagents from different lots.

It is not infectious and hazardous for the health (see Material Safety data Sheet – MSDS).

Do not eat, drink or smoke in the area where specimens and kit reagents are handled.

Read carefully the instructions notice before using this test.

Do not use beyond the expiration date which appears on the package label.

Do not use a test from a damaged protective wrapper.

#### LIMIT OF THE METHOD

Use only extracted DNA from **Whole blood collected in EDTA**.

Do not use contaminated DNA with extracted mucoproteins or hemoglobin: the latter inhibit the amplification reaction of nucleic acids and may cause invalid results .

The extreme sensitivity of gene amplification may cause false positives due to cross-contamination between samples and controls. Therefore, you should:

- physically separate all the products and reagents used for amplification reactions from those used for other reactions, as well as from post-amplification products;
- use tips with filters to prevent cross-contamination between samples;
- use disposable gloves and change them frequently;
- carefully open test tubes to prevent aerosol formation;
- close every test tube before opening another one.

The proper functioning of the amplification mix depends on the correct collection, correct transportation, correct storage and correct preparation of a biological sample.

As with any diagnostic device, the results obtained with this product must be interpreted taking in consideration all the clinical data and other laboratory tests done on the patient.

As with any diagnostic device, with this product there is a residual risk of obtaining invalid, false positives or false negatives results.

#### STORAGE AND STABILITY

Store the product MTHFR C677T at **-20°C**.

The **MTHFR C677T** kit is shipped on dry ice. The kit components should be frozen.

If one or more components are not frozen upon receipt or if the tubes have been compromised during transport, contact Clonit srl for assistance.

An intact and well stored product has a stability of 18 months from the date of production. Do not use beyond the expiration date which appears on the package label.

Repeat thawing and freezing of reagents (more than twice) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.

#### ANALYTICAL PROCEDURE

##### Human DNA Extraction

##### Manual Extraction

Ref. 51304/51306 - *QIAmp DNA mini kit* (QIAGEN).

Follow the instructions inside the kit *QIAmp DNA Mini Kit*. Elute the sample in 50 µl of buffer AE.

#### Automatic extraction

Ref. 62124 - *EZ1 DSP DNA Blood kit* on *EZ1 Advanced XL instrument*.

Follow the instructions inside the kit *EZ1 DSP DNA Blood kit*. Start from 200 µl of sample and elute it in 50 µl of buffer AE.

#### Manual extraction (SIEMENS)

Ref. 10629800 - *VERSANT® Sample Preparation 1.2 Reagents kit box 1*.

Ref. 10629801 - *VERSANT® Sample Preparation 1.2 Reagents kit box 2*.

Follow the instructions supplied by Siemens and elute it in 70 µl of Elution buffer. Transfer 55 µl of eluted sample to an appropriately size tube.

Sample can be stored at -20°C.

#### SOFTWARE SETTING:

##### *Lifetechnologies 7500 fast/StepOne plus*

Turn the instrument and the computer on and open the control software. Click on **"Advance Setup"**: by default the software will show the page **"experiment properties"**. Write in the **"experiment name"** the file name, choose the type of instrument (**ABI7500 o ABI7500fast**), the type of reaction (**Genotyping**), the type of used reagent (**Taqman®Reagents**) and the reaction time of analysis (**Standard ≈ 2 hours to complete a run**).

Open the page named **"page setup"**.

In the window **"Assign SNP assay to the selected wells"** open **"Create new SNP Assay"** and set:

**SNP Assay Name: MTHFR C677T:**

	Reporter	Quencer
<b>Allele 1 Name: MTHFR C677T WT</b>	FAM	NFQ
<b>Allele 2 Name: MTHFR C677T MUT</b>	VIC	NFQ

In the page **"plate setup"**, move on the area **"Assign Sample to the selected Wells"**: set the name of the analyzing samples, of positive controls and negative controls.

Choose an area of the plate where positive controls will be placed: select in the blank **"Assign SNP assay to the selected well"** and assign the MTHFR C677T SNP Assay. After set these tasks:

**"task Positive control Allele1/Allele1"** for MTHFR C677T Wild Type homozygote target;

**"task Positive control Allele2/Allele2"** for MTHFR C677T mutated homozygote target;

Choose an area of the plate where negative control will be placed: select **"Assign SNP assay to the selected well"** the "task Negative control" for MTHFR C677T SNP.

Select an area of the plate where samples will be placed: select the wells and set MTHFR C677T SNP. Link every well to a sample, through the window **"Assign samples to selected wells"**.

For each sample, select in the blank **"Assign SNP to selected wells"** the **"task Unknown (U)"** for the MTHFR C677T SNP.

Set ROX as passive reference, using it as normalizer of detected fluorescence.

Open **"Run Method"** (sheet **Graphic View**) and set the thermal cycling as follows:

cycles	Pre PCR Read	denaturation	Annealing extension	Post PCR Read
1	60°C 1min			
1		95°C 10min		
35		95°C 15sec	60°C 1min	
1				60°C 1min

In the window **"Reaction volume plate per well"** set a volume of 25 µl.

After having prepared the plate, and correctly inserted in the instrument, press the button **"Start Run"**.

#### ROTOR-GENE Q MDx

The experiments can be **set using the Quick Start Wizard or the Advanced Wizard**, which appears when the software is started.

Select the wizard **"Advanced"**. As a first step, select the model **"Two Step Reaction"** with a double click in the **"New Run"**.

In the next window, select the type of rotor installed on the instrument from the list that appears. Check the **"Locking Ring Attached"**, check the checkbox and then click **"Next"**.

Enter the name of the operator and the reaction volume of 25 µl, and then click **"Next"**.

In the next window click on **"edit profile"**. Set the following thermal cycle:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	60° C 1 min

Select the annealing / extension from the thermal profile and click on **"Acquiring A to cycling."**

In the next window, select yellow from the available channels and add it to acquiring channel along with the green channel and click **"OK"**. In the next window click on **"OK"** and then click **"Next"**.

click on "Edit Gain" button and set the following values for each channel:

Reporter	Gain
Green	8
Yellow	2

To begin the course, click on the button **"Start Run"**. You can save the model before you begin your run by clicking on **"Save Template"**.

After clicking on the button **"Start Run"** window appears **"Save As"**. The stroke can be saved in the desired position by the user.

Once the run started, the window **"Edit Samples"** allows you to set the name of samples and controls in the positions in which they were loaded on the instrument.

Select the locations where they were positioned the Wild Type and Mutant controls designate them as **MTHFR C677T WT positive CTR** and **MTHFR C677T Mutant Positive CTR**. Clicking on the box next to **"Type"** correspondent, in the dropdown menu **"Samples"** you can select the type of sample being analyzed. Select **"Positive Controls"**.

Select the location where you placed the Negative Control and name it as **Negative Control**. Clicking on the box next to **"Type"** correspondent, in the dropdown menu **"Samples"** you can select the type of sample being analyzed. Select **"Negative Controls"**

Select the location of each sample and enter the name or code of the patient. Clicking on the box next to **"Type"** correspondent, in the dropdown menu **"Samples"** you can select the type of sample being analyzed. Select **"Unknown"**

At the end of the operation click **"OK"** in the **"edit samples"** and wait until the end of the race for the analysis (see **"Interpretation of Results"**).

Analysis (see **"Interpretation of Results"**).

#### Versant kPCR AD or Stratagene MX3005P/MX3000P

Turn on the instrument, the computer and start the control software. Turn on the light at least 20 minutes before starting a new experiment. You can click on the lamp image for turning on the light from toolbar or you can select **"Lamp On"** from menù **"Instruments"**.

Verify the correct setting of fluorescent reporters gains: In the setting menu choose and then **"Filter set gain setting"**

Reporter	Gain
FAM	8
JOE/HEX	2
ROX	1

From the major screen of the software it will open a window **"New option – Select experiment/Project type"**: select **"Allelic Discrimination/SNPs Real Time"**. the software automatically open the window **plate setup**.

Choose a plate zone where you can put the positive control MTHFR C677T Wild Type and select in the bar on the right from the menu.

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
Pos. Control FAM	FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window **"well information"** here you can set the name of the calibrator: MTHFR-C677T WT

It's possible to set the name of the dye near the name of the analyzed target:

FAM	HEX
MTHFR C677T -WT	MTHFR C677T -MUT

Choose a plate zone where you can put the positive control MTHFR C677T mutated and select in the bar on the right from the menu.

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
Pos. Control HEX	FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window **"well information"** here you can set the name of the calibrator: MTHFR-C677T MUT

It's possible to set the name of the dye near the name of the analyzed target:

FAM	HEX
MTHFR C677T -WT	MTHFR C677T -MUT

Choose a plate zone where you can put the unknow samples and select in the bar on the right from the menu.

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
UnKnown	FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window **"well information"** here you can set the name of the sample: Every well **"Unknown"** should be singularly named clicking on single name it will open the window **"well information"** in this window it is possible to insert the sample name in analysis.

Choose a plate zone where you can put the negative control and select in the right bar on the right from the menu.

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
NTC	FAM/HEX/ROX	ROX	None

Select all the well corresponding to occupied places and clicking right key of the mouse selecting the voice **"well information"**: corresponding to **FAM** fluorophore set in the specific place (**Assay**): **MTHFR C677T WT**. Corresponding to fluorophore **HEX** set in the specific place (**Assay**): **MTHFR C677T MUT**.

Open page **"Thermal Profile setup"** and set the right thermal cycle:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	60° C 1 min

After preparing the plate and correctly insert it in the instument, push the button **"Start Run"**. the software will ask the file name for the saving data.

#### CFX 96 Real Time PCR

Turn the instrument and the computer on and start the control software. In the principal screen will appear the window **"Startup wizard"**: select **"CFX96"** and press **"ok"**. In the next window push **"create new"** and set the thermal protocol and the reaction volume (25µl):

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	60° C 1 min

Save the protocol and click the next button. The software will open in default the sheet **"plate"**. Click **"create new"**, select **"Fluorophores button"** to choose fluorophores (FAM and VIC). Select the locations where they were positioned the controls of known concentration and choose the **"Sample Type" Standards**. Click **"Load"** check boxes to load fluorophores and Type or select Target Name. Select the location where you placed the Negative Control. Choose the **"Sample Type" NTC**. Click **"Load"** check boxes to load fluorophores and Type or select Target Name. Select the location of each sample and enter the name or code of the patient. Choose the **"Sample Type" Unknown**. Click **"Load"** check boxes to load fluorophores and Type or select Target Name. Save the plate clicking the next button and start the experiment

#### PREPARATION OF THE REACTIONS:

Defrost a tube of **Amplification mMix**;

Defrost a tube of **MTHFR C677T probes Mix**;

Mix carefully by vortex 210µl of **Amplification mMix** and 126µl of **MTHFR C677T probes Mix** (the mix as produced is enough to prepare **16 reactions** of amplification: **2 positive controls, 1 negative control and 13 samples**).

Distribute, in the amplification plate, **20 µl of just reconstituted mix** in chosen positions and already setted on the instrument software.

Distribute, in the negative control position, **5µl** of solution taken by the **negative control** vial.

Distribute, in chosen position for each sample, **5µl** of **corresponding sample**.

Distribute, in chosen positions for the positive controls, **5µl** of **Positive Control Homozygote Wild-Type** and **Positive Control Homozygote Mutated**.

Seal up accurately the plate using an optichal adhesive film and verify that there aren't air bubbles in the mix, to avoid the amplification interferences.

Transfer the plate in the instrument and push the button **"Start Run"**.

## QUALITATIVE ANALYSIS


### Lifetechnologies 7500 Fast, StepOne Plus. Scatter Plot Analysis

At the end of the amplification reaction, the software automatically shows the obtained results in the "Allelic Discrimination Plot". The Lifetechnologies 7500Fast/StepOne Plus instrument automatically perform the genotyping of unknown samples by comparing controls versus homozygous wild-type and homozygous mutated contained in the kit. A proper results' analysis needs a correct settings of the instrumentation. For this aim, set:

Baseline fluorescence level from cycle 6;

	MTHFR- C677T WT - FAM Threshold	MTHFR- C677T MUT - VIC Threshold
7500 Fast	0.35	0.3
StepOne Plus	0.35	0.3

It is recommended to verify correct placement of each individual sample on the scatter plot. For writing the report containing all data obtained during the analysis, click the sheet "view table well".

The position of the negative control in other plot region could be a contamination of the reaction mix. Be sure that in the negative control not be increasing of specific fluorescence by examining targets (FAM and VIC).

### ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose it is necessary a correct setting of the software

Set the level of background fluorescence (Baseline) from cycle 6; Set the following threshold:

	MTHFR C677T -WT - FAM Threshold	MTHFR C677T -MUT - VIC Threshold
7500 Fast	0.25	0.25
StepOne Plus	0.25	0.25

Export data to Excel and set the formula for each sample and control:

**Allele2 Ct (MTHFR C677T -MUT) – Allele1 Ct (MTHFR C677T -WT)**

See paragraph "INTERPRETATION OF RESULTS"

### ROTOR - GENE Q MDx

#### Scatter Plot Analysis

At the end of the amplification reaction, click on **Analysis** from RotorGene Software. In the **Analysis window** select **Allelic Discrimination** sheet, click on **Cycling A green-Cycling A yellow** using **CTRL** button and click "show". The amplification plot will appear.

Set in the "CT calculation – Threshold" the following value:

	MTHFR C677T - WT - FAM Threshold	MTHFR C677T - MUT - VIC Threshold
RotorGene-Q	0.3	0.3

Click **genotypes** button and set:

	Reacting Channel	Reacting Channel
Wild type	cycling A (green)	
Heterozygous	cycling A (green)	cycling A (yellow)*
None		
Mutant		cycling A (yellow)*

In the **Analysis window** select "Scatter graph analysis" sheet, click on **Cycling A green-Cycling A yellow** and click "show". The scatter plot will appear with the 2 cloudes: Wild Type (bottom right) and Mutant (high left). Is it possible to print the final report.

### ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

At the end of the PCR run open the "Analysis" window. Select the "Quantification" sheet and click on "cycling A (green)". Select from the menu "Dynamic Tube" and subsequently "Slope correct". Set the correct setting of the threshold in the space provided "CT calculation – Threshold".

	MTHFR C677T -WT – Green (FAM) Threshold
RotorGene-Q	0.05

Open the "Analysis" window. Select the "Quantification" sheet and click on "cycling A (yellow)". Select from the menu "Dynamic Tube" and subsequently "Slope correct". Set the correct setting of the threshold in the space provided "CT calculation – Threshold".

	MTHFR C677T -MUT – Yellow (VIC) Threshold
RotorGene-Q	0.05

Export data to Excel, save the file as "Excel Analysis Sheet" and enter the following formula for each sample and control:

**Yellow Ct (MTHFR C677T -MUT) – Green Ct (MTHFR C677T -WT)**

See paragraph "INTERPRETATION OF RESULTS"

### Versant kPCR AD or Stratagene MX3005P/MX3000P

#### Scatter Plot Analysis

Click on button "Analysis" in the toolbar. The software immediately open sheet "Analysis Term Setting". Activate the keys FAM and HEX in below screen and select the samples and disposable controls. On the window "Analysis Term Setting" open the sheet "Results". Select in the right screen the voice "Amplification Plot" in the area "Area to Analysis". Set the values in the area "Threshold Fluorescence" the values:

	MTHFR C677T -WT - FAM Threshold	MTHFR C677T - MUT - VIC Threshold
Kpcr Versant	0.5	0.35

Select in the right screen the voice "Dual Color Scatter Plot" in the area "Area to Analysis". Select "Fluorescence" in the area "Display value for" and choose from menu "Rlast/Rfirst".

Select in the area "Allele Association":

A Allele: **MTHFR C677T Wild Type**

B Allele: **MTHFR C677T Mutated**

At the end it will be possible obtain the detailed account, clicking "Text Report" in the area "Area to Analysis".

### Only for Versant kPCR AD or Stratagene MX3005P/MX3000P

If the software can't discriminate the wild type and mutated positive controls click "Show all genotypes". It will be now possible to modify the genotyping windows, select the samples placed near the positive wild type and identify them as wild type. Select the samples placed near the positive mutated and identify them as mutated. Select the samples placed in an intermediate position between wild type positive control and mutated positive control and identify them as Heterozygote for the C677T Mutation.

### ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

Click "Analysis" in the toolbar.

Click the sheet "Results"; and choose the analysis "Amplification plot". Check the correct setting of the threshold in the window "Threshold fluorescence" and set the following values:

	MTHFR C677T -WT - FAM Threshold	MTHFR-MUT C677T - HEX Threshold
MX3000P – MX3005P	0.1	0.1

From the Text Report window you can export the results by clicking on the main menu: file, export  
Export data to Excel and set the formula for each sample and control:

**HEX Ct (MTHFR C677T -MUT) – FAM Ct (MTHFR C677T -WT)**

See paragraph "INTERPRETATION OF RESULTS"

### CFX96 real time PCR system

#### Scatter Plot Analysis

At the end of the PCR, select the "Allelic Discrimination" sheet.

On the bottom of the screen, set "Selected fluorophores": X = FAM and Y = VIC. Choose RFU from "Display Mode" and "Normalize data". Set the following Threshold:

	MTHFR C677T -WT - FAM Vertical Threshold	MTHFR C677T -MUT - VIC Horizontal Threshold
CFX 96	0.15	0.45

You can export the report pushing the paper block figure on the top of the screen.

### ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

At the end of the PCR, select the sheet "Quantitation". Analyzing a fluorophore at a time, open in the menu "Setting" the voice Baseline Threshold and set the following values of the Threshold in the "User Defined"

	MTHFR C677T -WT - FAM Threshold	MTHFR C677T -MUT - VIC Threshold
CFX 96	50	50

You can export the report by clicking on the figure notepad at the top of the screen.

Export data to Excel, save the file as "Excel Analysis Sheet" and enter the following formula for each sample and control:

**VIC Ct (MTHFR C677T -MUT) – FAM Ct (MTHFR C677T -WT)**

See paragraph "INTERPRETATION OF RESULTS"

### INTERPRETATION OF RESULTS

The use of positive and negative control in each amplification session allows to verify the correct functioning of the amplification mix and the absence of any contamination.

The instrument software is able to analyze the fluorescences that are emitted by the specific probe for MTHFR C677T Wild-Type (FAM) and by the specific probe for the MTHFR C677T mutated (VIC/HEX).

### Scatter Plot Analysis

A regular functioning of the amplification mix can be verified analyzing the correct position of positive controls and negative controls on the scatter plot.

- Positive control Homozygote Wild-Type: horizontal position on X axis (down on the right) Ct< 30
- positive Control homozygote Mutated: vertical position on Y axis (up on the left) Ct< 30
- negative Control: placed at the origin of cartesian plane (down on the left)

Genotyping tests (allelic discrimination) are endpoint experiments: fluorescence data are collected at the end of the reaction (Post PCR Read) and subtracted to initial read fluorescence (Pre PCR Read).

The software makes a scatter plot with obtained results: Y axis is the normalized fluorescence of Mutated Allele, while X axis shows the normalized fluorescence of Wild-Type Allele. The diagnosis obtained with the comparison between unknown samples and Homozygote Wild-Type and Homozygote Mutate, given by the kit. Selecting the controls, we obtain their disposition on the scatter plot, depending of their relation between fluorescence emitted by two probes FAM (Wild-Type) and VIC/HEX (Mutated).

### ΔCt Analysis

By analyzing ΔCt you can identify the correct genotype of the sample being analyzed by performing subtraction of the assigned Mutated allele Ct to assigned wild-type allele Ct. The genotype is determined by following the table below:

Genotype	ΔCt (Ct Mutated - Ct Wild Type)
Wild Type	ΔCt > 2
Mutated	ΔCt < -2
Heterozygous	-2 > ΔCt < 2

In the amplification of each positive control (Wild-Type and Mutated), the Ct values of the allele-specific probe are used to validate the assay. Make sure that the fluorescence emitted by the amplification allele has a correct Ct <27. If the result of the allele-specific amplification of each control has a Ct > 27 or undetermined the session can not be considered valid and must be repeated.

In the amplification of each sample, the Ct values of the allele-specific probe are used to validate the assay from the extraction process up to the stage of detection. Make sure that the fluorescence emitted by the allele-specific amplification of the sample, identified after analysis, has not a Ct > 27.

Sample	Wild Type Allele (FAM)	Mutated Allele (VIC)	Assay and Genotype
Wild Type	Ct < 27	Not relevant	Valid
Mutated	Not relevant	Ct < 27	Valid
Heterozygous	Ct < 27	Ct < 27	Valid

If a sample has a Ct > 27 means that there are problems in the extraction phase or in the amplification and therefore could be assigned to a sample genotype wrong. Repeat sample.

Sample	Wild Type Allele (FAM)	Mutated Allele (VIC)	Assay and Genotype
Wild Type	Ct > 27	Not relevant	Not Valid
Mutated	Not relevant	Ct > 27	Not Valid
Heterozygous	Ct > 27	Ct > 27	Not Valid

### PERFORMANCES

#### Clinic Sensitivity:

For the purposes of this evaluation it is considered as clinical sensitivity the skill of determining true heterozygote and homozygote mutated samples in the totality of screened samples. The test is performed following the method advices. Positive samples are confirmed with another CE disposable method.

Samples	Heterozygote		Homozygote Mutated	
	Pos	Neg	Pos	Neg
Donors EDTA blood	50	0	26	0

Clinical sensitivity is 100% for material extracted from EDTA blood.

#### Diagnostic Specificity:

For the purposes of this evaluation it is considered as diagnostic specificity the skill of the method of determining real negative samples. The diagnostic specificity of the system is valued analyzing human genomic samples tested and confirmed as negative with another disposable system.

Samples	N	Positive	Negative
Donors EDTA blood	41	0	41

#### Analytical specificity:

Test's specificity is guaranteed by the use of specific primers for determining MTHFR gene and of probes intentionally designed on C677T mutation.

The alignment of the choose regions for specific primers' hybridization with available sequences of present in database, demonstrated: their conservation and the complete specificity for the analyzed targets. Samples that are defined as positive for a determined genotype as much must be recognized by the amplification system.

#### INTERFERENCES:

Verify that in DNA extracted from the sample there aren't nucleoproteins and haemoglobin, in way of exclude possible inhibition of PCR reactions. The interference due to contaminants can be highlighted through the spectrophotometric analysis and obtained data report at 260 nm (maximum absorption of Nucleic Acids) and 280 nm (maximum absorption of Proteins). A pure DNA might have a rate of approximately 1.8.

#### QUALITY CONTROL

It is therefore recommended to insert as quality control of every extraction session, amplification and detection of a negative sample and of a positive sample which have already tested before or referential material with known concentration.

In accordance with the Clonit srl ISO EN 13485 Certified quality Management System, each lot of MTHFR C677T is tested against predetermined specification to ensure consistent product quality

#### BIBLIOGRAPHY

Soo-Sang Kang, Paul W. K. Wong, Armando Susmano, Judy Sora, Marija Norusis, and Neal Ruggie. Thermolabile Methylene-tetrahydrofolate Reductase: An Inherited Risk Factor for Coronary Artery Disease. Am. J. Hum. Genet. 48:536-545, 1991.

Magali Louis, Anne France Dekairelle and Jean-Luc Gala. Rapid combined genotyping of factor V, prothrombin and methylenetetrahydrofolate reductase single nucleotide polymorphisms using minor groove binding DNA oligonucleotides (MGB probes) and real-time polymerase chain reaction. Clin Chem Lab Med 2004;42(12):1364-1369.

#### TECHNICAL ASSISTANCE

For any question and support please contact our Technical support:

e-mail: [info@clonit.it](mailto:info@clonit.it)  
phone: +39 02 56814413

	In vitro diagnostic device
	Read the instruction's manual
	Range of temperature
	Use within (dd/mm/yyyy: year-month)
	Lot (xxxx)
	Code
	Manufacturer
	Contains sufficient for <n> tests

EDMA code: 16010490  
CND: W0106010499

The **MTHFR C677T** kit is CE marked diagnostic kit according to the European *in vitro* diagnostic directive 98/79/CE.

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## BIBLIOGRAFIA

Soo-Sang Kang, Paul W. K. Wong, Armando Susmano, Judy Sora, Marija Norusis, and Neal Ruggie. Thermolabile Methylenetetrahydrofolate Reductase: An Inherited Risk Factor for Coronary Artery Disease. Am. J. Hum. Genet. 48:536-545, 1991.

Magali Louis, Anne France Dekairrelle and Jean-Luc Gala. Rapid combined genotyping of factor V, prothrombin and methylenetetrahydrofolate reductase single nucleotide polymorphisms using minor groove binding DNA oligonucleotides (MGB probes) and real-time polymerase chain reaction. Clin Chem Lab Med 2004;42(12):1364–1369.

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### CLONIT S.r.l.