



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Elite™ Cell Cycle Assay Kit (Green Fluorescence)

CATALOG NUMBER: CA-C106, 100 assays

Description

The cell cycle has four sequential phases: G0/G1, S, G2, and M. During a cell's passage through cell cycle, its DNA is duplicated in S (synthesis) phase and distributed equally between two daughter cells in M (mitosis) phase. These two phases are separated by two gap phases: G0/G1 and G2. The two gap phases provide time for the cell to grow and double the mass of their proteins and organelles. They are also used by the cells to monitor internal and external conditions before proceeding with the next phase of cell cycle. The cell's passage through cell cycle is controlled by a host of different regulatory proteins.

This kit is designed to monitor cell cycle progression and proliferation by using our proprietary Nuclear Green™ LCS1 in permeabilized and fixed cells. The percentage of cells in each sample that are in G0/G1, S and G2/M phases, as well as the cells in the sub-G1 phase prior to apoptosis can be determined by flow cytometry. Cells stained with Nuclear Green™ LCS1 can be monitored with a flow cytometer at Ex/Em = 490 nm/520 nm (FL1 channel).

Features

Optimized for flow cytometry

Kit Components

- **Component A:** 200x Nuclear Green™ CCS1 1 vial (250 µl)
- **Component B:** Assay Buffer 1 bottle (50 ml)

Storage

Keep in freezer (-20 °C) and avoid exposure to light.

Assay Protocol

1. For each sample, prepare cells in 0.5 mL of warm medium or buffer of your choice at a density of 5×10^5 to 1×10^6 cells/ml.
Note: Each cell line should be evaluated on an individual basis to determine the optimal cell density for apoptosis induction.
2. Treat cells with test compound for a desired period to induce apoptosis or other cell cycle functions.
3. Add 2.5 µl of 200X Nuclear Green™ CCS1 (**Component A**), and incubate the cells in a 37 °C, 5% CO₂ incubator for 30 to 60 minutes.
Note 1: For adherent cells, gently lift the cells with 0.5 mM EDTA to keep the cells intact, and wash the cells once with serum-containing media prior to incubation with Nuclear Green™ CCS1.
Note 2: The appropriate incubation time depends on the individual cell type and cell concentration used. Optimize the incubation time for each experiment.
Note 3: It is not necessary to fix the cells before DNA staining since the Nuclear Green™ CCS1 is cell-permeable.
4. Optional: Centrifuge the cells at 1000 rpm for 4 minutes, then re-suspend the cells in 0.5 ml of assay buffer (**Component B**) or the buffer of your choice.
5. Monitor the fluorescence intensity by flow cytometry using the FL1 channel (Ex/Em = 490/525 nm). Gate on the cells of interest, excluding debris.



Data Analysis




Figure 1. DNA profile in growing and camptothecin treated Jurkat cells. Jurkat cells were treated without (red) or with 20 μ M camptothecin (blue) in a 37 °C, 5% CO₂ incubator for about 8 hours, and then dye loaded with Nuclear Green™ CCS1 for 60 minutes. The fluorescence intensity of Nuclear Green™ CCS1 was measured with a FACSCalibur (Becton Dickinson, San Jose, CA) flow cytometer using the FL1 channel. In growing Jurkat cells, nuclear stained with Nuclear Green™ CCS1 shows G1, S, and G2 phases (red). In camptothecin treated apoptotic cells (B), the fluorescence intensity of Nuclear Green™ CCS1 was decreased, and both S and G2 phases were diminished.