



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

(FOR RESEARCH USE ONLY. DO NOT USE IT IN CLINICAL DIAGNOSIS !)

FITC Labeling Kit

Catalog No: E-LK-F003

3 Reactions

This manual must be read attentively and completely before using this product.

If you have any problems, please contact our Technical Service Center for help.

Phone: 240-252-7368(USA) 240-252-7376(USA)

Email: techsupport@elabscience.com

Website: www.elabscience.com

Please kindly provide us the lot number (on the outside of the box) of the kit for more efficient service.

Introduction

The FITC Labeling Kit of Elabscience offering a collection of reagents required for FITC labeling is designed to label antibody with amidogen (NH₂-). The specially treated FITC in this kit can be used directly. The reagents are enough for approximately 3 labeling reactions each containing 0.1-2 mg of antibody or other protein. Each kit includes 6 Filtration tubes for desalination of antibody labeling without the need for dialysis. The whole procedure is simple and can be completed in 100min with proficient operation.

Product Features

1. **All-inclusive:** This kit provides all the reagents required for FITC labeling.
2. **Quick:** The whole procedure takes only 100 min.
3. **Convenient:** Desalination can be achieved with Filtration tube, dialysis or gel filtration is not necessary.
4. **Flexible:** The procedure can be easily adapted to both smaller and larger scales, with 0.1-2 mg of protein labeled each time.
5. **Perfect results:** This kit has been optimized to determine the optimum labeling ratio of FITC to antibody, lowering the possibility of protein inactivation resulted from insufficient labeling or excess FITC labeling.

Product component

Component	Amount
FITC	0.3 mg ×3 vials
Labeling Buffer	10 mL×1 vial
DMF	500 μL×1 vial
Filtration tube	0.5 mL×6

Materials required but not included in this kit

1. Adjustable high-precision transferpettor (10 μL, 50 μL, 200 μL, 1000 μL).
2. Incubator (37 °C).
3. Centrifuge (the centrifugal speed can be up to 12,000×g).

Storage

This kit can be stored at 2-8 °C for 1 year before opening.

Principle

Protein bonds covalently to FITC by thiourea connection as free amino groups of lysine residues in the protein molecule can react with nucleophilic FITC.

Calculation on the amount of Reactive FITC

The volume of FITC used in each reaction depends on the amount and concentration of the protein to be labeled. When labeling 2mg / ml of antibody (IgG, 150KD), taking antibodies and FITC with the mass ratio of 1 mg: 150 ug can achieve better marking effect.

Preparation before experiment

1. Read the manual carefully.
2. Calculate the volume of FITC to be added.
3. Bring the kit to room temperature for 20 min before experiment (**Note:** The unused FITC should be stored in the refrigerator).
4. Dissolve FITC: add 30 μ L of DMF to the vial of FITC, let it stand for 10 min until it dissolved fully. The concentration of FITC is 10 mg/mL now.

Assay procedure (we label 1 mg of protein in this assay):

1. Add 1 mg of protein sample and corresponding volume of Labeling Buffer to a Filtration tube to make the total volume is 0.5 mL. Centrifuge at 12,000 \times g for 10 min.
Note: ①The maximum volume of Filtration is 0.5 mL. ②The protein sample can be treated with centrifugal ultrafiltration first when at low concentration.
2. Add 15 μ L of FITC and appropriate volume of Labeling Buffer to the Filtration tube, making the final concentration of the protein solution is 2 mg/mL. Mix it thoroughly with a pipette and incubate the tube for 60 min at 37 $^{\circ}$ C.
3. Centrifuge at 12,000 \times g for 10 min.
4. Add appropriate Labeling Buffer to the Filtration tube to make the total volume is 0.5mL. Mix it thoroughly with a pipette and centrifuge at 12,000 \times g for 10 min. Repeat this step once again.
5. Add 0.2 mL of Labeling Buffer to the Filtration tube and mix it thoroughly with a pipette. Invert the filtration tube and put it into another centrifugal tube. Centrifuge at 6,000 \times g for 10 min.
6. Collect the solution in the centrifugal tube, namely antibody labeled by FITC.

Precautions

1. This kit can be also used to label antigen, HRP and polypeptides with amidogen (NH₂-). The labeling ratio depends on the amount of amidogen.
2. DMF should be preserved airtight in a dry place. Seal it with the parafilm immediately after use.
3. In the Step 5 above, Labeling Buffer is used to collect the labeled protein. You can also use other buffer or protective agents as you like.
4. This kit can be stored for 1 year before opening. Please use it within the expiration date.
5. The Filtration tube provided in this kit has a molecular weight cutoff (MWCO) of 10KD. So please be careful of the molecular weight of the antigen or polypeptide to be labeled.
6. In the Step 2 above, for other quality antibodies, the final concentration of antibody should be controlled to 2 mg/mL strictly, then calculate the volume of Reactive FITC required according to the quantity of the antibodies.
7. Ensure no free amino-group being introduced during the coupling process. (Tris, ammonia and sodium azide can react with the activated FITC, thus reducing the conjugation rate of protein- FITC).