SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ImmuGlo™

Anti-Exocrine Pancreas Antibody (ExPA) Test System

For *in vitro* diagnostic use

1194 40 Determinations

INTENDED USE

An indirect immunofluorescence test for the detection and semi-quantitation of exocrine pancreas antibodies in human serum to aid in the diagnosis of *Crohn's Disease (CD)*.

SUMMARY AND EXPLANATION

Crohn's Disease (CD) and ulcerative colitis (UC), collectively known as inflammatory bowel disease (IBD), are relatively common chronic inflammations of the gastrointestinal (GI) tract. Even though they are histopathologically and anatomically distinct, they have many common symptoms and are hard to differentiate. Because the treatment and prognosis are different for the two conditions, it is important to identify serological markers highly specific for CD and UC, respectively¹.

Exocrine pancreas autoantibodies (ExPA) have been suggested as a very specific Crohn's Disease serological marker. Clinical studies reported a prevalence of 30-50 % in CD patients²⁻¹¹. In spite of the relatively low sensitivity, ExPA detects a subpopulation of CD patients that are negative for other CD markers (ASCA IgG, ASCA IgA), thus being very useful in combination with other IBD tests.

ExPA are circulating antibodies that react with secretory granules in the cytoplasm of exocrine pancreas cells and their presence was not detected in normal individuals or in patients suffering from UC or other GI inflammations.

Significantly, ExPA were not only detected in CD patients but also in healthy first-degree relatives of CD patients that display a significantly increased risk of developing the disease⁴. Thus, ExPA testing may also have prognostic value for CD susceptibility.

PRINCIPLES OF PROCEDURE

Exocrine pancreas antibodies are detected by indirect immunofluorescence (IF). In this method, patients' sera are incubated on sections of monkey pancreas to allow binding of antibodies to the tissue substrate. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgG class are detected by incubation of these sections with a fluorescein-labeled conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters.

The presence of exocrine pancreas antibodies (ExPA) is demonstrated by a specific apple green reticulo-granular fluorescence pattern in the cytoplasm of the exocrine pancreas cells. The cytoplasm of the endocrine cells located in the islets of Langerhans does not present this staining pattern. The titer (the reciprocal of the highest dilution giving a positive reaction) is then determined by testing serial dilutions of patient sera¹².

EN REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious.⁷ Follow good laboratory practices in storing, dispensing and disposing of these materials.

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. NaN₃ is toxic if ingested. Report incidents immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

Materials provided

[REF] 1194 Exocrine Pancreas Antibody Test

Kits contains sufficient reagents to perform 40 determinations.

10 x	[SORB SLD 4]	4 well monkey pancreas substrate slide
1 x 0.5ml	[CONTROL + EXPA]*	Ready to use Positive Control . Contains human serum positive for exocrine pancreas antibodies.
1 x 0.5ml	[CONTROL -]*	Ready to use Negative Control . Contains human serum.
1 x 5ml	[IgG-CONJ FITC EB]*	Anti-human IgG FITC Conjugate containing Evan's Blue. Protect from light.
1 x 60 ml	[BUF]*	Buffered Diluent . Ready for use.
2 x vials	[BUF WASH]	Phosphate Buffered Saline. Powder Wash Buffer . Reconstitute to one liter each.
1 x 5ml	[MOUNTING MEDIUM]*	Mounting Medium . Do not freeze.
12 x	[COVER SLD]	Coverslips
1 x		Report form

Optional Components

1 x 5ml	[IgG-CONJ FITC]*	[REF] 2099X Anti-human IgG FITC Conjugate . Protect from light.
1 x 1ml	[EVANS]*	[REF] 2510 Evans Blue Counterstain

* Contains <0.1% NaN₃

Symbols used on labels

[LOT] Lot number

Catalog number

In vitro diagnostic use

Use by

Storage temperature

Read instructions before use

Number of tests

Manufacturer

Materials Required But Not Provided

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes
- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Absorbent paper towels
- Incubation chamber

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year.

PROCEDURE

Test Method

A. Screening:

- Step 1.** Dilute each patient serum **1:5** with the Buffered Diluent provided (0.1 ml serum + 0.4 ml diluent). Screening at more than one dilution helps to avoid a “prozone phenomenon.” For screening, **DO NOT dilute the Positive or Negative Controls**. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
- Step 2** Allow pouches containing substrate slides to equilibrate to room temperature for **10-15 minutes**. Carefully remove the slides from their pouch without touching the substrate.
- Step 3** Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
- Step 4** Using a micropipette or Pasteur pipette, apply **1 drop** (approximately 50 µl) of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
- Step 5** Using a micropipette or Pasteur pipette, apply **1 drop** of patient’s diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.

- Step 6** Place the lid on the incubation chamber and incubate slides **1 hour at room temperature**.
- Step 7** Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately **10 ml** of PBS using a pipette, or rinse slide in a beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash **10 minutes**. Repeat process with all remaining slides.
- Step 8** Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the **goat anti-human IgG FITC conjugate** dropper vial and gently squeeze to apply **1 drop** (approximately 50 µl) to each well. Repeat process with all remaining slides.
- Step 9** Replace the lid on the incubation chamber. Incubate **30 minutes** at room temperature.
- Step 10** Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for **10 minutes**. Repeat process with all remaining slides. **NOTE:** Improper washing may lead to increased background fluorescence.
- Step 11** Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **While slide is still wet mount the coverslip.** Place **3 drops** of the Mounting Medium evenly spaced on a coverslip and invert the slide onto the coverslip. To remove any air bubbles gently apply pressure along the edge of the coverslip. Avoid any movement of the coverslip. Repeat process with all remaining slides.
- Step 12** Examine for specific fluorescence under a fluorescence microscope at a magnification of **200x** or greater.

Slides may be read as soon as prepared. However, because of the presence of an anti-fading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed. Slides should be stored in the dark at 2°- 8°C.

B: End Point Determination (Titration)

Serum found positive in the screening test may be further tested to determine the titer by following **Steps 5 through 12**. Each test run should include the undiluted Negative Control and **undiluted, 1:2, 1:4, 1:8 and 1:16** dilutions of the Positive Control. Serial two-fold dilutions of the patient's serum should be made as well starting at **1:5**. The reciprocal of the highest dilution producing a positive reaction is the titer.

Preparation of Serial Dilutions

For each serum positive for exocrine pancreas antibodies upon screening, number four tubes 1 to 4 and add **0.4 ml** of Buffered Diluent to tube 1 and **0.2 ml** to subsequent tubes. Pipette **0.1 ml** of undiluted serum to tube 1 and mix thoroughly. Transfer **0.2 ml** from tube 1 to tube 2 and mix thoroughly. Continue transferring **0.2 ml** from one tube to the next after mixing to yield the expected dilutions as depicted in the following table:

EN Tubes	1	2	3	4
Serum	0.1 ml +			
Buffered Diluent	0.4 ml	0.2 ml ↗	0.2 ml ↗	0.2 ml ↗
Transfer		0.2 ml	0.2 ml	0.2 ml
Final dilution	1:5	1:10	1:20	1:40

Quality Control

Both a positive and negative control serum should be included with each test run. The negative control should show no significant fluorescence in the cytoplasm of the pancreas cells. The positive control should stain the cytoplasm of the exocrine pancreas cells but not the one of islet cells. Upon titration of the positive control, an endpoint titer of $\pm$ doubling dilution from the indicated titer should be obtained. If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control.
- Problems with the optical system of the fluorescence microscope. These may include: improper alignment, use of the bulb beyond the expected performance life, etc.
- Allowing the slide to dry during the procedure.
- Improper preparation of serial dilutions of control.


RESULTS

The results of the tests for islet cell antibodies should be reported as negative (<5) or positive with titer of 5 or higher.


Read for specific staining of the cytoplasm of exocrine pancreas cells (Figure 1). Various other tissue antibodies such as antinuclear antibodies (ANA), mitochondrial antibodies, islet cell (endocrine pancreas) antibodies, and smooth muscle antibodies may also be observed on pancreas sections. Sera giving any of these reactions should be reported as negative for antibodies to exocrine pancreas. Any sera giving nuclear staining reactions may be tested with the Antinuclear Antibody (ANA) Test (HEp-2 Cells) or the Antinuclear Antibody (ANA) Test (Mouse Liver Sections). Any sera giving smooth muscle or mitochondrial staining reactions may be tested with the Autoantibody Test System (Mouse Kidney/Stomach Sections). Any sera giving cytoplasmic islet cell staining may be additionally tested with Islet Cell Antibody (ICAb) Test System.

Figure 1. Florescence microscopy staining patterns obtained with a normal control serum (left) and with an ExPA positive CD patient sera.

ExPA negative staining reaction

ExPA positive staining pattern (CD)

LIMITATIONS OF THE PROCEDURE

Occasionally sera may exhibit strong staining for ANA or other autoantibodies. These may interfere with the ability to detect exocrine pancreas antibodies. In such cases, titrating the serum may permit the visualization of the exocrine pancreas antibodies. In other cases their titer may be lower than the ANA or other antibodies and hence may not be detected. The ImmuGlo™ exocrine pancreas antibody test should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only.

EXPECTED VALUES

The incidence of exocrine pancreas antibodies in IBD patients and normal individuals has been studied by various investigators. The incidence of ExPA in Crohn's disease was found to be between 30% and 50% (see Table 1)²⁻¹¹. The majority of clinical studies found ExPA to be highly specific for Crohn's disease, especially at higher titers. Few studies reported an ExPA association with particular subtypes of CD. For example, Koutrobakis et al. reported that 60% of the stenotic CD patients tested were ExPA positive⁵. The incidence of ExPA also varies with patient ethnicity and age. A recent study reported a higher incidence in Chinese CD patients (46%) compared to Australian Caucasian CD patients (22%)². When these results were sub-classified on age groups, the incidence of ExPA in Chinese CD patients under 35 was even higher (67%).

In an in-house small-scale study on 25 CD patients, 25 UC patients and 25 normal controls, we found that 56% of the CD patients were ExPA positive, while all UC patients and normal controls were ExPA negative. We tested the CD and UC patient sera for other IBD markers (pANCA, ASCA IgG, ASCA IgA) and we found that ExPA testing was able to detect a subpopulation of CD patients that were ASCA negative (Figure 2). Thus, when combining the ExPA results with ASCA testing, the sensitivity of the assay was increased from 52% to 76%, stressing the importance of testing multiple markers for an accurate distinction between UC and CD (Table 2).


Table 1. Incidence of exocrine pancreas antibodies in Crohn’s Disease patients.

Study	Number of CD Patients	% Positive
Lawrance et al. 2005	50 Chinese Hong Kong	46
	50 Australian Caucasian	22
Klebl et al. 2005	208	29
Joossens et al. 2004	169	32
Koutroubakis et al. 2005	77	42
Seibold et al. 1991	222	31
Stocker et al. 1987	59	39
Folwaczny et al. 1998	76	38
Seibold et al. 1997	212	27
IMMCO Diagnostics study 2006	25	56

Table 2. Sensitivity and specificity of various IBD markers in the IMMCO Diagnostics study.

IBD Marker	Sensitivity (%)	Specificity (%)
pANCA for UC	52	92
ExPA for CD	56	100
ASCA for CD	52	92
ExPA + ASCA for CD	76	92

Figure 2. Prevalence of IBD markers in Crohn’s disease patients (IMMCO Diagnostics study, n=25).

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
REFERENCES

1. Bossuyt X Serologic markers in inflammatory bowel disease. *Clinical Chemistry*; 2006, 52:171-181.
2. Lawrance IC, Hall A, Leong R, Pearce C, Murray K. A comparative study of goblet cell and pancreatic exocrine autoantibodies combined with ASCA and pANCA in Chinese and Caucasian patients with IBD. *Inflamm Bowel Dis*; 2005, 10:890-897.
3. Klebl FH, Bataille F, Huy C, Hofstadter F, Scholmerich J, Rogler G. Association of antibodies to exocrine pancreas with subtypes of Crohn's disease. *European J Gastroenterol & Hepatol*; 2005, 17:73-77.
4. Joossens S, Vermeire S, Van Steen K, Godefroidis G, Claessens G, Pierik M et al. Pancreatic autoantibodies in inflammatory bowel disease. *Inflamm Bowel Dis*; 2004, 10:771-777.
5. Koutrobakis IE, Drygiannakis D, Karmiris K, Drygiannakis I, Makreas S, Kouromalis EA. Pancreatic autoantibodies in Greek patients with inflammatory bowel disease. *Dig Dis Sci*; 2005, 50:2330-2334.
6. Torok HP, Folwaczny C. Pancreatic autoantibodies in Crohn's disease: a feasible diagnostic tool. *European J Gastroenterol & Hepatol*; 2005, 17:37-39.
7. Seibold F, Weber P, Jenss H, Wiedmann KH. Antibodies to a trypsin sensitive pancreatic antigen in chronic inflammatory bowel disease: specific markers for a subgroup of patients with Crohn's disease. *Gut*; 1991, 32:1192-1197.
8. Stocker W, Otte M, Ulrich S, Normann D, Finkbeiner H, Stocker K et al. Autoimmunity to pancreatic juice in Crohn's disease. Results of an autoantibody screening in patients with chronic inflammatory bowel disease. *Scand J Gastroenterol Suppl*; 1987, 139:41-52.
9. Stocker W, Otte M, Ulrich S, Normann D, Stocker K, Jantschek G. Autoantibodies against the exocrine pancreas and against intestinal goblet cells in the diagnosis of Crohn's disease and ulcerative colitis. *Dtsch Med Wochenschr*; 1984, 109:1963-1969.
10. Folwaczny C, Noehl N, Endres SP, Loeschke K, Fricke H. Antineutrophil and pancreatic antibodies in first-degree relatives of patients with inflammatory bowel disease. *Scand J Gastroenterol*; 1998, 33:523-528.
11. Seibold F, Mork H, Tanza S, Muller A, Holzhuter C, Weber P, Scheurlen M. Pancreatic autoantibodies in Crohn's disease: a family study. *Gut*; 1997, 40:557-558.
12. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York; 1987, 3rd Ed, 3-40.
13. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health. [HHS Pub. No. (CDC) 1993, 93-8395].

For technical assistance please contact:

IMMCO Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120

Telephone: (716) 691-0091

Fax: (716) 691-0466

Toll Free USA/Canada: 1-800-537-TEST

E-Mail: info@immcodiagnostics.com

or your local product distributor

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante Autorizzato/Representante
Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,
The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

www.emergogroup.com

MAR2011

Document No. PI4194