# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation

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[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)


[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 


**Immco Diagnostics, Inc.**  
60 Pineview Drive  
Buffalo, NY 14228-2120 USA  
Telephone: (716) 691-0091  
Fax: (716) 691-0466  
Toll Free USA/Canada: 1-800-537-TEST  
E-Mail: [info@immco.com](mailto:info@immco.com)

*or your local product distributor*

EMERGO Group, Inc.  
Molenstraat 15, 2513 BH, The Hague,  
The Netherlands  
[www.emergogroup.com](http://www.emergogroup.com)

EN

Product Codes #: 1602, 2201, 2202, 2203, 2204, 2205

#### INTENDED USE

For comparing the specificity and indirect immunofluorescent staining intensity of antinuclear antibodies (ANA) reactions on nuclear antigens of HEp-2 cell lines or rodent kidney substrate.

#### SUMMARY AND EXPLANATION

The detection of ANA performed by indirect immunofluorescence aids in the diagnosis of connective tissue diseases. Different staining patterns indicate the presence of specific antibodies associated with different diseases<sup>1</sup>. While the homogeneous staining pattern is seen in a variety of connective tissue diseases, the peripheral and speckled patterns are associated with SLE, speckled and/or nucleolar patterns are associated with *systemic sclerosis* and the centromere pattern is a marker for CREST syndrome, a benign form of *acrosclerosis*.

#### PRINCIPLES OF PROCEDURE

In the indirect immunofluorescence method patient and control sera are incubated on an appropriate substrate. Antibodies not bound to the antigen in the substrate are removed by rinsing. An incubation with fluorescein-labeled (FITC), anti-human IgG conjugate detects the binding to specific antibodies. Following a wash step the slides are ready to be cover-slipped. When observed under a fluorescence microscope equipped with appropriate filters, positive reactions appear as apple green fluorescence with a distinct staining pattern.

#### REAGENTS

**Description:** Human serum containing homogeneous or speckled or centromere or nucleolar or peripheral antinuclear antibodies and <0.1 % NaN<sub>3</sub> in a dropper vial.

**Volume:** Each control vial contains 0.5 ml.

**Storage:** Store at 2-8°C.

**Instructions:** Ready to use after equilibration to room temperature. Do not use if turbid or if precipitate is present.

**Precautions:** For *in vitro* Diagnostic Use. All human derived materials have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials.

**WARNING - Sodium azide (NaN<sub>3</sub>)** may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center. Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange components with those from sources other than the same Product Code from IMMCO DIAGNOSTICS. Do not use beyond expiration date.

## PROCEDURE

### Materials Provided

*Product Codes #: Panel Controls 1602 or individual Controls #2201, 2202, 2203, 2204, 2205*

1x0.5 ml ImmuGlo™ ANA Homogeneous Positive Control (#2201)

1x0.5 ml ImmuGlo™ ANA Speckled Positive Control, (#2202)

1x0.5 ml ImmuGlo™ ANA Centromere Positive Control (#2203)

1x0.5 ml ImmuGlo™ ANA Nucleolar Positive Control, (#2204)

1x0.5 ml ImmuGlo™ ANA Peripheral Positive Control (#2205)

### Material Required but not Provided

ImmuGlo™ IFA Tests (HEp-2 Cells, *Product Codes #1102 or 1103* and/or rodent kidney substrate *Product Codes #1107, 1107-1 and 1107-2*) and other materials required for performing the assay as specified in product insert.

### Optional Quality Control Reagent

ImmuGlo™ Optical Standard Slide, *Product Code 2550OS*.

### Test Method

Invert dropper vial and apply one drop (approximately 50 µl) of one of the five positive control sera to a well. Apply one drop of Negative Control to another well and one drop of patient's serum to additional wells. Follow steps described in package insert of test procedure.

### RESULTS

Examine substrate under a fluorescence microscope at 200x or greater magnification for nuclear fluorescence and staining patterns.

Compare patient's serum pattern and fluorescence intensities with the respective ANA Positive Controls and grade results as follows:

+4 Brilliant apple green fluorescence

+3 Bright apple green fluorescence

+2 Clearly distinguishable, positive fluorescence

+1 Lowest specific fluorescence that allows nuclei to be clearly differentiated from background

### REFERENCES

1. Reimer G, Cornell RC and Tan EM. The biochemical nature of nuclear antigens reactive with antinuclear antibodies. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 519-531, 1987.

### **Warranty Statement**

*These products are warranted to perform according to the labeling and procedure described herein. Results may be altered by any changes or modifications in the procedure. IMMCO Diagnostics Inc. disclaims any implied warranty, merchantability or fitness for any other purposes. In no event shall IMMCO Diagnostics Inc. be liable for any consequential damages arising out of the above express warranty.*

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