SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ImmuLisa™

PM/Scl Antibody ELISA

IVD For *in vitro* diagnostic use

PRODUCT INSERT

REF 5101 PM/Scl Antibody ELISA 96 Determinations

INTENDED USE

Enzyme-linked immunosorbent assay (ELISA) for the qualitative and semi-quantitative detection of IgG antibodies to the nucleolar antigen PM/Scl in human serum to aid in the diagnosis of polymyositis, scleroderma or an overlap of polymyositis and scleroderma in conjunction with other laboratory and clinical findings.

SUMMARY AND EXPLANATION

Anti-nuclear antibodies (ANA) are a characteristic feature in the serum of patients with connective tissue disorders (CTD) such as systemic lupus erythematosus (SLE), scleroderma (Scl), polymyositis (PM), dermatomyositis (DM), or PM/Scl overlap. ANA are visualized by various reaction patterns using indirect immunofluorescence methods. Of these, the most common include homogenous, speckled, centromere and nucleolar patterns. Each reaction pattern represents antibodies associated with a particular antigen(s) and hence a disease or disease subset. For example, the centromere reaction is associated with a subset of scleroderma. Similarly, the nucleolar is associated with scleroderma, polymyositis, dermatomyositis and overlap of these conditions. The nucleolar reaction patterns are of three types: homogenous, coarse speckled or clumpy. Each reaction pattern represents antibodies of a particular specificity. The homogenous nucleolar reaction pattern is associated with antibodies to PM/Scl. PM/Scl complex, also commonly referred to as human exosome, consists of a complex of nine proteins and some other proteins that associate with this complex. This core complex is most abundant in the nucleolus.¹⁻⁴ This complex is involved in ribosomal RNA processing and messenger RNA degradation. Antibodies to PM/Scl complex are present in approximately 30% of patients with PM/Scl overlap, 8% of patients with PM, 11% of patients with DM, and 2% of patients with scleroderma.

Polymyositis is a subset of idiopathic inflammatory myopathies (IBM) characterized by chronic muscle weakness and skeletal muscle inflammation. Other subsets of idiopathic inflammatory myopathy are dermatomyositis and inclusion body myositis. IBM is a rare disorder with a reported incidence of 0.5-

EN

1 per 100,000. PM and DM can occur either in childhood or at age 30-50 with female preponderance (F:M 2.5:1). Like many other autoimmune diseases, polymyositis has genetic link and is related to DRB1.

Autoantibodies are important in the diagnosis and understanding of inflammatory myopathies. The antibodies to PM/Scl are associated with PM/Scl overlap and to some extent in patients with PM, DM or Scl, and rarely in SLE.

Antibodies to PM/Scl can be detected by indirect immunofluorescence on HEp-2 cells with a specific reaction pattern, by immunodiffusion using rabbit thymus extract, by immunoprecipitation and by ELISA. The antigens to which PM/Scl antibodies react are to PM/Scl-75, one of the nine core complex proteins of the human exosome, as well as to an associated protein, PM/Scl-100. Recently, a specific epitope to which the PM/Scl antibodies react has been identified. The use of the peptide rather than the whole molecule may provide higher levels of sensitivity and specificity as compared to the full length molecules because of cross reactivity, antigenic epitope exposure and epitope density. This assay benefits from these findings and incorporates a suitable antigen molecule that presents with the best discrimination in detecting antibodies to PM/Scl.

PRINCIPLES OF PROCEDURE

The test is performed as a solid phase enzyme linked immunosorbent assay (ELISA). Individual Microwells are coated with PM/Scl antigen followed by blocking the unreacted sites to reduce non-specific binding. Control, calibrator and patient serum samples are incubated in the antigen coated wells allowing specific antibodies present in the serum to bind to the antigen. Unbound-antibody and other serum proteins are removed by washing the microwells. HRP labeled anti-human IgG conjugate is added to each well. The unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by spectrophotometer at 450 nm. Results are expressed in ELISA units per milliliter (EU/ml) and reported as positive or negative.

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

EN

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.⁵

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

Materials provided

ImmuLISA™ PM/Scl Antibody ELISA

REF 5101

Kits contains sufficient reagents to perform 96 determinations.

12 x 8	MICROPLATE PM/SCL	Microplate with individual breakaway microwells. Coated with PM/Scl antigen. Ready for use.
1 x 1.75 ml	CONTROL + PM/SCL	Ready to use Positive Control (red cap). Contains human serum positive for PM/Scl antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL -	Ready to use Negative Control (white cap). Contains human serum.
5 x 1.75 ml	CALIBRATOR A PM/SCL CALIBRATOR B PM/SCL CALIBRATOR C PM/SCL CALIBRATOR D PM/SCL CALIBRATOR E PM/SCL	Ready to use set of 5 Calibrators . Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing PM/Scl antibodies. Concentrations in EU/ml are printed on the labels.
1 x 15 ml	IgG-CONJ HRP	HRP goat anti-human IgG Conjugate . Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. Protect from light.
1 x 15 ml	STOP H ₂ SO ₄	Stop Solution. Ready for use.
2 x	BUF WASH	Powder Wash Buffer. Reconstitute to one liter each.
1 x		Protocol Sheets

Optional Components

1 x 60ml	BUF WASH	Liquid concentrated Wash Buffer. Reconstitute to one liter.
----------	------------------	--

Symbols used on labels

LOT	Lot number
REF	Catalog number
IVD	In vitro diagnostic use

EN

Use by

Storage temperature

Read instructions before use

Number of tests

Manufacturer

Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year.

PROCEDURE

Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- ***Good washing technique is critical.*** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. ***An automated microplate washer is recommended.***

EN

- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

Test Method

- Step 1** Let all reagents and specimens equilibrate to room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.

or

For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

Qualitative				Semi-Quantitative			
A	Blank	S5	Etc.	A	Blank	S1	Etc.
B	-Control	S6		B	-Control	S2	
C	+ Control	S7		C	+ Control	S3	
D	Cal D	S8		D	Cal A	S4	
E	S1	S9		E	Cal B	S5	
F	S2	S10		F	Cal C	S6	
G	S3	S11		G	Cal D	S7	
H	S4	S12		H	Cal E	S8	
	1	2	3		1	2	3

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500µl** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.
- Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.

EN

Step 14 Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.

Step 15 Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <10 EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

RESULTS

Calculations

The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

Abs. of Test Sample

----- X EU/ml of Calibrator D = EU/ml Test Sample

Abs. of Calibrator D

It is recommended that qualitative results be reported as “positive” or “negative.” Sample results greater than or equal to Calibrator D are considered positive.

2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in EU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as “positive,” “negative,” or “indeterminate” with EU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods.

Interpretation

Interpretation values were determined by testing 64 normal blood donors. The cutoff was established using mean of the normal subjects plus 2 SD and

assigned an arbitrary value of 20 EU/ml. IMMCO suggests use of the reference range below. Each laboratory should validate assay values for their own conditions.

PM/Scl Ab value	Interpretation
<20 EU/ml	Negative
20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	Positive

Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml values, multiply the units obtained by the dilution factor.

LIMITATIONS OF PROCEDURE

The assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only. The results obtained serve only as an aid in the diagnosis. Taken alone, these results should not be interpreted as diagnostic. This assay has not been validated in a pediatric population.

Specimens positive for antiphospholipid antibodies (APL) potentially cross-react with the PM/Scl antigen. It is recommended that additional laboratory testing of positive samples include assays for APL.

EXPECTED VALUES

Test results in a normal population are usually expected to be negative. If results are positive, additional testing may be warranted to substantiate the serum antibody findings.

Incidence of PM/Scl Antibodies⁴

Disease Group	n	n Positive	% Positive
PM/Scl Group			
Polymyositis/Sceroderma	40	22	55
Polymyositis	40	3	7.5
Scleroderma	205	27	13
Disease Controls			
Rheumatoid Arthritis	69	0	0
Systemic lupus erythematosus	114	3	2.5
Rheumatic disease controls	452	33	7

PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ PM/Scl ELISA for the detection of autoantibodies was evaluated by testing nucleolar positive and negative serum specimens from suspected connective tissue disorder subjects, as well as disease controls and "normal" human sera. These specimens were also tested on another ELISA test kit for the detection of PM/Scl antibodies. These results are summarized below.

A. Immulisa™ PM/Scl ELISA vs. other PM/Scl Antibody kit:

		Other PM/Scl ELISA		
IMMCO PM/SCL ELISA		Positive	Negative	Total
	Positive	20	6	26
	Negative	3	56	59
	Total	23	62	85
Positive Percent Agreement:		87.0%		
Negative Percent Agreement:		90.3%		
Overall Percent Agreement:		89.4%		

B. Cross Reactivity: A total of 80 potentially cross-reactive specimens from individuals with other autoimmune disorders or positive for other autoantibodies were tested for PM/Scl antibodies using the Immulisa™ PM/Scl ELISA.

Condition	n	Positive
ENA positive*	48	5
Celiac Disease	8	2
Autoimmune Vasculitis	8	0
Hashimoto's Disease	8	1
Total	72	8 (11%)

* Antibodies to Extractable Nuclear Antigens

Precision

Precision was tested with positive specimens selected throughout the range of the assay. Assay runs of 3 replicates of each specimen were conducted on 3 days. Repeatability was determined with 6 replicates of each specimen.

S #	Mean (EU/ml)	Total Imprecision		Between days		Within run (Repeatability)	
		SD (EU/ml)	CV%	SD (EU/ml)	CV%	SD (EU/ml)	CV%
1	14	1.01	7.0%	1.03	7.2%	1.40	9.6%
2	29	1.20	4.1%	1.28	4.4%	0.91	3.2%
3	58	2.45	4.2%	3.36	5.8%	0.84	1.4%
4	116	7.30	6.3%	9.07	7.7%	0.19	0.2%
5	231	19.29	8.4%	26.07	11.1%	2.02	0.9%

EN

Linearity

Multiple assay runs were performed using doubling dilutions of samples throughout the linear range of the assay. Average r^2 value was .992 for these runs with no run lower than 0.985. Assay runs are deemed acceptable if r^2 exceeds 0.95.

Interference

Interference was studied by mixing sera with known PM/Scl antibody levels with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342 $\mu\text{mol/L}$), and Rheumatoid Factor (100 EU/ml).

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IMMCO Diagnostics, Inc.
EMERGO Group, Inc.