SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ImmuLisa™ AGA Enhanced ELISA

Gliadin IgA or IgG Antibody ELISA

IVD For *in vitro* diagnostic use
CLIA Complexity: High
CDC Analyte Identification Code: 0528
CDC Test System Identification Code: 28513/28514

PRODUCT INSERT

| | | | |
|------------|-------|---------------|-------------------|
| REF | 5117A | AGA IgA ELISA | 96 Determinations |
| REF | 5117G | AGA IgG ELISA | 96 Determinations |

INTENDED USE

Enzyme linked immunosorbent assays (ELISA) for the detection and semi-quantitation of IgA or IgG anti-gliadin antibodies in human serum to aid in the diagnosis of patients with celiac disease and dermatitis herpetiformis.

SUMMARY AND EXPLANATION

Gluten sensitive enteropathy (GSE), i.e. celiac disease (CD) and dermatitis herpetiformis (DH), is a common clinically heterogenous gastrointestinal disorder with sensitivity to gluten, and can exhibit with non-classic or minimal symptoms.¹ There is some genetic component associated with GSE which is illustrated by showing, that approximately 5-10% of first degree relatives have symptomatic or asymptomatic CD. Children with short stature and patients with insulin dependent diabetes and other autoimmune disorders also demonstrate a greater likelihood of developing GSE. Strict avoidance of gluten in the diet is recommended to control the disease activity and early diagnosis in such patients may improve their overall prognosis.²

Published literature suggests the use of serological testing for patients suspected of GSE and to monitor dietary compliance.³⁻⁶ The European Society of Pediatric Gastroenterology and Nutrition (ESPGAN) has recommended the inclusion of various serological tests to reduce the number of intestinal biopsies needed to make a diagnosis.⁷ These include tests for anti-gliadin (AGA), anti-reticulín (ARA) and anti-endomysial (EMA) antibodies. Of these, AGA have been studied the most.⁸⁻²⁰ Using ELISA, both IgA and IgG class AGA are detected in sera of patients with GSE. Of these, IgG class AGA seem more sensitive but less specific indicators of GSE than IgA class AGA. IgA class AGA on the other hand are less sensitive but more specific for GSE.

PRINCIPLES OF PROCEDURE


Gliadin antigen is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrators and diluted patient sera are added to separate wells, allowing any gliadin antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human IgA or IgG conjugate is added to each well. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the appropriate class. After washing away any

For technical assistance please contact:

IMMCO Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immcodiagnostics.com

or your local product distributor

EMERGO Europe
Molenstraat 15, 2513 BH, The Hague,
The Netherlands
Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com

SEP2012
Document No. PI5117

unbound conjugate, specific enzyme substrate (TMB) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in ELISA units per milliliter (EU/ml).

REAGENTS

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.²¹

Stop Solution is a dilute sulfuric acid solution. Sulfuric acid (H₂SO₄) is poisonous and corrosive. Do not ingest and avoid contact with skin and eyes. Avoid exposure to bases, metals or other compounds that may react with acids.

TMB Enzyme Substrate contains an irritant that may be harmful if inhaled, ingested or absorbed through the skin. Do not ingest and avoid contact with skin and eyes.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

Materials provided

ImmuLisa™ AGA IgA ELISA

REF 5117A

ImmuLisa™ AGA IgG ELISA

REF 5117G

Kits contains sufficient reagents to perform 96 determinations.

| | | |
|--------------------|-----------------------|---|
| 12 x 8 | MICROPLATE AGA | Microplate with individual breakaway microwells. Coated with gliadin antigen. Ready for use. |
| 1 x 1.75 ml | CONTROL+ AGA-A | Ready to use Positive Control (red cap) for REF 5117A. Contains human serum positive for gliadin IgA antibodies. The expected concentration range in EU/ml is printed on the label. |
| 1 x 1.75 ml | CONTROL+ AGA-G | Ready to use Positive Control (red cap) for REF 5117G. Contains human serum positive for gliadin IgG antibodies. The expected concentration range in EU/ml is printed on the label. |
| 1 x 1.75 ml | CONTROL- | Ready to use Negative Control (white cap). Contains human serum. |

EFERENCES • ΠΑΡΑΠΟΜΠΕΣ • BIBLIOGRAFÍA • REFERENZEN • REFERENCES • BIBLIOGRAFIA • REFERÊNCIAS

- Pare P, Douville P, Coron D, et al. Adult celiac sprue: changes in the pattern of clinical recognition. *J Clin Gastroenterol*; 10:395-400, 1988.
- Hall MJ, Cooper BT, Rooney H et al. Coeliac disease and malignancy of the duodenum: diagnosis by endoscopy, successful treatment of malignancy and response to gluten free diet. *Gut* 32:90-92, 1991.
- Rautonen J, Rautonen N and Saufihti E. Antibodies to gliadin in children with celiac disease. *Acta Paediatr Scand*; 1991, 80:1200-1205.
- Bürgin-Wolff A, Lentze MJ. Relation of anti-gliadin antibodies to gluten sensitive enteropathy B. Diagnostic significance of antibodies against gliadin in celiac disease in children. In serologic diagnosis of celiac disease eds. Chorzelski TP, Beutner EH, Kumar V and Zalewski TP. CRC Press, Boca Raton FL, 1990 pp. 59-75.
- Lerner A, Kumar V, Iancu TC. Immunological diagnosis of childhood coeliac disease: comparison between antigliadin, antireticulin and antiendomysial antibodies. *Clin Exp Immunol*; 1994, 95:78-82.
- Patinen P, Björkstén F, Malmström M et al. Salivary and serum antigliadin antibodies in dermatitis herpetiformis. *Eur J Oral Sci*; 1995, 103.
- Walker-Smith JA, Guandalini S, Schmitz J, et al. Revised criteria for diagnosis of celiac disease. Report of working group of European Society of Pediatric Gastroenterology and Nutrition. *Arch Dis Child* 1990; 65:909-911.
- Savilahti E, Perkkio M, Kalimo K, et al. IgA antigliadin antibodies: a marker of mucosal damage in childhood celiac disease. *Lancet* 1973; i:320-322.
- Kumar V, Jain N, Lerner A, Beutner EH, Chorzelski TP, Lebenthal E. Comparative studies of different gliadin preparations in detecting antigliadin antibodies. *J Pediatr Gastroenterol Nutr* 1986; 5:730-734.
- Montgomery AMP, Goka AKJ, Kumar PJ, et al. Low gluten diet in the treatment of adult coeliac disease: effect on jejunal morphology and serum anti-gluten antibodies. *Gut* 1988; 29:1564-1568.
- Weiss JB, Austin RK, Schanfield MS, et al. Gluten-sensitive enteropathy. IgG heavy-chain (Gm) allotypes and the immune response to wheat gliadin. *J Clin Invest* 1983; 72:96-101.
- Mearin ML, Konincox CR, Biemond I, et al. Influence of genetic factors on the serum levels of antigliadin antibodies in celiac disease. *J Pediatr Gastroenterol Nutr* 1984; 3:373-377.
- Levenson SD, Austin RK, Dietler MD, et al. Specificity of antigliadin antibody in celiac disease. *Gastroenterology* 1985; 89:1-5.
- Lerner A, Lebenthal E. The controversy of antigliadin antibody (AGA) as a diagnostic tool in celiac disease. *J Pediatr Gastroenterol Nutr* 1991; 12:407-409.
- Lebenthal E, Heitlinger LA. Gliadin antibodies in celiac disease. *J Pediatr* 1983; 102:711-712.
- Kumar V, Lerner A, Jain N, Beutner EH. Are antigliadin antibodies specific for celiac disease? *J Pediatr Gastroenterol Nutr* 1984; 3:815.
- Tucker NT, Barghuthy FS, Prihoda TJ, Kumar V, Lerner A, Lebenthal E. Antigliadin antibodies detected by enzyme-linked immunosorbent assay as a marker of childhood celiac disease. *J Pediatr* 1988; 113:286-289.
- Unsworth DJ, Kieffer M, Holborow EJ, et al. IgA anti-gliadin antibodies in coeliac disease. *Clin Exp Immunol* 1981; 46:286-293.
- Bürgin-Wolff A, Bertele RM, Berger R et al. A reliable screening test for childhood celiac disease: fluorescent immunosorbent test for gliadin antibodies. A prospective multicenter study. *J Pediatr* 1983; 102:655-660.
- Kelly J, O'Farrelly C, Rees JPR, et al. Humoral response to alpha gliadin as serological screening test for coeliac disease. *Arch Dis Child* 1987; 62:469-473.
- Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health; 1993 (HHS Pub No [CDC] 93-8395).
- Unsworth DJ. Serological diagnosis of gluten sensitive enteropathy. *J Clin Pathol*, 1996, 49:704-711.

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EN

| | | |
|--------------------|---|--|
| 5 x 1.75 ml | CALIBRATOR A AGA-A CALIBRATOR B AGA-A CALIBRATOR C AGA-A CALIBRATOR D AGA-A CALIBRATOR E AGA-A | Ready to use set of 5 Calibrators for REF 5117A. Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing AGA IgA. Concentrations in EU/ml are printed on the labels. |
| 5 x 1.75 ml | CALIBRATOR A AGA-G CALIBRATOR B AGA-G CALIBRATOR C AGA-G CALIBRATOR D AGA-G CALIBRATOR E AGA-G | Ready to use set of 5 Calibrators for REF 5117G. Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing AGA IgG. Concentrations in EU/ml are printed on the labels. |
| 1 x 15 ml | IgA-CONJ HRP | HRP goat anti-human IgG Conjugate for REF 5117A. Ready for use. Color coded pink. |
| 1 x 15 ml | IgG-CONJ HRP | HRP goat anti-human IgG Conjugate for REF 5117G. Ready for use. Color coded pink. |
| 1 x 60 ml | DIL | Serum Diluent. Ready for use. Color coded blue. |
| 1 x 15 ml | SUBSTRATE TMB | TMB enzyme substrate. Ready for use. Protect from light. |
| 1 x 15 ml | STOP H2SO4 | Stop Solution. Ready for use. |
| 2 x vials | BUF WASH | Powder Wash Buffer. Reconstitute to one liter each. |
| 1 x | | Protocol Sheets |

Optional Components

| | | |
|------------------|-----------------|--|
| 1 x 60 ml | BUF WASH | Liquid concentrated Wash Buffer. Reconstitute to one liter. |
|------------------|-----------------|--|

Symbols used on labels

| | |
|---|------------------------------|
| LOT | Lot number |
| REF | Catalog number |
| IVD | In vitro diagnostic use |
|  | Use by |
|  | Storage temperature |
|  | Read instructions before use |
|  | Number of tests |
|  | Manufacturer |

Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels

EN

- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°-8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year.

PROCEDURE

Procedural Notes

- Carefully read the product insert before starting the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

Test Method

- Step 1** Let all reagents and specimens equilibrate to room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.
- or

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EN

For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

| Qualitative | | | | Semi-Quantitative | | | |
|-------------|-------|-----|------|-------------------|-------|----|------|
| A | Blank | S5 | Etc. | A | Blank | S1 | Etc. |
| B | - | S6 | | B | - | S2 | |
| C | + | S7 | | C | + | S3 | |
| D | Cal D | S8 | | D | Cal A | S4 | |
| E | S1 | S9 | | E | Cal B | S5 | |
| F | S2 | S10 | | F | Cal C | S6 | |
| G | S3 | S11 | | G | Cal D | S7 | |
| H | S4 | S12 | | H | Cal E | S8 | |
| | 1 | 2 | 3 | | 1 | 2 | 3 |

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500 µl** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.
- Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within 30 minutes of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630 nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3 . The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be

<10 EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining units. While performing Qualitative determinations, the optical density of the Calibrator D must be greater than that of the negative control and lesser than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

RESULTS

Calculations

The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{Units of Calibrator D} = \text{Units of Test Sample}$$

It is recommended that qualitative results be reported as “positive” or “negative.” Sample results greater than or equal to Calibrator D are considered positive. Qualitative positive results should be confirmed with a semi-quantitative assay run to determine units and assess results in comparison with applicable reference ranges as indicated in **Interpretation**.

2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in EU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as “positive,” “negative,” or “indeterminate” with EU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods, as indicated in **Limitations of Procedure**.

Interpretation

Interpretation values were determined by testing 64 normal blood donors and disease control specimens. The mean of these subjects plus 3SD was established as the assay cutoff and assigned an arbitrary value of 20 EU/ml. The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

| Gliadin Ab Value | Interpretation |
|------------------|----------------------------|
| <20 EU/ml | Negative |
| 20-25 EU/ml | Indeterminate (Borderline) |
| >25 EU/ml | Positive |

Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml values, multiply the units obtained by the dilution factor.

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EN

LIMITATIONS OF PROCEDURE

This test should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only.

The results obtained serve only as an aid in the diagnosis and should not be interpreted as diagnostic in themselves. In IgA deficient CD patients IgA-AGA may be negative.

EXPECTED VALUES

The expected values in a normal population are negative. However, it has been determined that some apparently healthy, asymptomatic individuals may test positive for IgA or IgG anti-gliadin antibodies.

Sets of clinical samples were tested on the Immulisa™ AGA ELISAs. Results demonstrating incidence in the populations for this study are provided below.

| Disease | n | Gliadin IgA | | Gliadin IgG | |
|-------------------------------|-----|-------------|-------|-------------|-------|
| | | n Pos | % Pos | n Pos | % Pos |
| Celiac disease | 306 | 133 | 43.5% | 191 | 62.4% |
| IgA Deficient Celiac disease | 64 | 0 | 0.0% | 63 | 98.4% |
| Pediatric Celiac disease | 123 | 31 | 25.2% | 70 | 56.9% |
| Dermatis herpetiformis | 30 | 11 | 36.7% | 16 | 53.3% |
| Rheumatoid Arthritis | 30 | 1 | 3.3% | 0 | 0.0% |
| Thyroiditis | 30 | 1 | 3.3% | 0 | 0.0% |
| Systemic lupus erythematosus | 30 | 2 | 6.7% | 3 | 10.0% |
| Vasculitis | 30 | 0 | 0.0% | 1 | 3.3% |
| Systemic sclerosis | 30 | 0 | 0.0% | 0 | 0.0% |
| Ulcerative Colitis | 15 | 0 | 0.0% | 1 | 6.7% |
| Crohn's disease | 15 | 1 | 6.7% | 2 | 13.3% |
| Pediatric Normal Human Sera | 130 | 0 | 0.0% | 4 | 3.1% |
| Adult Normal Human sera | 139 | 7 | 5.0% | 6 | 4.3% |
| Total NHS | 269 | 7 | 2.6% | 10 | 3.7% |
| Total NHS and Disease Control | 449 | 12 | 2.7% | 17 | 3.8% |

PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ Gliadin Antibody ELISA was evaluated by testing Celiac disease patients alongside disease controls and "normal" human sera. These specimens were also tested on commercially available test kits. Only specimens in the linear range of the assay were included in the method comparison. These results are summarized below.

A. Immulisa™ AGA ELISAs vs. other AGA ELISAs:

| | | Other AGA IgA Ab ELISA | | |
|---------------------------------------|-------|------------------------|------------------------|-------|
| | | Pos | Neg | Total |
| IMMCO AGA IgA Ab ELISA | Pos | 87 | 24 | 111 |
| | Neg | 1 | 350 | 351 |
| | Total | 88 | 374 | 462 |
| Pos Agreement | | 98.9% | (95% CI 92.9% - 99.9%) | |
| Neg Agreement | | 93.6% | (95% CI 90.5% - 95.8%) | |
| Overall % Agreement | | 94.6% | (95% CI 90.7% - 95.3%) | |

| | | Other Gliadin IgG Ab ELISA | | |
|---------------------------------------|-------|----------------------------|------------------------|-------|
| | | Pos | Neg | Total |
| IMMCO AGA IgG Ab ELISA | Pos | 251 | 6 | 257 |
| | Neg | 3 | 156 | 159 |
| | Total | 254 | 162 | 416 |
| Pos Agreement | | 98.8% | (95% CI 96.3% - 99.7%) | |
| Neg Agreement | | 96.3% | (95% CI 91.8% - 98.5%) | |
| Overall % Agreement | | 97.8% | (95% CI 95.9% - 98.9%) | |

B. Cross Reactivity: A total of 160 sera from individuals with other potentially cross-reactive autoimmune disorders were selected to test for gliadin antibodies using the Immulisa™ assay. The results appear below.

| Disease | n | Gliadin IgA | | Gliadin IgG | |
|------------------------------|-----|-------------|-------|-------------|-------|
| | | n Pos | % Pos | n Pos | % Pos |
| Rheumatoid Arthritis | 30 | 1 | 3.3% | 0 | 0.0% |
| Thyroiditis | 30 | 1 | 3.3% | 0 | 0.0% |
| Systemic lupus erythematosus | 30 | 2 | 6.7% | 3 | 10.0% |
| Vasculitis | 30 | 0 | 0.0% | 1 | 3.3% |
| Systemic sclerosis | 30 | 0 | 0.0% | 0 | 0.0% |
| Ulcerative Colitis | 15 | 0 | 0.0% | 1 | 6.7% |
| Crohn's disease | 15 | 1 | 6.7% | 2 | 13.3% |
| Total | 180 | 5 | 2.8% | 7 | 3.9% |

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Precision

Precision was tested with positive specimens selected throughout the range of the assay. Assay runs of 6 replicates of each specimen were conducted on 13 days. Repeatability was determined with 12 replicates of each specimen. Results are shown below.

| AGA IgA | S # | Mean (EU/ml) | Total Imprecision | | Between days | | Within run (Repeatability) | |
|------------|-----|-----------------|-------------------|------|---------------|------|-------------------------------|------|
| | | | SD (EU/ml) | CV% | SD (EU/ml) | CV% | SD (EU/ml) | CV% |
| | 1 | 10.4 | 0.536 | 5.2% | 0.559 | 5.4% | 0.251 | 2.4% |
| | 2 | 15.8 | 0.725 | 4.6% | 0.766 | 4.8% | 0.281 | 1.8% |
| | 3 | 24.5 | 1.872 | 7.6% | 1.949 | 8.0% | 1.266 | 5.1% |
| | 4 | 45.9 | 1.700 | 3.7% | 1.680 | 3.7% | 1.269 | 2.7% |
| | 5 | 84.4 | 2.485 | 2.9% | 2.591 | 3.1% | 1.291 | 1.5% |
| | 6 | 116.3 | 5.932 | 5.1% | 6.174 | 5.3% | 2.653 | 2.2% |

| AGA IgG | S # | Mean (EU/ml) | Total Imprecision | | Between days | | Within run (Repeatability) | |
|------------|-----|-----------------|-------------------|-------|---------------|-------|-------------------------------|------|
| | | | SD (EU/ml) | CV% | SD (EU/ml) | CV% | SD (EU/ml) | CV% |
| | 1 | 9.3 | 0.846 | 9.1% | 0.892 | 9.6% | 0.429 | 4.6% |
| | 2 | 15.4 | 1.108 | 7.2% | 1.131 | 7.4% | 0.316 | 1.9% |
| | 3 | 24.4 | 2.461 | 10.1% | 2.615 | 10.7% | 0.801 | 3.3% |
| | 4 | 37.1 | 4.429 | 11.9% | 4.362 | 11.6% | 0.696 | 2.1% |
| | 5 | 63.1 | 8.408 | 13.3% | 8.177 | 12.7% | 4.027 | 7.4% |
| | 6 | 127.9 | 11.132 | 8.7% | 11.288 | 8.7% | 5.649 | 4.7% |

Reproducibility

80 replicates of samples in the low negative range, near the cutoff, in the moderate positive range, and approximately +/-20% of the assay cutoffs for each AGA isotype were performed to determine qualitative reproducibility. All specimens produced 100% qualitative agreement, with the exception of one IgA sample with a value of 24.4 EU/ml (98% qualitative agreement).

Limit of Detection

The limit of detection (LoD) was determined based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples. LoD was determined to be 1.8 EU/ml for IgA and 3.7 EU/ml for IgG.

Linearity and Recovery

Linearity and recovery were tested by diluting positive specimens through the assay range in equidistant dilutions and comparing actual vs. expected results. The linear range of the assay was determined be 1.8 (LoD) – 160 EU/ml for IgA and 3.7(LoD) – 160 EU/ml for IgG. Results are summarized below.

| IgA | Test Range (EU/ml) | Slope (95% CI) | Y-intercept (95% CI) | R ² | % Recovery (Obtnd/Expctd) |
|-----|-----------------------|-------------------------|-------------------------|----------------|------------------------------|
| | 1.2 to 41.6 | 1.03 (.959 to 1.109) | 0.4 (-1.4 to 2.2) | 0.995 | 95% to 110% |
| | 2.3 to 101.6 | 1.05 (.929 to 1.168) | 2.2 (-4.6 to 9.0) | 0.987 | 109% to 115% |
| | 4.7 to 273.4 | 0.95 (.887 to 1.014) | 3.0 (-7.8 to 13.8) | 0.996 | 93% to 114% |

EN

| IgG | Test Range (EU/ml) | Slope (95% CI) | Y-intercept (95% CI) | R ² | % Recovery (Obtnd/Expctd) |
|-----|--------------------|-------------------------|--------------------------|----------------|---------------------------|
| | 2.7 to 47.0 | 0.97 (.891 to 1.040) | 0.89 (-1.3 to 3.1) | 0.994 | 94% to 117% |
| | 12.8 to 116.8 | 0.95 (.957 to 1.046) | -0.24 (-7.5 to 7.0) | 0.991 | 90% to 96% |
| | 96.8 to 227.0 | 0.96 (.776 to 1.141) | -6.69 (-39.8 to 26.4) | 0.965 | 87% to 92% |

Interference

Interference was studied by mixing sera with known gliadin antibody levels for each isotype with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342 µmol/L), Rheumatoid Factor (100 EU/ml), Triglycerides (37 mmol/L), and Cholesterol (13 mmol/L).

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