**SZABO
SCANDIC**

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic

Optical Standard (OS) Slide

IVD

For Laboratory Use

PRODUCT INSERT

REF 2550 OS Optical Standard Slide 8 Wells

INTENDED USE

The Optical Standard (OS) Slide is intended to aid laboratories in the calibration and optimization of microscopes equipped for immunofluorescence studies.

SUMMARY AND EXPLANATION

Immunofluorescence (IF) microscope methods are simple, sensitive and specific techniques used in a wide variety of laboratory tests. Interpretation of IF reactions, however, are highly influenced by various factors, the most significant being the microscope itself. Microscope performance variance may result from different conditions that involve alignment of the optical system, selection of filter(s), age/type of light source and the objective selected. Any of these factors may significantly affect the observed fluorescence intensity emitted by a specimen. To minimize the effect of such variations on the reading and interpretation of IF reactions, NCCLS recommends the use of an optical standard slide.¹

Methods were proposed by T. Ooishi and N. Uyesaka to ameliorate disparities in microscope performance through the use of latex microspheres labeled with varying concentrations of FITC.² This fluorescein bead method offers applications in calibration, however, the substrate has proven difficult to read and interpret, especially by technicians used to semi-quantitation of antibody levels on tissues, cells and other biological substrates. The fluorescence from fluorophores attached to analytes can differ greatly from the same fluorophore immobilized on a microsphere.³ For these reasons Immco™ has established an optical standard based on a biological substrate.

An earlier version of the Immco™ OS Slide utilized microspheres of graduated intensity. However, studies of various OS Slide formats have concluded that calibration and standardization of optical systems for IF applications in the clinical laboratory is ideally performed using biological cells bound and labeled with graduated titers of antibodies. These studies, facilitated by image analysis software available from the National Institutes of Health,⁴ prompted a change in Immco™ OS Slide format.

The new and improved ImmuGlo™ Optical Standard Slide offers the laboratory a practical and innovative approach to optimize and monitor the performance of the fluorescent microscope using cells commonly employed in the clinical laboratory. Utilization of this method along with controls of known titer can significantly improve the reliability and inter-scope correlation of IF results.

PRINCIPLE OF PROCEDURE

The ImmuGlo™ Optical Standard Slide contains human polymorphonuclear leukocytes. They are fixed onto the glass in eight wells and exhibit graduated apple-green fluorescence intensities from four plus (++++) to negative (-) when observed under a fluorescent light microscope with fluorescein filter settings. This allows the microscopist to determine desired end point. Variation from this end point reflects sub-optimal performance of the fluorescent light microscope, prompting a corrective action.

REAGENTS

Storage and Cleaning

Keep slide in its cassette in the dark at 2-8°C. Before use, let it equilibrate for 15 minutes at room temperature. If the slide requires cleaning, carefully hold it at the edge and wipe both sides of slide with alcohol swab. Apply minimal pressure to the coverglass to prevent any lateral movement.

Materials Provided

ImmunoGlo™ Optical Standard slide [REF] 2550 OS

8-well microscope slide containing fixed human polymorphonuclear leukocytes of graduated fluorescence intensities.

Materials required but not provided

Microscope equipped with fluorescent light source and excitation filters.

PROCEDURE

STEP 1 Switch on fluorescent light microscope and let light source equilibrate for 15 minutes.

STEP 2 Place OS slide on microscope stage and read apple-green fluorescence intensity of neutrophil nuclei at 400x in each well. Please note normal light refraction when focusing on nuclei.

STEP 3 Determine end point and compare results with the following figure:


Figure 1. OS Slide Layout


Figure 2. Fluorescent Light Microscope

EXPECTED VALUES

Neutrophils in all wells should display apple-green nuclear reactions of varying intensities. Well number 1 should read four plus (+++). Intensity readings should gradually decrease as well numbers increase, reading negative (-) in well number 7. Endpoint (+/-) should be obtained in well number 6. Deviation of a microscope system by more than one well from the expected endpoint indicates inadequate performance of the fluorescent light microscope and/or interpretation subjectivity. Remedial actions are suggested below:

- Check and adjust optical alignment per microscope manufacturer's instructions
- Clean objective and eyepiece lenses
- Check filters. Excitation and barrier filters may not be appropriate for the fluorochrome used
- Check the microscope light source (bulb). Replace the bulb and recalibrate the microscope
- Consult with microscope manufacturer if adjustments are unsuccessful and problems persist

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REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Quality assurance for the indirect immunofluorescence test for autoantibodies to nuclear antigen (IF-ANA); approved guideline. I/LA2-A vol. 16, No.11, 1996.
2. Ooishi T et al. A new standard fluorescent microsphere for quantitative flow cytometry. J Immunolog Meth. 84, 143-154, 1985.
3. Gaigalas AK et al. The development of fluorescence intensity standards. J Res Natl Inst Stand Technol. 106, 381-389, 2001.
4. Refer to <http://rsb.info.nih.gov/ij/>
5. Wick et al. Immunofluorescence. University of Innsbruck, Austria, 1978.
6. Stout GW et al. Automation of an indirect fluorescent antibody test for syphilis. Ann N Y Acad Sci. 177:453-458, 1971.
7. Ten Veen JH et al. An approach toward automation of antinuclear antibody (ANA) determination. Ann N Y Acad Sci. 177:459-66, 1971.
8. Camargo ME et al. A microscopic immunofluorescence technique with soluble protein antigens fixed to cellulose particles. Int Arch Allergy Appl Immunol. 39 (2-3): 292-300, 1970.
9. Rostami R et al. Quantitative studies of immunofluorescent staining. Int Arch Allergy Immunol. 98:200-204, 1992.
10. Haiijman JJ et al. Microfluoremetric immunoassays with antigens bound to sepharose beads. VI international conference on immunofluorescence and related staining techniques. Ann NY Acad Sci. 254:137-148, 1975.

for technical assistance please contact:

IMMCO Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120 USA
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immco.com
or your local product distributor

EMERGO EUROPE
Prinsesegracht 20
2514 AP The Hague
The Netherlands