**SZABO  
SCANDIC**

Part of Europa Biosite

## Produktinformation

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[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](http://linkedin.com/company/szaboscandic)

## ImmunoLISA™

# ENA Antibody ELISA

### RNP Antibody ELISA

### Sm Antibody ELISA

**[IVD]** For *in vitro* diagnostic prescription use only

#### PRODUCT INSERT

<b>REF</b>	5126	RNP Antibody ELISA	96 Determinations
<b>REF</b>	5127	Sm Antibody ELISA	96 Determinations

#### INTENDED USE

An enzyme linked immunoassays (ELISA) for the detection and semi-quantitation of Sm or RNP antibodies in human serum as an aid in diagnosis of Systemic Lupus Erythematosus (SLE) and associated connective tissue disorders in conjunction with other laboratory and clinical findings.

#### SUMMARY AND EXPLANATION

ENA are soluble ribonucleoprotein (snurps) complexes. Autoantibodies directed against various ENA have proven to be of value in the diagnosis and monitoring of various systemic connective tissue diseases. Anti-Sm antibodies are disease specific and thus are a marker for systemic lupus erythematosus (SLE). Antibodies to Sm occur in approximately 30-40% of SLE patients. They are rare in other systemic connective tissue diseases and if present, indicate either overlap of disease or patients that have not yet fulfilled the ARA criteria for SLE.<sup>1-6</sup> Other antibodies such as those directed against SS-A (Ro), SS-B (La) and RNP are not disease specific.

Antibodies to RNP occur in 35-45% of SLE patients and in over 95% of patients with mixed connective tissue disease (MCTD). Occasionally they are also found in scleroderma, rheumatoid arthritis and drug induced LE.<sup>1-6</sup>

Patients with anti-RNP antibodies have a lower incidence of renal disease as compared to patients with anti-Sm antibodies. Antibodies to SS-A (Ro) and SS-B (La) occur in approximately 30-40% and 10-15% of SLE patients and 60-70% and 40-60% of patients with Sjögren's syndrome respectively.

These antibodies can be detected by various methods. ELISA methodology has many advantages over older technologies such as immunodiffusion: assay performance times are reduced, individual subjectivity in reading results is eliminated, quantitation is achieved without serum titration, there is potential for automation and ELISAs offer greater sensitivity.<sup>10</sup>

#### PRINCIPLES OF PROCEDURE

The test is performed as a solid phase immunoassay. Microwells are coated with purified Sm or RNP antigen followed by a blocking step to reduce non-specific protein binding during the assay run. Controls, calibrators and patient sera are incubated in the antigen coated wells to allow specific antibodies present in the serum to bind to the antigen. Unbound antibodies and other serum proteins are removed by washing the microwells. Bound antibodies are detected by adding an enzyme labeled anti-human IgG conjugate to the microwells. Unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in ELISA units per milliliter (EU/ml) and reported as positive or negative.

#### REAGENTS

##### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

## EN

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

## Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.<sup>11</sup>

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

## Materials provided

Kits contains sufficient reagents to perform 96 determinations.

ImmunoLISA™ RNP Antibody ELISA

[REF] 5126

12 x 8 [MICROPLATE|RNP]  
Microplate with individual breakaway microwells. Coated with RNP antigens. Ready for use.

1 x 1.75 ml [CONTROL+|RNP]  
Ready to use Positive Control (red cap). Contains human serum positive for RNP antibodies. The expected concentration range in EU/ml is printed on the label.

1 x 1.75 ml [CONTROL-]  
Ready to use Negative Control (white cap). Contains human serum.

5 x 1.75 ml [CALIBRATOR|A|RNP]  
[CALIBRATOR|B|RNP]  
[CALIBRATOR|C|RNP]  
[CALIBRATOR|D|RNP]  
[CALIBRATOR|E|RNP]  
Ready to use set of 5 Calibrators. Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing RNP antibodies. Concentrations in EU/ml are printed on the labels.

1 x 15 ml [IgG-CONJ|HRP]  
HRP goat anti-human IgG Conjugate. Ready for use. Color coded pink.

1 x 60 ml [DIL]  
Serum Diluent. Ready for use. Color coded purple.

1 x 15 ml [SUBSTRATE|TMB]  
TMB enzyme substrate. Ready for use. Protect from light.

1 x 15 ml [STOP|H2SO4]  
Stop Solution\*. Ready for use.

2 x vials [BUF|WASH]  
Powder Wash Buffer. Reconstitute to one liter each.

1 x Protocol Sheets

ImmunoLISA™ Sm Antibody ELISA

[REF] 5127

12 x 8 [MICROPLATE|SM]  
Microplate with individual breakaway microwells. Coated with Sm antigen. Ready for use.

1 x 1.75 ml [CONTROL+|SM]  
Ready to use Positive Control (red cap). Contains human serum positive for Sm antibodies. The expected concentration range in EU/ml is printed on the label.

1 x 1.75 ml [CONTROL-]  
Ready to use Negative Control (white cap). Contains human serum.

5 x 1.75 ml [CALIBRATOR|A|SM]  
[CALIBRATOR|B|SM]  
[CALIBRATOR|C|SM]  
[CALIBRATOR|D|SM]  
Ready to use set of 5 Calibrators. Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing Sm antibodies. Concentrations

EN

	CALIBRATOR ESM	in EU/ml are printed on the labels.
1 x 15 ml	IgG-CONJ HRP	HRP goat anti-human <b>IgG Conjugate</b> . Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. <b>Protect from light.</b>
1 x 15 ml	STOP H2SO4	Stop Solution*. Ready for use.
2 x	BUF WASH	Powder Wash Buffer. <b>Reconstitute to one liter each.</b>
1 x		Protocol Sheets

### Optional Components

1 x 60ml	BUF WASH	Liquid concentrated Wash Buffer. <b>Reconstitute to one liter.</b>
----------	----------	--

### Symbols used on labels

LOT	Lot number
REF	Catalog number
IVD	In vitro diagnostic use
	Use by
	Storage temperature
	Consult instructions for use
	Number of tests
	Manufacturer
	Date of Manufacture

\*Danger; Causes severe skin burns and eye damage. Causes severe eye damage. Wear protective gloves, protective clothing, eye protection, face protection. Wash exposed skin thoroughly after handling. If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician

### Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

### SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year.

### PROCEDURE

#### Procedural Notes

- Carefully read the product insert before starting the assay.

## EN

- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

**Test Method**

**Step 1** Let all reagents and specimens equilibrate to room temperature.

**Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.

**Step 3** For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.

or

For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

Qualitative

A	Blank	S5	Etc.
B	-Control	S6	
C	+ Control	S7	
D	Cal D	S8	
E	S1	S9	
F	S2	S10	
G	S3	S11	
H	S4	S12	

1      2      3

Semi-Quantitative

A	Blank	S1	Etc.
B	-Control	S2	
C	+ Control	S3	
D	Cal A	S4	
E	Cal B	S5	
F	Cal C	S6	
G	Cal D	S7	
H	Cal E	S8	

1      2      3

**Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500µl** of Serum Diluent.

**Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.

**Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.

**Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.

**Step 7** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.

**Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.

**Step 9** Pipette **100 µl** of Conjugate into microwells.

**Step 10** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.

**Step 11** Wash all microwells as in Step 8.

**Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.

**Step 13** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.

EN

- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

### Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <10 EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

## RESULTS

### Calculations

The concentrations of the patient samples can be determined by either of two methods:

#### 1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

It is recommended that qualitative results be reported as "positive" or "negative." Sample results greater than or equal to Calibrator D are considered positive.

#### 2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in EU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as "positive," "negative," or "indeterminate" with EU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods for detection of antibodies to ENA.

### Interpretation

Interpretation values were determined by testing 80 normal blood donors. The cutoff was established using mean of the normal subjects plus 2.5 SD for the RNP Ab ELISA and 3.5 SD for the Sm Ab ELISA. The cutoff was assigned an arbitrary value of 20 EU/ml. IMMCO suggests use of the reference range below. Each laboratory should validate assay values for their own conditions.

anti-RNP Ab value	anti-Sm Ab value	Interpretation
<20 EU/ml	<20 EU/ml	Negative
20-25 EU/ml	20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	>25 EU/ml	Positive

### Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml values, multiply the units obtained by the dilution factor.

## LIMITATIONS OF PROCEDURE

The assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only. The results obtained serve only as an aid in the

EN

diagnosis. Taken alone, these results should not be interpreted as diagnostic. This assay has not been validated in a pediatric population.

## EXPECTED VALUES

Test results in a normal population are usually expected to be negative. The incidence of ENA antibodies in various systemic connective tissue diseases is summarized in the following table:

### Diagnostic Significance of Antibodies to RNP and Sm Antigens

Antibody	Disease Association (% Incidence)
RNP	SLE (20-30%) MCTD (95-100%)
Sm	SLE (10-40%)

SLE = systemic lupus erythematosus

MCTD = mixed connective tissue disease

Note: The frequency of each antibody specificity in a disease represents compilation from the literature.<sup>3</sup> The incidence varies depending on the patient population.

## PERFORMANCE CHARACTERISTICS

The utility of the ImmuLisa™ RNP and Sm Antibody ELISAs were evaluated by testing well-characterized ENA antibody positive specimens alongside disease controls and "normal" human sera. These specimens were also tested on commercially available test kits. Only specimens in the linear range of the assay were included in the method comparison. These results are summarized below.

### A. Method Comparison: ImmuLisa™ RNP Antibody ELISA vs. other RNP Antibody kit:

Other RNP Antibody ELISA			
	Positive	Negative	Total
IMMCO	Positive	43	3
RNP	Negative	2	95
Ab ELISA	Total	45	98
			143

Positive Percent Agreement: 95.6% (95% CI 83.6% to 99.2%)

Negative Percent Agreement: 96.9% (95% CI 90.7% to 99.2%)

Relative Agreement: 96.5% (95% CI 91.6% to 98.7%)

### B. Method Comparison: ImmuLisa™ Sm Antibody ELISA vs. other Sm Antibody kit:

Other Sm Ab ELISA			
	Positive	Negative	Total
IMMCO	Positive	53	12
Sm Ab	Negative	2	167
ELISA	Total	55	179
			234

Positive Percent Agreement: 96.4% (95% CI 86.4% to 99.4%)

Negative Percent Agreement: 93.3% (95% CI 88.3% to 96.3%)

Relative Agreement: 94.0% (95% CI 90.0% to 96.6%)

### C. Cross Reactivity: A total of 80 potentially co-incident antibody positive or cross-reactive specimens were selected from individuals positive for other ENA antibodies or suffering from other autoimmune disorders were tested for RNP and Sm antibodies using the ImmuLisa™ assays.

Condition	n	RNP Ab Pos n (%)	Condition	n	Sm Ab Pos n (%)
Associated ENA Ab Pos					
Jo-1	8	0	Jo-1	8	0
La	8	0	La	8	0
Pm-Scl	8	1	Pm-Scl	8	0
Ro	8	0	Ro	8	0
Scl-70	8	4	Scl-70	8	1

EN				RNP		
Sm	8	1			8	3
Total	48	6 (12.5)			48	4 (12)
Potentially cross-reactive autoimmune diseases						
Celiac Disease	8	0	Celiac Disease	8	0	
Hashimoto's	8	0	Hashimoto's	8	0	
Rheumatoid Arthritis	8	0	Rheumatoid Arthritis	8	0	
Vasculitis	8	0	Vasculitis	8	0	
Total	32	0 (0)		32	0 (0)	

## Precision

Precision was tested with positive specimens selected throughout the range of the assay. Assay runs of 4 replicates of each specimen were conducted on 4 days. Repeatability was determined with 8 replicates of each specimen.

Kit	S #	Mean (EU/ml)	Total Imprecision		Between days		Within run (Repeatability)	
			SD		SD		SD	
			(EU/ml)	CV%	(EU/ml)	CV%	(EU/ml)	CV%
RNP Ab Assay	1	8.4	0.419	5.0%	0.353	4.2%	0.553	6.5%
	2	17.7	1.721	9.7%	1.670	9.8%	0.385	2.0%
	3	23.7	2.289	9.7%	2.193	9.5%	2.184	8.9%
	4	61.6	5.103	8.3%	5.431	9.0%	2.581	4.0%
	5	95.6	4.630	4.8%	4.995	5.2%	3.591	3.8%
	6	198.1	8.584	4.3%	7.720	3.8%	7.781	4.0%
Sm Ab Assay	1	2.3	0.235	10.2%	0.273	11.9%	0.150	6.5%
	2	14.1	0.624	4.4%	0.464	3.2%	0.361	2.7%
	3	19.7	1.136	5.8%	1.373	7.0%	0.444	2.2%
	4	30.0	1.644	5.5%	1.844	6.1%	1.043	3.5%
	5	66.7	3.791	5.7%	3.742	5.7%	2.432	3.5%
	6	120.8	6.133	5.1%	6.701	5.6%	4.687	3.8%

## Reproducibility

8 replicates of samples in the negative range, ~20% below cutoff, ~20% above cutoff and in the moderate positive range of the assay were performed to determine intra-assay qualitative reproducibility. 4 replicates of the same samples were tested in 4 runs to determine inter-assay reproducibility. Assay results for these specimens produced 100% qualitative agreement.

## Limit of Detection

The limit of detection (LoD) was determined based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples. LoD for RNP Ab was 3.0 EU/ml. LoD for Sm Ab was 4.2 EU/ml.

## Linearity and Recovery

Studies were performed using equidistant dilution series of positive samples with values throughout the calibrator range to determine linear range of the assay. The linear range of the assays were determined be 3.0 (LoD) – 400 EU/ml for RNP Ab and 4.2 (LoD) – 255 EU/ml for Sm Ab. Results are summarized below.

Test Range (EU/ml)	Slope (95% CI)	Y-intercept (95% CI)	R <sup>2</sup>	% recovery (obtained/expected)
RNP Ab				
4.4 to 136	0.98 (0.95 to 1.00)	2.1 (1.8 to 3.0)	0.9979	88.9 to 103.3
263 to 611.1	1.27 (1.10 to 1.51)	-15.6 (-16.1 to -14.2)	0.9817	96.1 to 139.8
7.2 to 70.3	0.97 (0.88 to 1.14)	0.5 (0.3 to 0.7)	0.9958	97.3 to 106.9
Sm				
5.2 to 75.8	1.01 (0.93 to 1.51)	-0.5 (-0.8 to 0.1)	0.9987	97.3 to 105.7

**EN**

4.8 to 93.9	1.00 (0.90 to 1.51)	1.1 (0.3 to 1.8)	0.9963	90.9 to 102
5.3 to 317.9	0.72 (0.61 to 1.51)	-2.2 (-4.0 to -0.2)	0.9753	101.9 to 159.6

**Interference**

Interference was studied by mixing sera with known RNP and Sm antibody levels with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342 µmol/L), and Rheumatoid Factor (100 EU/ml).

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## REFERENCES

1. Tan EM. Autoantibodies to nuclear antigens (ANA): their immunobiology and medicine. *Adv Immunol* 33: 167-240, 1982.
2. Tan EM. Special antibodies for the study of systemic lupus erythematosus. *Arth Rheum* 25: 753-756, 1982.
3. Tan EM, Chan E, Sullivan KF and Rubin RL. Antinuclear antibodies (ANAs): diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity. *Clin Immunol Immunopathol* 47: 121-141, 1988.
4. Reichlin M. Current perspectives on serological reactions in SLE patients. *Clin Exp Immunol* 44: 1-10, 1981.
5. Reichlin M and Harley JB. Antibodies to extractable nuclear antigens: clinical significance. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 555-563, 1987.
6. Kumar V, Beutner EH and Chorzelski TP. Autoimmunity and the skin. In "Concepts in Immunopathology", Vol 1, Cruse JM and Lewis RE, Eds, Karger, Basel, 318-353, 1985.
7. McCarty GA. Autoantibodies and their relation to rheumatic diseases. *Medical Clinics of North America* 70: 237-261, 1986.
8. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. *Arth Rheum* 29: 457-460, 1986.
9. Williamson GG and Boyle JA. Antigenic relatedness of small ribonucleoprotein particles. *Biochem Biophys Acta* 798:149-155, 1984.
10. Maddison PJ, Skinner RP, Vlachoyiannopoulos P, Brennan DM and Hough D. Antibodies to nRNP, Ro (SS-A) and La (SS-B) detected by ELISA: their specificity and inter-relations in connective tissue disease sera. *Clin Exp Immunol* 62: 337- 245, 1985.
11. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, (HHS Pub No [CDC] 88-8395), 1988.

*For technical assistance please contact:*

**■ IMMCO Diagnostics, Inc.**  
60 Pineview Drive  
Buffalo, NY 14228-2120 USA  
Telephone: (716) 691-0091  
Fax: (716) 691-0466  
Toll Free USA/Canada: 1-800-537-TEST  
E-Mail: [info@immco.com](mailto:info@immco.com)

*or your local product distributor*

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands