SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Immulisa™

ENA Antibody ELISA

SS-A (Ro) Antibody ELISA

SS-B (La) Antibody ELISA

[IVD] For *in vitro* diagnostic prescription use only

CLIA Complexity: High

PRODUCT INSERT

[REF] 5128 SS-A (Ro) Antibody ELISA 96 Determinations

[REF] 5129 SS-B (La) Antibody ELISA 96 Determinations

INTENDED USE

Immuli

sa Enhanced™ SS-A (Ro) Antibody ELISA: Enzyme linked immunoassay (ELISA) for the qualitative and semi-quantitative detection of SS-A (Ro) (52 kD and 60 kD) IgG antibodies in human serum as an aid in diagnosis of Systemic Lupus Erythematosus (SLE) and Sjögren's Syndrome in conjunction with clinical findings and other laboratory tests.

Immuli

sa Enhanced™ SS-B (La) Antibody ELISA: Enzyme linked immunoassay (ELISA) the qualitative and semi-quantitative detection of SS-B (La) IgG antibodies in human serum as an aid in diagnosis of Systemic Lupus Erythematosus (SLE) and Sjögren's Syndrome in conjunction with clinical findings and other laboratory tests.

SUMMARY AND EXPLANATION

ENA are soluble ribonucleoprotein (snurps) complexes. Autoantibodies directed against various ENA have proven to be of value in the diagnosis and monitoring of various systemic connective tissue diseases.¹⁻¹⁰ Antibodies to SS-A (Ro) and SS-B (La) occur in approximately 30-50% and 10-15% of SLE patients and 60-95% and 40-90% of patients with Sjögren's syndrome respectively. Antibodies to SS-B (La) occur frequently in association with SS-A (Ro) antibodies.¹⁸⁻²¹

These antibodies can be detected by various methods. ELISA methodology has many advantages over older technologies such as immunodiffusion: assay performance times are reduced, individual subjectivity in reading results is eliminated, quantitation is achieved without serum titration, there is potential for automation and ELISAs offer greater sensitivity.²²

PRINCIPLES OF PROCEDURE

The test is performed as a solid phase immunoassay. Microwells are coated with recombinant SS-A (Ro) 52 kD and 60 kD antigen or SS-B (La) antigen followed by a blocking step to reduce non-specific protein binding during the assay run. Controls, calibrators and patient sera are incubated in the antigen coated wells to allow specific antibodies present in the serum to bind to the

EN

antibodies and other serum proteins are removed by washing the microwells. Bound antibodies are detected by adding an enzyme labeled anti-human IgG conjugate to the microwells. Unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in ELISA units per milliliter (EU/ml) and reported as positive or negative.

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.²³

Stop Solution is a dilute sulfuric acid solution. Sulfuric acid (H₂SO₄) is poisonous and corrosive. Do not ingest and avoid contact with skin and eyes. Avoid exposure to bases, metals or other compounds that may react with acids.

TMB Enzyme Substrate contains an irritant that may be harmful if inhaled, ingested or absorbed through the skin. Do not ingest and avoid contact with skin and eyes.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

Caution: United States federal law restricts this device to sale by or on the order of a licensed practitioner.

Materials provided

Kits contains sufficient reagents to perform 96 determinations.

EN

ImmuLisa™ SS-A (Ro) (52 kD and 60 kD) Antibody ELISA REF 5128

12 x 8	MICROPLATE SSA	Microplate with individual breakaway microwells. Coated with recombinant Ro 52 kD and 60 kD antigen. Ready for use.
1 x 1.75 ml	CONTROL + SSA	Ready to use Positive Control (red cap). Contains human serum positive for SS-A (Ro) (52 kD and 60 kD) antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL -	Ready to use Negative Control (white cap). Contains human serum.
5 x 1.75 ml	CALIBRATOR A SSA CALIBRATOR B SSA CALIBRATOR C SSA CALIBRATOR D SSA CALIBRATOR E SSA	Ready to use set of 5 Calibrators . Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing SS-A (Ro) 52kD and 60kD antibodies. Concentrations in EU/ml are printed on the labels.
1 x 15 ml	IgG-CONJ HRP	HRP goat anti-human IgG Conjugate . Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. Protect from light .
1 x 15 ml	STOP H2SO4	Stop Solution*. Ready for use.
2 x	BUF WASH	Powder Wash Buffer. Reconstitute to one liter each .
1 x		Protocol Sheets

ImmuLisa™ SS-B (La) Antibody ELISA REF 5129

12 x 8	MICROPLATE SSB	Microplate with individual breakaway microwells. Coated with recombinant La antigen. Ready for use.
1 x 1.75 ml	CONTROL + SSB	Ready to use Positive Control (red cap). Contains human serum positive for SS-B (La) antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL -	Ready to use Negative Control (white cap). Contains human serum.
5 x 1.75 ml	CALIBRATOR A SSB CALIBRATOR B SSB CALIBRATOR C SSB CALIBRATOR D SSB CALIBRATOR E SSB	Ready to use set of 5 Calibrators . Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing SS-B (La) antibodies. Concentrations in EU/ml are printed on the labels.
1 x 15 ml	IgG-CONJ HRP	HRP goat anti-human IgG Conjugate . Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. Protect from light .
1 x 15 ml	STOP H2SO4	Stop Solution*. Ready for use.
2 x	BUF WASH	Powder Wash Buffer. Reconstitute to one liter each .

Optional Components

1 x 60ml

BUF	WASH
-----	------

Liquid concentrated Wash Buffer. **Reconstitute to one liter.**

Symbols used on labels

Lot number

Date of Manufacture

Catalog number

*Danger; Causes severe skin burns and eye damage. Causes severe eye damage.

In vitro diagnostic use

Wear protective gloves, protective clothing, eye protection, face protection.

Use by

Wash exposed skin thoroughly after handling. If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Storage temperature

Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

Read instructions before use

Number of tests

Manufacturer

Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year.

EN PROCEDURE

Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

Test Method

Step 1 Let all reagents and specimens equilibrate to room temperature.

Step 2 Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.

Step 3 For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.

or

For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

Qualitative			
A	Blank	S5	Etc.
B	-Control	S6	
C	+ Control	S7	
D	Cal D	S8	
E	S1	S9	
F	S2	S10	
G	S3	S11	
H	S4	S12	
	1	2	3

Semi-Quantitative			
A	Blank	S1	Etc.
B	-Control	S2	
C	+ Control	S3	
D	Cal A	S4	
E	Cal B	S5	
F	Cal C	S6	
G	Cal D	S7	
H	Cal E	S8	
	1	2	3

EN

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500ul** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.
- Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3 . The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <10 EU/ml. The optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. The positive control must give values in the range stated on the vial.

EN RESULTS

Calculations

The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

It is recommended that qualitative results be reported as “positive” or “negative.” Sample results greater than or equal to Calibrator D are considered positive.

2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in EU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as “positive,” “negative,” or “indeterminate” with EU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods for detection of antibodies to ENA.

Interpretation

Interpretation values were determined by testing normal blood donors and disease control sera. The cutoff was established using mean of the normal subjects plus 4.5 SD for the SS-A (Ro) (52 kD and 60 kD) Ab ELISA and 3 SD for the SS-B (La) Ab ELISA. The cutoff was assigned an arbitrary value of 20 EU/ml. Negative, Indeterminate and Positive assay ranges are provided below.

SS-A (Ro) (52 kD

and 60 kD) Ab value

SS-B (La) Ab value

Interpretation

<20 EU/ml

<20 EU/ml

Negative

20-25 EU/ml

20-25 EU/ml

Indeterminate (Borderline)

>25 EU/ml

>25 EU/ml

Positive

Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. The measuring range of these assays are from the EU/ml defined in the Limit of Detection (see below) to Calibrator A (160 EU/ml) for semi-quantitative assays. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve

EN

when retested. For determining EU/ml values, multiply the units obtained by the dilution factor.

LIMITATIONS OF PROCEDURE

The assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only. The results obtained serve only as an aid in the diagnosis. Taken alone, these results should not be interpreted as diagnostic. This assay has not been validated in a pediatric population.

EXPECTED VALUES

Test results in a normal population are usually expected to be negative. The prevalence of SS-A (Ro) and SS-B (La) antibodies in various systemic connective tissue diseases is summarized in the following table:

Diagnostic Significance of Antibodies to SS-A (Ro) and SS-B (La) Antigens

Antibody	Disease Association (% Incidence)
SS-A (Ro)	SLE (30-50%) Sjögren's syndrome (60-95%)
SS-B (La)	SLE (10-15%) Sjögren's syndrome (15-60%)

Note: The frequency of each antibody specificity in a disease represents compilation from the literature.^{3,10-17} The prevalence varies depending on the patient population.

Sets of clinical samples were tested on the Immulisa™ SS-A (Ro) (52 kD and 60 kD) and SS-B (La) Antibody ELISAs. Results demonstrating prevalence in the populations for this study are provided below.

Condition	SS-A (Ro) Ab			SS-B (La) Ab		
	n	n pos	% pos	n	n pos	% pos
SLE	160	69	41.3%	137	19	13.9%
SLE with APS	15	3	20.0%	14	3	21.4%
Sjögren's	104	81	77.9%	68	41	60.3%
Myositis	91	31	34.1%	48	0	0.0%
Systemic sclerosis	154	30	19.5%	70	3	4.3%
Healthy subjects	142	4	2.8%	133	3	2.3%

SLE=systemic lupus erythematosus, APS=antiphospholipid syndrome

EN

Clinical sensitivity and specificity for the populations tested are summarized below:

Analyte	Disease	Sensitivity	Specificity
SS-A (Ro) Ab	SLE	39.4%	91.5%
SS-A (Ro) Ab	Sjögren's	77.9%	91.5%
SS-B (La) Ab	SLE	14.6%	97.7%
SS-B (La) Ab	Sjögren's	60.3%	97.7%

PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ SS-A (Ro) (52 kD and 60 kD) and SS-B (La) Antibody ELISAs were evaluated by testing specimens from subjects with SLE, Sjögren's, systemic sclerosis, myositis and a variety of autoimmune and infectious disease control populations. These specimens were also tested on commercially available test kits. Only specimens in the linear range of the assay were included in the method comparison. These results are summarized below.

- A. Method Comparison: Immulisa™ SS-A (Ro) (52 kD and 60 kD) Antibody ELISA vs. other SS-A (Ro) (52 kD and 60 kD) Antibody kit:

IMMCO SS-A (Ro) Ab ELISA		Other SS-A (Ro) Antibody ELISA		
		Positive	Negative	Total
	Positive	138	24	162
	Negative	9	452	461
	Total	147	476	623

Positive Percent Agreement: 93.9% (95% CI 88.4% - 97.0%)

Negative Percent Agreement: 95.0% (95% CI 92.5% - 96.7%)

Total Agreement: 94.7% (95% CI 92.6% - 96.3%)

- B. Method Comparison: Immulisa™ SS-B (La) Antibody ELISA vs. other SS-B (La) Antibody kit:

IMMCO SS-B (La) Ab ELISA		Other SS-B (La) Antibody ELISA		
		Positive	Negative	Total
	Positive	62	8	70
	Negative	5	498	503
	Total	67	506	573

Positive Percent Agreement: 92.5% (95% CI 82.7% - 97.2%)

Negative Percent Agreement: 98.4% (95% CI 96.8% - 99.3%)

Total Agreement: 97.7% (95% CI 96.1% - 98.7%)

EN

C. Cross Reactivity: Numerous potentially cross-reactive autoimmune and infectious disease sera were tested for SS-A (Ro) (52 kD and 60 kD) and SS-B (La) antibodies using the Immulisa™ assays. Studies demonstrated minimal cross-reactivity for these specimens as indicated below.

Condition	SS-A (Ro) Ab		SS-B (La) Ab	
	n	n pos (%)	n	n pos (%)
IBD ¹	40	0	40	0
Myositis ²			48	0
MCTD	17	0	10	0
Rheumatoid arthritis	33	1	33	1
Other AI controls ³	127	0	123	0
Infectious disease ⁴	110	3	110	3
Other disease controls ⁵	20	1	20	1
Total	347	5 (1.4%)	384	5 (1.3%)

1. Ulcerative colitis (20), Crohn's disease (20)
2. Myositis: polymyositis (34), dermatomyositis (14)
3. Other autoimmune disease controls: celiac disease (26 SS-A, 25 SS-B), dermatitis herpetiformis (10), Wegener's (8), Churg-Strauss (10), Hashimoto's (8), Graves (10), antiphospholipid syndrome (15 SS-A, 12 SS-B), Crohn's disease (20), ulcerative colitis (20)
4. Infectious disease controls: syphilis (30), hepatitis (30), toxoplasmosis (10), cytomegalovirus (10), herpes (20), rubella (10)
5. Other disease controls: osteoarthritis (10), thrombocytopenia (10)

Precision

Precision was tested with positive specimens selected throughout the range of the assay. Assay runs of 6 replicates of each specimen were conducted over multiple days. Repeatability was determined with 12 replicates of each specimen.

Kit	S #	Mean (EU/ml)	Total Imprecision		Between days		Within run (Repeatability)	
			SD (EU/ml)	CV%	SD (EU/ml)	CV%	SD (EU/ml)	CV%
SS-A (Ro) Assay	1	9.6	1.2	12.4%	1.2	12.3%	1.2	12.3%
	2	17.2	1.5	8.9%	1.6	9.1%	1.2	7.0%
	3	20.3	1.2	6.1%	1.2	6.0%	1.4	6.6%
	4	23.2	1.8	7.7%	1.8	7.9%	1.6	6.8%
	5	85.2	6.0	7.0%	6.0	7.1%	5.5	6.3%
	6	116.4	7.1	6.1%	7.2	6.2%	5.4	4.5%
	7	137.2	13.2	9.6%	13.4	9.9%	8.2	5.7%

SS-B (La) Assay	1	11.9	1.0	8.1%	0.9	7.9%	0.8	7.1%
	2	15.2	1.2	8.0%	1.2	8.0%	1.1	7.7%
	3	19.8	1.7	8.4%	1.7	8.6%	1.1	5.7%
	4	25.2	1.9	7.6%	1.9	7.4%	1.7	6.9%
	5	51.4	4.7	9.1%	4.8	9.4%	3.6	6.9%
	6	80.4	7.7	9.5%	7.9	9.9%	5.8	7.0%
	7	152.1	10.1	6.7%	10.6	7.0%	6.7	4.4%

Reproducibility

Intra-assay qualitative reproducibility was tested with 90 replicates of samples in the negative range, ~20% below cutoff, at ~cutoff, ~20% above cutoff and in the moderate positive range of the assays. For SS-A/Ro (52 kD and 60 kD) Ab the cutoff +/- 20% specimens produced 97.8% qualitative agreement, the cutoff specimen produced 51.1% negative and 48.9% positive results and other specimens produced 100% qualitative agreement. For SS-B/La Ab the cutoff specimen produced 53.3% negative and 46.7% positive results and other specimens produced 100% qualitative agreement.

Inter-lot reproducibility was tested with positive specimens selected throughout the range of the assay. Each specimen was tested in 3 runs of 3 replicates over 3 days on each of 3 lots. Assay results for these specimens produced 100% qualitative agreement.

Limit of Detection

The limit of detection (LoD) was determined based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples. LoD for SS-A (Ro) (52 kD and 60 kD) Ab was 2.0 EU/ml. LoD for SS-B (La) was 1.6 EU/ml.

Linearity and Recovery

Linearity and recovery were tested by diluting positive specimens through the assay range in equidistant dilutions and comparing actual vs. expected results. The linear ranges of the assays were determined be 2.0 (LoD) – 160 EU/ml for SS-A (Ro) (52 kD and 60 kD) Ab and 1.6 (LoD) – 160 EU/ml for SS-B (La) Ab. Results are summarized below:

Test Range (EU/ml)	Slope (95% CI)	Y-intercept (95% CI)	R ²	% recovery (obtained/expected)
SS-A (Ro) Ab				
2.6 to 35.8	0.96 (0.89 to 1.03)	0.76 (-0.7 to 2.2)	0.9953	88% to 104%
3.0 to 175.1	0.98 (0.93 to 1.03)	2.91 (2.76 to 3.06)	0.9919	86% to 108%
18.2 to 156.0	0.99 (0.92 to 1.07)	-1.32 (-8.5 to 5.8)	0.9948	94% to 110%
30.3 to 248.5	0.98 (0.81 to 1.14)	-7.7 (-35.1 to -19.7)	0.9724	100% to 121%
1.7 to 21.0	1.02 (0.96 to 1.09)	-0.59 (-1.43 to 0.25)	0.9961	98% to 115%

SS-B (La) Ab				
2.8 to 22.2	0.98 (0.89 to 1.07)	0.26 (-1.0 to 1.6)	0.9914	92% to 106%
8.7 to 81.8	0.88 (0.76 to 1.01)	2.61 (-4.2 to 9.4)	0.981	84% to 116%
8.7 to 159.3	0.90 (0.79 to 1.01)	3.92 (-6.6 to 14.5)	0.9894	94% to 117%
6.4 to 197.4	1.00 (0.90 to 1.11)	-5.80 (-18.64 to 7.05)	0.9894	100% to 119%
1.1 to 13.7	1.07 (0.93 to 1.22)	-0.66 (-1.9 to 0.58)	0.9815	92% to 122%

EN

Interference

Interference was studied by mixing sera with known SS-A (Ro) (52 kD and 60 kD) and SS-B (La) antibody levels with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342 $\mu\text{mol/L}$), and Rheumatoid Factor (100 EU/ml).

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REFERENCES

1. Tan EM. Autoantibodies to nuclear antigens (ANA): their immunobiology and medicine. *Adv Immunol* 33: 167-240, 1982.
2. Tan EM. Special antibodies for the study of systemic lupus erythematosus. *Arth Rheum* 25: 753-756, 1982.
3. Tan EM, Chan E, Sullivan KF and Rubin RL. Antinuclear antibodies (ANAs): diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity. *Clin Immunol Immunopathol* 47: 121-141, 1988.
4. Reichlin M. Current perspectives on serological reactions in SLE patients. *Clin Exp Immunol* 44: 1-10, 1981.
5. Reichlin M and Harley JB. Antibodies to extractable nuclear antigens: clinical significance. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 555-563, 1987.
6. Kumar V, Beutner EH and Chorzelski TP. Autoimmunity and the skin. In "Concepts in Immunopathology", Vol 1, Cruse JM and Lewis RE, Eds, Karger, Basel, 318-353, 1985.
7. McCarty GA. Autoantibodies and their relation to rheumatic diseases. *Medical Clinics of North America* 70: 237-261, 1986.
8. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. *Arth Rheum* 29: 457-460, 1986.
9. Williamson GG and Boyle JA. Antigenic relatedness of small ribonucleoprotein particles. *Biochem Biophys Acta* 798:149-155, 1984.
10. Tan EM. Antinuclear antibodies: diagnostic markers for autoimmune diseases and probes for cell biology. *Adv Immunol.* 1989;44:93-151.
11. von Mühlen CA, Tan EM. Autoantibodies in the diagnosis of systemic rheumatic diseases. *Semin Arthritis Rheum.* 1995;24:323-58.
12. Meilof JF, Smeenk RJ. Autoantibodies and their target antigens in Sjögren's syndrome. *Neth J Med.* 1992;40:140-7.
13. Scofield RH, Farris AD, Horsfall AC, Harley JB. Fine specificity of the autoimmune response to the Ro/SSA and La/SSB ribonucleoproteins. *Arthritis Rheum.* 1999;42:199-209.
14. Meyer O. [Anti-SSA/Ro and anti-SSB/La antibodies. What's new?]. *Ann Med Interne (Paris).* 2002;153:520-9.
15. Parker JC, Burlingame RW, Bunn CC. Prevalence of antibodies to Ro-52 in a serologically defined population of patients with systemic sclerosis. *Journal of Autoimmune Diseases* 2009, Mar 6; 6:2
16. Schulte-Pelkum J, Fritzier M, Mahler M. Latest update on the Ro/SS-A autoantibody system. *Autoimmun Rev.* 2009;8:632-7. Epub 2009 Feb 12.
17. Routsias JG, Tzioufas AG. Sjögren's syndrome--study of autoantigens and autoantibodies. *Clin Rev Allergy Immunol.* 2007;32:238-51.
18. Clark G, Reichlin M and Tomasi TB. Characterization of a soluble cytoplasmic antigen reactive with sera from patients with systemic lupus erythematosus. *J Immunol*1102:117-124, 1969.
19. Alspaugh MA and Maddison P. Resolution of the identity of certain antigen antibody systems in systemic lupus erythematosus and Sjogren's syndrome: an interlaboratory collaboration. *Arth Rheum* 22: 796-798, 1979.
20. Sontheimer RD, Thomas JR and Gilliam JN. Subacute cutaneous lupus erythematosus subset. *Arch Dermatol* 115: 1409-1415, 1979.
21. Sontheimer RD, Maddison PJ and Reichlin M. Serologic and HLA associations in subacute cutaneous lupus erythematosus, a clinical subset of lupus erythematosus. *Ann Intern Med* 97:664-671, 1982.
22. Maddison PJ, Skinner RP, Vlachoyiannopoulos P, Brennard DM and Hough D. Antibodies to nRNP, Ro (SS-A) and La (SS-B) detected by ELISA: their specificity and inter-relations in connective tissue disease sera. *Clin Exp Immunol* 62: 337- 245, 1985.
23. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health. 2007; (HHS Pub No [CDC] 93-8395).

For technical assistance please contact:

IMMCO Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120 USA

Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immco.com

or your local product distributor

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands