# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation

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[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# ImmLISA™ Enhanced RF ELISAs

## Detection of IgA, IgG or IgM Rheumatoid Factor (RF)

**IVD** For *in vitro* diagnostic use

### PRODUCT INSERT

REF	5138A	RF IgA ELISA	96 Determinations
REF	5138G	RF IgG ELISA	96 Determinations
REF	5138M	RF IgM ELISA	96 Determinations

### INTENDED USE

ImmLISA Enhanced™ RF IgA Antibody ELISA: Enzyme linked immunoassay (ELISA) for the qualitative or semi-quantitative detection of Rheumatoid Factor IgA antibodies in human serum to aid in the diagnosis of rheumatoid arthritis (RA) in conjunction with other laboratory tests and clinical findings.

ImmLISA Enhanced™ RF IgG Antibody ELISA: Enzyme linked immunoassay (ELISA) for the qualitative or semi-quantitative detection of Rheumatoid Factor IgG antibodies in human serum to aid in the diagnosis of rheumatoid arthritis (RA) in conjunction with other laboratory tests and clinical findings.

ImmLISA Enhanced™ RF IgM Antibody ELISA: Enzyme linked immunoassay (ELISA) for the qualitative or semi-quantitative detection of Rheumatoid Factor IgM antibodies in human serum to aid in the diagnosis of rheumatoid arthritis (RA) in conjunction with other laboratory tests and clinical findings.

### SUMMARY AND EXPLANATION

The measurement of RF is important in the diagnosis and prognosis of rheumatoid arthritis, as high titers of RF occur in sera of patients who tend to develop extra-articular complications.<sup>1,2</sup> The majority of routine laboratory tests measure IgM RF by its ability to agglutinate sheep red blood cells, latex or similar particles coated with IgG.<sup>1-4</sup> However, studies suggest that RF of other immunoglobulin isotypes are also present in rheumatoid arthritis.<sup>5-9</sup>

RF is present in 70-90% of patients with rheumatoid arthritis and is included in the classification criteria.<sup>10</sup> According to the revised ARA criteria, a low positive RF result contributes 2 out of 6 points and a high positive RF contributes 3 out of 6 points for classification of definite RA.

RF as detected by agglutination is of the IgM isotype. Other methods such as ELISA have improved specificity, sensitivity and reliability over existing and routinely used agglutination methods.<sup>3-5</sup> ELISA methods can detect RF of various immunoglobulin isotypes. Such a distinction is not possible with traditional agglutination tests. Significant levels of an isotype of RF commonly found in rheumatic diseases and elevations of various isotypes together is almost diagnostic for rheumatoid arthritis.<sup>8</sup> A combination of IgM, IgG and IgA RF results increases sensitivity for RA as compared to IgM RF results alone.<sup>9</sup>

## PRINCIPLES OF PROCEDURE

The test is performed as a solid phase immunoassay. Microwells are coated with purified rabbit IgG followed by a blocking step to reduce non-specific binding during the assay run. Controls, calibrators and patient sera are incubated in the antigen coated wells to allow specific antibodies present in the serum to bind to the antigen. Unbound antibodies and other serum proteins are removed by washing the microwells. Bound antibodies are detected by adding an enzyme labeled anti-human IgA, IgG or IgM conjugate to the microwells. Unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in ELISA Units per milliliter (EU/ml) for RF IgA and RF IgG and International Units per milliliter (IU/ml) for RF IgM and reported as positive or negative.

## REAGENTS

### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

Open kit studies have been performed with materials/reagents and demonstrate stability through a minimum of 30 days.

### Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.<sup>11</sup>

Stop Solution is a dilute sulfuric acid solution. Sulfuric acid (H<sub>2</sub>SO<sub>4</sub>) is poisonous and corrosive. Do not ingest and avoid contact with skin and eyes. Avoid exposure to bases, metals or other compounds that may react with acids.

TMB Enzyme Substrate contains an irritant that may be harmful if inhaled, ingested or absorbed through the skin. Do not ingest and avoid contact with skin and eyes.

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

Caution: United States federal law restricts this device to sale by or on the order of a licensed practitioner.

**Materials provided**

ImmLiSa™ RF IgA ELISA  
 ImmLiSa™ RF IgG ELISA  
 ImmLiSa™ RF IgM ELISA


REF	5138A
REF	5138G
REF	5138M

Kits contains sufficient reagents to perform 96 determinations.

12 x 8	MICROPLATE RF	Microplate with individual breakaway microwells. Coated with purified rabbit IgG antigen. Ready for use.
1 x 1.75 ml	CONTROL + RFA	Ready to use <b>Positive Control</b> (red cap). Contains human serum positive for RF IgA. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL + RFG	Ready to use <b>Positive Control</b> (red cap). Contains human serum positive for RF IgG. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL + RFM	Ready to use <b>Positive Control</b> (red cap). Contains human serum positive for RF IgM. The expected concentration range in IU/ml is printed on the label.
1 x 1.75 ml	CONTROL -	Ready to use <b>Negative Control</b> (white cap). Contains human serum.
5 x 1.75 ml	CALIBRATOR A RFA	Ready to use set of 5 <b>Calibrators</b> . Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing RF IgA antibodies. Concentrations in EU/ml are printed on the labels.
	CALIBRATOR B RFA	
	CALIBRATOR C RFA	
	CALIBRATOR D RFA	
	CALIBRATOR E RFA	
5 x 1.75 ml	CALIBRATOR A RFG	Ready to use set of 5 <b>Calibrators</b> . Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing RF IgG antibodies. Concentrations in EU/ml are printed on the labels.
	CALIBRATOR B RFG	
	CALIBRATOR C RFG	
	CALIBRATOR D RFG	
	CALIBRATOR E RFG	
5 x 1.75 ml	CALIBRATOR A RFM	Ready to use set of 5 <b>Calibrators</b> . Calibrator A (green cap) 80 IU/ml, Calibrator B (violet cap) 40 IU/ml, Calibrator C (blue cap) 20 IU/ml, Calibrator D (yellow cap) 10 IU/ml, and Calibrator E (orange cap) 1 IU/ml. Derived from human serum containing RF IgM antibodies. Concentrations in IU/ml are printed on the labels.
	CALIBRATOR B RFM	
	CALIBRATOR C RFM	
	CALIBRATOR D RFM	
	CALIBRATOR E RFM	
1 x 15 ml	IgA-CONJ HRP	HRP goat anti-human <b>IgA Conjugate</b> . Ready for use. Color coded pink.
1 x 15 ml	IgG-CONJ HRP	HRP rabbit anti-human <b>IgG Conjugate</b> . Ready for use. Color coded pink.
1 x 15 ml	IgM-CONJ HRP	HRP goat anti-human <b>IgM Conjugate</b> . Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. <b>Protect from light</b> .

EN

1 x 15 ml  Stop Solution. Ready for use.

2 x vials  Powder Wash Buffer. **Reconstitute to one liter each.**

1 x Protocol Sheets

### Optional Components

1 x 60ml  Liquid concentrated Wash Buffer. **Reconstitute to one liter.**

### Symbols used on labels

-  Lot number
-  Catalog number
-  In vitro diagnostic use
-  Use by
-  Storage temperature
-  Read instructions before use
-  Number of tests
-  Manufacturer

### Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

### SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. It is recommended that frozen specimens be tested within one year. Avoid repeated freezing and thawing of samples.

### PROCEDURE

#### Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.

- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

**Test Method**

- Step 1** Let all reagents and specimens equilibrate to room temperature.
- Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.
- Step 3** For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.
- or**
- For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

	Qualitative				Semi-Quantitative		
A	Blank	S5	Etc.	A	Blank	S1	Etc.
B	-	S6		B	-	S2	
C	+	S7		C	+	S3	
D	Cal D	S8		D	Cal A	S4	
E	S1	S9		E	Cal B	S5	
F	S2	S10		F	Cal C	S6	
G	S3	S11		G	Cal D	S7	
H	S4	S12		H	Cal E	S8	
	1	2	3		1	2	3

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500ul** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.
- Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** (± 5 min) at room temperature.

- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

### Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <10 EU/ml for IgA and IgG and less than 5 IU/ml for IgM. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

## RESULTS

### Calculations

The concentrations of the patient samples can be determined by either of two methods:

#### 1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{Unit of Calibrator D} = \text{Unit Test Sample}^*$$

It is recommended that qualitative results be reported as "positive" or "negative." Sample results greater than or equal to Calibrator D are considered positive.

Note: Units for RF IgA and IgG are EU/ml. Units for IgM are IU/ml

#### 2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in U/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as “positive,” “negative,” or “indeterminate” with EU/ml or IU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods.

### Interpretation

Interpretation values were determined by testing normal blood donors and non-rheumatoid arthritis disease control specimens. The mean of the normal subjects plus 2SD was established as the assay cutoff for IgA and IgG and assigned an arbitrary value of 20 EU/ml. The mean of the normal subjects plus 2.5SD was established as the assay cutoff for IgM and assigned a value of 10 IU/ml. Negative, Indeterminate and Positive assay ranges are provided below.

RF IgA or IgG value	RF IgM value	Interpretation
<20 EU/ml	<10 IU/ml	Negative
20-25 EU/ml	10-12.5 IU/ml	Indeterminate (Borderline)
>25 EU/ml	>12.5 IU/ml	Positive

### Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A.

### LIMITATIONS OF PROCEDURE

The assay should not be performed on grossly hemolyzed, microbially contaminated, icteric or lipemic samples. This method should be used for testing human serum samples only. The results obtained serve only as an aid in the diagnosis. Taken alone, these results should not be interpreted as diagnostic.

### EXPECTED VALUES

ELISA testing for RF has been found to be considerably more sensitive than conventional agglutination tests. Approximately 90% of patients with rheumatoid arthritis have elevated RF by the ELISA system compared to approximately 70% by agglutination.<sup>7</sup> Most patients with rheumatoid arthritis have elevation of 2 or more isotypes of RF as compared to other rheumatic diseases such as systemic lupus erythematosus. In addition, rheumatoid patients have higher levels of RF as compared to other rheumatic conditions.<sup>9,12</sup> The following tables summarize prevalence of RF reported in literature.

#### Prevalence of RF Isotypes in Rheumatoid Arthritis<sup>8,9</sup>

RF Isotype	~Prevalence (%)
IgM	66
IgG	44
IgA	51
IgM, IgG and/or IgA	70-90

#### Prevalence of RF in Patients with Other Rheumatic Diseases<sup>6,10,12-22</sup>

Condition	~ Prevalence (%)
Sjögren's syndrome	60%
Systemic lupus erythematosus	35%
Systemic sclerosis	25%
Mixed connective tissue disease	15%
Juvenile rheumatoid arthritis	15%
Polymyositis/dermatomyositis	10%

Sets of clinical samples were tested on the Immulisa™ RF ELISAs. Results demonstrating incidence in the populations for this study are provided below.

Clinical Populations	RF IgA ELISA			RF IgG ELISA			RF IgM ELISA		
	#	#	% Pos	#	#	% Pos	#	#	% Pos
	Tested	Pos	% Pos	Tested	Pos	% Pos	Tested	Pos	% Pos
Rheumatoid Arth.	249	131	52.6%	249	126	50.6%	249	176	70.7%
Juvenile Arth	10	0	0.0%	10	1	10.0%	10	1	10.0%
Osteoarthritis	30	0	0.0%	30	0	0.0%	30	1	3.3%
Psoriatic Arth	33	2	6.1%	33	3	9.1%	33	2	6.1%
Spondyloarthritis	33	1	3.0%	33	0	0.0%	33	1	0.0%
Churg-Strauss	10	1	10.0%	10	1	10.0%	10	1	10.0%
Mixed CT Disease	10	1	10.0%	10	2	20.0%	10	4	40.0%
Sjögren's	20	9	45.0%	20	2	10.0%	20	11	55.0%
Systemic Lupus Eryth	30	11	36.7%	30	6	20.0%	30	5	16.7%
Systemic Sclerosis	20	3	15.0%	20	4	20.0%	20	5	25.0%
Wegener's	8	1	12.5%	8	1	12.5%	8	1	12.5%
Misc AI*	36	0	0.0%	36	2	5.6%	36	1	2.8%
Misc ID**	70	4	5.7%	70	1	1.4%	70	6	8.6%
Healthy normals	122	5	4.1%	120	7	5.8%	144	5	3.5%

\* includes celiac (8), Graves (10), Hashimoto (8), ulcerative colitis (10)

\*\* includes syphilis (10), rubella (10), mononucleosis (5), Lyme (5), HCV (10), CMV (10), HSV (10), toxoplasmosis (10)

The clinical sensitivity and specificity for RA using these assays is summarized below\*:

RF IgA		RF IgG		RF IgM	
Sensitivity	Specificity	Sensitivity	Specificity	Sensitivity	Specificity
52.6%	89.4%	50.6%	92.9%	70.7%	87.4%

\* Excludes NHS.

## PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ RF ELISAs were evaluated by testing well-characterized rheumatoid arthritis serum specimens alongside disease controls and "normal" human sera. These specimens were also tested on commercially available ELISA kits. These results are summarized below.

A. Method Comparison: Immulisa™ RF ELISAs vs. other RF antibody ELISAs.

		Other RF IgA ELISA		
		Pos	Neg	Total
IMMCO RF	Pos	113	9	122
IgA ELISA	Neg	5	250	255
	Total	118	259	377
	Pos % Agreement		95.8%	(95% CI 89.9 - 98.4)
	Neg % Agreement		96.5%	(95% CI 93.3 - 98.3)
	Overall % Agreement		96.3%	(95% CI 93.7 - 97.9)

Disease Associated (RA) Specimens: 230  
 Disease Control (non-RA) Specimens:

Other RF IgG ELISA

		Pos	Neg	Total
IMMCO RF IgG ELISA	Pos	88	24	112
	Neg	7	271	278
Total		95	295	390

Pos % Agreement 92.6% (95% CI 84.9 - 96.7)  
 Neg % Agreement 91.9% (95% CI 88.0 - 94.6)  
 Overall % Agreement 92.1% (95% CI 88.8 - 94.5)

Disease Associated (RA) Specimens: 154  
 Disease Control (non-RA) Specimens: 236

Other RF IgM ELISA

		Pos	Neg	Total
IMMCO RF IgM ELISA	Pos	124	16	140
	Neg	4	235	239
Total		128	251	379

Pos % Agreement 96.9% (95% CI 91.7 - 99.0)  
 Neg % Agreement 93.6% (95% CI 89.7 - 96.2)  
 Overall % Agreement 94.7% (95% CI 91.8 - 96.7)

Disease Associated (RA) Specimens: 139  
 Disease Control (non-RA) Specimens: 239

B. Cross Reactivity: A total of 106 potentially cross-reactive specimens from individuals with autoimmune disorders and infectious diseases not reported to demonstrate high incidence of autoantibodies for RF were tested using the ImmuLisa™ RF ELISAs.

Clinical Populations	RF IgA ELISA			RF IgG ELISA			RF IgM ELISA		
	# Tested	# Pos	% Pos	# Tested	# Pos	% Pos	# Tested	# Pos	% Pos
Celiac Disease	8	0	0.0%	8	0	0.0%	8	0	0.0%
Graves	10	0	0.0%	10	1	10.0%	10	1	10.0%
Hashimoto	8	0	0.0%	8	0	0.0%	8	0	0.0%
Ulcerative Colitis	10	0	0.0%	10	1	10.0%	10	0	0.0%
Syphilis	10	0	0.0%	10	0	0.0%	10	0	0.0%
Rubella	10	0	0.0%	10	0	0.0%	10	1	10.0%
Mononucleosis	5	1	20.0%	5	0	0.0%	5	0	0.0%
Lyme	5	0	0.0%	5	0	0.0%	5	0	0.0%
HCV	10	1	10.0%	10	1	10.0%	10	2	20.0%
CMV	10	0	0.0%	10	0	0.0%	10	0	0.0%
HSV	10	1	10.0%	10	0	0.0%	10	2	20.0%
Toxoplasmosis	10	1	10.0%	10	0	0.0%	10	1	10.0%

Total	106	4	3.8%	106	3	2.8%	106	7	6.6%
-------	-----	---	------	-----	---	------	-----	---	------

### Precision

Precision was tested with 7 positive specimens selected throughout the range of the assay for each isotype. Assay runs of 6 replicates of each specimen were conducted over approximately 21 days. Repeatability was determined with 12 replicates of each specimen.

Sample	Mean (units/ml)	Total Imprecision		Between days		Within run (Repeatability)	
		SD (units/ml)	CV%	SD (units/ml)	CV%	SD (units/ml)	CV%
IgA 1	9.2	0.54	5.9%	0.37	4.0%	0.39	4.2%
IgA 2	14.9	0.81	5.5%	0.57	3.8%	0.58	3.9%
IgA 3	21.2	1.20	5.7%	0.71	3.3%	0.97	4.6%
IgA 4	26.4	1.42	5.4%	1.22	4.6%	0.73	2.8%
IgA 5	74.6	5.00	6.7%	3.86	5.2%	3.18	4.3%
IgA 6	115.6	5.28	4.6%	4.09	3.5%	3.34	2.9%
IgA 7	158.1	6.36	4.0%	3.89	2.5%	5.03	3.2%
IgG 1	8.8	0.76	8.7%	0.43	4.9%	0.63	7.2%
IgG 2	16.4	1.04	6.3%	0.71	4.3%	0.76	4.6%
IgG 3	19.4	1.06	5.5%	0.6	3.1%	0.88	4.5%
IgG 4	22.7	1.08	4.8%	0.54	2.4%	0.94	4.1%
IgG 5	45.9	2.16	4.7%	1.11	2.4%	1.86	4.0%
IgG 6	79.7	4.92	6.2%	4.15	5.2%	2.65	3.3%
IgG 7	133.0	6.36	4.8%	4.14	3.1%	4.83	3.6%
IgM 1	4.8	0.17	3.6%	0.13	2.7%	0.11	2.4%
IgM 2	8.5	0.33	3.9%	0.11	1.3%	0.31	3.7%
IgM 3	10.3	0.79	7.7%	0.48	4.7%	0.63	6.1%
IgM 4	11.9	0.70	5.9%	0.54	4.5%	0.44	3.7%
IgM 5	25.6	1.08	4.2%	0.81	3.2%	0.71	2.8%
IgM 6	45.5	1.87	4.1%	1.46	3.2%	1.17	2.6%
IgM 7	80.2	3.64	4.5%	2.99	3.7%	2.08	2.6%

### Qualitative Reproducibility

Ninety replicates of samples in the negative range, ~20% below cutoff, ~20% above cutoff and in the moderate positive range of the assay were performed to determine qualitative reproducibility for each isotype. These samples were tested in multiple runs. Assay results for each specimen produced 100% qualitative agreement.

### Limit of Detection

The limits of detection (LoD) for RF antibody using these assays were determined to be 3.7 EU/ml for IgA, 2.2 EU/ml for IgG and 1.3 IU/ml for IgM based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples.

### Linearity

Studies were performed using equidistant dilution series of positive samples with values throughout the calibrator range to determine linear range of the assay. The linear range of RF IgA was determined to be 3.3 – 160 EU/ml. The linear range of RF IgG was determined to be 2.6 – 160 EU/ml. The linear range of RF IgM was determined to be 1.2 – 80 IU/ml. Results are summarized below.

Sample	Test Range (units/ml)	Slope (95% CI)	Y-Intercept (95% CI)	R <sup>2</sup>	% Recovery
IgA 1	4.8 - 13.7	0.96 (0.86 to 1.06)	0.83 (-0.12 to 1.77)	1.00	91.2 to 100.0
IgA 2	10.5 - 64.9	0.96 (0.84 to 1.09)	3.66 (-1.16 to 8.48)	0.98	83.4 to 101.3
IgA 3	17.3 - 155.0	1.00 (0.88 to 1.14)	2.55 (-9.4 to 14.5)	0.99	88.0 to 104.1
IgG 1	3.6 - 13.8	1.00 (0.92 to 1.08)	0.33 (-0.40 to 1.07)	1.00	92.2 - 100.9
IgG 2	7.0 - 68.7	1.01 (0.89 to 1.12)	1.81 (-2.85 to 6.47)	0.99	87.5 - 100.1
IgG 3	16.1 - 154.7	1.07 (0.95 to 1.20)	-1.59 (-13.3 to 10.2)	1.00	88.2 - 105.6
IgM 1	2.0 - 8.1	0.97 (0.82 to 1.26)	-0.50 (-2.60 to 1.25)	0.99	83.0 - 102.3
IgM 2	8.1 - 61.4	0.96 (0.91 to 1.01)	1.43 (-0.56 to 3.43)	1.00	87.6 - 103.1
IgM 3	13.0 - 73.3	1.00 (0.84 to 1.15)	2.24 (-4.88 to 9.36)	1.00	89.4 - 100.0

### Interference

Interference was studied by mixing sera with known RF isotype levels with potentially interfering serum samples and studying deviation from expected results. The following substances were tested at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342 µmol/L), Cholesterol (13 mmol/L) and Triglycerides (37 mmol/L). It was noted that bilirubin, cholesterol and triglyceride spiked specimens demonstrated up to +10% greater difference than expected coefficient of variance percent for these assays. Please note sample requirements listed in Limitations of Procedure.

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## REFERENCES

1. Jónsson T et al. Is measurement of rheumatoid factor isotypes clinically useful? *Ann Rheum Dis.* 1993; 52: 161-164.
2. van Zeben D et al. Clinical significance of rheumatoid factors in early rheumatoid arthritis: results of a follow up study. *Ann Rheum Dis.* 1992; 51:1029-1035.
3. Shmerling RH et al. The rheumatoid factor: an analysis of clinical utility. *Am J Med.* 1991; 91:528-534.
4. Wolfe F et al. The latex test revisited-rheumatoid factor testing in 8,287 rheumatic disease patients. *Arthr Rheum.* 1991; 34:951-960.
5. Bampton JLM et al. Measurement of rheumatoid factors by enzyme linked immunosorbent assay (ELISA) and comparison with other methods. *Ann Rheum Dis.* 1985; 44: 13-19.
6. Markusse HM et al. Rheumatoid factor isotypes in serum and salivary fluid of patients with primary Sjögren's Syndrome. *Clin Immun Immunopath.* 1993; 66:26-32.
7. Jónsson P. Studies on the clinical significance of rheumatoid factor isotypes. Thesis University of Iceland, Reykjavik Iceland, 1993.
8. Jónsson T et al. Combined elevation of IgM and IgA rheumatoid factor has high diagnostic specificity for rheumatoid arthritis. *Rheumatol Int.* 1998; 18: 119-22.
9. Vallbracht I et al. Diagnostic and clinical value of anti-cyclic Citrullinated peptide antibodies compared with rheumatoid factor isotypes in rheumatoid arthritis. *Ann Rheum Dis.* 2004; 63: 1079-1084.
10. Aletaha D et al. Rheumatoid arthritis classification criteria: An American College of Rheumatology/European League Against Rheumatism collaborative initiative. *Arthritis Rheum* 2010; 62: 2569-2581.
11. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health. 2007; (HHS Pub No [CDC] 93-8395).
12. Sheldon J. Laboratory testing in autoimmune rheumatic diseases. *Best Pract Res Clin Rheumatol.* 2004; 18:249-69.
13. Atkinson JC et al. Serum anti-SS-B/La and IgA rheumatoid factor are markers of salivary gland disease activity in primary Sjögren's syndrome. *Arth Rheum.* 1992; 35:1368-1372.
14. Anderson LG et al. The spectrum of benign to malignant lymphoproliferation in Sjögren's Syndrome. *Clin Exp Immunol.* 1972; 10:199-221.
15. Saulsbury FT. Heavy and light chain composition of serum IgA and IgA rheumatoid factor in Henoch-Schönlein purpura. *Arthr Rheum.* 1992; 35:1377-1380.
16. Saulsbury FT. IgA rheumatoid factor in Henoch Schönlein purpura. *J Pediatr.* 1986; 108:71-76.
17. Bas S et al. Anti-cyclic citrullinated peptide antibodies, IgM and IgA rheumatoid factors in the diagnosis and prognosis of rheumatoid arthritis *Rheumatology.* 2003; 42:677-680.
18. Westwood OMR et al. Rheumatoid factors: what's new? *Rheumatology.* 2006; 45:379-385.
19. Ingegnoli F et al. Rheumatoid factors: clinical applications. *Dis Markers.* 2013; 35:727-34.
20. Barcelos F et al. Anti-cyclic citrullinated peptide antibodies and rheumatoid factor in Sjögren's syndrome. *Acta Reumatol Port.* 2009; 34:608-12.
21. Huo AP et al. Predictive and prognostic value of antinuclear antibodies and rheumatoid factor in primary Sjogren's syndrome. *Int J Rheum Dis.* 2010; 13:39-47.
22. Mimura Y et al. Rheumatoid factor isotypes in mixed connective tissue disease. *Clin Rheumatol.* 2006; 25:572-4.

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**  
**60 Pineview Drive**  
**Buffalo, NY 14228-2120**  
**Telephone: (716) 691-0091**  
**Fax: (716) 691-0466**  
**Toll Free USA/Canada: 1-800-537-TEST**  
**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante Autorizzato/Representante  
Autorizado/Représentant Autorisé  
EMERGO Group, Inc.  
Molenstraat 15, 2513 BH, The Hague,  
The Netherlands  
Tel (+31) 345 8570, Fax (+31) 346 7299  
[www.emergogroup.com](http://www.emergogroup.com)

December 2015  
Document No. PI5138