# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation

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[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# ImmuLisa™

## PR3 Antibody ELISA

### Proteinase 3 (PR3) Antibody ELISA

**IVD** For *in vitro* diagnostic use

#### PRODUCT INSERT

**REF** 5162 PR3 Antibody ELISA96 Determinations

#### INTENDED USE

An enzyme linked immunosorbent assay (ELISA) for the detection and semi-quantitation of antibodies to proteinase 3 (PR3) in human serum to aid in the diagnosis of anti-neutrophil cytoplasmic antibody (ANCA) associated vasculitides, including Wegener's granulomatosis, in conjunction with other laboratory and clinical findings.

#### SUMMARY AND EXPLANATION

Anti-neutrophil cytoplasmic antibodies (ANCA) are a group of antibodies directed against proteins in the granules of neutrophils. The presence of ANCA in patients with vasculitis was first observed in 1982 by Davies.<sup>1</sup> These antibodies can be detected by indirect immunofluorescence on ethanol fixed neutrophils, producing a characteristic cytoplasmic staining pattern (cANCA) or by ELISA. cANCA reactions target Proteinase 3 (PR3) and to some extent minor antigens such as cathepsin G and elastase. PR3 is a neutral serine proteinase localized in the azurophilic granules of the neutrophils with a MW of 28kD.<sup>2</sup> Antibodies against the PR3 antigen serve as a marker for Wegener's Granulomatosis (WG), a systemic necrotizing vasculitis.<sup>3</sup> WG is systemic small vessel vasculitis characterized by necrotizing granulomatous inflammation of the upper and lower respiratory tract with pauci-immune necrotizing crescentic glomerulonephritis. Early diagnosis is of paramount importance, since Wegener's granulomatosis often has a progressive course rapidly increasing extent and severity of organ involvement. PR3 antibodies have clinical importance in the diagnosis and monitoring of disease activity, and predicting relapse. Several studies have established a direct correlation between PR3 antibody levels and the active phase of WG. The concentration of serum PR3 antibodies rise dramatically during disease exacerbations, and relapses are usually accompanied by significant increases in antibody levels.<sup>6,7</sup>

#### PRINCIPLES OF PROCEDURE

The test is performed as a solid phase immunoassay. Microwells are coated with purified PR3 antigen followed by a blocking step to reduce non-specific binding during the assay run. Controls, calibrators and patient sera are incubated in the antigen coated wells to allow specific antibodies present in the serum to bind to the PR3 antigen. Unbound antibodies and other serum proteins are removed by washing the microwells. Bound antibodies are detected by adding an enzyme labeled anti-human IgG conjugate to the microwells. Unbound conjugate is removed by washing. Specific enzyme substrate (TMB) is then added to the wells and the presence of antibodies is detected by a color change produced by the conversion of TMB substrate to a colored reaction product. The reaction is stopped and the intensity of the color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 450 nm. Results are expressed in International Units per milliliter (IU/ml) and reported as positive or negative.

#### REAGENTS

Store all reagents at 2-8°C. **Do not freeze.** Reagents are stable until the expiration date when stored and handled as directed.

Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use.

EN

Reconstitute the wash buffer to 1 liter with distilled or deionized water. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date.

Coated microwell strips are for one time use only. Unused microwell strips should be carefully resealed in the pouch containing desiccants to prevent condensation and stored at 2-8°C.

**Precautions**

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials.<sup>10</sup>

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

**Materials provided**

Immulin™ PR3 Antibody ELISA REF 5162

Kits contains sufficient reagents to perform 96 determinations.

12 x 8	MICROPLATE PR3	Microplate with individual breakaway microwells. Coated with purified PR3 antigen. Ready for use.
1 x 1.75 ml	CONTROL + PR3	Ready to use <b>Positive Control</b> (red cap). Contains human serum positive for PR3 antibodies. The expected concentration range in IU/ml is printed on the label.
1 x 1.75 ml	CONTROL -	Ready to use <b>Negative Control</b> (white cap). Contains human serum.
5 x 1.75 ml	CALIBRATOR A PR3	Ready to use set of 5 <b>Calibrators</b> . Calibrator A (green cap) 156 IU/ml, Calibrator B (violet cap) 62.5 IU/ml, Calibrator C (blue cap) 25 IU/ml, Calibrator D (yellow cap) 10 IU/ml, and Calibrator E (orange cap) 1 IU/ml. Derived from human serum containing PR3 antibodies. Concentrations in IU/ml are printed on the labels.
	CALIBRATOR B PR3	
	CALIBRATOR C PR3	
	CALIBRATOR D PR3	
	CALIBRATOR E PR3	
1 x 15 ml	IgG-CONJ HRP	HRP goat anti-human <b>IgG Conjugate</b> . Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. <b>Protect from light.</b>
1 x 15 ml	STOP H2SO4	Stop Solution*. Ready for use.
2 x	BUF WASH	Powder Wash Buffer. <b>Reconstitute to one liter each.</b>
1 x		Protocol Sheets

**Optional Components**

1 x 60ml BUF|WASH Liquid concentrated Wash Buffer. **Reconstitute to one liter.**

**Symbols used on labels**

- LOT Lot number
- REF Catalog number
- IVD In vitro diagnostic use
-  Use by
-  Storage temperature
-  Consult instructions for use
-  Number of tests
-  Manufacturer
-  Date of Manufacture
-  \*Danger; Causes severe skin burns and eye damage. Causes severe eye damage. Wear protective gloves, protective clothing, eye protection, face protection. Wash exposed skin thoroughly after

handling. If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

### Materials Required But Not Provided

- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm
- Automatic microplate washer capable of dispensing 200 µl

### SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. It is recommended that frozen specimens be tested within one year. Avoid repeated freezing and thawing of samples.

### PROCEDURE

#### Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

#### Test Method

**Step 1** Let all reagents and specimens equilibrate to room temperature.

**Step 2** Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.

**Step 3** For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.

**or**

For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

Qualitative			
A	Blank	S5	Etc.
B	-	S6	
C	+	S7	

Semi-Quantitative			
A	Blank	S1	Etc.
B	-	S2	
C	+	S3	

EN

D	Cal D	S8	
E	S1	S9	
F	S2	S10	
G	S3	S11	
H	S4	S12	
	1	2	3

D	Cal A	S4	
E	Cal B	S5	
F	Cal C	S6	
G	Cal D	S7	
H	Cal E	S8	
	1	2	3

- Step 4** Prepare a **1:101** dilution of the patient samples by mixing **5 µl** of the patient sera with **500ul** of Serum Diluent.
- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.  
**Note:** Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer’s instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** ( $\pm 5$  min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

### Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be  $<0.3$ . The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be  $<10$  IU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining IU/ml. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

## RESULTS

### Calculations

The concentrations of the patient samples can be determined by either of two methods:

#### 1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{IU/ml of Calibrator D} = \text{IU/ml Test Sample}$$

It is recommended that qualitative results be reported as “positive” or “negative.” Sample results greater than or equal to Calibrator D are considered positive.

#### 2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in IU/ml on the X-axis against the absorbance on the Y-axis and draw a

EN

point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as “positive,” “negative,” or “indeterminate” with IU/ml unit values. Indeterminate/borderline results should be retested and evaluated along with other laboratory methods, such as assays for detection of ANCA by IFA.

### Interpretation

Interpretation values were determined by testing 136 normal blood donors and non-autoimmune vasculitides disease control specimens. The mean of the normal subjects plus 3 SD was established as the assay cutoff and assigned a value of 10 IU/ml derived from the CDC international standard for PR3. MPO disease controls were not included in the cut-off determination of the PR3 assay. IMMCO suggests use of the reference range below. Each laboratory should validate assay values for their own conditions.

<b>anti-PR3 Ab value</b>	<b>Interpretation</b>
<10 IU/ml	Negative
10-12.5 IU/ml	Indeterminate (Borderline)
>12.5 IU/ml	Positive

### Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining IU/ml values, multiply the units obtained by the dilution factor.

### LIMITATIONS OF PROCEDURE

The assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only. The results obtained serve only as an aid in the diagnosis. Taken alone, these results should not be interpreted as diagnostic. This assay has not been validated in a pediatric population. MPO disease controls were not included in the cut-off determination of the PR3 assay

### EXPECTED VALUES


The expected values in a normal population are negative. However, 2-4% of apparently healthy, asymptomatic individuals may test positive for PR3 antibodies. In contrast, some patients with active disease may have undetectable levels of these antibodies. Immunosuppressive therapy, initiation or alteration in treatment should not be started on the basis of just positive PR3 antibody results, but rather on careful clinical observations.

The following table depicts the frequency of PR3 and MPO specific ANCA in sera from 112 ANCA associated vasculitides patients.<sup>11</sup> The incidence of PR3 antibodies varies depending upon the patient population. A compilation of incidence taken from the literature appears below.

Incidence of anti-PR3 and anti-MPO in ANCA associated vasculitides

Antibody association	Wegener's granulomatosis	Microscopic polyangiitis	Churg-Strauss syndrome
ANCA positive by IFA	78%	59%	67%
anti-PR3 positive	90%	0%	10%
anti-MPO positive	0%	62%	17%

The clinical utility of ANCA detection with cANCA specificity with positive anti-PR3 antibody specificity in various vasculitic disorders is illustrated below.<sup>5</sup>

Sets of clinical samples were tested on the ImmuLisa™ PR3 Antibody ELISA and another commercially available MPO antibody ELISA. Results demonstrating incidence in the populations for this study are provided below.

Patient Group	IMMCO			Other PR3 Ab ELISA		
	n	n Pos	% Pos	n	n Pos	% Pos
Disease Associated						
Wegener's granulomatosis	59	57	96.6%	59	56	94.9%
Glomerulonephritis	29	2	6.9%	29	2	6.9%
Undifferentiated ANCA positive	11	0	0.0%	11	0	0.0%
Disease Control						
Non-ANCA associated vasculitis	24	1	4.2%			
Inflammatory bowel disease	16	0	0.0%			
Systemic lupus erythematosus	32	1	3.1%	32	0	0.0%
Rheumatoid arthritis	8	0	0.0%	8	0	0.0%
Other autoimmune disease	16	0	0.0%	16	0	0.0%
Healthy normals	80	0	0.0%	80	0	0.0%

## PERFORMANCE CHARACTERISTICS

The utility of the ImmuLisa™ PR3 Antibody ELISA was evaluated by testing well-characterized PR3 antibody positive serum specimens from cANCA antibody positive subjects alongside disease controls and “normal” human sera. These specimens were also tested on commercially available ELISA kits. These results are summarized below.

### A. Method Comparison: ImmuLisa™ PR3 Antibody ELISA vs. other PR3 antibody ELISA:

		Other PR3 ELISA		
		Positive	Negative	Total
IMMCO	Positive	41	4	45
PR3	Negative	2	58	60
ELISA	Total	43	62	105

Positive Percent Agreement: 95.3% (95% CI 82.9% to 99.2%)

Negative Percent Agreement: 93.5% (95% CI 83.5% to 97.9%)

Overall Percent Agreement: 94.3% (95% CI 87.5% to 97.7%)

B. Cross Reactivity: A total of 136 potentially cross-reactive specimens from individuals with other autoimmune disorders or positive for other autoantibodies were tested for PR3 antibodies using the ImmuLisa™ PR3 Antibody ELISA.

Condition	n	Positive n
Glomerulonephritis	29	2
Undifferentiated pANCA positive	11	0
Non-ANCA associated vasculitis	24	1
Inflammatory bowel disease	16	0

EN		
Systemic lupus erythematosus	32	0
Rheumatoid Arthritis	8	0
Other autoimmune disease	16	0
<b>Total</b>	<b>136</b>	<b>3 (2.2%)</b>

### Precision

Precision was tested with 6 positive specimens selected throughout the range of the assay. Assay runs of three replicates of each specimen were conducted on three days. Repeatability was determined with 6 replicates of each specimen.

Sample	Mean (IU/ml)	Total Imprecision		Between days		Within run (Repeatability)	
		SD (IU/ml)	CV%	SD (IU/ml)	CV%	SD (IU/ml)	CV%
1	8.7	0.809	9.3%	0.809	9.5%	0.758	8.4%
2	12.2	0.812	6.7%	0.974	8.0%	0.547	4.5%
3	27.2	2.211	8.1%	2.183	8.3%	1.218	4.2%
4	57.0	3.017	5.3%	2.872	5.0%	3.458	6.1%
5	112.4	4.211	3.7%	5.059	4.5%	2.821	2.5%
6	139.3	6.199	4.4%	6.654	4.8%	6.043	4.3%

### Reproducibility

Six replicates of samples in the negative range, ~10% below cutoff, ~20% above cutoff and in the moderate positive range of the assay were performed to determine intra-assay qualitative reproducibility. Three replicates of the same samples were tested in three runs to determine inter-assay reproducibility. Assay results for these specimens produced 100% qualitative agreement.

### Linearity

Studies were performed using equidistant dilution series of positive samples with values throughout the calibrator range to determine linear range of the assay. The linear range of the assay was determined to be from 2.2 – 156 IU/ml. Results are summarized below. It is suggested that samples producing results greater than the top calibrator be diluted and retested.

Sample	Test Range	Slope (95% CI)	Y-Intercept (95% CI)	R <sup>2</sup>	% Recovery
1	3.7 to 77.2	1.07 (0.94 to 1.20)	-0.031 (-0.134 to 0.071)	0.986	89.8 to 111.9
2	3.8 to 158.2	1.012 (0.934 to 1.090)	-0.053 (-0.153 to 0.047)	0.994	88.9 to 102.2
3	3.4 to 43.4	1.023 (0.929 to 1.118)	0.009 (-0.045 to 0.063)	0.991	100 to 112.3

### Limit of Detection

The limit of detection (LoD) for PR3 antibody using this assay was determined to be 2.2 IU/ml based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples.

### Interference

Interference was studied by mixing sera with known PR3 antibody levels with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (2 g/L), Bilirubin (342 µmol/L), and Rheumatoid Factor (100 EU/ml).

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## REFERENCES

1. Davies DJ, Moran JE, Niall JF et al. Segmental necrotising glomerulonephritis with antineutrophil antibody: possible arbovirus etiology? *Br Med J*; 1982, 285:606.
2. Kao RC, Wehner NG, Stubitz KM et. al. Proteinase 3: a distinct human polymorphonuclear leucocyte proteinase that produces emphysema in hamster. *J Clin Invest*; 1988, 82:1963-1973.
3. van der Woude FJ, Rasmussen N, Lobatto S et. al. Autoantibodies against neutrophils and monocytes: tool for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet*; 1985, 1:425-429.
4. Cohen Tervaert JW, van der Woude FJ, Fauci AS et. al. Association between active Wegener's granulomatosis and anti-cytoplasmic antibodies. (ACPA). *Arch Int Med*; 1989, 149: 2461-2465.
5. Gross WL, Schmitt WH, Csernok E. ANCA and associated diseases: immunodiagnostic and pathogenetic aspects. *Clin Exp Imm*; 1993, 91:1-12.
6. Nölle B, Specks U, Lüdemann J et. al. Anticytoplasmic autoantibodies, their immunodiagnostic value in Wegener's granulomatosis. *Ann Int Med*; 1989, 111:28-40.
7. Cohen Tervaert J W, Huitema MG, Hene RJ et. al. Prevention of relapses in Wegener's granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre; *Lancet*; 1990, 336:709-711.
8. Kallenberg CG, Mulder AH, Cohen Tervaert JW. Antineutrophil cytoplasmic antibodies: a still-growing class of autoantibodies in inflammatory disorders. *Am J Med*; 1992, 93:675-682.
9. Hagen EC, Ballieux BE, van Es LA et al. Anti-neutrophil cytoplasmic autoantibodies: a review of the antigens involved, the assays, and the clinical and possible pathogenetic consequences. *Blood*; 1993, 81:1996-2002.
10. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health. 2007; (HHS Pub No [CDC] 93-8395).
11. Stoffel MP, Csernok C, Herzberg T et. al. Anti-neutrophil cytoplasmic antibodies (ANCA) directed against bactericidal/permeability increasing protein (BPI): a new seromarker for inflammatory bowel disease and associated disorders. *Clin Exp Immunol*; 1996, 104:54-59.
12. Jennette JC, Falk RJ. Diagnosis and management of glomerulonephritis presenting as acute renal failure. In. *Renal failures and transplantation, Medical clinics of North America* Eds. Mandal AK, and Hebert LA. Saunders Co; 1990, 74:893-908.

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**  
60 Pineview Drive  
Buffalo, NY 14228-2120 USA  
Telephone: (716) 691-0091  
Fax: (716) 691-0466  
Toll Free USA/Canada: 1-800-537-TEST  
E-Mail: [info@immco.com](mailto:info@immco.com)

*or your local product distributor*

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands