SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ImmcoStripe™ ANA Advanced LIA

Line immunoassay (LIA) for the detection of anti-nuclear antibodies

IVD

PRODUCT INSERT

REF 6011 Anti-Nuclear (ANA-Advanced) Antibody LIA 20 Determinations

INTENDED USE

A line immunoassay for the detection and identification of antibodies to nuclear and cytoplasmic antigens associated with systemic and certain organ specific autoimmune diseases.

SUMMARY AND EXPLANATION

Laboratory methods to detect levels of antinuclear antibodies (ANA) in sera are commonly used to aid in the diagnosis of connective tissue disorders including systemic lupus erythematosus (SLE), mixed connective tissue disease (MCTD), Sjögren's syndrome and scleroderma.¹⁻⁵ ANA occur in about 95% of SLE patients as well as in patients with other connective tissue diseases. ANA may also occur in other disorders such as chronic active hepatitis and primary biliary cirrhosis.⁶⁻⁸

ANA are not disease specific; however, identification of particular ANA specificities may be used to support diagnosis of a particular disease or disease subset. ANA specificities can be determined by gel precipitation, hemagglutination, IFA, ELISA, western blot and other methods.

Antibodies to U-1 RNP are characteristic of MCTD patients. Antibodies to Sm antigen occur exclusively in SLE. SS-A (Ro) antibodies are detected in Sjögren's syndrome and SLE.¹ Antibodies to Jo-1 (histidyl-tRNA synthetase) occur in patients with myositis.^{1,3} SS-B (La) is primarily reported in SLE and Sjögren's syndrome. Antibodies to both the 52 kD SS-A (Ro) and the 48 kD SS-B (La) are found in virtually all patients with neonatal lupus erythematosus and complete heart block, and the presence of both in SLE indicates a later age of disease onset and a lower frequency of lupus nephritis.⁹ Scl-70 antibodies react with human topoisomerase and are found in patients with diffuse scleroderma. Antibodies to nucleosome remodeling and genome maintenance associated proteins such as Mi2 and Ku respectively are an indication of SLE. Autoantibodies against nuclear or cytoplasmic antigens directed against ribonucleoproteins involved in protein synthesis ("anti-synthetases") such as Jo1, PL-7, PL-12, EJ and OJ or translational transport (SRP54) are detectable in about a fifth of polymyositis (PM) / dermatomyositis (DM) patients. The ANA Advanced LIA panel is able to detect antibodies to most of the above described antigens as well as AMA-M2, a major target of anti-mitochondrial autoantibodies, a characteristic feature of PBC (primary biliary cirrhosis). Anti-DFS70 antibodies have been associated with normal population. The ImmcoStripe™ ANA-Advanced LIA panel is capable of simplifying the detection of auto-antibodies against 22 key antigens that may also be observed by IFA screening on HEP-2 cell substrates. Following antigens are immobilized on the test strips: PM-Scl100, PM-Scl75, SSA/Ro52, SSA/Ro60, Jo1, Ribo-P, Nucleosomes, DNA, Histones, Sm, U1SnRNP68, U1SnRNP A, U1SnRNP C, SSB/La, Scl70, CENP-B, PCNA, Mi2, Ku, SRP54, AMA-M2 and DFS70 (LEDGF). In addition to the antigen lines, a cut-off, serum and conjugate control line is provided in the bottom of the strip.

PRINCIPLES OF PROCEDURE

To perform the test, strips are incubated with diluted patient serum. In positive sera antibodies specifically bind to one or more of the test lines on the strip. The strips are washed according to the protocol, and then the pre-diluted, ready-to-use conjugate is added to the test strips. After incubation and wash steps, the ready-to-use substrate is added to the strips. During a 10 minute incubation, conjugate and substrate binding produces visible blue/purple lines for serum, conjugate and cut-off control lines. If the sample is positive for any of the antigen coated test lines, it will show a reaction more intense than the cut-off line. Reactions are read visually and reported as positive, negative or equivocal (comparable to cut-off line).

EN

Materials Provided

ANA Advanced LIA **REF** 6011

Kits contain sufficient reagents to perform 20 determinations.

20 x **STRIP**ANA Adv**LIA**

1 x 120 µl **CONTROL**+**ANA**

1 x 30 ml **CONJ**LIA

1 x 30 ml **SUBSTRATE**

1 x 40 ml **SAMPLE**DIL

1 x 50 ml **BUF**WASH**LIA**

2 x

1 x

Line Immunoassay Test Strips, containing antinuclear antigen coated test lines and control lines. Ready for use.

Positive Control (red cap). Contains human serum positive for ANA.

IgG Conjugate.

Enzyme Substrate (amber bottle). Ready for use. **Protect from light**.

Sample Diluent

Wash Buffer Concentrate. **Reconstitute to one liter** with deionized or distilled water or as needed proportionally.

LIA 10 well Assay Trays

Report/Scoring Sheet

Symbols used on labels

Lot number

Catalog number

In vitro diagnostic use

Use by

Storage temperature

Consult instructions for use

Number of tests

Manufacturer

Date of Manufacture

Materials Required But Not Provided

- Clean 1000 ml graduated cylinder
- Non-serrated forceps (Filter forceps)
- Rocker or rotating platform shaker
- Absorbent paper or paper towels
- Deionized or distilled water
- Squeeze bottles to hold diluted wash buffer or distilled water
- Pipettes capable of delivering 10 to 1000 µl
- Disposable pipette tips
- Timer

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C; **do not freeze**.

All reagents must be brought to room temperature (18-25°C) and mixed thoroughly prior to use. Do not use if reagent is not clear or if insoluble precipitate is present. The reagents are stable until the indicated expiration date when stored at 2-8°C and protected from contamination, or as stated below after opening and/or reconstitution.

- Antigen coated test strips **STRIP**ANA Adv**LIA** are ready for use. Please allow the test strip bag to reach the room temperature before opening to avoid condensation and associated deterioration. Please re-pack unused test strips and store at 2-8°C in dark and dry conditions.
- Sample Diluent **SAMPLE**DIL is ready to use. After opening, Sample Diluent **SAMPLE**DIL is stable for at least 8 weeks when stored properly and protected from microbial or chemical contamination.

EN

- Reconstitute 1 part [BUF|WASH|LIA] into 19 parts of distilled or deionized water to produce 1 liter of Wash Buffer. Wash Buffer is stable for at least 8 weeks after reconstitution when stored properly and protected from microbial contamination.
- [CONJ|LIA] and [SUBSTRATE] are stable for at least 8 weeks after opening when stored properly and protected from microbial contamination. [SUBSTRATE] is light sensitive and must be stored in the provided amber colored bottle.

Antigen strips can only be used once. Do not interchange components of different lots. Do not use reagents beyond expiration date indicated on labels.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of the above materials.¹⁰

WARNING – Proclin 300 is a preservative. Upon disposal of liquids containing Proclin 300, flush with large volumes of water to dilute the components below active levels.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Specimens with gross hemolysis, elevated lipids or microbial contamination may interfere with the performance of the test and therefore must not be used. Store specimens at 2°-8°C for no longer than one week. For longer storage, serum specimens should be frozen. It is recommended that frozen specimens be tested within one year. Avoid repeated freezing and thawing of samples.

PROCEDURE

Procedural Notes

- Read Product Insert carefully before starting with the assay.
- Let serum specimens and test reagents equilibrate to room temperature for approximately 30 minutes prior to starting the test procedure. Return all unused specimens and reagents to the refrigerator promptly after use.
- Proper washing technique is critical to the satisfactory performance of the assay.
- Handle test strips with clean forceps or gloves only. Avoid touching the white antigen coated areas.
- The test lines are placed above the cut-off, serum and conjugate control lines as described in the schematic (Figure 1). Serum and conjugate control lines appear on the same piece of nitrocellulose at the bottom of the strip.
- Assign specimen identification numbers to the respective strips on the Report Form. Each strip has the strip number and lot number printed on the bottom for traceability.
- Complete all other relevant information on the Report Form prior to starting the assay.

Test Method

- Step 1** Using gloves or blunt forceps, peel off the required number of strips. Care should be taken not to touch the antigen coated areas with bare hands or pointed forceps.
- Step 2** Place required number of [STRIP|ANA Adv|LIA] labeled side up into individual wells of the assay tray.
- Step 3** Pipet 1.5 ml of [SAMPLE|DIL] into each well; make sure that the strips are completely submerged under the liquid.
- Step 4** Incubate the strips in [SAMPLE|DIL] for at least 10 minutes on a rocker. The blue color in coated antigen and control locations starts to disappear as the membrane is soaked.
- Step 5** Pipet 15 µl of serum or positive control sample into the appropriate wells to obtain a 1:101 dilution. Incubate 60 minutes (±5 min.) at room temperature on a rocker or rotating shaker.

Step 6 WASH: Aspirate sample solution into waste container. Thoroughly wash strips with the **reconstituted** Wash Buffer by squirting approximately 2ml of solution directly onto strips. Wash strips with gentle agitation for 5 minutes and aspirate solution into waste container. Repeat the wash two more times.

Caution: Complete washing of the strips between incubations is crucial to obtain valid results. Improper washing will result in high background staining. Do not allow the strips to become dry at any step during the assay.

Step 7 Pipet 1.0 ml of **CONJ|LIA** into each well. Incubate 30 minutes (± 5 min) at room temperature on rocker or rotating shaker.

Step 8 Repeat Step 6.

Step 9 Pipet 1.0 ml **SUBSTRATE|TMB** into each well and incubate with gentle shaking 10 minutes at room temperature in reduced light. The serum and conjugate control lines develop intense color after incubation in substrate. The cut-off control line develops into a blank to faintly colored line after the incubation.

Step 10 To stop the reaction, rinse strips 2x with distilled/deionized water by squirting approximately 2ml of water directly onto strips followed by aspiration. Do not soak/wash for more than 10 minutes as this may result in decreased sensitivity of the developed colored lines.

Step 11 Using blunt forceps remove strips from assay tray and place them gently onto absorbent paper and allow them to dry. Let the strips dry before analysis or affixing them on the report/scoring sheet.

Quality Control

Procedural Controls: Each strip has three procedural controls for the addition of serum and conjugate and a cut-off line for determining the weak or negative reactions.

Positive and Negative controls are available as optional components and may be run for additional quality control.

Individual labs are expected to optimize the substrate development time by ± 4 minutes based on the blot processor or manual methodology. It is recommended that the cut-off line should be faintly visible to the eye, post incubation with substrate.

Interpretation

The test strips contain control lines at the bottom and test line above the controls. The bottom end of the test strip (near serial number) has three control lines: the cut-off line, the serum control line and the conjugate control line from top to bottom. The cut-off allows the technician to determine the test result as positive, negative or indeterminate (+/-). The two procedure control lines ensure the addition of specimen, conjugate and substrate.

Compare the reaction of the test lines with those of the controls. Use of a magnifying glass can assist in observation of weak reactions.

- As labeled in Figure 1, the serum and conjugate control lines should be clearly positive indicating a successful experiment. The cut-off is a faint line with variation in intensity based on the experimental conditions. The schematic in Figure 1 shows an example test line. In the ANA Advanced LIA (Figure 2) there are 22 test lines and 3 control lines. The test line development depends on the sample. Positive reactions can occur in varying intensities from weak to strong. **Weak reactions should be compared with intensity of the provided cut-off line within the strip.** Reactions that are distinctly darker or denser than the intensity of the cut-off line should be considered positive.
- Strips may show a homogeneous or discolored background due to various interfering factors in lipemic or hemolytic sera. This effect can also be seen if the test strips are not sufficiently blocked or accidentally allowed to dry up during the assay.
- In case of weak positive and negative reactions, the reacted line intensity should be compared to cut-off line to determine the result as negative (weaker intensity than the cut-off line) or equivocal (+/-; indistinguishable from cut-off line). Visualization of weak reactions is improved when the strips are completely dry. Equivocal samples should be confirmed by an alternate method.


Figure 1: Schematic of reacted LIA strips with one test line.

EN

- Dried strips can be assembled in the provided report/scoring sheet. The plastic protective flap is permanently affixed to the report sheet on the left edge. Carefully peel the plastic flap in the right to left direction like a page of the book. Place the reacted strips on the adhesive tape in the respective slot and cover the plastic flap back in place. The protective plastic flap is designed to be reusable for multiple sessions of experiments and the strips can be assembled in respective slots. The technician can use the form to record the lot numbers of used reagents, specimen number and results/comments.

Figure 2: Schematic of Report Sheet

LIMITATIONS OF PROCEDURE

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Do not store specimen at 2-8°C more than a week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

The ANA Advanced LIA should only be used as an aid to diagnosis. Positive results may be found in other autoimmune conditions or certain infectious diseases. Hence, results should be evaluated and interpreted by a medical authority in light of the patient's clinical history and other laboratory findings.

EXPECTED VALUES

At 1:40 titer, elevated levels of anti-nuclear and cytoplasmic antibodies may be present in up to 30% individuals in a "normal" population as determined by IFA methodology. At a titer of 1:320 positivity rate in normal population has been reported to drop to approximately 3%. Incidence of autoantibody positivity varies to a great extent based on the specific antigens, studies and selected cohort of specimens. The expected values in a normal population are negative on LIA. LIA panels are designed to provide optimal sensitivity and specificity. It is recommended to confirm the positive result by an alternate methodology.

PERFORMANCE CHARACTERISTICS

Antigen Specificities

The ANA Advanced LIA is able to detect autoantibodies to the following antigens: PM-Scl100, PM-Scl75, Ro-52, Ro-60, Jo-1, Ribo-P, Nucleosomes, DNA, Histones, Sm, U1RNP68, U1RNP A, U1RNP C, SS-B/La, Scl-70, Cenp-B, PCNA, Mi-2, Ku, SRP-54 and AMA-M2 and DFS70 (LEDGF). Refer to Figure 2 for individual antigen and control line positions.

Cross-reactivity

A panel of potentially cross-reactive autoimmune disease sera from conditions not associated with ANA was tested using the ANA-Advanced LIA test. 1 out of 550 determinations demonstrated a positive reaction indicating a specificity of 99.8% in this population.

Interference

Interference was studied by mixing sera with known levels of autoantibodies for each analyte with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: Hemoglobin (5 mg/ml), Bilirubin (0.4 mg/ml), Rheumatoid Factor (200 EU equivalent) and Triglycerides (25 mg/ml). Interference studies have been done according to the CLSI guidelines (publication EP7-A2).

Reproducibility

Assays of samples in the negative range, equivocal and positive range were performed to determine qualitative reproducibility from run to run and operator to operator. Results produced 100% qualitative agreement.

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
1. Kumar V, Beutner EH and Chorzelski TP. Autoimmunity of the skin. In "Concepts in Immunopathology", Vol 1, Cruse JM and Lewis RE Jr, Eds, Karger, Basel. 1985, 318-53.
2. Reimer G, Cornell RC and Tan EM. The biochemical nature of nuclear antigens reactive with antinuclear antibodies. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP, Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed. 1987, 519-531.
3. Tan EM, Chan EKL, Sullivan KF and Rubin RL. Antinuclear antibodies (ANAs): Diagnostically specific immune markers and clues toward the understanding of systemic autoimmunity. Clin Immunol Immunopathol; 1988, 47:121-141.
4. Von Muhlen CA, Tan EM. Autoantibodies in the diagnosis of systemic rheumatic diseases. Semin Arthritis Rheum; 1985, 21:323-58.
5. Tan, EM. Antinuclear Antibodies: Diagnostic Markers for Autoimmune Diseases and Probes for Cell Biology. "Advances in Immunology" Vol. 44; 1987, pp 93-151.
6. Manns M, Gerken G, Kyriatsoulis A and Meyer zum Büschenfeld KH. Significant autoimmune markers of autoimmune liver disorders: Current status. J Clin Lab Anal; 1987, 1:362-370.
7. Mackay IR. Autoimmunity and the liver. Clin Aspects Immunity. 2: 8- 17, 1988.
8. McMillan SA, Alderdice JM, McKee CM et al. Diversity of autoantibodies in patients with anti-mitochondrial antibody and their diagnostic value. J Clin Path 4: 232-236, 1987.
9. Ben-Chetrit, E, Fox RI, and Tan EM. Dissociation of Immune Responses to the SS-A(Ro) 52kD and 60kD polypeptides in Systemic Lupus Erythematosus and Sjögren's Syndrome. Arth. Rheum. 33, 349-355, 1990.
10. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health. 2007; (HHS Pub No [CDC] 93-8395).

For technical assistance please contact:

 **IMMCO Diagnostics, Inc.**
60 Pineview Drive
Buffalo, NY 14228-2120 USA
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immco.com

or your local product distributor

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Rev. OCT2017

Document No. PI6011