SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

FUNGAL IMMUNODIFFUSION SYSTEM

REF Catalog No. 100100

IVD For *in vitro* diagnostic use only

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[REF] Catalog No. 100100 [IVD] For *in vitro* diagnostic use only

INTENDED USE

The Meridian Bioscience Fungal Immunodiffusion Reagents are standardized, purified preparations for the *in vitro* determination of precipitating antibodies to four systemic fungal pathogens: *Histoplasma capsulatum* (*H. cap.*), *Blastomyces dermatitidis* (*B. dermat.*), *Coccidioides immitis* (*C. imm.*) and *Aspergillus fumigatus* (*A. fum.*). These reagents are optimized for use in the Ouchterlony double diffusion method.

EXPLANATION

H. cap., *B. dermat.*, *C. imm.* and *A. fum.* are causative agents of deep-seated mycoses. These fungi present a diagnostic challenge to the microbiology laboratory and the physician.

Radiographically, lesions produced by the systemic fungi may be difficult to distinguish from tuberculous lesions, neoplasms, cancerous tissue and each other. Symptoms are often unremarkable and may mimic various pneumonias, sarcoidosis, cancer and other maladies. Culturally and histologically the organisms may be difficult to demonstrate, even after repeated attempts. With the exception of *A. fum.* growth is very slow, requiring 2-6 weeks (3,4,5).

Frequently serology offers the only evidence available to guide treatment, suggest prognosis or lead to the selection of more definitive diagnostic techniques such as intensive culture or biopsy. In addition, quantitative serology, such as complement fixation testing, can provide important information on the effects of chemotherapy (11).


The fungal immunodiffusion test is a simple and reliable tool for the evaluation of suspected mycoses. Anticomplementary sera may be tested using this technique. The test also provides specificity data on reactions obtained by the complement fixation method. No expensive equipment is required and the technique is simple enough to be performed by any laboratory, thus providing an excellent screening method.

BIOLOGICAL PRINCIPLES

The Meridian Bioscience Fungal Immunodiffusion System is based on the principle of double diffusion as described by Oudin and Ouchterlony. An antibody and its homologous soluble antigen are placed in separate wells cut in an agarose diffusion medium and allowed to diffuse outward. Between the two wells, a concentration gradient of each of the reaction components is established ranging from antigen excess closest to the antibody well to antibody excess closest to the antigen well. A visible line of precipitate forms at the point of equivalence (8,9,10).

Antigens or antibodies may be tested for "identity" by placing a test well of the substance in question adjacent to the wells of a known system. If the antigen-antibody complexes are identical, the precipitin lines form an unbroken line of identity with the known system. Partial and nonidentity reactions are also possible. (See Figure 1)


Figure 1.

2

A partial identity reaction occurs when certain components of the antigens (or antibodies) are identical and others are not. The "spur" represents the components which are unrelated. A nonidentity reaction will occur when the antigen-antibody complexes are different. The resulting "X" or crossed reaction indicates that two unrelated complexes are present. The well patterns of the Meridian immunodiffusion plate are arranged to provide each patient test well with a reference known system so that identity reactions are readily apparent. (see Figure 2)

Figure 2. Arrangement of Sera and Antigens in Immunodiffusion Well Pattern.

Positive control serum wells 1 & 4; Patient sera, wells 2, 3, 5 & 6; Antigen, well 7.

COMPONENTS

A. **Antigens.** The following antigens are available in a 0.2ml size as a part of the Fungal Immunodiffusion System (Catalog #100100) or separately in a 1.0ml size. Coccidioidin antigen is also available in a 5.0ml size.

1. **Histoplasmin.** A purified mycelial cultural filtrate from the growth of *Histoplasma capsulatum* containing the "H" and "M" antigens. The preparation contains 100 units/ml* of each antigen. Lipids, media components, and most other antigenic components have been removed. (Note: the "H" band appears closest to the antibody well while the "M" band is closest to the antigen well.)

2. **Blastomycin.** A purified extract of the yeast phase growth of *Blastomyces dermatitidis* containing the "A" antigen at a concentration of 100 units/ml*.

3. **Coccidioidin.** A purified cultural filtrate containing the "F" antigen of *Coccidioides immitis* at a concentration of 100 units/ml*.

4. **Aspergillin.** A purified, alcohol-precipitated, carbohydrate preparation from mycelial phase cultures of *Aspergillus fumigatus*, *A. niger* and *A. flavus* containing approximately 450 micrograms/ml of carbohydrate. Nomenclature has not been designated for these antigens at this time.

These antigens contain thimerosal at a final concentration of 0.01% as a preservative. Care should be exercised, however, to avoid microbial contamination.

*These units are internal Meridian standards designed only to insure lot to lot reproducibility of antigen concentration.

B. **Control Sera.** The following lyophilized control sera are provided in 0.4ml size as a part of the Fungal Immunodiffusion System or are available separately in the 1.0ml size. All control sera are raised in hyper-immune (not infected) goats against purified antigen preparations.

1. **Histoplasmosis positive control.** Contains specific antibodies directed against the "H" and "M" antigens of *H. capsulatum*. (Note: the "H" band appears near the antibody well, while the "M" band is closest to the antigen well.)

2. **Blastomycosis positive control.** Contains specific antibodies directed against the "A" antigen of *B. dermatitidis*.

3. **Coccidioidomycosis positive control.** Contains specific antibodies directed against the "F" antigen of *C. immitis*.

4. **Aspergilliosis positive control.** Contains antibodies directed against *A. fumigatus*, *A. niger* and *A. flavus*. Standard nomenclature has not as yet been designated for these antigens. Three bands are apparent against the Meridian aspergillin.

These positive control sera are preserved with 0.1% sodium azide. Nevertheless, care should be exercised to avoid microbial contamination.

3

The control sera are preserved with sodium azide. Accumulation of this chemical in metal plumbing fixtures may represent an explosion hazard. It is therefore recommended that excess control serum simply be discarded in an appropriate waste receptacle or flushed down a drain with copious amounts of water.

RISK AND SAFETY PHRASES

CONTROL SERA: HARMIFUL - SODIUM AZIDE

RISK PHRASES:

- 32 Harmful if swallowed.
- 32 Contact with acids liberates very toxic gas.

SPECIMEN COLLECTION

For optimal results, sterile serum is obtained from the blood of a patient. If a delay is encountered in specimen processing, refrigeration for up to 72 hours is permissible. Specimens may be stored up to six months at -20°C with no loss of activity, provided they are not repeatedly thawed and refrozen. Specimens in transit between laboratories should be maintained at 2°-8°C for optimal results. Specimens may be preserved with 0.01% thimerosal or 0.1% sodium azide if necessary.

PERFORMING THE TEST

This test should be performed by qualified personnel per local regulatory requirements.


A. Reconstitution of positive control sera.

1. Using a 1.0ml pipet or tuberculin syringe add 0.38ml (for 0.4ml vials) or 0.95ml (for 1.0ml vials) of reagent quality water to the control serum vials.
2. Gently "flick" the bottom of the vials to mix.
3. Allow the vials to stand at room temperature for about 1 hour, occasionally inverting to mix. (Do not shake)
4. If, upon its first use, the serum appears inactive, it may be that sufficient standing time was not allowed. In this case repeat the test.
5. Aliquot and freeze any control serum which will not be consumed within one month. (see SHELF LIFE AND STORAGE)

B. Test Procedure

1. Please refer to figures 2, 3 and 4 for proper designation of the test well patterns.
2. Using the provided capillary pipets (Fungal Immunodiffusion System only), fill the control serum wells of an ID plate with the appropriate control sera as follows:
 - a. *Histoplasma* control serum - Series I, wells 1 & 4.
 - b. *Blastomyces* control serum - Series II, wells 1 & 4.
 - c. *Coccidioides* control serum - Series III, wells 1 & 4.
 - d. *Aspergillus* control serum - Series IV, wells 1 & 4.

Figure 3.

C. Immunodiffusion plates. Nine ID plates are provided in the Fungal Immunodiffusion System. When four sera per plate are tested, sufficient materials are provided to test 36 patient sera. Additional plates may be obtained by ordering Meridian Catalog #101009. The immunodiffusion plates consist of 0.9% agarose in glycine/phosphate buffer at pH 7.4 +/- 0.2. Note: recent advances in plate formulation have eliminated the need for Immunodiffusion Band Intensifying Fluid.

NOTE: Due to the standardization necessary to the production of high quality fungal serology reagents, performance of antigens, control sera and plates with materials other than those produced by Meridian Bioscience, Inc. cannot be guaranteed. The user assumes full responsibility for any modifications to the procedures published herein.

MATERIALS AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

1. Reagent quality water.
2. Moist chamber. Any convenient container such as a petri dish, plastic box or glass jar with a tight fitting cover containing moist filter paper or paper toweling is satisfactory provided the ID plates are stationary, level and remain hydrated during incubation.
3. Reading light. Darkfield plate readers are commercially available, however a satisfactory system can be devised using a high-intensity lamp in which the front of the bulb shield-reflector is covered with black construction paper with a small hole (1-2 cm in diameter) cut for illumination. Alternatively, the plate may be held at a 45° angle to almost any bright indirect light source for adequate visualization.

SHELF LIFE AND STORAGE

Fungal antigens are best stored at refrigerator temperature (2°-8°C). Repeated freezing and thawing is detrimental to these antigens. When stored at 2°-8°C the antigens are stable until the expiration as indicated on each bottle.

Positive control sera are stable in the lyophilized state until the expiration as indicated on each bottle when stored at 2°-8°C. Once reconstituted it is suggested that they be aliquoted and frozen if they will not be consumed within one month. At -20°C a minimum nine month life may be expected. Repeated freezing and thawing should be avoided. When the positive control sera are in use, the period at room temperature should be kept as short as possible.

Immunodiffusion plates are individually packaged and if unopened and stored at 2°-8°C, are stable until the expiration as indicated on each plate. If plate packages are opened and the plates must be stored prior to use, they must be kept in a tightly sealed, humid container to avoid drying. A resealable plastic container with a small moist sponge or toweling is ideal.

USER QUALITY CONTROL

User quality control is performed on each test run by the inclusion of a positive control serum for each antigen. A negative control may usually be found among the various patients tested; however, Meridian Catalog #102901 Negative Control Serum may be used if desired.

The results obtained with the controls should be recorded in an appropriate log book to maintain high quality testing and comply with regulatory requirements.

PRECAUTIONS

All components of the Meridian Fungal Immunodiffusion System are for *in vitro* diagnostic use only. Occasionally, bands due to C reactive protein may be evident in the serum of patients when reacted with *Aspergillus* sp. antigen. While this reaction will not form an identity line with the supplied control serum, interpretation may be difficult due to the diffuse nature of some bands. Flooding the plate with 1.0M sodium citrate and allowing it to remain at room temperature will eliminate this reaction. When handling blood specimens, adequate measures should be taken to prevent the dissemination of etiologic agents potentially present in the specimen.

- Record the name, date, and/or lab number of the first patient on line 2 of the left hand column (Figure 4) of the recording form and, using a capillary tube, fill the well 2 of each of series I, II, III and IV with the patient's serum. Using a fresh capillary tube for each patient, repeat this procedure for additional sera in wells 3, 5, and 6 of each series.
- Using a fresh capillary tube, fill the center well of series I with Histoplasma antigen. Repeat this process for Blastomycosis antigen (series II), Coccioidin antigen (series III), and Aspergillus antigen (series IV).
- Number and date the ID plate and the analysis form.
- Place the ID plate in a moist chamber and incubate at room temperature for 24 hours. Note: bands may be evident in as little as six hours for some antigens. It will still be necessary to confirm the results of the other antigens later.
- After 24 hours read and record the ID bands on the analysis form (Figure 4). (see **Reading the Test**) An interim report should be issued at this point if no identity reactions are observed. Positive results should be reported immediately.

Figure 4.

Well Number	SERIES I (Histoplasma) Reading	SERIES II (Blastomycosis) Reading	SERIES III (Coccioidin) Reading	SERIES IV (Aspergillus) Reading
1. Control Serum	Histo Control	Blasto Control	Cocci Control	Asper Control
2. CARTER, R. 3287	H and W BANDS	NEGATIVE	NEGATIVE	NEGATIVE
3. LESTER, C. 2519	NEGATIVE	NEGATIVE	NEGATIVE	NEGATIVE
4. Control Serum	Histo Control	Blasto Control	Cocci Control	Asper Control
5. BARRIS, T. 1562	NEGATIVE	NEGATIVE	NEGATIVE	3 - BANDS
6. HASKAMP, D. 1154	NEGATIVE	NEGATIVE	1 - BAND	NEGATIVE
7. Antigen	Histoplasmin	Blastomycitin	Coccioidin	Aspergillin

After 24 hour incubation, arrow bands observed on the well diagram
 NOTE: The correct bands have been drawn on the well diagram above. If these bands are not observed on the plates after 24 hour incubation, the test is not valid and must be repeated.

- An additional 48 hour incubation period is recommended to confirm a negative result with any reagent. A final report is made at the conclusion of this period.

READING THE TEST

Any bright indirect light source may be used to aid in seeing the precipitin bands. A dark background is preferred. For instance, the plate may be held next to the shade of a high intensity lamp whose light beam is directed straight downward, over a black counter-top.

A bright direct light source may also be used if the plate is held in the light beam at a 45 degree angle against a dark background.

Light boxes, providing indirect light against a dark background, prepared specifically for reading immunodiffusion plates are also commercially available.

Particular attention should be paid to the orientation of bands produced by the patient serum in relation to the control bands. A smooth junction of the bands is indicative of an identity reaction. If the control band is seen to bend toward a position in front of the patient well, it is indicative of patient antibody at a low titer.

Should no reaction be observed between the positive control and the antigen the test must be repeated.

Occasionally non-identity reactions may be noted in the aspergillin series (Series IV). While these reactions will not form an identity reaction with the anti-Aspergillus sp. control serum, their diffuse nature may cause some confusion. Such reactions may be tested by flooding with 1 M sodium citrate for 5 minutes. Reactions which disappear after this treatment are due to C reactive protein and may be disregarded.

INTERPRETATION

A band of identity with a known positive control is indicative of patient antibody against the antigen in question. Partial identity reactions are regarded as positive for antibody against the antigen only if no other identity reaction is present on the plate. Non-identity reactions are regarded as a negative test (11).

REPORTING OF TEST RESULTS

If identity or partial identity reactions are observed the test may be reported as "Positive for antibodies to (name of antigen)" or "Identity reaction with (name of antigen)." Non-identity reactions and sera which are unreactive may be reported as "Negative for antibodies to (name of antigen) by immunodiffusion."

EXPECTED VALUES

In general, an identity reaction against a given antigen is indicative of active or recent past infection. Partial identity reactions are also important indicators of probable disease, especially when no identity reaction is observed among other agents. Non-identity reactions may be apparent when the disease state is caused by a mycotic agent other than the one tested. A non-identity reaction, however, is a negative test for the organism in question.

If activity against any of the fungal agents is observed, a vigorous attempt should be made to demonstrate the organism culturally for confirmation.

Histoplasmosis

The clinically significant antigens for *H. cap.* are designated "H" and "M" antigens. (Note: the "H" band appears closest to the antibody well while the "M" band is closest to the antigen well.) The "M" band has been found in about 63% of patients with active histoplasmosis, while the "H" band is apparent in only 27%. The "H" band is rarely observed in the absence of antibodies against the "M" antigen. The "M" band may appear in patients who have recently recovered from histoplasmosis as well as in the serum of previous histoplasmosis patients who have been recently skin tested. Since the test may be negative in as many as 10% of culturally demonstrated cases, the absence of either an "H" or an "M" band does not rule out histoplasmosis (4, 5).

Blastomycosis

The ID test for antibodies directed against *B. dermatitidis* is positive in approximately 80% of the culturally provable cases, rendering a negative test of little value. Care must be taken to demonstrate an identity reaction with the "A" antigen since cross-reactions from other fungi are common (particularly *H. cap.*) (11).

Coccioidin mycosis

Antibodies directed against the "F" antigen are indicative of active or recent past (up to 1 year) infection. The ID test is usually positive within 4 weeks after infection and remains positive throughout the clinically active disease. Latex and complement fixation testing may provide important additional information regarding patient status (4).

Cross-reactions may be seen in patients harboring other systemic fungi (especially *H. cap.*) so care must be exercised when reading for identity reactions (7, 11, 13).

Aspergilliosis

Less is known about the serological-clinical relationship of aspergilliosis than the other systemic mycoses. Demonstration of one or more precipitin bands against the *Aspergillus* antigen is indicative of infection, colonization or allergy to an *Aspergillus* sp. In general, the greater the number of precipitin bands observed, the higher the probability of infection caused by a member of this genus. Approximately 20% of patients with allergic bronchopulmonary aspergilliosis will form no precipitin bands with *Aspergillus* antigen. A negative test does not exclude the possibility of aspergilliosis (6).

LIMITATIONS OF THE TEST

The high rate of negative serologic tests observed among culturally demonstrable cases limits the predictive value of a negative test.

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BIBLIOGRAPHY

1. Bauman, D. S. and C. D. Smith. 1975. Comparison of immunodiffusion and complement fixation tests in the diagnosis of histoplasmosis. *J. Clin. Microbiol.* 2:77-80.
2. Crowle, A. J. 1973. Immunodiffusion. Second Edition. Academic Press. New York. Chapters 1, 2, and 4.
3. Heiner, D. C. 1968. Diagnosis of histoplasmosis using precipitin reactions in agar gel. *Pediatrics* 22:616-627.
4. Huppert, M. and J. W. Bailey. 1965. The use of immunodiffusion tests in coccidioidomycosis. *Amer. J. Clin. Path.* 44:364-368.
5. Kaufman, L. 1973. Value of immunodiffusion tests in the diagnosis of systemic mycotic diseases. *Ann. Clin. Lab. Sci.* 3:141-46.
6. Kaufman, L. 1975. Practical applications of immunologic procedures for the diagnosis of opportunistic fungus infections. In: *Proceedings of the Second International Conference on Opportunistic Fungal Infections*. C. C. Thomas Publishing, Springfield, pp. 47-57.
7. Kaufman, L. and M. J. Clark. 1974. Value of the concomitant use of complement fixation and immunodiffusion tests in the diagnosis of coccidioidomycosis. *Appl. Microbiol.* 28:641-43.
8. Ouchterlony, O. 1949. Antigen-antibody reactions in gels and the practical application of this phenomenon in the laboratory diagnosis of diphtheria. *Med. Diss.*, Stockholm.
9. Ouchterlony, O. 1968. Handbuch der Immunodiffusion und Immunoelektrophoresis. Ann Arbor Publishers, Inc. Ann Arbor, gelose. *Premiere parsa.* Ann Inst. Pasteur 75:30-52.
10. Oudin, J. 1946. L'analyse immunochimique qualitative. *Méthode par diffusion des antigènes au sein de l'immunosérum précipitant gelose.* *Premiere parsa.* Ann Inst. Pasteur 75:30-52.
11. Palmer, D. F., et al. 1977. Serodiagnosis of Mycotic Diseases. C. C. Thomas Publishing, Springfield, pp.78.
12. Smith, C. E. 1943. Coccidioidomycosis. *Med. Clin. North Amer.* 27:790-807.
13. Smith, C. E., M. T. Saito, R. R. Beard, K. M. Kepp, R. W. Clark, and B. V. Eddie. 1950. Serological tests in the diagnosis of coccidioidomycosis. *Am. J. Hyg.* 52:1-2.
14. Tosh, F. E. 1975. Problems in the serological diagnosis of opportunistic fungus diseases. In: *Proceedings of the Second International Conference on Opportunistic Fungal Infections*. C. C. Thomas Publishing, Springfield, pp. 36-43.

USA/Corporate Office
 Manufactured by 3471 River Hills Drive
 Cincinnati, Ohio 45244
 Telephone: (513) 271-3700
 Orders/Customer Service: 1-800-543-1980
 Technical Support Center: 1-800-343-3858
 Information Fax: (513) 272-5432
 Ordering Fax: (513) 271-0124

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Tel.: +33 (0) 493187210
 Fax: +33 (0) 493187211
 e-mail: mdeur.info@fnac.net

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product

	Manufactured by		For Performance Evaluation Only
	Authorized Representative		Temperature Limitations
	Catalog Number		Use By/Expiration Dating Information
	In vitro Diagnostics		Sufficient for 'X' Tests
	Batch Code/Lot Information		CE Symbol

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Halderheiweg 6
 5282 SN Boxtel
 Holland

Tel.: +31 (411) 62 11 66
 Fax: +31 (411) 62 48 41
 e-mail: meridian@wxs.nl

For technical assistance, please call our Technical Support Center at 1-800-343-3858 between the hours of 8 AM and 6 PM, Eastern Standard Time. To place an order, please call our Customer Service Department at 1-800-543-1980.

SN 10400

Printed in the U.S.A. on recycled paper
 Copyright © Meridian Bioscience, Inc.

Rev. 11/03