SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

External Control Material for use with Alethia Malaria DNA Amplification Assays

REF 479970

IVD

Rx Only

INTENDED USE

The Alethia Malaria External Control Kit contains Positive and Negative Control Reagents for use with the Alethia Malaria or Alethia Malaria PLUS DNA Amplification Assays. External controls are used as part of a routine quality control program to aid the user in detection of unexpected conditions that may lead to test errors.

SUMMARY AND EXPLANATION OF THE TEST

Quality control testing is performed to detect factors such as reagent deterioration, adverse environmental or test conditions, or variance in operator performance that can cause test errors. External control materials are not built into the test system, but are prepared, amplified and detected in the same manner as patient specimens.

BIOLOGICAL PRINCIPLES

The Alethia Malaria and Alethia Malaria PLUS DNA Amplification test kits should be tested with known positive and negative external control reagents on a periodic basis to verify that the kits are performing correctly. The frequency with which a laboratory performs external controls will be affected by:

- The proficiency level of the laboratory,
- The laboratory's own internal requirements,
- The requirements of the laboratory's accrediting agencies,
- The number of new and different operators performing the test,
- Whether a new kit lot is being added to testing,
- Whether the kit lot is from a different shipment, and
- Whether deviations from the manufacturer's stated storage or handling conditions have occurred.

When unacceptable quality control test results are obtained, all test results should be considered invalid. QC test failures are an indication that either reagents, test environment or operator performance have changed.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Alethia Malaria Positive Control:** Tris-buffered solution with plasmid containing DNA inserts (*Plasmodium* sp. and human mitochondrial DNA inserts) and azide (0.09%) as a preservative.
- Alethia Malaria Negative Control:** Tris-buffered solution with plasmid containing human mitochondrial DNA inserts) and azide (0.09%) as a preservative.

MATERIALS PROVIDED SEPARATELY

- Alethia Malaria PLUS DNA Amplification Test Kit**, Meridian Bioscience, Inc. Catalog Number 481125
- Alethia Malaria DNA Amplification Test Kit**, Meridian Bioscience, Inc. Catalog Number: 480925

MATERIALS NOT PROVIDED

- Disposable latex gloves, powder free
- DNase/RNase free, aerosol resistant pipette tips

EQUIPMENT NOT PROVIDED

- Interval timer
- Vortex Mixer (optional)
- Micropipette capable of dispensing 50 µL
- Micropipette capable of dispensing 250 µL (catalog number 481125 only)
- Alethia Incubator/Reader™, Meridian Bioscience, Inc. Catalog Number: 610189

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- This is a quality control reagent and is used to evaluate the performance of the Alethia Malaria and Alethia Malaria PLUS DNA Amplification Assays. It is not directly used to test patient samples.
- Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
- Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
- Quality Control Programs for Molecular Testing Laboratories should be employed.
- The Alethia Malaria Test Devices contain lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
- The Alethia Malaria Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
- Dispose of used Alethia Test Devices, **M-prep** Columns, and tubes immediately after processing, leaving the device latch securely in place. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit components at the temperature indicated on the label.

REAGENT PREPARATION

- Bring Alethia Malaria External Positive and Negative Controls, and all Alethia Malaria kit components to room temperature (19-30 C) before use. Incorrect results may be obtained if control materials and components are not brought to room temperature prior to use.
- Use one Alethia Malaria Test Device for each Positive Control and Negative Control to be tested.

For controls to be processed by Alethia Malaria PLUS: Refer to the Alethia Malaria Package Insert (SN11035) for information regarding column preparation, orientation and use. Use 1 **M-prep** column per sample. Remove the column top cap, followed by the bottom twist-off tip, and place the column tip into an ST tube to drain. Drained **M-prep** columns should be used within 1 hour.

QUALITY CONTROL SAMPLE PREPARATION

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

NOTE: Ensure that the Alethia Incubator/Reader is powered and the required performance verifications have been completed prior to initiation of Quality Control Sample Preparation. Refer to the Alethia Incubator/Reader Operator's Manual for further information regarding instrument set-up and operation.

- Preparation of controls using Alethia Malaria (Catalog number 480925)**
Note: The Positive and Negative Control require separate Buffer tubes and SMP PREP IV.
 - Mix Positive Control by inversion 5 times or by vortexing 10 seconds. Using a micropipette, add 50 µL of appropriate Positive or Negative Control to one Buffer I tube. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Incubate the sample for 2 minutes.
 - Mix by inversion 5 times or by vortexing for approximately 10 seconds. Using a micropipette, immediately transfer 50 µL of prepared sample into SMP PREP IV. Mix by inversion 5 times or vortex for 10 seconds.
 - Gently squeeze the SMP PREP IV and slowly collect 5 to 10 drops into a clean Tube I. Proceed to Quality Control Test Procedure.
- Preparation using Alethia Malaria PLUS (Catalog number 481125)**
Note: The Positive and Negative Control require separate Buffer I tube and M-prep Columns.
 - Mix Positive Control by inversion 5 times or by vortexing 10 seconds. Using a micropipette, add 50 µL of appropriate Positive or Negative Control to one Buffer I tube. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Incubate the sample for 2 minutes.
 - Mix by inversion 5 times or by vortexing for approximately 10 seconds. Using a micropipette, immediately transfer 250 µL of the prepared sample to the top of an appropriately labeled and prepared **M-prep** Column. Wait approximately 2 minutes, or until the sample has been absorbed by the column and flow stops.
 - Using a micropipette, add 250 µL of **M-prep** Buffer II to the top of the **M-prep** Column. Discard the pipette tip. The column will have a red appearance after the addition of **M-prep** Buffer II. Wait approximately 2 minutes, or until the red-colored buffer is absorbed by the column and flow stops.
 - Remove the last drop of liquid from the column tip with the ST Tube. Discard the tube.
 - Place a clean ST Tube under the **M-prep** Column. Using a micropipette, add 250 µL of **M-prep** Buffer III to the top of the **M-prep** Column. Discard the pipette tip. Wait approximately 2 minutes or until flow stops.
 - Remove the last drop of liquid from the column tip with the ST Tube. Label the tube and proceed to the Quality Control Test Procedure.

QUALITY CONTROL TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single Alethia Incubator/Reader run.

- Remove 1 Alethia Malaria Test Device from its protective pouch for the Negative Control. Carefully open the device, holding the chambers such that the lyophilized reagents will not fall out upon opening. Place the device on a flat surface or in a rack that can accommodate the device.
- Using a micropipette, transfer 50 µL of the prepared Negative Control to both the TEST (Left/White Bead) and CONTROL (Right/Yellow Bead) chambers of the Alethia Malaria Test. Take care to not introduce air to the reaction mixture. Do not mix reactions with pipette.
- Close the Alethia Test Device and fasten the latches securely.
- Repeat Test Procedure Steps 1-3 for the prepared Positive Control sample to be tested.
- Tap device(s) on the bench top or mix to remove air bubbles. Carefully examine the Test Device(s) for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
- Insert the Alethia Test Devices into the Alethia Incubator/Reader and initiate run **using the Malaria Program**. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat control testing as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat control testing as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia Incubator/Reader Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

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Meridian Bioscience, Inc.
 Corporate Office
 3471 River Hills Drive
 Cincinnati, Ohio 45244 USA
 Telephone: 513.271.3700
 Orders/Customer Service:
 800.543.1980
 Technical Support Center:
 800.343.3858
 Information Fax: 513.272.9432
 Ordering Fax: 513.271.0124

Manufactured By

Meridian Bioscience Europe S. r.l.
 Via dell'Industria, 7
 20020 Villa Cortese, Milano
 ITALY
 Tel: +39 0331 43 36 36
 Fax: +39 0333 43 36 16
 Email: info@meridianbioscience.eu
 WEB: www.meridianbioscience.com/au

Authorized Representative

Meridian Bioscience Europe s.a./n.v.
 2 Avenue du Japon - 1420 Braine l'Alleud
 BELGIUM
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 Email: info@meridianbioscience.eu

Meridian Bioscience Europe France
 34 rue de Ponthieu - 75008 Paris
 FRANCE
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 Email: info@meridianbioscience.eu

Meridian Bioscience Europe b.v.
 Postbus 301 - 5460 AH Veghel
 NETHERLANDS
 Tel: +31 (0) 41 62 11 66
 Fax: +31 (0) 41 62 48 41
 Email: info@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum
CE	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung
	Female / Femminile / De sexe féminin / Hembra / Frau
STERILE R	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Sterilisation par irradiation aux rayons gamma / Esterilizado por irradiación gamma / Sterilization durch Gammastrahlen
STERILE EO	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Sterilisation par oxyde d'éthylène / Esterilizado por óxido de etileno / Sterilization durch Ethylenoxid
	Contains sufficient for n tests / Contenuto sufficiente per n saggi / Contenu suffisant pour n tests / Contenido suficiente para n ensayos / Inhalt ausreichend für n Prüfungen
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät
MEDIA	Media / Terreno di trasporto / Milieux / Medio / Medium
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum
	LASER RADIATION: Avoid Exposure to beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Éviter toute exposition au faisceau / Radiación Laser: Evite Exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement Laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr
CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
RoHS	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Peligrosas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Caution, consult accompanying documents / Attenzione, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
	Recycle - do not dispose of as general waste / Riciclare - non eliminare come rifiuto generico / Recycle - ne pas jeter dans une poubelle / Recicle - no desheche como basura general / Recycling - dieses Produkt nicht über den Hausmüll entsorgen
HT TUBE	Heat Treatment Tube / Provetta per il trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhrchen zur Hitzebehandlung
	For IVD Performance Evaluation Only / Solamente per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	HOT SURFACE: Keep hands away from hot surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
IPX-0	CAUTION: Protect from water / ATTENZIONE: Protegger dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer
MIN OIL	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl
	Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
ST TUBE	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhrchen mit Schnappverschluss
COL	Preparation Column / Colonna di preparazione del campione / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufarbeitung
BUF SMP	Sample Buffer / Soluzione tampone per il campione / Tampon de l'échantillon / Tampón de muestra / Probenpuffer
SMP PREP	Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenvorbereitung
R Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig
TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß

SN11038

REV. 10/18

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.