**SZABO
SCANDIC**

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic

FOR EXPORT USE ONLY. NOT FOR SALE IN THE USA.

Malaria and Malaria PLUS DNA Amplification Assays

DNA Amplification Assays for the Detection of *Plasmodium* sp.

REF

480925, 481125

IVD

INTENDED USE

The Alethia Malaria and Alethia Malaria PLUS DNA amplification assays, performed on the Alethia Incubator/Reader™, are qualitative in vitro diagnostic tests for the direct detection of Plasmodium sp. DNA in human venous EDTA whole blood specimens from individuals with signs and symptoms of malarial infection. Results from Alethia Malaria assays are intended to be used as an aid in the diagnosis of human malaria infection.

Alethia Malaria assays utilize loop-mediated isothermal DNA amplification (LAMP) technology to detect Plasmodium sp. DNA by targeting segments of the Plasmodium genome. Alethia Malaria assays do not distinguish between Plasmodium species.

Alethia Malaria and Alethia Malaria PLUS are intended for use in hospital, reference or state laboratory settings. The devices are not intended for nonlaboratory point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

The Alethia Malaria and Alethia Malaria PLUS molecular assays are based on loop-mediated amplification (LAMP)^{1,2} technology. The assays target a region of the Plasmodium genome that is conserved across *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae*, and *Plasmodium knowlesi*. The Alethia Malaria assays target a 214 base pair (bp) sequence of a Plasmodium sp. mitochondrial DNA noncoding region.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian Alethia Incubator/Reader Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The Alethia Malaria kit includes Alethia Malaria Test Devices, Sample Preparation Apparatus IV (SMP PREP IV), Alethia Buffer I, and tubes. Whole blood specimens are first treated with Buffer I; cells are disrupted and nucleic acids are released. The lysate is added to SMP PREP IV and filtered into a clean Tube I by gently squeezing. Sample effluent collected after filtration contains nucleic acids for use with the Alethia Malaria Test Device.

The Alethia Malaria PLUS kit includes Alethia Malaria Test Devices, **M-prep** Columns, Alethia Buffer I, **M-prep** Buffer II, **M-prep** Buffer III, and ST Tubes. Whole blood specimens are processed using size exclusion chromatography by gravity flow to separate and purify nucleic acids. **M-prep** Columns provide the matrix for purification of nucleic acids while the included buffers facilitate cell lysis and sample elution. Whole blood specimens are first treated with Buffer I; cells are disrupted and nucleic acids are released. Sample material is applied to the **M-prep** Column and allowed to move through the column by gravity, with the largest molecules moving through the column matrix first, or fastest. **M-prep** Buffer II is applied to the **M-prep** Column, facilitating movement of the sample through the column matrix. **M-prep** Buffer II contains a non-reactive colorant that provides a visual indicator of sample progression. **M-prep** Buffer III is applied to the **M-prep** Column after Buffer II has entered the gel matrix. Sample effluent collected after the addition of **M-prep** Buffer III contains nucleic acids for use with the Alethia Malaria Test Device. Neither Alethia Malaria nor Alethia Malaria PLUS require specialized laboratory equipment.

The Alethia Malaria Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with Plasmodium sp.-specific primers and a CONTROL chamber with human mitochondrial DNA-specific primers. Human mitochondrial DNA in whole blood samples and the human mitochondrial DNA-specific primers in the Test Device CONTROL chambers function as the Internal Control for the assay. During specimen preparation, human mitochondrial DNA is extracted with the *Plasmodium* sp. mitochondrial DNA to allow for parallel processing of target DNA and Control DNA through amplification and detection. The Internal Control monitors DNA extraction, amplification inhibition, assay reagent performance and sample processing effectiveness. The Control target must be amplified and detected in the final reaction or the test is considered invalid and results are not reported.

The Alethia Incubator/Reader monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start ($Signal_{start}, S_1$) and at the assay Run End ($Signal_{final}, S_2$). The Alethia Incubator/Reader calculates the change in light transmission between Run End and Run Start ($S_2 - S_1$) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S_2/S_1 ratios less than 70% are reported as 'POSITIVE'; TEST chamber S_2/S_1 ratios greater than or equal to 70% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_2/S_1 ratios less than 85% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S_2/S_1 ratios greater than or equal to 85% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

Malaria is a public health challenge world-wide, with an estimated 3.2 billion people at risk for infection in 97 countries.³ The World Health Organization (WHO) estimates 198 million cases occurred in 2013 that led to approximately 584,000 deaths.³ Malaria is caused by the intracellular parasite *Plasmodium* and is transmitted to humans by female *Anopheles* mosquitoes.³ Once humans are infected, *Plasmodium* exhibit multiple life stage morphologies for initial growth in liver cells, followed by invasion of red blood cells.⁴ *Plasmodium* replicate within red blood cells, causing cell rupture and clinical symptoms of the disease.⁴

There are five *Plasmodium* species known to cause malaria in humans (*Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium malariae*, *Plasmodium ovale*, and *Plasmodium knowlesi*), with the majority of infections caused by *P. falciparum* and *P. vivax*. Malaria treatment is dependent on the *Plasmodium* species, clinical status of the patient, and antimicrobial susceptibility.⁵ Infections caused by *P. falciparum* and *P. knowlesi* are typically more severe and require rapid treatment.⁵

The WHO recommends diagnostic testing for all individuals with suspected malaria.³ Thin and thick smear microscopy is currently the gold standard diagnostic method. However, accuracy and precision of microscopy is highly dependent on the technician skill (limit of detection ranging from 4-100 parasites/ μ L or higher), with false negative and false positive rates between 10-25% and approximately 2-3 fold variability in quantitation.^{6,7,8} In addition, errors in the *Plasmodium* species identification are common in routine microscopy examination.⁸ Immunoassay-based rapid diagnostic tests (RDTs) are also an available diagnostic method. Quality assessment testing performed by the WHO demonstrated that at low parasite densities (200 parasites/ μ L), only 76% and 42% of available RDTs could detect *P. falciparum* and *P. vivax*, respectively.³

The Alethia Malaria assays provide rapid results with higher sensitivity than microscopy and RDTs.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

For Alethia Malaria, Catalog Number 480925

1. **Alethia Malaria Test Device:** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either *Plasmodium* sp.-specific primers (TEST Chamber) or human mitochondrial DNA-specific primers (CONTROL Chamber).
2. **Alethia Buffer I:** Lysis solution containing 0.2N sodium hydroxide.
3. **Alethia Sample Preparation Apparatus IV (SMP PREP IV):** Tris-buffer solution containing 0.09% azide as a preservative.
4. **Alethia Tube I:** 1.5 mL tubes

For Alethia Malaria PLUS, Catalog Number 481125

1. **Alethia Malaria Test Device:** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either *Plasmodium* sp.-specific primers (TEST Chamber) or human mitochondrial DNA-specific primers (CONTROL Chamber).
2. **Alethia Buffer I:** Lysis solution containing 0.2N sodium hydroxide.
3. **M-prep Buffer II:** Solution containing phenol red and 0.09% azide as a preservative.
4. **M-prep Buffer III:** Tris-buffered solution with 0.09% azide as a preservative.
5. **M-prep Column:** 5.0 cm long column with twist-off tip and plug seal cap containing chromatography resin in a tris-buffered solution containing sodium azide (0.09%) as a preservative.
6. **ST Tubes:** 2.0 mL screw-top microcentrifuge tubes.

MATERIALS PROVIDED SEPARATELY

1. Alethia Malaria External Control Kit, Meridian Bioscience, Inc. Catalog Number: 479970

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips
3. Blood collection tubes with EDTA anticoagulant

EQUIPMENT NOT PROVIDED

1. Interval Timer
2. Vortex mixer (optional)
3. Micropipette capable of dispensing 50 μ L
4. Micropipette capable of dispensing 250 μ L (Catalog Number 481125 only)
5. Alethia Incubator/ReaderTM, Meridian Bioscience, Inc. Catalog Number: 610189

PRECAUTIONS

1. All reagents are for *in vitro* diagnostic use only.
2. Do not interchange kit reagents and Test Devices between lots. Individual lots of Tube I and ST Tubes are interchangeable within each kit type provided they are within the assigned expiration date when used.
3. The Alethia Incubator/Reader may produce incorrect results if the Malaria assay program is not used for testing.
4. Follow Biosafety Level 2 and Good Laboratory practices during testing.⁹ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
6. Quality Control Programs for Molecular Testing Laboratories, including proper use and care of equipment, should be employed.¹⁰
7. The Alethia Malaria Test Devices contain lyophilized reagents. The protective pouches should not be opened until ready to perform the assay.
8. The Alethia Malaria Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
9. Dispose of used Alethia Malaria Test Devices, **M-prep** Columns, and tubes immediately after processing. Leave the Test Device latch securely in place. DO NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD AND PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com, for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit components at the temperature indicated on the label.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (19-30 C) before use. Ensure Buffer I has been warmed to room temperature completely prior to use and no precipitate is visible. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.


COLUMN PREPARATION (For Alethia Malaria PLUS Catalog Number 481125 only)

Visually inspect each **M-Prep** Column and ensure the column bed and filters were not disturbed during transportation. Columns with disturbed resin beds or improperly placed filters may not perform correctly.

Column Set-up:

1. Use 1 **M-prep** Column for each sample to be tested. Remove the top cap and the bottom twist-off tip.
2. Place the column tip into an ST Tube. The ST Tube should be held at a slight angle (approximately 15 degrees) from the **M-prep** Column.
NOTE: Do not seat the column directly in the collection tube as this may create a seal and cause improper flow through the column. Refer to **Figure 1** for proper column and collection tube placement.
3. Let the **M-prep** Column drain into the ST Tube. Drained **M-prep** Columns should be used within 1 hour.
4. Proceed with Sample Preparation for **Alethia Malaria PLUS**.

FIGURE 1: M-prep Column and ST Tube placement

SPECIMEN COLLECTION AND PREPARATION

Sample type: Human venous whole blood samples with EDTA as a preservative.

Sample Collection: Venous whole blood samples should be collected into a specimen tube containing EDTA according with institutional guidelines for collection of clinical specimens for malaria infection. Ensure the blood collection tube is inverted immediately after collection at least 8 times, or according to the manufacturer's instructions.

EDTA venous whole blood samples may be stored at 2-30 C after collection and during transportation to the laboratory. Samples should be tested as soon as possible, but may be stored for up to 7 days at room temperature (19-30 C) or up to 14 days refrigerated (2-8 C) prior to testing. Samples that will not be tested within this time frame should be frozen immediately at -20 C for up to 30 days until tested. Samples may be frozen and thawed 2 times after storage at -20 C prior to testing with the Alethia Malaria assays.

SPECIMEN PREPARATION

NOTE: Ensure that the Alethia Incubator/Reader is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the Alethia Incubator/Reader Operator's Manual for further information regarding instrument set-up and operation.

NOTE: Ensure specimens are at room temperature (19-30 C) before specimen preparation.

Specimens to be tested with the Alethia Malaria assay may be prepared using one of the following methods:

1. Alethia Malaria Specimen Preparation (Catalog Number 480925)

- a. Invert EDTA whole blood specimen 2-3 times to mix.
- b. Add 50 μ L of the collected venous whole blood sample (with EDTA) to one tube of Buffer I. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Hold the sample for 2 minutes.
- c. Mix by inversion 5 times or by vortexing for approximately 10 seconds and immediately transfer 50 μ L of lysate into SMP PREP IV. Mix by inversion 5 times or by vortexing for 10 seconds.
- d. Gently squeeze the SMP PREP IV and slowly collect 5 to 10 drops into a clean Tube I. Visually verify that the eluate is tinted red to reddish brown. Label the tube with the specimen identification and proceed to the Test Procedure.

2. Alethia Malaria PLUS Specimen Preparation (Catalog Number 481125)

NOTE: Sample elution steps with **M-prep** Columns should take no longer than 30 minutes. Samples that take longer than 30 minutes to elute should be discarded and re-tested with the original patient sample.

- a. Invert EDTA whole blood specimen 2-3 times to mix.
- b. Add 50 μ L of the collected venous whole blood sample (with EDTA) to one tube of Buffer I. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Hold the sample for 2 minutes.
- c. Mix by inversion 5 times or by vortexing for approximately 10 seconds and, using a micropipette, immediately transfer 250 μ L of the prepared sample to the top of an appropriately labeled and prepared **M-prep** Column. Wait approximately 2 minutes, or until the sample has been absorbed by the column and flow stops. This step should take no longer than 30 minutes.
- d. Using a micropipette, add 250 μ L of **M-prep** Buffer II to the top of the **M-prep** Column. Discard the pipette tip. The column will have a red appearance after the addition of **M-prep** Buffer II. Wait approximately 2 minutes, or until the red-colored buffer is absorbed by the column and flow stops. This step should take no longer than 30 minutes.

- e. Remove the last drop of liquid from the column tip with the ST Tube. Discard the tube.
- f. Place a clean ST Tube under the **M-prep** Column. Using a micropipette, add 250 µL of **M-prep** Buffer III to the top of the **M-prep** Column. Discard the pipette tip. Wait approximately 2 minutes or until flow stops. This step should take no longer than 30 minutes.
- g. Remove the last drop of liquid from the column tip with the ST Tube. Visually verify that the eluate is tinted red to reddish brown. Label the tube with sample identification information and proceed to the Test Procedure.

TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single Alethia Incubator/Reader run.

1. Remove 1 Alethia Malaria Test Device from its protective pouch per sample. Carefully open the device, holding the chambers such that the lyophilized reagents will not fall out upon opening. Place the device on a flat surface or in a rack that can accommodate the device.
2. Using a micropipette, transfer 50 µL of the sample to both the TEST (Left/White Bead) and CONTROL (Right/Yellow Bead) chambers of the Alethia Malaria Test. Take care to not introduce air to the reaction mixture. Do not mix reactions with pipette.
3. Close the Alethia Test Device and fasten the latches securely.
4. Tap device(s) on the bench top or mix to remove air bubbles. Carefully examine the Test Device(s) for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
5. Repeat Test Procedure Steps for all samples to be tested.
6. Insert the Alethia Test Devices into the Alethia Incubator/Reader and initiate run **using the Malaria Program**. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains <i>Plasmodium</i> sp. target DNA.
	NEGATIVE	No <i>Plasmodium</i> sp. DNA detected.
	INVALID	No reportable result. Repeat the test using the original sample. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia Incubator/Reader Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Each device contains an internal control that controls for amplification inhibition, assay reagents, DNA preparation, and sample processing effectiveness. Human mitochondrial DNA, which serves as the internal control DNA, is isolated from the whole blood sample and processed through all steps of the procedure. Primers for amplification of internal control DNA are present in the Control Chamber of the Alethia Test Device.
2. Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
3. Alethia Malaria External Positive Control Reagent is supplied separately (Catalog 479970). Alternatively, previously characterized clinical or contrived *Plasmodium* sp. positive blood samples can be used as an external positive control. A qualified negative human whole blood sample with EDTA may be used as an external negative control. It is recommended that the reactivity of each new lot and each new shipment of Alethia Malaria and Alethia Malaria PLUS be verified on receipt and before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia Malaria and Alethia Malaria PLUS test kit should not be used in patient testing if the external controls do not produce the correct results.
4. A separate Test Device must be used for each external control reagent.

EXPECTED VALUES

The incidence of *Plasmodium* sp. as detected by the Alethia Malaria assay during the 2015 clinical study was 68.1% (147/216) using the Alethia Malaria sample preparation method (Catalog number 480925) and 70.6% (149/211) using the Alethia Malaria PLUS Method (Catalog number 481125). The overall incidence of each *Plasmodium* species is provided below.

Prevalence by <i>Plasmodium</i> species				
Species	Alethia Malaria (n=216)		Alethia Malaria PLUS (n=211)	
	Total Positive	Prevalence	Total Positive	Prevalence
<i>P. falciparum</i>	137	63.4%	134	63.5%
<i>P. vivax</i>	0	0%	0	0%
<i>P. ovale</i>	1	0.5%	1	0.5%
<i>P. malariae</i>	1	0.5%	1	0.5%
<i>P. knowlesi</i>	0	0%	0	0%
Unknown*	8	3.7%	13	6.2%

*Samples were negative by microscopy and therefore species could not be determined.

LIMITATIONS OF THE PROCEDURE

1. This product can only be used with the Alethia Incubator/Reader instrument.
2. Alethia Malaria and Alethia Malaria PLUS do not distinguish between Plasmodium species.
3. Hematocrit (>57.5%) and hemoglobin (>30 g/dL) above normal physiological levels may produce invalid results with the Alethia Malaria assays.
4. Performance of the Alethia Malaria assays has not been established for packed red blood cell specimens.
5. Performance of Plasmodium knowlesi with the Alethia Malaria assay was established using purified genomic DNA only; whole organism testing was not performed.
6. Alethia Malaria and Alethia Malaria PLUS are qualitative assays and do not provide quantitative values or information about organism load.
7. The detection of nucleic acids is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
8. Organism nucleic acid may persist *in vivo*, independent of organism viability. Alethia Malaria and Alethia Malaria PLUS do not distinguish between viable and nonviable organisms.
9. As with all molecular based diagnostic tests, (A) False negative results may occur from the presence of inhibitors, technical error, sample mix-up or low numbers of organisms in the clinical specimen; (B) False positive results may occur from the presence of cross-contamination by target organisms, their nucleic acids or amplified product, and from non-specific signals.

SPECIFIC PERFORMANCE CHARACTERISTICS

The performance characteristics of the Alethia Malaria and Alethia Malaria PLUS DNA Amplification Assays were established in clinical studies conducted in November 2015 in Senegal, Africa. Performance characteristics of the assay were compared to the reference method, thin and thick film microscopy. Plasmodium species identification for all positive samples was determined by microscopy.

Samples included prospective and retrospective venous whole blood specimens with EDTA from patients with signs and symptoms of Plasmodium infection. A total of 216 eligible, de-identified whole blood specimens were evaluated. 66 specimens were tested prospectively and 150 specimens were prospectively-collected and stored frozen prior to testing by the Alethia Malaria assay (retrospective). Each specimen was prepared by both Alethia Malaria and Alethia Malaria PLUS sample preparation methods. An additional five samples were considered ineligible for the Alethia Malaria PLUS assay due to process error. Therefore, there are 211 eligible samples for Alethia Malaria PLUS. All samples were tested prospectively by microscopy.

Analysis of assay performance with prospective and retrospective specimens indicated there is no difference in performance for fresh and frozen specimens with the Alethia Malaria and Alethia Malaria PLUS assays. The tables below summarize the performance of the Alethia Malaria and Alethia Malaria PLUS DNA Amplification Assays for both prospective and retrospective samples combined.

Alethia Malaria Performance Compared to Microscopy

	Microscopy			Alethia	Performance				
	Pos	Neg	Total		INV ^a			95% CI	
Alethia Malaria	Pos	139	8	147	0(1)	Sensitivity	100.0%	139/139	97.3-100.0%
	Neg	0	67	67	2(2)	Specificity	89.3%	67/75	80.3-94.5%
	Total	139	75	214	2				

^aInitial invalid results are reported within the parentheses. The final number of invalid samples remaining after repeat testing is shown before the parenthesis.

Alethia Malaria PLUS Performance Compared to Microscopy

	Microscopy			Alethia	Performance				
	Pos	Neg	Total		INV			95% CI	
Alethia Malaria PLUS	Pos	136	13	149	0	Sensitivity	100.0%	136/136	97.3-100.0%
	Neg	0	62	62	0	Specificity	82.7%	62/75	72.6-89.6%
	Total	136	75	211	0				

Samples were collected from patients ranging in age from 4 to 77 years. 15 were between 1 – 12 years of age; 75 were between 13 – 21 years of age; 125 were ≥22 years old; the age of one patient was not defined. There was no noted performance difference based on age. The study population included 84 (38.9%) female and 131 (60.6%) male patients; the gender of one patient was not defined. There is no expectation that assay performance is influenced by gender.

ANALYTICAL SENSITIVITY

Limit of Detection (LoD) was determined using *Plasmodium falciparum* and *Plasmodium vivax*. The LoD was confirmed using 20 replicates and a stated probability (eg. 95%, where 19/20 replicates are positive) of obtaining positive responses.

	Alethia Malaria	Alethia Malaria PLUS
Plasmodium Species	Parasites/ μ L	Parasites/ μ L
<i>P. falciparum</i> (3D7)	2	0.25
<i>P. vivax</i> (India VII)	0.125	0.063

ASSAY REACTIVITY

Quantified stocks of *P. malariae*, *P. ovale*, and *P. knowlesi* DNA were diluted in negative whole blood to approximately three times the limit of detect for each Alethia Malaria and Alethia Malaria PLUS sample preparation method (approximately 72 copies/ μ L or 6 parasites/ μ L for Alethia Malaria; approximately 9 copies/ μ L or 0.75 parasites/ μ L for Alethia Malaria PLUS). All species reacted at the concentrations tested.

REPRODUCIBILITY

Reproducibility studies were carried out by three internal laboratories. Blind-coded panels of 10 samples were supplied to participating laboratories. The panels included contrived *Plasmodium falciparum* (strain 3D7) samples manufactured as moderate positive samples, low positive samples, or high negative samples. The panel also included one negative whole blood sample. Testing was performed by at least 2 different operators at each laboratory on the same day (intra-assay variability) for five days (inter-assay variability). Three lots each of Alethia Malaria or Alethia Malaria PLUS and 8 Alethia Incubator/Reader instruments were used in this study. External Positive and Negative Controls were tested with each panel; a qualified negative blood sample was used as the Negative Control.

Reproducibility Study Summary: Alethia Malaria Sample Preparation Method								
Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
Moderate Positive (8 parasites/ μ L)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive (2 parasites/ μ L)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
High Negative (0.458 parasites/ μ L)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Negative	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
<hr/>								
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

Reproducibility Study Summary: Alethia Malaria PLUS Sample Preparation Method								
Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
Moderate Positive (1 parasite/ μ L)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive (0.5 parasites/ μ L)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
High Negative (0.029 parasites/ μ L)	29/30	96.7%	28/30	93.3%	30/30	100.0%	87/90	96.7%
Negative	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
<hr/>								
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

CROSSREACTIVITY STUDIES

Crossreactivity studies employed positive (*P. falciparum* 3D7) and negative whole blood specimens inoculated with bacterial or fungal organisms to a minimum concentration of 1.0×10^8 CFU/mL, virus at a minimum of 1.0×10^5 TCID₅₀/mL, or protozoa to a minimum concentration of 1.0×10^5 organisms/mL. Where whole organisms were not available, 1.0×10^6 copies/mL for genomic DNA was tested. None of the following organisms or their genetic material reacted with the Alethia Malaria assays:

Babesia microtei, *Borrelia burgdorferi*, *Haemophilus influenza*, *Klebsiella pneumoniae*, *Leishmania donovani*, *Leptospira interrogans*, *Mycobacterium smegmatis*, *Salmonella typhi*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Streptococcus pneumoniae*, *Treponema pallidum*, *Trypanosoma cruzi*, *Trypanosoma rangeli*, *Vibrio parahaemolyticus*, *Cytomegalovirus (CMV)*, *Epstein-Barr virus (EBV)*, *Hepatitis B virus*, *Hepatitis C virus*, *Herpes simplex virus 1 (HSV 1)*, *HIV-1*, *Human Papilloma Virus (HPV)*, and *Rubella virus*.

The following organisms or their DNA could not be obtained for testing and were evaluated through in silico analysis. None of the following organisms are expected to react with the Alethia Malaria assays based on this analysis:

Anaplasma phagocytophilum, *Clostridium botulinum*, *Orientia tsutsugamushi*, *Rickettsia conorii*, *Rickettsia rickettsii*, *Rickettsia typhi*, Chikungunya virus, Dengue virus (types 1-4), West Nile virus, and Yellow Fever virus.

Human genomic DNA was nonreactive at 1.0×10^6 copies/mL.

TESTS FOR INTERFERING SUBSTANCES

Interference testing was performed in the presence of worst-case concentrations of chemical and biological substances introduced directly into contrived low positive (*P. falciparum* 3D7) and negative whole blood samples. No interference was observed with the following substances for either Alethia Malaria sample preparation methods:

Acetaminophen (1 mg/mL), Amoxicillin (0.1 mg/mL), Artesunate (1 mg/mL), Aspirin (1 mg/mL), Atovaquone (0.1 mg/mL), Cephalexin (0.1 mg/mL), Chloroquine (1 mg/mL), Ciprofloxacin (0.1 mg/mL), Clindamycin (1 mg/mL), Doxycycline hydrochloride (1 mg/mL), Erythromycin (0.1 mg/mL), Hydroxychloroquine sulfate (1 mg/mL), Ibuprofen (1 mg/mL), Lumefantrine (1 mg/mL), Mefloquine (1 mg/mL), Primaquine phosphate (1 mg/mL), Proguanil hydrochloride (1 mg/mL), Pyrimethamine (1 mg/mL), Quinine sulfate (1 mg/mL), Sodium Citrate (0.11 M), elevated Bilirubin (>0.15 mg/mL), elevated leukocytes (buffy coat) >10% v/v, serum albumin (>0.03 g/mL), and triglycerides (>9.9 mg/mL).

Hematocrit (>57.5%) and hemoglobin (>30 g/dL) above normal physiological levels may produce invalid results with the Alethia Malaria assays.

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REFERENCES

1. Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. Mol Cell Probes 2002;16:223-29.
2. Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. J Biochem Biophys 2004;59:145-47.
3. World Malaria Report 2014. Global Malaria Programme, World Health Organization. http://www.who.int/malaria/publications/world_malaria_report/en/.
4. Cowman AF and Crabb BS. Invasion of Red Blood Cells by Malaria Parasites. Cell 2006; 124: 755-766.
5. Treatment of Malaria (Guidelines for Physicians). Centers for Disease Control and Prevention, July 2013. www.cdc.gov/malaria/resources/pdf/clinicalguidance.pdf.
6. O'Meara WP, Barcus M, Wongsrichanalai C, Muth S, Maguire JD, Jordan RG, Prescott WR, and McKenzie FE. Reader Technique as a Source of Variability in Determining Malaria Parasite Density by Microscopy. Malaria J 2006;5:118.
7. McKenzie FE, Sirichaisinthop J, Miller RS, Gasser RA, and Wongsrichanalai C. Dependence of Malaria Detection and Species Diagnosis by Microscopy on Parasite Density. Am J Trop Med Hyg 2003;69(4):372-376.
8. Wongsrichanalai C, Barcus MJ, Muth S, Utamihardja A, and Wernsdorfer WH. A Review of Malaria Diagnostic Tools: Microscopy and Rapid Diagnostic Test (RDT). Am J Trop Med Hyg 2007;77(6):119-127.
9. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
10. CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.

Manufactured By

Meridian Bioscience, Inc.
Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S. r. L
Via dell'Industria, 7
20020 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.com/eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus: 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.