SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

PREMIER™

CRYPTOCOCCAL ANTIGEN

REF Catalogue Number 602096

IVD In Vitro Diagnostic Medical Device

EIA for the Detection of Cryptococcal Antigen in Serum and CSF

INTENDED USE

The Premier Cryptococcal Antigen enzyme immunoassay (EIA) is a screening or a semi-quantitative test system for the detection of capsular polysaccharide antigens of *Cryptococcus neoformans* in serum and cerebrospinal fluid (CSF).

SUMMARY AND EXPLANATION OF THE TEST

Cryptococcosis is typically initiated as a subclinical pulmonary infection caused by the encapsulated soil yeast *C. neoformans*.¹ Although *C. neoformans* can infect any organ, the most serious complications occur when the organism invades the central nervous system of the immunocompromised patient.¹⁻³ Four major serotypes (A-D) exist within the capsular polysaccharide antigen, with 80-85% of isolates being either A or B.⁴

Detection of cryptococcal antigens in serum or CSF is used successfully in the diagnosis of meningitis associated with infection by *C. neoformans*.^{1-3,7,8} When both serum and CSF are evaluated the sensitivity of antigen detection is greater than 90%.^{1,2,8} Antigen detection is more sensitive than the India ink mount^{5,8} and is a valuable aid in establishing a diagnosis when the culture is negative.⁶

Determining the relative changes in antigen levels has prognostic value and is used to monitor the efficacy of drug therapy.^{1,2,5,6}

The Premier Cryptococcal Antigen EIA detects and provides a reliable semi-quantitative determination of all cryptococcal antigen serotypes present in specimens.

PRINCIPLE OF THE TEST

The Premier Cryptococcal Antigen EIA utilizes anticryptococcal polyclonal antibodies adsorbed to microwells and a detection system based upon a monoclonal-peroxidase conjugate. If cryptococcal antigens are present in the sample, a complex is formed between the antigens, enzyme conjugate, and the adsorbed antibody. After washing to remove unbound conjugate, a substrate solution is added. Color develops in the presence of bound enzyme. The amount of color is related to the quantity of cryptococcal antigen present.

REAGENTS/MATERIALS PROVIDED

- 1. Premier Enzyme Conjugate** – Cryptococcal antigen specific monoclonal antibody conjugated to horseradish peroxidase in buffered protein solution containing 0.02% Thimerosal as a preservative.
- 2. Premier 20X Wash Buffer II** – Concentrated wash buffer containing 0.2% Thimerosal as a preservative.
- 3. Premier Substrate II** – Buffered solution containing urea peroxide and tetramethylbenzidine.
- 4. Premier Stop Solution II** – 2N Sulfuric acid. **CAUTION:** Avoid contact with skin. Flush with water if contact occurs.
- 5. Premier Positive Control** – Cryptococcal antigen in buffered protein solution containing 0.01% Thimerosal and 0.16% Gentamicin as preservatives.
- 6. Premier Sample Diluent** – Buffered protein solution containing 0.02% Thimerosal as a preservative.
- 7. Premier Antibody Coated Microwells** – Breakaway plastic microwell strips coated with polyclonal antibodies specific for cryptococcal antigen.
- 8. Microwell Strip Holder**

The maximum number of tests obtained from this kit is listed on the outer box.

MATERIALS NOT PROVIDED

1. Pipettes capable of delivering 50 µL, 100 µL and 200 µL.
2. Test tubes (12 x 75 mm) for dilution of sample.
3. Distilled or deionized water.
4. EIA plate reader capable of reading absorbance at 450 nm or 450/630.
5. Squirr bottle or automatic EIA plate washer.
6. Timer.
7. Graduated cylinder for making 1X Wash Buffer.

PRECAUTIONS

1. All reagents are for *in vitro* diagnostic use only.
2. Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards. Wear disposable gloves while handling patient specimens and performing the test procedure.
3. All reagents should be gently mixed and at room temperature before use.
4. Do not interchange reagents from different kit lot numbers or use expired reagents.
5. Hold reagent vials vertically to ensure proper drop size and delivery.
6. Return caps to correct vials.
7. The Positive Control contains inactivated cryptococcal antigen. This reagent, used wash buffer, used microwells, etc., should be disposed of as potentially hazardous material. The microwell plate holder should be disinfected and retained for future assays.
8. Avoid skin contact with Stop Solution (2N Sulfuric Acid). Flush with water immediately if contact occurs.

9. **Unused microwells must be placed back inside the foil ziplock pouch and sealed immediately. It is important to protect strips from moisture.**
10. Avoid splashing when dispensing diluted sample into microwells by placing the pipette tip about halfway down the microwell and dispensing slowly down the side of the microwell.
11. Perform microwell washing as directed in the assay procedure. **Inadequate washing may cause elevated background.**
12. Incubation times and reagent concentrations have been optimized for sensitivity and specificity. Avoid any deviations from set incubation times or dilutions.
13. Do not store specimens in a frost-free type freezer. Freezing and thawing of a specimen repeatedly can affect test results.

RISK AND SAFETY PHRASES

Sulfuric Acid

Stop Solution II - Xi Irritant

R PHRASES

R:36/38 Irritating to eyes and skin

Methanol

Substrate II - F – Highly Flammable, T – Toxic

R PHRASES

R:11 Highly flammable.
R:66 Repeated exposure may cause skin dryness or cracking.
R:20/21/22 Harmful by inhalation, in contact with skin and if swallowed.
R:39/23/24/25 Toxic: danger of very serious irreversible effects through inhalation, in contact with skin and if swallowed.

S PHRASES

S:45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
S:36 Wear suitable protective clothing.
S:37 Wear suitable gloves.

Thimerosal

20X Wash Buffer II - Xn Harmful

R PHRASES

R:33 Danger of cumulative effects.
R:20/21/22 Harmful by inhalation, in contact with skin and if swallowed.

S PHRASES

S:36 Wear suitable protective clothing.
S:37 Wear suitable gloves.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use.

REAGENT PREPARATION

1. Bring the entire kit including microwell pouch to room temperature before use.
2. Prepare enough 1X Wash Buffer for use by diluting 1 part 20X Wash Buffer with 19 parts purified water. The 1X Wash Buffer can be stored at room temperature for up to one month.

SPECIMEN PREPARATION

1. **Cerebrospinal Fluid**
 - a. Collect specimen aseptically following accepted procedures. (**See PRECAUTION #2**).
 - b. Centrifuge at 1000 xg for 15 minutes to ensure the removal of all white cells and particulate matter.
 - c. Carefully aspirate the CSF into a sterile container and seal.
 - d. Specimen may be tested immediately, stored up to 3 days at 2-8 C, or preserved by freezing. Do not store in a frost-free freezer.
2. **Serum**
 - a. Collect whole blood aseptically according to in-house procedures. (**See PRECAUTION #2**). The specimen must not contain anticoagulants as this may invalidate the test.
 - b. Permit blood to clot for 10 minutes or more at room temperature in a collection tube.
 - c. Centrifuge at 1000 xg for 10 minutes.
 - d. Carefully aspirate the serum into a sterile container and seal.
 - e. Specimen may be tested immediately, stored up to 3 days at 2-8 C, or preserved by freezing. Do not store in a frost-free freezer.

NOTE: Pretreatment of specimens is not recommended for the Premier Cryptococcal Antigen EIA. Studies indicate that pronase or heat treatment of specimens (56 C, 15 minutes) does not appear to adversely effect screening assay results. Occasionally, pronase or heat treatment alters EIA Titers. Therefore, if pronase treatment is used, it is essential that pronase be used for **all** serial specimens from an individual patient which are utilized to monitor therapy.

Specimens to be used for future comparisons should be aliquotted and frozen, as storage longer than 72 hours at 2-8 C may affect titers.

SCREENING ASSAY PROCEDURE

This test should be performed by qualified personnel per local regulatory requirements.

The screening assay procedure can be performed directly with serum or CSF prepared as above. The screening assay is intended solely for differentiation of positive and negative specimens. Values obtained with the screening assay cannot be used for calculation of EIA Titers. For accurate determination, the semi-quantitative assay should be performed.

- Snap off a sufficient number of microwells for samples and two controls. Insert microwells into the holder and record sample positions. **(See PRECAUTION #9).**
- Add 50 µL of each specimen, Sample Diluent (negative control), and 1 drop of Positive Control to the bottom of separate microwells. Mix by gently shaking. **Caution should be taken not to cross contaminate samples, as this could cause erroneous results.**
- Incubate at room temperature (20-30 C) for 10 minutes.
- Remove samples from each microwell by aspiration or by inverting the microwells and firmly tapping on an absorbent pad. Fill each microwell with 1X Wash Buffer (approximately 300 µL). Invert the microwells and firmly tap on an absorbent pad. Repeat this wash procedure three times for a total of four washes. **(See PRECAUTION #11).**
- Add 1 drop of Enzyme Conjugate to each microwell. Mix by gently shaking. Incubate at room temperature (20-30 C) for 10 minutes.
- Remove solution from each microwell by aspiration or by inverting the microwells and firmly tapping on an absorbent pad. Fill each microwell with 1X Wash Buffer (approximately 300 µL). Invert the microwells and firmly tap on an absorbent pad. Repeat this wash procedure three times for a total of four washes. **(See PRECAUTION #11).**
- Add 2 drops of Substrate Solution to each microwell. Mix by gently shaking 15-20 seconds.
- Incubate at room temperature (20-30 C) for 10 minutes.
- Add 2 drops of Stop Solution to each microwell. Mix by gently shaking and wait 2 minutes before reading.
- Observe Reactions:
 - Visual determination – Read within 15 minutes after adding Stop Solution.
 - Spectrophotometric Determination – Zero EIA reader on air. Wipe the underside of microwells with a lint free tissue. Read absorbance at 450 nm or 450/630 nm within 15 minutes after adding Stop Solution.
- Disinfect and retain microwell holder. Discard the used assay materials as biohazardous waste.

QUALITY CONTROL

At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing, or leakage.

The positive and negative controls must be used with each batch of specimens to provide quality assurance of the reagents.

The negative control (Sample Diluent) absorbance should be > 0.000 and < 0.100 at 450 nm or > 0.000 and < 0.070 at 450/630 nm. The Positive Control absorbance should be > 0.700 at either 450 nm or 450/630 nm.

If absorbance values do not correspond to visual readings, wipe the underside of the microwell, check the alignment of the microwells within the plate holder as well as the alignment of the plate holder, and reread.

It is suggested that the results of each quality control check be recorded in an appropriate log book to facilitate monitoring testing procedures and to comply with regulatory agencies.

If the expected values are not observed, please contact Meridian Technical Services Department at 1-800-343-3858.

INTERPRETATION OF SCREENING ASSAY

- Visual Reading
 - Negative = Colorless
 - Positive = Definite yellow color
 - Use an EIA reader if there is uncertainty in determining a definite yellow color.
- Spectrophotometric Single Wavelength (450 nm)
 - Negative = $OD_{450} < 0.100$
 - Indeterminant = $OD_{450} \geq 0.100$ and < 0.150
 - Positive = $OD_{450} \geq 0.150$
- Spectrophotometric Dual Wavelength (450/630 nm)
 - Negative = $OD_{450/630} < 0.070$
 - Indeterminant = $OD_{450/630} \geq 0.070$ and < 0.100
 - Positive = $OD_{450/630} \geq 0.100$

Indeterminant results should be repeated. If the repeated results are still indeterminate, a second specimen should be obtained and tested. Extremely strong reactions may yield a purple precipitate within a few minutes of stopping; these reactions should be reported as positive.

SEMI-QUANTITATIVE ASSAY PROCEDURE

This test should be performed by qualified personnel per local regulatory requirements.

The semi-quantitative assay procedure can be performed on any positive CSF or serum specimen. This procedure is intended to yield values necessary for calculation of an accurate EIA Titer.

If drug therapy is to be monitored using sequential specimens, the preferred method is to repeat each specimen at the same time as the new specimen is titered. It is not recommended that EIA titers from different assay runs be compared directly to monitor therapy. To prevent multiple freeze-thawing of specimens, several aliquots of each specimen should be frozen.

A. Dilutions

The following dilution series is recommended when performing the semi-quantitative procedure. Combine 100 µL of specimen and 100 µL of Sample Diluent in a test tube, labeled #1, and vortex 1-3 seconds. From this mixture prepare a series of five-fold dilutions as follows:

- Place 200 µL of Sample Diluent in each of 4 test tubes labeled 2-5.

- Using a clean pipet, transfer 50 µL of the diluted patient specimen from tube #1 into tube #2 and mix well.
- Transfer 50 µL from tube #2 into tube #3 and mix well. Continue this dilution procedure through tube #5. Further 1:5 dilutions may be made, if necessary. Final dilutions are: 1:2, 1:10, 1:50, 1:250, 1:1250

B. Assay Procedure

For increased precision, all reagents must be delivered using a pipetter. Use the dropper bottle to deliver the required number of drops (see below) into a clean test tube. Transfer the required volumes of reagents to the microwells using a pipetter.

NOTE: Do not remove dropper tips from reagent bottles; contamination may compromise activity.

Number of Drops Required

	Enzyme Conjugate	Substrate	Stop Solution
1 Titer (7 wells)	7	22	12
2 Titers (12 wells)	11	35	19
3 Titers (17 wells)	15	49	27
4 Titers (22 wells)	19	62	30
10 titers (52 wells)	42	143	79

- Snap off a sufficient number of microwells for diluted samples and two controls. Insert into the microwell holder and record sample positions. **(See PRECAUTION #9).**
- Pipet 50 µL of each diluted sample, Positive Control and Sample Diluent (negative control) into the bottom of separate microwells. **Caution should be taken not to cross contaminate samples, as this could cause erroneous results.**
- Incubate at room temperature (20-30 C) for 10 minutes.
- Remove solution from each microwell by aspiration or by inverting microwells and firmly tapping on an absorbent pad. Fill each microwell with 1X Wash Buffer (approximately 300 µL). Invert microwells and firmly tap on an absorbent pad. Repeat this wash procedure three times for a total of four washes. **(See PRECAUTION #11).**
- Pipet 50 µL of Enzyme Conjugate into each microwell, including control microwells. Mix by gently shaking. Incubate at room temperature (20-30 C) for 10 minutes.
- Remove solution from each microwell by aspiration or by inverting the microwells and firmly tapping on an absorbent pad. Fill each microwell with 1X Wash Buffer (approximately 300 µL). Invert the microwells and firmly tap on an absorbent pad. Repeat this wash procedure three times for a total of four washes. **(See PRECAUTION #11).**
- Pipet 100 µL of Substrate Solution into each microwell. Mix by gently shaking 15-30 seconds.
- Incubate at room temperature (20-30 C) for 10 minutes.
- Pipet 100 µL of Stop Solution into each microwell. Mix by gently shaking and wait 2 minutes before reading.
- Spectrophotometric Determination – Zero EIA reader on air. Wipe the underside of wells with a lint free tissue. Read absorbance at 450 nm or 450/630 nm within 15 minutes after adding Stop Solution.
- Disinfect and retain microwell holder. Discard the used assay materials as biohazardous waste.

INTERPRETATION OF SEMI-QUANTITATIVE ASSAY

NOTE: Quality Control procedures listed after the Screening Assay should be followed when performing the semi-quantitative assay.

When using a single wavelength reader to obtain results for the semi-quantitative assay, the absorbance of the negative control (Sample Diluent) must be subtracted from all other absorbance values before calculation of the EIA Titer.

Dual wavelength absorbance values need not be adjusted.

A. Selections of Absorbance Values for Calculations

- Use only absorbance values which are from the semi-quantitative assay to calculate the EIA titer.
- Select the highest absorbance value within the acceptable range of 0.100 and 1.500 to calculate the EIA titer. If the next greater dilution also yields an absorbance within the acceptable range, this value should also be selected.
- If the 1:2 dilution of the specimen yields an absorbance < 0.100, repeat the assay on the undiluted specimen.
- Using the selected absorbance, calculate the EIA titer as outlined in B below. If two absorbances are selected, calculate each titer separately and average the results.

NOTE: Appearance of a purple precipitate does not invalidate a positive result, but readings from such wells cannot be used for calculation of EIA Titer.

B. Calculation

- Calculate the EIA titer for a particular specimen with the following formula:

$$\text{Absorbance Value} \times \text{Multiplication Factor} = \text{EIA titer}$$

Tube/Well Number	Multiplication Factors*				
	1	2	3	4	5
Dilution	1:2	1:10	1:50	1:250	1:1250
Multiplication Factor	20	100	500	2,500	12,500

*Use a Multiplication Factor of 10 to calculate an EIA titer for undiluted specimens.

2. Sample Calculations:
For example, if a patient specimen had an absorbance of 1.2 at a 1:50 dilution, the EIA Titer would be: $1.2 \times 500 = 600$ which is reported as 1:600


C. Reporting EIA Titers

Numbers derived from the EIA calculations may be reported to physicians as the EIA Titer. A higher titer reflects a greater concentration of antigen. However, do not attempt to compare EIA Titer with latex titer. Due to serotype variations and differing test technologies, EIA Titers are not numerically equivalent to latex titers. Decreases in titers of less than four-fold are not considered clinically significant.⁶

EXPECTED VALUES

Cryptococcal antigen in CSF or serum of the untreated patient indicates an active infection.¹⁻³ The Premier Cryptococcal Antigen EIA should have both diagnostic and prognostic value since antigen concentrations rise during disease progression and fall during successful treatment.^{1,2,5,6} Failure of titers to decline may indicate inadequate therapy. However, in some treated patients, antigen levels may remain stationary for extended periods of time when no viable organisms can be demonstrated.⁸

Comparison of EIA and Latex Titers on an Individual Patient

The primary use for semi-quantitative values derived with the latex or the EIA is to follow the course of disease and treatment within an individual patient. These data suggest that the Premier Cryptococcal Antigen EIA will be a valuable method for monitoring the efficacy of drug therapies.

AIDS patients may have very high EIA Titers. The dilution series recommended in the procedure may not be adequate to calculate an EIA Titer. Further 1:5 dilutions may be made as needed.

LIMITATIONS OF THE PROCEDURE

The Premier Cryptococcal Antigen EIA is intended for use with CSF or serum specimens. Use of this assay is not recommended with other specimens.

A negative result does not preclude diagnosis of cryptococcosis, particularly if only a single specimen has been tested and the patient shows symptoms consistent with cryptococcosis.

A positive result implies the presence of cryptococcal antigen; however, all test results should be reviewed in light of other clinical data by the physician. *Trichosporon beigeli* yielded a false positive in a single case with latex⁹ and a similar cross-reaction is possible with the EIA.

SPECIFIC PERFORMANCE CHARACTERISTICS

The Premier Cryptococcal Antigen EIA was evaluated at two major clinical laboratories (279 specimens) in the United States and at Meridian Bioscience, Inc., (88 specimens, including one indeterminate specimen. Insufficient quantity precluded retesting of the indeterminate, and this specimen is not listed below). The CALAS® Latex Agglutination System was used as the reference assay. The combined test results are summarized in the table below.

		Latex	
		+	-
Premier EIA	+	108	6*
	-	0	252

Sensitivity 100%
Specificity 98%
Agreement 98%

*In four of six discrepancies listed above, previous or subsequent specimens confirmed that the patients had cryptococcosis.

The EIA appeared to be more sensitive by detecting low concentrations of cryptococcal antigen at early or late stages of disease. Limits of detection for serotypes A and B have been demonstrated to be 0.63 ng/mL for EIA and 7.6 ng/mL for the latex assay.¹⁰

Patients with AIDS have emerged as a major group susceptible to cryptococcosis. Cryptococcosis is the fourth most common infection complicating AIDS.² A comparison of clinical diagnosis and Premier Cryptococcal Antigen results in an AIDS patient population is listed below.

		Clinical Diagnosis Cryptococcosis	
		+	-
Premier EIA	+	44	2
	-	0	25

ASSAY SPECIFICITY

The following list represents organisms that were negative when tested in the EIA:

Microorganisms tested for Cross-reactivity

Bacteria	Fungal Antigen Preparation	Virus	Yeast
<i>Escherichia coli</i> K1	<i>Aspergillus flavus</i>	Coxsackie A-16, B-5	<i>Candida albicans</i>
<i>Haemophilus influenzae</i>	<i>Aspergillus fumigatus</i>	Echovirus II	<i>Rhodotorula rubra</i>
<i>Neisseria meningitidis</i> 301, 362	<i>Aspergillus niger</i>	Poliovirus I, II, III	<i>Torulopsis glabrata</i>
<i>Pseudomonas aeruginosa</i>	<i>Blastomyces dermatitidis</i>		
<i>Staphylococcus aureus</i>	<i>Coccidioides immitis</i>		
<i>Staphylococcus aureus</i> Cowan strain 1	<i>Histoplasma capsulatum</i>		
<i>Staphylococcus epidermidis</i>			
<i>Streptococcus agalactiae</i>			
<i>Streptococcus pneumoniae</i>			

The only organism found to cross-react in the EIA was *T. beigeli*. This organism is not likely to be found in typical specimens being tested for cryptococcal antigen.⁹

ASSAY PRECISION

Intra-assay Variability – 16 replicates of 3 known positive sera were tested in one assay to determine the intra-assay reproducibility.

Sample	Mean Absorbance	% C.V.
1	0.67	2.1
2	1.74	2.2
3	2.70	3.0

Coefficients of Variation for the same 3 samples using the screening assay with dropper-tips (16 replicates each) were 12, 3.5, and 3.0% respectively.

Inter-assay Variability – 16 replicates of 3 known positive sera were tested in multiple assays by more than one operator.

Sample	Mean Absorbance	% C.V.
4	0.27	11.8
5	0.73	11.5
6	1.76	7.2

Coefficients of Variation for the same 3 samples using the screening assay with dropper-tips (16 replicates each) were 30, 5.8, and 18.9% respectively.

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BIBLIOGRAPHY

- Perfect, J.R. 1989. Cryptococcosis. Infectious Disease Clinics of North America. 3:77-102.
- Patterson, T.F. and V.T. Andriole. 1989. Current Concepts in Cryptococcosis. Eur. J. Clin. Microbiol. Infect. Dis. 8:457-465.
- Rinaldi, M.G. 1982. Cryptococcosis. Lab. Medicine 13:11-19.
- Murphy, J.M. 1989. Cryptococcosis, in *Immunology of the Fungal Diseases*. Chapter 4, pp. 98-138. CRC Press., Boca Raton, FL.
- Bennett, J.E., H.F. Hasenclever and B.S. Tynes. 1964. Detection of cryptococcal polysaccharide in serum and spinal fluid: value in diagnosis and prognosis. Trans. Assoc. Am. Physicians 77:145-150.
- Diamond, D. and E. Bennett. 1974. Prognostic factors in cryptococcal meningitis. Ann. Int. med. 80:176-181.
- Goodman, J.S., L. Kaufman and M.G. Koening. 1971. Diagnosis of Cryptococcal meningitis: Value of immunologic detection of cryptococcal antigen. New Eng. J. of Med. 285:434-436.
- Kaufman, L. and S. Blumer. 1977. Cryptococcosis: The Awakening Giant. Proc. of the Fourth International Conference on the Mycoses. Pan American Health Organization Scientific Publication #356. pp. 176-182.
- McManus, E.J. and J.M. Jones. 1985. Detection of a *Trichosporon beigelii* antigen cross-reactive with *Cryptococcus neoformans* capsular polysaccharide in serum from a patient with disseminated *Trichosporon* infection. J. Clin. Microbiol. 21:681-685.
- Data on file (presented at ICAAC, October 1990, reprints available. Manuscript to be submitted to *Journal of Clinical Microbiology*, 1991).

SN11125

Rev. 12/07

Manufactured By

Meridian Bioscience, Inc.
USA/Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244
Telephone: (513) 271-3700
Orders/Customer Service:
(800) 343-1980
Technical Support Center:
(800) 343-3858
Information Fax:(513) 272-5432
Ordering Fax:(513) 271-0124

EC REP

Meridian Bioscience Europe
Via dell'Industria, 7
20020 Villa Cortese (MI)
Italy
Authorized Representative
Tel.: +39 0331 433636
Fax: +39 0331 433616
e-mail: info@mdeur.com

Meridian Bioscience Europe France
Le Quadra
455, Promenade des Anglais
06299 Nice Cedex-3
France
Tel.: +33 (0) 49 3187210
Fax: +33 (0) 49 3187211
e-mail: info@meridianbioscience.fr

Meridian Bioscience Europe s.a. / n.v.
Rue de l'Industrie 7
1400 Nivelles
Belgium
Tel.: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
e-mail: info@mdeur.be

Meridian Bioscience Europe b.v.
Halderheiweg 6
5282 SN Boxtel
The Netherlands
Tel.: +31 (411) 621166
Fax: +31 (411) 624841
e-mail: meridian.info@planet.nl

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:
Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	SUBS	Substrate / Substrato / Substrat / Substrato / Substrat
	Batch Code / Codice del lotto / Code du lot / Codigo de lote / Chargenbezeichnung	SOLN STOP	Stopping Solution / Soluzione di Stop / Solución d'arrêt / Solution de parada / Stoplösung
	In Vitro Diagnostic Medical Device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugue enzymatique / Conjugado enzimático / Enzymkonjugat
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
		CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
		CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
	Catalogue number / Numero de catalogo / Référence du catalogue / Número de catalogo / Bestellnummer	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campioni / Diluente echantillons / Diluyente de muestra / Probenverdünnungspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller	EC REP	Authorized Representative in the European Community / Mandatario nella Comunità Europea / Mandatario dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	Contains sufficient for n tests / Contenuto sufficiente per n assg / Contenu suffisant pour n tests / Contenido suficiente para n ensayos / Inhalt ausreichend für n Prüfungen		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	REAG	Reagent / Reagent / Réactifs / Reactivos / Reagenzien

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.