# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation

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
[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Premier™ Adenoclone®-Type 40/41

EIA for the Detection of Adenovirus Serotypes 40 and 41 in Human Fecal Samples

**REF** Catalog No. 696006

**IVD** For *in vitro* diagnostic use only

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## Premier Adenoclone®-Type 40/41

EIA for the detection of Adenovirus Serotypes 40 and 41 in Human Fecal Samples

Catalog No. 696006

### INTENDED USE

Premier Adenoclone®-Type 40/41 is an Enzyme Immuno-assay (EIA) for the qualitative detection of enteric adenovirus serotypes 40 and 41 in human fecal specimens.

### EXPLANATION

Human adenoviruses consist of 41 known types, classified primarily by serology and DNA electrophoretotyping. The viruses are nonenveloped, DNA-containing icosahedrons composed of 252 capsomers. The virion contains a number of proteins associated with genus and type specificity. Those associated with genus (group) specificity are the alpha portion of the hexon, the penton base, the delta portion of the fiber and polypeptide IX. Those associated with species (type) specificity are the epsilon portion of the hexon and the gamma portion of the fiber. Those associated with both genus and type specificity are the major core protein and polypeptide IX (11). The monoclonal antibodies in Premier Adenoclone®-Type 40/41 specific EIA are directed against the type-specific portion of the hexon of adenovirus types 40 and 41. The type-specific monoclonal antibodies in the EIA will not react with any of the other adenovirus serotypes (1).

The enteric adenovirus types 40 and 41 are the only adenoviruses which have been consistently associated with gastroenteritis in infants and young children, and may be the second major cause of gastroenteritis after rotaviruses (2-4, 8, 13-17). Until recently, these viruses could be detected only by electron microscopy (EM). These adenovirus types can now be detected by cultivation in Graham 293 cells, a human embryonic kidney cell line transformed by adenovirus type 5 DNA, but cannot readily be cultivated in all lines traditionally used for adenovirus isolation, e.g. HEP-2 cells (5). Identification of adenovirus in tissue culture isolation is followed by the appearance of characteristic cytopathic effect (CPE). Generally CPE occurs within 3 to 7 days but may require up to 28 days (7,11). Once the virus is isolated, it must be further characterized by analysis of restriction endonuclease digests of the virus DNA followed by gel electrophoresis (6).

The Premier Adenoclone®-Type 40/41 test kit for direct stool specimens is a rapid and specific test for adenovirus 40 and 41 that eliminates the need for tissue culture and DNA electrophoretotyping.

### BIOLOGICAL PRINCIPLES

Premier Adenoclone®-Type 40/41 utilizes monoclonal antibodies in a solid phase sandwich type assay. Breakaway microwells are coated with a monoclonal antibody directed against the group specific antigen for all known human adenoviruses. An aliquot of fecal suspension is added to the microwell and incubated simultaneously with anti-adenovirus types 40 and 41 monoclonal antibodies, conjugated to horseradish peroxidase, resulting in the adenovirus antigen being sandwiched between the solid phase and enzyme conjugate. After a 60 minute incubation at room temperature, the sample well is washed with distilled or deionized water to remove unbound specimen and enzyme labelled antibodies. Enzyme substrate A (urea peroxide) and substrate B (tetramethylbenzidine) are added to the wells and incubated for 10 minutes at room temperature. Any bound enzyme-conjugate in the wells converts the colorless substrate to a blue color.

### MATERIALS PROVIDED

Premier Adenoclone®-Type 40/41 contains 48 tests.

1. 48 Microwells — Coated with Group Specific Anti-Adenovirus Murine Monoclonal Antibody.
2. 50 ml Sample Diluent / Negative Control — Buffered saline with 0.02% thimerosal as preservative.
3. 2.0 ml Positive Control — Inactivated Adenovirus (AD41) in buffer with 0.02% thimerosal as preservative.
4. 6 ml Enzyme Conjugate — Horseradish Peroxidase conjugated to anti-adenovirus types 40 and 41 murine monoclonal antibodies in a buffered protein solution with gentamicin and 0.02% thimerosal as preservatives.
5. 6 ml Part A Substrate — Substrate buffer contains urea peroxide.
6. 6 ml Part B Substrate — Substrate buffer contains tetramethylbenzidine.
7. 6 ml Stop Solution — Contains 1 N sulfuric acid. **Caution:** Avoid contact with skin. If skin contact occurs, flush with water.
8. 1 Microwell Holder
9. 48 Sample Transfer Pipettes

### MATERIALS NOT PROVIDED

1. 12 x 75 mm test tubes, test tube rack.
2. Deionized or distilled water.
3. Absorbent paper.
4. Precision micropipette tips to deliver 100  $\mu$ l and 1000  $\mu$ l (optional).
5. Waste container with a 1:10 dilution of household bleach. For autoclaving, use an iodophor disinfectant.
6. Microwell plate reader capable of reading absorbance at 450 nm (optional).
7. Device for dispensing wash solution such as multi-channel pipette, syringe with manifold, wash bottle, etc.
8. Timer (minimum 1 hour).

## PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Do not mouth pipette samples or reagents. Avoid contact with broken skin or mucous membranes.
3. Do not smoke, eat or drink in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling samples and wash hands after assay is complete. Patient specimens, assay controls and all materials coming into contact with them should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories," 1988.
5. Avoid skin contact with stop solution (1N sulfuric acid). It may cause irritation and burns. Flush with water if contact occurs.
6. Dispose of all materials used to perform the test by autoclaving at 121°C for at least one hour. Liquid waste may be disposed by mixing with a 1:10 dilution of household bleach for a minimum of 30 minutes. **CAUTION:** Liquid waste containing stop solution must be neutralized before addition of household bleach.
7. Avoid splashing or generation of aerosols.
8. Do not use **Premier Adenoclone®-Type 40/41** reagents beyond the kit expiration date. Each reagent has been optimized for maximum performance. Dilution or adulteration of these reagents may result in a loss of sensitivity. Incubation time and temperatures other than those specified may give erroneous results.
9. Do not interchange or mix different lots of **Premier Adenoclone®-Type 40/41** reagents.
10. Avoid microbial contamination of reagents or incorrect results may occur. Contamination of samples could cause erroneous results.
11. Use separate pipettes or pipette tips for each sample, control and reagent.
12. DO NOT REUSE MICROWELLS.

## RISK AND SAFETY PHRASES

**PART B SUBSTRATE:** TOXIC — METHANOL

RISK PHRASES:

23/25 Toxic by inhalation and if swallowed.

SAFETY PHRASES:

- 45 In case of accident or if you feel unwell, seek medical advice immediately (show label where possible).

## STORAGE INSTRUCTIONS

Store kit reagents at 2°-8°C. Bring kit reagents to room temperature (20°-30°C) before use and promptly return to 2°-8°C after use. Keep microwells in pouch until the pouch reaches room temperature. Return all unused microwells to their original foil pouch. The microwells are stable for the shelf life of the kit after opening.

## INDICATION OF INSTABILITY OR DETERIORATION

The following conditions may indicate reagent deterioration:

1. Any evidence of microbial contamination or heavy precipitation.
2. Any blue color in the substrate solutions before addition to microwells.
3. A negative control value greater than 0.150 absorbance units at 450 nm may indicate deterioration of reagents.
4. A positive control value of less than 0.3 absorbance units may indicate deterioration of reagents.

If any of the above conditions are observed, contact our Technical Support Center at 1-800-343-3858.

## SPECIMEN HANDLING

Stool specimens or rectal swabs should be collected in containers that do not contain media, preservatives, animal serum or detergent as any of these additives may interfere with the **Premier Adenoclone®-Type 40/41** assay. Caution should be taken when using rectal swabs to insure that sufficient fecal sample (30 to 50 mg of raw stool) is obtained.

Diluted samples may be stored at 2°-8°C for 3 days without interfering with the assay performance. For long term storage of undiluted specimens, -20°C or colder is recommended. Repeated freezing and thawing of samples is not recommended and may cause erroneous results. Do not store in self-defrosting freezers.

## REAGENT PREPARATION

1. Bring all reagents to room temperature (20°-30°C) before use.
2. Return all reagents to 2°-8°C immediately after use.
3. Do not allow microwells to dry between steps.
4. Prepare decontaminating vessel for discarding reagents and materials.
5. Reproducibility in any EIA assay is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the EIA test procedure.

## SPECIMEN PREPARATION


Add 1 ml of sample diluent to properly marked tube, using a transfer pipette or precision pipette. Add samples by one of the following:

1. Solid stool — press sample into transfer pipette to first mark.
2. Liquid stool — aspirate sample into transfer pipette to first mark.
3. Rectal swab — swirl swab in one ml of sample diluent to release fecal material. Firmly press swab against side of tube to remove liquid. Mix thoroughly. *Do not dilute further.*

Use the sample transfer pipette to add two drops (100 µl) of diluted stool specimen to each well.

## Sample Transfer Pipette

## PROCEDURE

*This test should be performed by qualified personnel per local regulatory requirements.*

1. Snap off a sufficient number of wells for samples and the controls and insert into the microwell holder. Record sample position.
2. Add 2 drops (100  $\mu$ l) of diluted fecal sample, negative control (sample diluent) or positive control to the bottom of separate microwells.
3. Add 2 drops (100  $\mu$ l) of enzyme conjugate to each microwell. Mix by gently swirling on table top.
4. Incubate at room temperature for 60  $\pm$  5 minutes.
5. Pour the liquid out of the wells into a discard vessel. Tap the microwell holder upside down vigorously against absorbent paper to ensure complete removal of liquid from the wells.
6. Fill all the wells to overflowing with deionized water and pour the liquid out as in Step 5.
7. Repeat the washing procedure (Steps 5 & 6) four more times for a total of 5 washes.
8. Add 2 drops (100  $\mu$ l) of Substrate A (urea peroxide) solution, to each microwell.
9. Add 2 drops (100  $\mu$ l) of Substrate B (TMB) solution to each microwell.
10. Incubate at room temperature for 10 minutes.
11. Discard the used wells and absorbent paper into biohazardous waste after reading and recording results.

Visual determinations can be made after incubating for 10 minutes in step 10, and should be read within 10 minutes. Samples with blue color greater than the negative control are positive. Samples showing equal or less color than negative control are negative.

### Optional - Spectrophotometric Procedure

Spectrophotometric determinations can be made by adding 2 drops (100  $\mu$ l) of stop solution (1 N sulfuric acid) to each microwell after the incubation in step 10. Read the absorbance of each microwell at 450 nm using a >600 nm reference filter (optional) against an air blank within 60 minutes.

## QUALITY CONTROL

**Premier Adenoclone<sup>®</sup>-Type 40/41** is supplied with a positive and negative (sample diluent) control which are to be used every time an assay is performed. These reagents are used to insure that the kit reagents are working properly. The positive control well should be clearly blue with an absorbance at 450 nm of greater than 0.30. The negative control must have no visible color, and an absorbance at 450 nm of less than 0.150. If the positive and negative controls do not perform properly, contact our Technical Support Center at 1-800-343-3858.

## INTERPRETATION OF RESULTS

1. Visual Determination:  
Any sample with blue color more intense than that of the negative control is considered positive. Any sample with color equal to or less intense than the negative control is considered negative.
2. Spectrophotometric Determination:  
Specimens with absorbance units ( $A_{450}$ ) equal to or greater than 0.150 are considered positive. Specimens with absorbance less than 0.150 are considered negative. Occasionally, a discrepant result may occur between visual and spectrophotometric determinations with samples containing low amounts of antigen. Spectrophotometric determination, being an objective method, is slightly more accurate.

NOTE: A precipitate may form in high positive samples. This will not affect the results.

## LIMITATIONS OF THE PROCEDURE

**Premier Adenoclone<sup>®</sup>-Type 40/41** is highly specific and sensitive for adenovirus types 40 and 41 antigen. Failure to detect adenovirus 40/41 may be a result of factors such as collection of specimen at an improper time in the disease, improper sampling or handling of the specimen.

A negative result for a stool sample may not preclude the presence of an infection with non-enteric adenovirus in other body fluids. If respiratory or ophthalmic infection is suspected specimens from the site of infection should be cultured and confirmed with Meridian Bioscience's **Premier Adenoclone<sup>®</sup>** (Catalog #696007).

All positive results must be interpreted with caution since adenovirus is capable of latency and recrudescence (12). Asymptomatic shedding may occur for months after infection (13). Enteric adenoviruses may be found in the stools of asymptomatic children (13). Test results should be interpreted in conjunction with information available from epidemiological studies or clinical evaluation of patient or other diagnostic procedures (12).

While the relationship between adenovirus 40/41 and gastroenteritis is well established, co-infection with other viral and bacterial pathogens (including non-enteric adeno-viruses) is possible. Therefore, bacteriological tests should be performed in parallel with this test to rule out bacteriological etiology. Meridian Bioscience's **Premier Rotaclone<sup>®</sup>** rotavirus diagnostic kit (Catalog #696004) may be performed in parallel to rule out rotaviral gastroenteritis.

False positives could occur with high levels of *S. aureus* possessing protein A such as Cowan strain. *Staphylococcal enterocolitis* is an uncommon disease generally described in adult patients and is extremely rare in infants and children (9,10).

## EXPECTED VALUES

Rate of positivity may vary depending on geographic location, method of specimen collection, handling and transportation, test employed and general health environment of patient population under study.

The frequency of adenovirus infection will vary with the clinical syndrome and the age of the individual. It has been reported that approximately 5 to 10% of viral gastroenteritis in infants and children under 2 years is caused by adenovirus types 40 and 41 (2,4,6,8,13).

Performance characteristics have not yet been fully determined for direct testing of neonatal stool samples. Lack of false positive results should be determined in each individual laboratory before neonatal stools are tested directly.

## PERFORMANCE CHARACTERISTICS

Premier Adenoclone®-Type 40/41 for enteric adenovirus antigen test was tested in collaboration with institutions in Arizona, Massachusetts, Ontario and Manitoba. The results were compared with initial Electron Microscopy followed by DNA restriction endonuclease analysis of adenovirus positive samples. All of the samples were tested by Premier Adenoclone®-Type 40/41 EIA and scored visually then read spectrophotometrically. There was 100% correlation between visual and spectrophotometric determinations.

		EM/DNA Restriction Analysis	
		+	-
Premier Adenoclone®	+	58	3
-Type 40/41 EIA	-	1	169
Sensitivity		98%	
Specificity		98%	
Total Agreement		98%	

## LIMITS OF DETECTION

Virus particle counts from known adenovirus 40 and 41 stool samples were determined by electron microscopy. Serial dilutions were tested by Premier Adenoclone®-Type 40/41. These results show that adenovirus particle counts as low as 10<sup>4</sup> per ml can be detected by Premier Adenoclone®-Type 40/41.

Virus Particle Count	Ad 40 (A <sub>450</sub> )	Ad 41 (A <sub>450</sub> )
2 x 10 <sup>6</sup>	>2.0	—
1 x 10 <sup>6</sup>	>2.0	>2.0
5 x 10 <sup>5</sup>	>2.0	>2.0
2.5 x 10 <sup>5</sup>	>2.0	>2.0
1.25 x 10 <sup>5</sup>	1.510	>2.0
6.25 x 10 <sup>4</sup>	1.001	>2.0
3.1 x 10 <sup>4</sup>	0.595	1.445
1.6 x 10 <sup>4</sup>	0.307	0.781
8 x 10 <sup>3</sup>	0.167	0.425
4 x 10 <sup>3</sup>	0.108	0.219
2 x 10 <sup>3</sup>	—	0.142

## CROSS REACTIVITY

1. Non-enteric Adenovirus Subgroups: Adenoviral 3 + CPE culture supernatants were tested by Premier Adenoclone®-Type 40/41. The data shows no cross reactivity with the Premier Adenoclone®-Type 40/41.

PLEASE NOTE: CPE culture supernates are not approved samples for use with this product. The information derived from tests with culture filtrates that is cited here is provided for information purposes only. Do not use culture supernates to verify specificity or sensitivity.

Adenoviruses		Premier Adenoclone® -Type 40/41	Premier Adenoclone® Group specific*
Subgroup	Type	(A <sub>450</sub> )	(A <sub>450</sub> )
A	12	0.043	>3.0
A	18	0.042	>3.0
A	31	0.037	>3.0
B	3	0.050	>3.0
B	7	0.061	>3.0
B	21	0.027	>3.0
C	1	0.042	>3.0
C	2	0.041	>3.0
C	6	0.039	>3.0
D	22	0.030	>3.0
D	37	0.033	>3.0
D	39	0.035	>3.0
E	4	0.047	>3.0
Negative control		0.049	0.060
Positive control		1.071	0.580

\*The Adenovirus subgroups were tested by Meridian Bioscience's Premier Adenoclone® (Catalog #696007) group specific EIA, which detects all 41 types of human adenovirus to demonstrate the presence of adenovirus.

II. Microorganisms and Viruses: Common intestinal pathogens and other organisms occasionally present in feces were tested in Premier Adenoclone® -Type 40/41 and showed no cross reactivity in the test.

**Microorganisms**

<i>Acholeplasma laidlawii</i>	<i>Enterobacter aerogenes</i>	<i>Mycoplasma hominus</i>	<i>Pseudomonas aeruginosa</i>
<i>Campylobacter lariidis</i>	<i>Enterobacter cloacae</i>	<i>Mycoplasma hyorhinis</i>	<i>Salmonella dysenteria</i>
<i>Campylobacter fetus</i>	<i>Escherichia coli</i>	<i>Mycoplasma orale</i>	<i>Serratia marcescens</i>
<i>Candida albicans</i>	<i>Klebsiella pneumoniae</i>	<i>Mycoplasma pneumoniae</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i>	<i>Mycoplasma arginini</i>	<i>Proteus mirabilis</i>	<i>Staphylococcus aureus*</i> (Cowan)

\* False positives could occur with high levels of *S. aureus*, possessing protein A such as Cowan strain. *Staphylococcal enterocolitis* is an uncommon disease generally described in adult patients and is extremely rare in infants and children (9, 10).

**Viruses**

Astrovirus	Echovirus 8	Papovavirus SV-40	Respiratory Syncytial virus
Calicivirus	Echovirus 11	Parainfluenza-1	Rhinovirus 3
Chlamydia trachomatis (serotype A)	Echovirus 16	Parainfluenza-2	Rhinovirus 17
Coxsackievirus A-9	Echovirus 20	Parainfluenza-3	Rhinovirus 25
Coxsackievirus B-3	Echovirus 22	Poliovirus I	Rhinovirus 43
Coxsackievirus B-4	Human Coronavirus	Poliovirus II	Rhinovirus 62
Coxsackievirus B-5	Influenza B	Poliovirus III	Rotavirus
Cytomegalovirus	Measles	Reovirus 1	Simian foamyvirus (type 2 & 8)
Echovirus 4	Myxo/Paramyxoviruses SV-5	Reovirus 2	Small Round virus (SRV)
	Norwalk virus	Reovirus 3	Varicella zoster virus

A urine sample was tested in the Premier Adenoclone®-Type 40/41 and found to be negative.

**ASSAY PRECISION**

**Inter-assay variation**

The inter-assay precision was determined by diluting known adenovirus 40 and adenovirus 41 positive samples to obtain high, medium, and low absorbance values. These diluted samples were tested in 21 different assays on different days by several technicians and the mean and coefficient of variation determined.

Adenovirus Type 40		High	Medium	Low	Negative
	Mean	1.456	0.856	0.262	0.042
	S.D.	0.174	0.088	0.025	0.017
	CV	12.0%	10.3%	9.5%	40.5%

Adenovirus Type 41		High	Medium	Low
	Mean	0.892	0.565	0.193
	S.D.	0.072	0.068	0.019
	CV	8.1%	12.0%	9.8%

**Intra-assay variation**

The intra-assay precision was determined by diluting known adenovirus 40 and adenovirus 41 positive samples to obtain high, medium, and low absorbance values. These diluted samples were tested 21 times within a single assay and the mean and coefficient of variation determined.

Adenovirus Type 40		High	Medium	Low	Negative
	Mean	1.522	0.894	0.193	0.064
	S.D.	0.038	0.038	0.011	0.022
	CV	3.8%	4.3%	5.7%	34.3%

Adenovirus Type 41		High	Medium	Low
	Mean	0.924	0.647	0.156
	S.D.	0.086	0.039	0.012
	CV	9.3%	6.0%	7.7%

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## BIBLIOGRAPHY

1. Herrmann, J.E., Perron-Henry, D.M., Stobbs-Walro, D., and Blacklow, N.R. Preparation and characterization of monoclonal antibodies to enteric adenovirus types 40 and 41. *Arch. Virol.* 94:259-265, 1987.
2. Uhnno, I., Wadell, G., Svensson, L., and Johansson, M. Two new serotypes of adenovirus from human infantile diarrhea. *Dev. Biol. Stand.* 55:311-318, 1983.
3. Brandt, C.D., Kim, H.W., Yoiken, R.H., Kapikian, A.Z., Arrobio, J.O., Rodriguez, W.J., Wyatt, R.G., Chanock, R.M., and Parrott, R.H. Comparative epidemiology of two rotavirus serotypes and other viral agents associated with pediatric gastroenteritis. *Amer. J. Epidemiol.* 110: 243-254, 1979.
4. Uhnno, I., Wadell, G., Svensson, L., Johansson, M.E. Importance of enteric adenoviruses 40 and 41 in acute gastroenteritis in infants and young children. *J. Clin. Microbiol.* 20:365-372, 1984.
5. Takiff, H.F., Straus S.E., and Garon, C.F. Propagation and *in vitro* studies of previously non-cultivable enteric adenoviruses in 293 cells. *Lancet* ii:832-834, 1981.
6. Wadell, G. Molecular epidemiology of human adenoviruses. In: Doerfler W., ed. *The molecular biology of adenoviruses 2. 30 years of adenovirus research 1953-1983. Current topics in microbiology and immunology.* Berlin: Springer-Verlag, 110:191-220, 1984.
7. Herrmann, J.E., Personal communication, 1987, University of Massachusetts Medical School, Worcester, Massachusetts.
8. Brandt, C.D., Kim, H.W., Rodriguez, W.J., Arrobio, J.O., Jeffries, B.C., Stallings, E.P., Lewis, C., Miles, A.J., Chanock, R.M., Kapikian, A.K., Parrott, R.H. Pediatric viral gastroenteritis during eight years of study. *J. Clin. Microbiol.* 18:71-78, 1983.
9. Gutman, L.T., Idriss, Z.H., Gehlbach, S., Blackman, L. Neonatal staphylococcal enterocolitis; association with indwelling feeding catheters and *S. aureus* colonization. *J. Pediatr.*, 88:836-839, 1976.
10. Hay, P. and McKenkie, P. Side effects of oxytetracycline therapy. *Lancet* i:945-950, 1954.
11. Horwitz, M.S. Adenoviruses and Their Replication. In B.N. Fields ed. *Virology*, Raven Press, New York, pp. 433-438.
12. Horwitz, M.S. Adenoviruses and Their Replication. In B.N. Fields ed. *Virology*, Raven Press, New York, pp. 480.
13. Cukor, G., and Blacklow, N.R. Human viral gastroenteritis. *Microbiological Reviews.* 48:157-179, 1984.
14. Chiba, S., Wakamura, I., Urasawa, S., Nakata, S., Taniguchi, K., Fujinaga, K., Nakao, T. Outbreak of infantile gastroenteritis due to type 40 adenovirus. *Lancet* ii:954-957, 1983.
15. Johansson, M.E., Uhnno, I., Kidd, A.H., Madeley, C.R., Wadell, G. Direct identification of enteric adenovirus, a candidate new serotype, associated with infantile gastroenteritis. *J. Clin. Micro.* 12:95-100, 1980.
16. Wigand, R., Baumeister, H.G., Maass, G., Kuhn, J., Nammer, H.J. Isolation and identification of enteric adenoviruses. *J. Med. Virol.* 11:233-240, 1983.
17. de Jong, J.C., Wigand, R., Kidd, A.H., Wadell, G., Kapsenberg, J.G., Muzerie, C.J., Wermenbol, A.G., Firtzlaff, R.G. Candidate adenoviruses 40 and 41: fastidious adenoviruses from human infant stool. *J. Med. Virol.* 11:215-231, 1983.

**USA/Corporate Office**  
 Manufactured by 3471 River Hills Drive  
 Cincinnati, Ohio 45244  
 Telephone: (513) 271-3700  
 Orders/Customer Service: 1-800-543-1980  
 Technical Support Center: 1-800-343-3858  
 Information Fax: (513) 272-5432  
 Ordering Fax: (513) 271-0124

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product

	Manufactured by		For Performance Evaluation Only
	Authorized Representative		Temperature Limitations
	Catalog Number		Use By/Expiration Dating Information
	In vitro Diagnostics		Sufficient for "X" Tests
	Batch Code/Lot Information		CE Symbol

Authorized Representative  
**Meridian Bioscience Europe**  
 Via dell'Industria, 7  
 20020 Villa Cortese (MI)  
 Italy  
 Tel: +39 0331 433636  
 Fax: +39 0331 433616  
 e-mail: info@mdeur.com

**Meridian Bioscience Europe s.a./n.v.**  
 Rue de l'Industrie 7  
 1400 Nivelles  
 Belgium  
 Tel.: +32 (67) 895959  
 Fax: +32 (67) 895958  
 e-mail: info@mdeur.be

**Meridian Bioscience Europe France**  
 455, Promenade des Anglais  
 06299 Nice Cedex-3  
 France  
 Tel.: +33 (0) 493187210  
 Fax: +33 (0) 493187211  
 e-mail: mdeur.info@fnac.net

**Meridian Bioscience Europe b.v.**  
 Halderheiweg 6  
 5282 SN Boxtel  
 Holland  
 Tel.: +31 (411) 62 11 66  
 Fax: +31 (411) 62 48 41  
 e-mail: meridian@wxs.nl

**OTHER MERIDIAN PRODUCTS**

MERIDIAN PRODUCT

Merifluor® Cryptosporidium/Giardia  
 Premier™ C. difficile Toxin A  
 Premier™ H. pylori  
 Premier™ EHEC  
 Premier™ Rotacclone®  
 Premier™ Human Lyme  
 ImmunoCard® Mycoplasma  
 Para-Pak® C&S

CATALOG NUMBER

250050  
 601096  
 606096  
 608096  
 696004  
 696016  
 709030  
 900612

For technical assistance, please call our Technical Support Center at 1-800-343-3858 between the hours of 8 AM and 6 PM, Eastern Standard Time. To place an order, please call our Customer Service Department at 1-800-543-1980.

SN 11244

Printed in the U.S.A. on recycled paper  
 Copyright © Meridian Bioscience, Inc.

Rev. 5/03