



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

THE ABSOLUTE SOLUTION FOR *C. DIFFICILE* LABORATORY DIAGNOSIS



ImmunoCard® *C. difficile* GDH is a single test, rapid enzyme immunoassay for the detection of *C. difficile* common antigen, glutamate dehydrogenase. This screening test allows your lab to identify suspected cases of *C. difficile* disease in a rapid and cost-effective manner for the lab.

HIGHLY ACCURATE

- Rapid, cost-effective identification of suspected CDI specimens
- Early and confident elimination of negative results
 - Sensitivity: 97.6%
- High Negative Predictive Value
 - NPV: 99.6%

RAPID RESULTS

- Test outcomes within 30 minutes

ABSOLUTE PORTFOLIO

- GDH, toxin, and molecular testing options for *C. difficile* testing
- Algorithmic program empowers your lab with molecular confirmation in one hour




Meridian Bioscience offers the Absolute portfolio of diagnostic options for *C. difficile*




In a test tube, add 200 μ L of Sample Diluent using the dropper assembly.




Add 25 μ L of liquid stool (first mark on the transfer pipette) or 2 mm diameter of solid stool to the test tube. Vortex mixture.




Hold Enzyme Conjugate vertically and add 3 drops to the test tube. Vortex. Let the tube stand at 20–26 C for 15 minutes.




Using a transfer pipette, add 150 μ L (second mark on the pipette) of diluted specimen to the right and left sample (lower) ports. Incubate 5 minutes at 20–26 C.



Hold Wash Reagent vertically and add 3 drops to each reaction (upper) port.



Once Wash Reagent is absorbed, hold Substrate Reagent vertically and add 3 drops to both reaction (upper) ports. Incubate for 5 minutes at 20–26 C. Visually read the results within 30 seconds of the end of incubation.



This illustration is representative of the current Package Insert at the time of publication. Please refer to the most current version of the Package Insert for complete instructions.



ImmunoCard®
***C. difficile* GDH**

Ordering Information

Assays	Kit Size	Catalog #
ImmunoCard® <i>C. difficile</i> GDH	50 tests	716050
Premier® <i>C. difficile</i> GDH	96 microwells	611096
<i>illumigene</i> ® <i>C. difficile</i>	50 tests	280050
Premier® Toxins A&B	96 microwells	616096
ImmunoCard® Toxins A&B	50 tests	712050



For more information, contact a specialist at info@meridianbioscience.eu



Meridian
Bioscience® Europe
Inspired Science. Trusted Solutions.™

MERIDIAN BIOSCIENCE EUROPE

Italy – Milan Tel: +39 0331 43 36 36
France – Paris Tel: +33 1 42 56 04 40
Belgium – Nivelles Tel: +32 67 89 59 59
The Netherlands – Boxtel Tel: +31 411 62 11 66
Web: www.meridianbioscience.eu
Email: info@meridianbioscience.eu

For other countries, please visit: www.meridianbioscience.eu/distributors

WORLD HEADQUARTERS

Meridian Bioscience, Inc.
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Web: www.meridianbioscience.com

AUSTRALIA

Bioline (Aust) Pty Ltd
Tel: 1 800 993 120
Email: info.aust@bioline.com