SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

immunocard STAT!®

Rota-Adeno-Noro2

EXTERNAL CONTROL SET

External Control Set for use with Immunocard STAT! ROTA-ADENO-NORO2

REF 750301

IVD

INTENDED USE

Immunocard STAT! ROTA-ADENO-NORO2 External Control Set is a set including Positive and Negative Controls to be used with the Immunocard STAT! ROTA-ADENO-NORO2 kit (Catalog# 750320), as part of a routine quality control program.

The ROTA-ADENO Positive Control contains purified viral proteins from real samples and purified antigens, diluted in a buffered solution. The ROTA-ADENO Negative Control is a buffered solution stabilized by means of proteins. It contains sodium azide as preservative. These controls are intended to be used with **Immunocard STAT! ROTA-ADENO-NORO2** for the Rota-Adeno portion of the assay.

The NORO2 Positive Control contains virus-like particles (non-infectious) belonging to genotypes GI.1 (Norwalk) and GII.4 (MD145), diluted in a buffered solution. The NORO2 Negative Control is a buffered solution stabilized by means of proteins. It contains sodium azide as preservative. These controls are intended to be used with **Immunocard STAT! ROTA-ADENO-NORO2** for the Norovirus portion of the assay.

The External Controls are supplied ready to use.

SUMMARY AND EXPLANATION OF THE TEST

Quality control testing is performed to detect factors such as reagent deterioration, adverse environmental or test conditions, or variance in operator performance that can cause test errors. External control reagents, such as Immunocard STAT! ROTA-ADENO-NORO2 External Control Set, are reagents that are not built into the test system but are sampled and tested in the same manner as patient specimens.

BIOLOGICAL PRINCIPLES

Immunocard STAT! ROTA-ADENO-NORO2 External Control Set is used to create a known Positive and Negative Control panel. The frequency with which a laboratory performs external controls will be affected by:

1. The proficiency level of the laboratory,
2. The laboratory's own internal requirements,
3. The requirements of the laboratory's accrediting agencies,
4. The number of new and different operators performing the test,
5. Whether a new kit lot is being added to testing,
6. Whether the kit lot is from a different shipment,
7. Whether deviations from the manufacturer's stated storage or handling conditions have occurred.

When unacceptable quality control test results are obtained, all test results should be considered invalid. QC test failures can mean that one of the reagents or components is no longer reactive at the time of use, the test was not performed correctly or that reagents or sample were not added.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this kit is listed on the outer box.

1. **Immunocard STAT! ROTA-ADENO Positive Control:** buffered solution containing purified viral proteins from real samples and purified antigens.
2. **Immunocard STAT! ROTA-ADENO Negative Control:** buffered solution containing sodium azide as preservative.
3. **Immunocard STAT! NORO2 Positive Control:** buffered solution containing virus-like particles (non-infectious) belonging to genotypes GI.1 (Norwalk) and GII.4 (MD145).
4. **Immunocard STAT! NORO2 Negative Control:** buffered solution containing sodium azide as preservative.

The stability of the control set is 24 months if stored at 2-8 C. Meet always the Good Laboratory Practices when manipulating controls. The ROTA-ADENO Positive Control has been evaluated against the most risk infectious agents (Hepatitis A&B, HIV) leading to a negative result. The NORO2 Positive Control does not come from human culture cells and therefore is free from most risk infectious agents. The Negative Controls do not contain any transmissible infectious agent. The reagents are ready to use.

MATERIALS NOT PROVIDED

1. **Immunocard STAT! ROTA-ADENO-NORO2** kit (Catalog# 750320)
2. Vortex
3. Interval Timer
4. Micropipette (p200)

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. This is a quality control reagent set and is used to evaluate the performance of **Immunocard STAT! ROTA-ADENO-NORO2**. It is not directly used to test patient samples.
3. Store these reagents at 2-8 C when not in use. Do not freeze.
4. These reagents can be used with any **Immunocard STAT! ROTA-ADENO-NORO2** kit lot number.

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

QUALITY CONTROL TEST PROCEDURE

The test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Let the Positive and Negative Controls reach room temperature (19-27 C) for about 30 minutes before performing a test.
2. Use one **Immunocard STAT! ROTA-ADENO-NORO2** Test Device for each Positive and Negative Control reagent to be tested for Rotavirus-Adenovirus and/or Norovirus. Remove each Test Device from its foil pouch. Label each Test Device with the name of the Control to be tested.
3. Using a new pipette tip, transfer 170 µL of the ROTA-ADENO Positive or Negative Control to the sample area of reaction strip for the detection of Rotavirus-Adenovirus located in the left side of an appropriately labelled Test Device.
4. Using a new pipette tip, transfer 170 µL of the NORO2 Positive or Negative Control to the sample area of reaction strip for the detection of Norovirus located in the right side of an appropriately labelled Test Device.
5. Incubate the test for 15 minutes at room temperature.
6. At the end of testing, read results as described in the section INTERPRETATION OF RESULTS of **Immunocard STAT! ROTA-ADENO-NORO2** insert.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure.

LIMITATIONS OF THE PROCEDURE

The **Immunocard STAT! ROTA-ADENO-NORO2 External Control Set** is validated for use only with the **Immunocard STAT! ROTA-ADENO-NORO2** Rapid Assay.

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Meridian Bioscience Europe S.r.l