SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ImmunoCard STAT!® CGE

(The rapid discriminator of *Entamoeba histolytica*¹)

DEFINITIVE ANSWERS, CONFIDENT RESULTS™

- Easy, effective and efficient detection of *Cryptosporidium parvum*, *Giardia lamblia* & *Entamoeba histolytica*
- Room Temperature storage and test results within 10 minutes
- Expedite the initiation of treatment and pathogen control measures

ImmunoCard
STAT! CGE

 **Meridian**
Bioscience® Europe
Inspired Science. Trusted Solutions.®

TEST PROCEDURE

- Using the dropper with the Sample Diluent/Negative Control, add 1ml of Sample Diluent to a test tube.
- Add a small portion of approximately 5-6mm size (50mg) or 100µl stool using appropriate applicator or pipette. Vortex 15 seconds.
- Wait 3-5 minutes until the solid particles settle at the bottom of the tube or centrifuge.
- Add 125µl of the prepared specimen to the sample port of the device and read the results at 10 minutes.

This illustration is representative of the current Package Insert at the time of publication. Please refer to the most current version of the Package Insert for complete instructions.

RESULTS INTERPRETATION

Negative **Positives**

- Control
- Entamoeba
- Giardia
- Cryptosporidium

Ordering Information

Product	Quantity	Catalogue #
ImmunoCard STAT!® CGE	20 Tests	751420
ImmunoCard STAT!® CGE Positive Control	3 Tests	751403

For more information contact a specialist at info@meridianbioscience.eu
Visit us on the web at www.meridianbioscience.eu

REFERENCES

¹WHO recommends that *E. histolytica*/*E. dispar* is differentiated and such patients should not be treated on the basis of microscopy findings alone. **Wkly. Epidemiol. Rec.** 1997, 72, 97-99

MERIDIAN BIOSCIENCE EUROPE

Italy - Milan Tel: +39 0331 43 36 36
 France - Paris Tel: +33 (0)1 42 56 04 40
 Belgium - Nivelles Tel: +32 (0)67 89 59 59
 The Netherlands - Boxtel Tel: +31 (0)411 62 11 66
 United Kingdom - London Tel: +44 (0)20 8453 7970
 Web: www.meridianbioscience.eu
 Email: info@meridianbioscience.eu

WORLD HEADQUARTERS

Meridian Bioscience, Inc.
 3471 River Hills Drive
 Cincinnati, Ohio 45244 USA
 Web: www.meridianbioscience.com

AUSTRALIA

Bioline (Aust) Pty Ltd
 Tel: +61 (0)2 9209 4180
 Web: www.meridianbioscience.com
 Email: info.aust@bioline.com

For other countries, please visit www.meridianbioscience.eu/distributors