SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

NovaLisa[®]

Avidity Measles Virus IgG

ELISA Supplementary Instruction

Only for in-vitro diagnostic use

English	2
Deutsch	6
Français	10
Italiano	14
Español	18
Português.....	22
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía / Bibliografia	26
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos.....	27
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procedure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	28

ENGLISH

1. INTRODUCTION

The presence of IgG antibodies to Measles virus indicates the occurrence of the infection but does not distinguish between recent and past infection. Specific IgM antibodies are first detected approximately in ten days and peak at about four weeks post infection. They may persist for several months after acute infections. Based on the evidence that antibody avidity gradually increases after exposure to an immunogen, avidity of IgG antibodies can be used as a marker for distinguishing recent primary from long-term infections. Avidity describes the binding strength of a specific antibody to its antigen. Low-avidity IgG antibodies indicate a primary infection, whereas the presence of IgG antibodies with high avidity points to persistency or reactivation of infection.

2. INTENDED USE

The Avidity Measles Virus IgG ELISA is intended to indicate the Measles-specific IgG avidity in human serum or plasma (citrate, heparin) to differentiate between acute and past infection.

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microplates are coated with specific antigens to bind corresponding antibodies of the sample (dual pipetting). After washing the wells to remove all unbound sample material, one well is incubated with avidity reagent and the corresponding well with washing buffer. The avidity reagent removes the low-avidity antibodies from the antigens whereas the high-avidity ones are still bound to the specific antigens. After second washing step to remove the rest of avidity reagent and low-avidity antibodies, a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a third washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Avidity Reagent:** 1 bottle containing 15 ml of an urea solution, coloured blue, ready to use; black cap.
- **Avidity Measles Virus IgG Control Low:** 1 vial containing 2 ml control (human serum or plasma); coloured yellow; ready to use; blue cap.
- **Avidity Measles Virus IgG Control High:** 1 vial containing 2 ml control (human serum or plasma); coloured yellow; ready to use; red cap.

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 instruction for use Avidity Measles Virus IgG ELISA (Product Number: AMEA7330)
- 1 Instruction for use Measles Virus IgG ELISA (Product Number: MEAG0330)
- 1 empty labelled bottle (white with white cap) for ready to use Washing Buffer

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Avidity Reagent

If crystals have formed in the reagent warm up to 37 °C e.g. in a water bath and mix gently until they disappear.

6.2. Washing Buffer

It is recommended to fill 15 ml ready to use Washing Buffer into supplied bottle (s. 4.2) to use it in step Step 5 of the test preparation.

Note: Ready to use Washing Buffer is stable for 5 days at room temperature (20...25 °C).

7. SPECIMEN COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) **Measles Virus IgG positive** samples with this assay.

Note: For samples with high absorbance values (OD > 2.000) appropriate higher dilutions should be used.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Diluent. Dispense 10 µl sample and 1 ml IgG Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three to five and the volume of Washing Buffer from 300 µl to 350 µl to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

For avidity determination dual pipetting of standards/controls and diluted samples is needed.

1. Dispense 100 µl standards/controls and diluted samples into their respective wells. Leave wells A1/A2 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µl of Avidity Reagent in wells B1, C1, D1, E1 etc, except for the Substrate Blank well A1.
Dispense 100 µl of Washing Buffer in wells B2, C2, D2, E2 etc, except for the Substrate Blank well A2.
6. **Incubate for exactly 10 min at 37 ± 1 °C.**
7. Repeat step 4.
8. Dispense 100 µl Conjugate into all wells except in the blank wells (A1/A2).
9. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
10. Repeat step 4.
11. Dispense 100 µl TMB Substrate Solution into all wells.
12. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
13. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate, thereby a colour change from blue to yellow occurs.
14. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.2. Measurement

Adjust the ELISA microwell plate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA microwell plate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate Blank:** Absorbance value **< 0.100**
- **Control Low:** Avidity (%): **< 45 %**
- **Control High:** Avidity (%): **> 55 %**

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

For each patient sample or control calculate the ratio between the absorbance of the well dispensed with Avidity Reagent and the absorbance of the well dispensed with Washing Buffer multiplied by 100:

$$\frac{\text{Absorbance (sample or control) Avidity Reagent}}{\text{Absorbance (sample or control) Washing Buffer (diluted 1+19)}} \times 100 = \text{Avidity (\%)}$$

Note: For samples with high absorbance values (OD > 2.000) appropriate higher dilutions should be used.

9.3. Interpretation of Results

Result	Avidity	Interpretation
Low-avidity IgG	< 45 %	An avidity index of less than 45 % indicates a primary infection acquired within the past 2 months.
Equivocal	45 – 55 %	No clinical interpretation can be deduced from an equivocal result. It is recommended to take a second sample within an appropriate period of time (e.g. 2 weeks) and repeat testing. If the result of the repeated test is still equivocal, precise statements regarding the time of infection cannot be made.
High-avidity IgG	> 55 %	The presence of high-avidity IgG indicates a past infection or reinfection.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.
In immunocompromised patients and newborns serological data only have restricted value.
A result of high avidity can not exclude the possibility of a recent infection.

9.3.1. Antibody Isotypes and State of Infection

IgG	IgM	IgG-Avidity	Probable result
+	-	low	Vague, further investigation necessary
+	-	high	Indicatives of a past infection
+	+	low	Suggests a current or very recent infection
+	+	high	Suggests a past infection with persisting IgM or reactivation of infection

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.
For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Diagnostic Performance

The evaluation of the diagnostic performance of the Avidity Measles Virus IgG test was performed in comparison to well defined samples. The resulting relative agreement was 100 %.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of the in vitro diagnostic medical devices is intended by the manufacturer to secure suitability, performances and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

13. ORDERING INFORMATION

Prod. No.: AMEA7330 Avidity Measles Virus IgG (48 Determinations)

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Condorelli, Francesca; Ziegler, Thedi (1993): Dot immunobinding assay for simultaneous detection of specific immunoglobulin G antibodies to measles virus, mumps virus, and rubella virus. In *Journal of Clinical Microbiology* 31 (3), pp. 717–719.

Garg, R. K. (2002): Subacute sclerosing panencephalitis. In *Postgraduate Medical Journal* 78 (916), pp. 63–70. DOI: 10.1136/pmj.78.916.63.

Mercader, Sara; Garcia, Philip; Bellini, William J. (2012): Measles virus IgG avidity assay for use in classification of measles vaccine failure in measles elimination settings. In *Clinical and vaccine immunology : CVI* 19 (11), pp. 1810–1817. DOI: 10.1128/CVI.00406-12.

Messling, Veronika von; Milosevic, Dragana; Cattaneo, Roberto (2004): Tropism illuminated: lymphocyte-based pathways blazed by lethal morbillivirus through the host immune system. In *Proceedings of the National Academy of Sciences of the United States of America* 101 (39), pp. 14216–14221. DOI: 10.1073/pnas.0403597101.


Moss, William J.; Clements, C. John; Halsey, Neal A. (2003): Immunization of children at risk of infection with human immunodeficiency virus. In *Bulletin of the World Health Organization* 81 (1), pp. 61–70.

Pinsky, Norman A.; Huddleston, Jeanne M.; Jacobson, Robert M.; Wollan, Peter C.; Poland, Gregory A. (2003): Effect of multiple freeze-thaw cycles on detection of measles, mumps, and rubella virus antibodies. In *Clinical and diagnostic laboratory immunology* 10 (1), pp. 19–21. DOI: 10.1128/CDLI.10.1.19-21.2003.

Robert Koch Institut (RKI) (2015): Empfehlungen der Ständigen Impfkommission (STIKO) am Robert Koch-Institut/Stand: August 2015. In *Epidemiologisches Bulletin* 34, pp. 327–362. DOI: 10.17886/EpiBull-2015-001.2.

World Health Organization (WHO) (2009): The Immunological Basis for Immunization Series. Module 7: Measles.

SYMBOLS KEY / SYMBOLSCHLÜSSEL / EXPLICATION DES SYMBOLES / LEGENDA / SIMBOLOS / TABELA DE SIMBOLOS

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Avidity Measles Virus IgG

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all samples and standards/controls on the plate layout supplied in the kit.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank A1/A2	Avidity Control Low B1	Avidity Control Low B2	Avidity Control High C1	Avidity Control High C2	Sample (diluted 1+100) e.g. D1	Sample (diluted 1+100) e. g. D2
Avidity Control Low	-	100 µl	100 µl	-	-	-	-
Avidity Control High	-	-	-	100 µl	100 µl	-	-
Sample (diluted 1+100)	-	-	-	-	-	100 µl	100 µl
Cover wells with foil Incubate for 1 h at 37 °C Wash each well three times with 300 µl of Washing Buffer							
Avidity Reagent	-	100 µl	-	100 µl	-	100 µl	-
Washing Buffer	-	-	100 µl	-	100 µl	-	100 µl
Incubate for exactly 10 min at 37°C Wash each well three times with 300 µl of Washing Buffer							
Conjugate	-	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
Incubate for 30 min at room temperature (20...25 °C) Do not expose to direct sunlight Wash each well three times with 300 µl of Washing Buffer							
TMB Substrate solution	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
Incubate for exactly 15 min at room temperature (20...25 °C) in the dark							
Stop Solution	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
Photometric measurement at 450 nm (reference wavelength: 620 nm)							

NovaTec Immundiagnostica GmbH

Waldstraße 23 A6
63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760
Email: info@NovaTec-ID.com
Internet: www.NovaTec-ID.com

Fax: +49 (0) 6074-487629

AMEA7330engl,dt,fr,it,es,port-09012018-JH-gr.Lot095