**SZABO  
SCANDIC**

Part of Europa Biosite

## Produktinformation

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
[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](http://linkedin.com/company/szaboscandic)

**NovaLisa™**

# Candida albicans

## IgG - ELISA

CE

Enzyme immunoassay for the qualitative determination of IgG-class antibodies against Candida albicans in human serum or plasma

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For further languages please contact our authorized distributors.

Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía	Page / Seite / Page / Pagina / Página	30
Symbols Key / Symbolschlüssel / Explication des symboles / Legenda / Símbolos	Page / Seite / Page / Pagina / Página	31
Summary of Test Procedure/ Kurzanleitung Testdurchführung/ Résumé de la procedure de test/ Schema della procedura/ Resumen de la técnica	Page / Seite / Page/ / Pagina / Página	32

---

Product Number: CANG0060 (96 Determinations)

---

## 1. INTRODUCTION

Candida albicans is an ubiquitous present yeast like fungus belonging to the family of ascomycota. It is a facultative pathogen belonging to the normal microbial flora of skin and mucosa. Beside the yeast form which appears mainly with superficial infections, so called pseudomycelia are another morphological form of ascomycota. Pseudomycelia prevail with invasive mycosis. Although not pathogenic in healthy humans the fungus may be opportunistic in those suffering from a variety of disorders, and in those treated intensively with broad-spectrum antibiotics or immunosuppressive measures. Candidiasis is caused to about 90% by C. albicans. It is an acute, superficial or subacute, invasive infection. Acute infections with the fungus may produce lesions in the mouth, vagina, skin and nails. Subacute infections may affect bronchi, lungs, heart or meningeal. In immunosuppressant patients with cellular immunodeficiency, e.g., AIDS patients, C. albicans may lead to severe necroses of infected tissues.

Species	Symptoms	Complications	Mechanism of Infection
Candida albicans	Depending upon the location of C. albicans infectious the kind of lesion varies. Oral infections for example show typical creamy white patches of thrush whereas in bronchial candidiasis cough is the most characteristic symptom	Sepsis	Endogene: Usually candidiasis is due to increased susceptibility to a member of the normal human flora. Exogene: Transmission between individuals: Infection of newborns from the birth canal; transmission by non-sterile hospital equipment

### Detection:

- Microscopy: Gram stain, isolation of the organism in cultures
- Serology: Detection of specific antibodies by ELISA.

## 2. INTENDED USE

The NovaTec Candida albicans IgG ELISA is intended for the qualitative determination of IgG class antibodies against Candida albicans in human serum or plasma (citrate).

## 3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of IgG-class antibodies against Candida albicans is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiter strip wells are precoated with Candida albicans antigens to bind corresponding antibodies of the specimen. After washing the wells to remove all unbound sample material horseradish peroxidase (HRP) labelled anti-human IgG conjugate is added. This conjugate binds to the captured Candida albicans-specific antibodies. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is proportional to the amount of Candida albicans -specific IgG antibodies in the specimen. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450 nm is read using an ELISA microwell plate reader.

## 4. MATERIALS

### 4.1. Reagents supplied

- **Candida albicans Coated Wells (IgG):** 12 breakapart 8-well snap-off strips coated with Candida albicans antigen; in resealable aluminium foil.
- **IgG Sample Diluent \*\*\*:** 1 bottle containing 100 ml of buffer for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.2 mol/l; ready to use; red cap.
- **Washing Solution (20x conc.)\*:** 1 bottle containing 50 ml of a 20-fold concentrated buffer (pH 7.2 ± 0.2) for washing the wells; white cap.
- **Candida albicans anti-IgG Conjugate\*\*:** 1 bottle containing 20 ml of peroxidase labelled rabbit antibody to human IgG; coloured blue, ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB); ready to use; yellow cap.
- **Candida albicans IgG Positive Control\*\*\*:** 1 bottle containing 2 ml; coloured yellow; ready to use; red cap.
- **Candida albicans IgG Cut-off Control\*\*\*:** 1 bottle containing 3 ml; coloured yellow; ready to use; green cap.
- **Candida albicans IgG Negative Control\*\*\*\*:** 1 bottle containing 2 ml; coloured yellow; ready to use; blue cap.

\* contains 0.1 % Bronidox L after dilution

\*\* contains 0.2 % Bronidox L

\*\*\* contains 0.1 % Kathon

## **4.2. Materials supplied**

- 1 Strip holder
- 1 Cover foil
- 1 Test protocol
- 1 distribution and identification plan

## **4.3. Materials and Equipment needed**

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620nm
- Incubator 37°C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Disposable tubes
- Timer

## **5. STABILITY AND STORAGE**

---

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

## **6. REAGENT PREPARATION**

---

*It is very important to bring all reagents, samples and controls to room temperature (20...25°C) before starting the test run!*

### **6.1. Coated snap-off strips**

The ready to use breakapart snap-off strips are coated with Candida albicans antigen. Store at 2...8°C. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

### **6.2. Candida albicans anti-IgG Conjugate**

The bottle contains 20 ml of a solution with anti-human-IgG horseradish peroxidase, buffer, stabilizers, preservatives and an inert blue dye. The solution is ready to use. Store at 2...8°C. *After first opening stability expiry date when stored at 2...8°C.*

### **6.3. Controls**

The bottles labelled with Positive, Cut-off and Negative Control contain a ready to use control solution. It contains 0.1% Kathon and has to be stored at 2...8°C. *After first opening stability expiry date when stored at 2...8°C.*

### **6.4. IgG Sample Diluent**

The bottle contains 100 ml phosphatebuffer, stabilizers, preservatives and an inert yellow dye. It is used for the dilution of the patient specimen. This ready to use solution has to be stored at 2...8°C. *After first opening stability expiry date when stored at 2...8°C.*

### **6.5. Washing Solution (20xconc.)**

The bottle contains 50 ml of a concentrated buffer, detergents and preservatives. Dilute washing solution 1+19; e.g. 10 ml washing solution + 190 ml fresh and germ free redistilled water. The diluted buffer will keep for 5 days if stored at room temperature. *Crystals in the solution disappear by warming up to 37 °C in a water bath. After first opening the concentrate is stable until the expiry date*

### **6.6. TMB Substrate Solution**

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C, away from the light. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away. After first opening stability expiry date when stored at 2...8°C.*

### **6.7. Stop Solution**

The bottle contains 15 ml 0.2 M sulphuric acid solution (R 36/38, S 26). This ready to use solution has to be stored at 2...8°C. *After first opening stability until expiry date..*

## **7. SPECIMEN COLLECTION AND PREPARATION**

---

Use human serum or plasma (citrate) samples with this assay. If the assay is performed within 5 days after sample collection, the specimens should be kept at 2...8°C; otherwise they should be aliquoted and stored deep-frozen (-20 to -70°C). If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing.*

Heat inactivation of samples is not recommended.

### **7.1. Sample Dilution**

Before assaying, all samples should be diluted 1+100 with IgG Sample Diluent. Dispense 10µl sample and 1ml IgG Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

## **8. ASSAY PROCEDURE**

---

### **8.1. Test Preparation**

Please read the test protocol carefully before performing the assay. Result reliability depends on strict adherence to the test protocol as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend to increase the washing steps from three to five and the volume of washing solution from 300 $\mu$ l to 350 $\mu$ l to avoid washing effects. Prior to commencing the assay, the distribution and identification plan for all specimens and controls should be carefully established on the result sheet supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Please allocate at least:

- |         |              |                             |
|---------|--------------|-----------------------------|
| 1 well  | (e.g. A1)    | for the substrate blank,    |
| 1 well  | (e.g. B1)    | for the negative control,   |
| 2 wells | (e.g. C1+D1) | for the cut-off control and |
| 1 well  | (e.g. E1)    | for the positive control.   |

*It is recommended to determine controls and patient samples in duplicate, if necessary.*

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each control and sample.

Adjust the incubator to 37° ± 1°C.

1. Dispense 100 $\mu$ l controls and diluted samples into their respective wells. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37±1°C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 $\mu$ l of Washing Solution. Avoid overflows from the reaction wells. The soak time between each wash cycle should be >5sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!

*Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.*

5. Dispense 100 $\mu$ l Candida albicans anti-IgG Conjugate into all wells except for the blank well (e.g. A1). Cover with foil.
6. **Incubate for 30 min at room temperature. Do not expose to direct sunlight.**
7. Repeat step 4.
8. Dispense 100 $\mu$ l TMB Substrate Solution into all wells
9. **Incubate for exactly 15 min at room temperature in the dark.**
10. Dispense 100 $\mu$ l Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution.

*Any blue colour developed during the incubation turns into yellow.*

*Note: Highly positive patient samples can cause dark precipitates of the chromogen! These precipitates have an influence when reading the optical density. Predilution of the sample with physiological sodium chloride solution, for example 1+1, is recommended. Then dilute the sample 1+100 with dilution buffer and multiply the results in NTU by 2.*

11. Measure the absorbance of the specimen at 450/620nm within 30 min after addition of the Stop Solution.

### **8.2. Measurement**

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

*If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!*

**Measure the absorbance** of all wells at **450 nm** and record the absorbance values for each control and patient sample in the distribution and identification plan.

*Dual wavelength reading using 620 nm as reference wavelength is recommended.*

Where applicable calculate the **mean absorbance values** of all duplicates.

## **9. RESULTS**

---

### **9.1. Run Validation Criteria**

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate blank** in A1: Absorbance value **< 0.100**.
- **Negative control** in B1: Absorbance value **< 0.200 and < cut-off**
- **Cut-off control** in C1 and D1: Absorbance value **0.150 – 1.30**.
- **Positive control** in E1: Absorbance value **> cut-off**.

If these criteria are not met, the test is not valid and must be repeated.

## 9.2. Calculation of Results

The mean absorbance value of the cut-off control has to be calculated.

The cut-off is calculated by multiplying the cut-off control value with 2.5.

Example: Absorbance value Cut-off control 0.39 + absorbance value Cut-off control 0.37 = 0.76 / 2 = 0.38

Cut-off control value = 0.38

Calculated Cut-off = 0.38 x 2.5 = 0.95

## 9.3. Interpretation of Results

Samples are considered NEGATIVE if the absorbance value is lower than the calculated cut-off.

Samples with an absorbance value up to 20% above the calculated cut-off should **not** be considered as clearly positive  
→ grey zone

It is recommended to repeat the test again 2 - 4 weeks later with a fresh sample. If results in the second test are again in the grey zone the sample has to be considered NEGATIVE.

Samples are considered POSITIVE if the absorbance value is more than 20% higher than the calculated cut-off.

### 9.3.1. Results in NovaTec Units

Patient (mean) absorbance value x 10 = [NovaTec Units = NTU]  
Cut-off control value

Example:  $\frac{1.216 \times 10}{0.38} = 32$  NTU (NovaTec Units)

Calculated Cut-off: 25 NTU

Grey zone: 25 - 30 NTU

Negative: <25 NTU

Positive: >30 NTU

The borderline value of *Candida albicans* IgG was established in a way that most of tested serum samples from unselected blood donors showed negative or borderline results. Therefore positive IgG-results indicate an acute infection and not the remaining IgG antibodies from a previous infection. If already moderately elevated antibody titers should be detected and interpreted as positive results, a borderline region of 10-25 NTU can be also applied. In this connection, however, it has to be unconditionally observed that titers which are due to seroprevalence are no longer completely blanked out. Therefore a serological positive result has to be looked upon then very critically, especially if the clinical picture and/or results of other detection methods are not available.

## 10. SPECIFIC PERFORMANCE CHARACTERISTICS

### 10.1. Precision

<u>Interassay</u>	<u>n</u>	<u>Mean</u>	<u>Cv (%)</u>
Pos. Serum	20	18	5.1
<b>Intraassay</b>			
<u>Intraassay</u>	<u>n</u>	<u>Mean</u>	<u>Cv (%)</u>
Pos. Serum	8	1.0	4.3

### 10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is >95 %.

### 10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 93 %.

### 10.4. Interferences

Interferences with hemolytic, lipemic or icteric sera are not observed up to a concentration of 10 mg/ml hemoglobin, 5 mg/ml triglycerides and 0.2 mg/ml bilirubin.

**Note:** The results refer to the groups of samples investigated; these are not guaranteed specifications.

## 11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values. Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

In immunocompromised patients and newborns serological data only have restricted value.

## **12. PRECAUTIONS AND WARNINGS**

---

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of the in vitro diagnostic medical devices is intended by the manufacturer to secure suitability, performances and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- The NovaLisa™ ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

**WARNING:** In the used concentration Bronidox L has hardly any toxicological risk upon contact with skin and mucous membranes!

**WARNING:** Sulphuric acid irritates eyes and skin. Keep out of the reach of children. Upon contact with the eyes, rinse thoroughly with water and consult a doctor!

### **12.1. Disposal Considerations**

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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## **BIBLIOGRAPHY / LITERATUR / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFÍA**

---

- IRURETAGOYENA, R.J., REGÚEZ, P., QUINDÓS, G. AND PONTÓN, J. (2000) ANTIBODIES TO CANDIDA ALBICANS GERM TUBES IN TWO INTENSIVE CARE PATIENTS WITH INVASIVE CANDIDIASIS. REV. IBEROAM. MICOL. 17, 93-96
- ELLEPOLA, A.N.B. AND MORRISON, C.J. (2005) LABORATORY DIAGNOSIS OF INVASIVE CANDIDIASIS. J. MICROBIOL. 43 (S), 65-84
- QUINDÓS, G., MORAGUES, M.D. AND PONTÓN, J. (2004) IS THERE A ROLE FOR ANTIBODY TESTING IN THE DIAGNOSIS OF INVASIVE CANDIDIASIS? REV. IBEROAM. MICOL. 21, 10-14
- ARAJ, G.F., HOPFER, R.L., CHESTNUT, S., FAINSTEIN, V. AND BODEY, G.P. (1982) DIAGNOSTIC VALUE OF THE ENZYME-LINKED IMMUNOSORBENT ASSAY FOR DETECTION OF CANDIDA ALBICANS CYTOPLASMIC ANTIGEN IN SERA OF CANCER PATIENTS. J. CLIN. MICROBIOL. 16, 46-52
- AKPAN, A. AND MORGAN, R. (2002) ORAL CANDIDIASIS. POSTGRAD. MED. J. 78, 455-459
- MEUNIER-CARPENTIER, F., KIEHN, T.E. AND ARMSTRONG, D. (19981) FUNGEMIA IN THE IMMUNOCOMPROMISED HOST. AM. J. MED. 71, 363-370
- WITKIN, S.S., YU, I.R. AND LEDGER, W.J. (1983) INHIBITION OF CANDIDA ALBICANS- INDUCED LYMPHOCYTE PROLIFERATION BY LYMPHOCYTES AND SERA FROM WOMEN WITH RECURRENT VAGINITIS. AM. J. OBSTET. GYNECOL. 147, 809-811

Symbols Key/ Symbolschlüssel/ Explication des symboles / Legenda / Símbolos	
	Manufactured by / Hergestellt von/ Fabriqué par/ Prodotto da/ Fabricado por
	In Vitro Diagnostic Medical Device/ In Vitro Diagnosticum/ Dispositif médical de diagnostic <i>in vitro</i> / Diganostico <i>in vitro</i> / Producto para diagnóstico In vitro
	Lot Number/ Chargenbezeichnung/ Numéro de lot/ Lotto/ Número de lote
	Expiration Date/ Verfallsdatum/ Date de péremption/ Scadenza/ Fecha de caducidad
	Storage Temperature/ Lagertemperatur/ Température de conservation/ Temperatura di conservazione / Temperatura de almacenamiento
	CE Mark/ CE-Zeichen/ Marquage CE / Marchio CE/ Marca CE
	Catalogue Number/ Katalog Nummer/ Référence du catalogue/ Numero di codice/ Número de Catálogo
	Consult Instructions for Use/ Gebrauchsanweisung beachten/ Consulter la notice d'utilisation/ Consultare le istruzioni/ Consulte las Instrucciones de Uso
	Microplate/ Mikrotiterplatte/ Microplaquette/ Micropiastra/ Microplaca
	Conjugate/ Konjugat/ Conjugué/ Conjugato/ Conjugado
	Control serum, negative/ Kontrollserum, negative/ Sérum de contrôle négatif/ siero di controllo, negativo /Suero control negativo
	Control serum, positive/ Kontrollserum, positiv/ Sérum de contrôle positif/ siero di controllo, positivo/ Suero de control positivo
	Cut off control serum/ Cut off Kontrollserum/ Sérum de contrôle du cut-off/ siero di controllo, cut-off/ Suero control Cut-off
	Sample diluent buffer IgG/ IgG-Probenverdünnungspuffer/ Tampon diluant pour échantillon IgG/ soluzione tampone per i campioni IgG/ solución tampón para muestras IgG
	Stop solution/ Stopplösung/ Solution d'arrêt/Soluzione bloccante
	TMB Substrate solution/ TMB-Substratlösung/ Substrat TMB/ soluzione substrato TMB/ solución substrato TMB
	Washing solution 20x concentrated/ Waschlösung 20x konzentriert/ Solution de lavage concentré 20 x/ soluzione di lavaggio concentrazione x20/ solución de lavado concentrado x20
	Contains sufficient for "n" tests/ Ausreichend für "n" Tests/ Contenu suffisant pour "n" tests/ Contenuto sufficiente per "n" saggi/ Contenido suficiente para "n" tests

# SCHEME OF THE ASSAY

Candida albicans IgG-ELISA

## Test preparation

Prepare reagents and samples as described.  
Establish the distribution and identification plan for all specimens and controls on the result sheet supplied in the kit.  
Select the required number of microtiter strips or wells and insert them into the holder.

## Assay procedure

	Substrate blank (e.g. A1)	Negative control	Positive control	Cut-off control	Sample (diluted 1+100)
Negative control	-	100µl	-	-	-
Positive control	-	-	100µl	-	-
Cut-off control	-	-	-	100µl	-
Sample (diluted 1+100)	-	-	-	-	100µl
Cover wells with foil supplied in the kit					
<b>Incubate for 1 h at 37°C</b>					
Wash each well three times with 300µl of washing solution					
Conjugate	-	100µl	100µl	100µl	100µl
Cover wells with foil supplied in the kit					
<b>Incubate for 30 min at room temperature</b>					
Wash each well three times with 300µl of washing solution					
TMB Substrate	100µl	100µl	100µl	100µl	100µl
<b>Incubate for exactly 15 min at room temperature in the dark</b>					
Stop Solution	100µl	100µl	100µl	100µl	100µl
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

**NovaTec Immundiagnostica GmbH**

**Technologie & Waldpark**

Waldstr. 23 A6  
D-63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760  
Fax: +49 (0) 6074-487629  
Email : info@NovaTec-ID.com  
Internet: www.NovaTec-ID.com

CANG0060engl,dt,fr,it,es07042009-AD