SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Giardia lamblia antigen

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	7
Español	12
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía / Bibliografia	18
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	18
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	19
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procedure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	20

Product Number: GIA0160S (96 Determinations)

ENGLISH

1. INTRODUCTION

Giardia lamblia is one of the most common human intestinal protozoan pathogens worldwide. The incidence strongly depends on the geographic region and reaches 2-7 % in central Europe and exceeds 50 % in tropical countries.

The life cycle of *Giardia lamblia* is characterized by two stages: the trophozoite and the cyst stage. The trophozoite is the motile dividing stage and inhabits the upper small intestine. Ascending infections of the gallbladder may also occur. The cyst is the infective form of the parasite. It develops in the intestine and is excreted with the faeces. Cysts are transmitted via contaminated food or drinking water but also from person to person. The clinical picture of a *Giardia lamblia* infection ranges from the asymptomatic carrier state to acute diarrhea which is often accompanied by abdominal pain and flatulence. Chronic giardiasis can cause severe malabsorption syndrome.

Giardiasis is usually diagnosed by microscopic detection of trophozoites and/or cysts in faecal smears after commonly used staining techniques or direct immune fluorescence. These methods are time-consuming, require trained personnel and can only detect parasites with intact morphology. Immunologic methods like enzyme immunoassays detecting *Giardia lamblia* antigens may overcome these problems.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Giardia lamblia</i>	Giardiasis (Lambliasis)	Diarrhoea, nausea, weight loss, malabsorption, abdominal cramps, flatulence and anaemia	Transmission through faecal-oral ingestion of cysts. Adhesive disk allows organism to attach to the surface of intestinal cells, blocking the absorption of fats and vitamins

2. INTENDED USE

The *Giardia lamblia* antigen ELISA is intended for the qualitative determination of *Giardia lamblia* antigens in faeces.

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of *Giardia lamblia* antigens is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microplates are coated with polyclonal anti-*Giardia lamblia* antibodies. Diluted patient specimen (stool), and ready to use controls are added to these wells and during the first incubation, *Giardia* antigens present in the stool supernatant are captured. After washing the wells to remove all unbound sample and control a horseradish peroxidase labelled anti-*Giardia* antibody is added which binds to the captured antigen. After incubation and washing off unbound enzyme the antigen-conjugate-complex is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity is proportional to the amount of *Giardia lamblia* antigens in the patient specimen. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Coated microplate:** 12 breakapart 8-well snap-off strips coated with polyclonal anti-*Giardia lamblia* antibodies; vacuum sealed, in resealable aluminium foil.
- **Sample Diluent:** 1 bottle containing 100 ml of ready to use buffer for sample dilution, coloured yellow; ready to use; black cap.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.25 mol/l; ready to use; yellow cap.
- **Washing Buffer (10x conc.):** 1 bottle containing 100 ml of a 10-fold concentrated buffer for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 12 ml of peroxidase labelled anti-*Giardia lamblia* antibodies, coloured green; ready to use; brown cap; < 2 % NMP (for hazard and precautionary statements see 12.1).
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-Tetramethylbenzidine (TMB); ready to use; blue cap.
- **Positive Control:** 1 bottle containing 2 ml, coloured blue; ready to use; red cap.
- **Negative Control:** 1 bottle containing 2 ml, coloured blue; ready to use; green cap.

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620 nm
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C. After first opening the reagents are stable for another 2 months if stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run! Mix gently without causing foam!

6.1. Coated Microplate

The breakapart snap-off strips are coated with polyclonal antibodies against *Giardia lamblia*. The strips are vacuum sealed. Allow the sealed plate to reach room temperature before opening. Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (10x conc.)

Dilute Washing buffer 1+9 with distilled water; e.g. 10 ml Washing buffer + 90 ml distilled water. The diluted buffer is stable for 30 days if stored at 2...8 °C.

6.3. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

The test is intended for the detection of *Giardia lamblia* in diluted stool specimen. Either fresh or frozen samples may be used in this test. If samples are stored frozen, thaw sample quickly, warm it to room temperature and mix thawed samples well before dilution. If the assay is performed within 72 hours after sample collection, the specimen should be kept at 2...8°C; otherwise they should be aliquoted and stored deep-frozen (-20 °C). Avoid repeated freezing and thawing.

7.1. Sample Dilution

Before assaying, all fresh or frozen samples should be diluted 1+10 with Sample Diluent. Pipette 1000 µl sample diluent into a clean tube. Using a disposable stirring rod transfer about 100 mg (about 2-3 mm diameter) of faeces if solid or pipette 100 µl if liquid into the tube with sample diluent and suspend thoroughly, e.g. on a vortex. Allow floating particles to sediment for 10 min at most. If a stool suspension sediments longer than 10 min the sample should be mixed again immediately before starting the assay.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and controls should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each control and sample.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each control and sample.

1. Dispense 75 µl ready to use controls and **100 µl** diluted stool samples (**supernatant**) into their respective wells. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 30 min at room temperature (20...25 °C).**
4. When incubation has been completed, remove the foil, decant the content of the wells and wash each well five times with 300 µl of Washing Buffer. Avoid overflows from the reaction wells. The soak time between each wash cycle should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.
5. Dispense 2 drops (or 75 µl) of Conjugate into all wells except for the blank well (e.g. A1). Cover with foil.
6. **Incubate for 30 min at room temperature (20...25 °C).**
7. Repeat step 4.
8. Dispense 2 drops (or 75 µl) of TMB Substrate Solution into all wells.
9. **Incubate for 10 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 2 drops (or 75 µl) of Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance of the specimen at 450/620 nm within 30 min after addition of the Stop Solution.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the Substrate Blank, subtract the absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each control and patient sample in the distribution and identification plan.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Negative Control:** mean absorbance value ≤ 0.200
- **Positive Control:** mean absorbance value ≥ 0.800

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The cut-off is calculated by addition of 0.20 absorbance units to the measured absorption of the mean value of the two negative control determinations.

Example: $0.12 \text{ OD Negative Control} + 0.14 \text{ OD Negative Control} = 0.26 \div 2 = 0.13$
Cut-off = absorbance mean value of the Negative Control + 0.20
Cut-off = $0.13 + 0.20 = 0.33$

9.3. Interpretation of Results

Samples are considered **positive** if the absorbance value is equal to or higher than the Cut-off value.

Samples are considered **negative** if the absorbance value is lower than the Cut-off value.

10. SPECIFIC PERFORMANCE CHARACTERISTICS

10.1. Precision

Intraassay	n	Mean (E)	CV (%)
#1	12	1.438	6.86
#2	12	1.007	6.41
#3	12	0.673	6.36
#4	12	0.243	5.04

Interassay	n	Mean (E)	CV (%)
#1	11	1.488	4.77
#2	11	1.034	4.64
#3	11	0.656	7.77
#4	11	0.346	10.9

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. 369 known negative samples were tested with this assay in comparison to IFT resulting in a specificity of 99.5% (367/369).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. 40 known positive samples were tested with this assay in comparison to IFT resulting in a sensitivity of 97.5% (39/40).

10.4. Lower Detection Limit

The lower detection limit of the Giardia lamblia antigen assay was determined by titration of faecal samples spiked with Giardia lamblia cysts and trophozoites from culture. The lower detection limit for Giardia lamblia was determined 5.0 x 10³ cysts and 2.5 x 10⁴ trophozoites per ml of diluted faecal sample.

10.5. Cross Reactivity

Faecal samples positive for one of the following intestinal parasites and other pathogens resp., did not show any cross reactions in this assay:

Adenovirus, Ancylostoma duodenale, Ascaris lumbricoides, Blastocystis hominis, Cryptosporidium parvum, Entamoeba dispar, Entamoeba histolytica, Rotavirus.

Negative stool samples have been spiked with >10⁸ colony forming units of the following microorganisms and tested negative with the Giardia lamblia antigen ELISA.

<i>Aeromonas hydrophila</i>	ATCC 7966
<i>Bacillus cereus</i>	ATCC 11778
<i>Bacillus subtilis</i>	ATCC 6633
<i>Bacteroides fragilis</i>	ATCC 25285
<i>Candida albicans</i>	ATCC 10231
<i>Campylobacter coli</i>	ATCC 33559
<i>Campylobacter jejuni</i>	ATCC 33291
<i>Citrobacter freundii</i>	ATCC 8090
<i>Clostridium sordellii</i>	ATCC 9714
<i>Enterobacter aerogenes</i>	ATCC 13048
<i>Enterobacter cloacae</i>	ATCC 13047
<i>Enterococcus faecalis</i>	ATCC 29212
<i>Escherichia coli</i>	ATCC 25922
<i>Klebsiella pneumoniae</i>	ATCC 13883
<i>Peptostreptococcus anaerobius</i>	ATCC 27337
<i>Proteus vulgaris</i>	ATCC 8427
<i>Pseudomonas aeruginosa</i>	ATCC 10145
<i>Salmonella enterica</i> Serovar <i>enteritidis</i>	ATCC 13076
<i>Salmonella enterica</i> Serovar <i>typhimurium</i>	ATCC 14028
<i>Shigella flexneri</i>	ATCC 12022
<i>Shigella sonnei</i>	ATCC 25931
<i>Staphylococcus aureus</i>	ATCC 25923
<i>Staphylococcus epidermidis</i>	ATCC 12228
<i>Vibrio parahaemolyticus</i>	ATCC 17802

Note: The results refer to the groups of samples investigated; these are not guaranteed specifications.

11. LIMITATIONS OF THE PROCEDURE

There is no correlation between the measured absorbance and the seriousness of the infection. It is also not allowed to correlate absorbance of the sample with that of the positive control. Cross-contamination of reagents and samples can produce false results. Incorrect dilutions, not sufficiently homogenized samples and samples, which stayed for sedimentation for more than 10 minutes can cause false results.

A negative ELISA result does not exclude a *Giardia lamblia* infection, because the excretion of cysts is periodic. Thus at least one further stool specimen of the regarding person should be demanded in case of a negative test result but clinical suspect.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

12. PRECAUTIONS AND WARNINGS

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of the in vitro diagnostic medical devices is intended by the manufacturer to secure suitability, performances and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- Some of the reagents contain small amounts of Thimerosal (< 1% w/v) and Kathon (1.0% v/v) as preservatives. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components as potentially hazardous.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

12.1. Safety note for NMP-containing reagents

The Conjugate contains NMP.

Therefore, the following hazard and precautionary statements apply.

Danger

H360D

May damage the unborn child.

P280

Wear protective gloves, protective clothing, eye protection.

P308 + P313

IF exposed or concerned: Get medical advice/ attention.

Further information can be found in the safety data sheet.

12.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

13. ORDERING INFORMATION

Prod. No.: GIA0160S Giardia lamblia antigen ELISA (96 Determinations)

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Cardona, Guillermo A.; Carabin, Hélène; Goni, Pilar; Arriola, Larraitz; Robinson, Guy; Fernández-Crespo, Juan C. et al. (2011): Identification and molecular characterization of Cryptosporidium and Giardia in children and cattle populations from the province of Alava, North of Spain. In *The Science of the total environment* 412-413, pp. 101–108. DOI: 10.1016/j.scitotenv.2011.09.076.

Faubert, Gaétan (2000): Immune Response to Giardia duodenalis. In *Clinical Microbiology Reviews* 13 (1), pp. 35–54. DOI: 10.1128/CMR.13.1.35-54.2000.

Janoff, Edward N.; Craft, J. Carl; Pickering, Larry K.; Novotny, Thomas; Blaser, Martin J.; Knisley, Cathy V.; Reller, L. Barth (1989): Diagnosis of Giardia lamblia infections by detection of parasite-specific antigens. In *Journal of Clinical Microbiology* 27 (3), pp. 431–435.


Manual of clinical microbiology (1995). 6. ed., [1. Dr.]. Edited by Patrick R. Murray. Washington, DC: ASM Press.

Wolfe, Martin S. (1992): Giardiasis. In *Clinical Microbiology Reviews* 5 (1), pp. 93–100.

ABBREVIATIONS / ABKÜRZUNGEN / ABRÉVIATIONS / ABBREVIAZIONI / ABREVIACIONES / ABREVIATURAS

NMP	N-Methyl-2-pyrrolidone
------------	------------------------

SYMBOLS KEY / SYMBOLSCHLÜSSEL / EXPLICATION DES SYMBOLES / LEGENDA / SIMBOLOS / TABELA DE SIMBOLOS

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Giardia lamblia antigen ELISA

Test Preparation

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all specimens and controls on the plate layout supplied in the kit.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank (A1)	Negative Control	Positive Control	Sample (diluted 1+10)
Negative control	-	75 µl	-	-
Positive control	-	-	75 µl	-
Sample (diluted 1+10)	-	-	-	100 µl
Cover wells with foil supplied in the kit Incubate for 30 min at room temperature (20...25 °C) Wash each well five times with 300 µl of Washing Buffer				
Conjugate	-	2 drops (or 75 µl)	2 drops (or 75 µl)	2 drops (or 75 µl)
Cover wells with foil supplied in the kit Incubate for 30 min at room temperature (20...25 °C) Do not expose to direct sunlight Wash each well five times with 300 µl of Washing Buffer				
TMB Substrate solution	2 drops (or 75 µl)	2 drops (or 75 µl)	2 drops (or 75 µl)	2 drops (or 75 µl)
Incubate for 10 min at room temperature (20...25 °C) in the dark				
Stop Solution	2 drops (or 75 µl)	2 drops (or 75 µl)	2 drops (or 75 µl)	2 drops (or 75 µl)
Photometric measurement at 450 nm (reference wavelength: 620 nm)				

NovaTec Immundiagnostica GmbH

Waldstraße 23 A6
 63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760 Fax: +49 (0) 6074-487629
 Email: info@NovaTec-ID.com
 Internet: www.NovaTec-ID.com

GIA0160Sengl,dt,es-17072018-JH-gr.Lot052