# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation

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[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

**NovaLisa®**

# Strongyloides

**ELISA**

CE

**Only for in-vitro diagnostic use**

English .....	2
Deutsch .....	6
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía .....	13
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas .....	13
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos .....	14
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procédure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste .....	15

---

Product Number:      STRO0690 (96 Determinations)

---

## ENGLISH

### 1. INTRODUCTION

---

*Strongyloides* is a genus containing some 50 species of obligate gastrointestinal parasites of vertebrates. *Strongyloides stercoralis* is the scientific name of a human parasitic roundworm causing the disease of strongyloidiasis. Its common name is pinworm in the UK and threadworm in the US. The *Strongyloides stercoralis* nematode can parasitize humans. The adult parasitic stage lives in tunnels in the mucosa of the small intestine.

*S. stercoralis* can be found in areas with tropical and subtropical climates but cases also occur in temperate area, more frequently in rural areas. *S. stercoralis* has a very low prevalence in societies where fecal contamination of soil or water is rare. Many people infected are usually asymptomatic at first. Symptoms include dermatitis: swelling, itching, larva currens, and mild hemorrhage at the site where the skin has been penetrated. If the parasite reaches the lungs, the chest may feel as if it is burning, and wheezing and coughing may result, along with pneumonia-like symptoms (Löffler's syndrome). The intestines could eventually be invaded, leading to burning pain, tissue damage, sepsis, and ulcers. In severe cases, edema may result in obstruction of the intestinal tract, as well as loss of peristaltic contractions. *Strongyloides* infection in immunocompromised individuals (particularly following the administration of steroids, for example following transplant surgery) can result in disseminated strongyloidiasis, in which worms move beyond the confines of the gut into other organs. This is fatal unless anti-*Strongyloides* therapy is given.

Locating juvenile larvae, either rhabditiform or filariform, in recent stool samples will confirm the presence of this parasite. Other techniques used include direct fecal smears, culturing fecal samples on agar plates, serodiagnosis through ELISA, and duodenal fumigation.

### 2. INTENDED USE

---

The Strongyloides ELISA is intended for the qualitative determination of antibodies against Strongyloides in human serum or plasma (citrate or heparin).

### 3. PRINCIPLE OF THE ASSAY

---

The qualitative immunoenzymatic determination of IgG and IgM antibodies against Strongyloides is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiter strip wells are precoated with recombinant Strongyloides antigens to bind corresponding antibodies of the specimen. After washing the wells to remove all unbound sample material horseradish peroxidase (HRP) labelled anti-human IgG and IgM conjugate is added. This conjugate binds to the captured Strongyloides -specific antibodies. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is proportional to the amount of Strongyloides -specific antibodies in the specimen. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450 nm is read using an ELISA microwell plate reader.

### 4. MATERIALS

---

#### 4.1. Reagents supplied

- **Strongyloides Coated Wells:** 12 breakapart 8-well snap-off strips coated with recombinant Strongyloides antigens in resealable aluminium foil.
- **Sample Diluent:** 1 bottle containing 100 ml of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap. < 0.1 % MIT, < 0.1% CMIT, < 0.1% NaN<sub>3</sub>
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.2 mol/l; ready to use; red cap.
- **Washing Solution (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated phosphate buffer (0.2 M); pH 7.2 ± 0.2; for washing the wells; white cap. < 1% Ethanol, < 0.5% Bromonitrodioxane
- **Strongyloides Conjugate:** 1 bottle containing 20 ml of peroxidase labelled antibody to human IgG/IgM; in phosphate buffer (10 mM); coloured blue, ready to use; black cap. < 1% Ethanol, < 0.5% Bromonitrodioxane
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB) (< 0.04%); ready to use; yellow cap. < 0.0001% CMIT, < 0.0001% MIT, < 0.01% H<sub>2</sub>O<sub>2</sub>
- **Strongyloides Positive Control:** 1 bottle containing 2 ml control (human serum or plasma); coloured yellow; ready to use; red cap. < 0.1% Bromonitrodioxan, < 0.1% MIT
- **Strongyloides Cut-off Control:** 1 bottle containing 3 ml control (human serum or plasma); coloured yellow; ready to use; green cap. < 0.1% Bromonitrodioxan, < 0.1% MIT
- **Strongyloides Negative Control:** 1 bottle containing 2 ml control (human serum or plasma); coloured yellow; ready to use; blue cap. < 0.1 % MIT, < 0.1% CMIT, < 0.1% NaN<sub>3</sub>

#### 4.2. Materials supplied

- 1 Strip holder
- 1 Cover foil
- 1 Test protocol
- 1 Distribution and identification plan

### 4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Disposable tubes

## 5. STABILITY AND STORAGE

---

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

## 6. REAGENT PREPARATION

---

It is very important to bring all reagents, samples and standards/controls to room temperature (20...25 °C) before starting the test run!

### 6.1. Coated snap-off strips

The ready to use break-apart snap-off strips are coated with recombinant Strongyloides antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

### 6.2. Washing Solution (20x conc.)

Dilute Washing Solution 1 + 19; e. g. 10 ml Washing Solution + 190 ml fresh and germ free redistilled water. The diluted buffer is stable for 5 days at room temperature. Crystals in the concentrate disappear by warming up to 37 °C in a water bath.

### 6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

## 7. SPECIMEN COLLECTION AND PREPARATION

---

Use human serum or plasma (citrate or heparin). If the assay is performed within 5 days after sample collection, the specimens should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

### 7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with Sample Diluent. Dispense 10 µl sample and 1 ml Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

## 8. ASSAY PROCEDURE

---

### 8.1. Test Preparation

Please read the test protocol carefully **before** performing the assay. Result reliability depends on strict adherence to the test protocol as described. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three to five and the volume of Washing Solution from 300 µl to 350 µl to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all specimens and standards/controls should be carefully established on the result sheet supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Please allocate at least:

1 well	(e. g. A1)	for the Substrate Blank,
1 well	(e. g. B1)	for the Negative Control,
2 wells	(e. g. C1+D1)	for the Cut-off Control and
1 well	(e. g. E1)	for the Positive Control

It is recommended to determine standards/controls and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µl standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl of Washing Solution. Avoid overflows from the reaction wells. The soak time between each wash cycle should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µl Strongyloides Conjugate into all wells except for the Substrate Blank well (e. g. A1). Cover with foil.
6. **Incubate for 30 min at room temperature.** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µl TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.**
10. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate.  
Any blue colour developed during the incubation turns into yellow.
11. Measure the absorbance of the specimen at 450/620 nm within 30 min after addition of the Stop Solution.

## 8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **Substrate Blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the Substrate Blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

**Measure the absorbance** of all wells at **450 nm** and record the absorbance values for each standard/control and patient sample in the distribution and identification plan.

Dual wavelength reading using 620 nm as reference wavelength is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

## 9. RESULTS

### 9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate Blank** in A1: Absorbance value **< 0.100**
- **Negative Control** in B1: Absorbance value **< 0.200 and < Cut-off**
- **Cut-off Control** in C1 and D1: Absorbance value **0.150 – 1.300**
- **Positive Control** in E1: Absorbance value **> Cut-off**

If these criteria are not met, the test is not valid and must be repeated.

### 9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43

Cut-off = 0.43

#### 9.2.1. Results in Units [NTU]

$\frac{\text{Patient (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec Units} = \text{NTU}]$

Example:  $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU (Units)}$

### 9.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as <b>negative</b> .
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

In immunocompromised patients and newborns serological data only have restricted value.

## 10. SPECIFIC PERFORMANCE CHARACTERISTICS

### 10.1. Precision

<u>Interassay</u>	<u>n</u>	<u>Mean (NTU)</u>	<u>Cv (%)</u>
Positive Sample	24	22,11	5,52
Positive Sample	24	22,47	6,28
Negative Sample	24	6,09	8,52

<b>Intraassay</b>	<b>n</b>	<b>Mean (OD)</b>	<b>Cv (%)</b>
Positive Sample	24	0,604	5,74
Positive Sample	24	0,710	9,68
Negative Sample	24	0,323	4,61

### 11. 10.2. DIAGNOSTIC SPECIFICITY

---

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 95,8%.

### 12. 10.3. DIAGNOSTIC SENSITIVITY

---

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 87,9%

#### 12.1. Cross Reactivity

Cross reactions with antibodies against other worms and parasites cannot be excluded.

Note: The results refer to the groups of samples investigated; these are not guaranteed specifications.

### 13. LIMITATIONS OF THE PROCEDURE

---

Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values.

### 14. PRECAUTIONS AND WARNINGS

---

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of the in vitro diagnostic medical devices is intended by the manufacturer to secure suitability, performances and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately to the bottom of wells.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.
- The concentrations of the hazardous materials mentioned in point 4.1. are very low. Therefore there is hardly any toxicological risk. Nevertheless rinse with plenty of water upon contact with eyes, skin or mucous membranes and consult a doctor in case of irritations. All solutions should be handled with adequate care.

#### 14.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Strongyloides-ELISA

## Test Preparation

Prepare reagents and samples as described.  
 Establish the distribution and identification plan for all specimens and controls on the result sheet supplied in the kit.  
 Select the required number of microtiter strips or wells and insert them into the holder.

## Assay Procedure

	Substrate Blank (e. g. A1)	Negative Control	Cut-off Control	Positive Control	Sample (diluted 1+100)
Negative Control	-	100 µl	-	-	-
Cut-off Control	-	-	100 µl	-	-
Positive Control	-	-	-	100 µl	-
Sample (diluted 1+100)	-	-	-	-	100 µl
Cover wells with foil supplied in the kit <b>Incubate for 1 h at 37 °C</b> Wash each well three times with 300 µl of Washing Solution					
Conjugate	-	100 µl	100 µl	100 µl	100 µl
Cover wells with foil supplied in the kit <b>Incubate for 30 min at room temperature</b> Wash each well three times with 300 µl of Washing Solution					
TMB Substrate	100 µl	100 µl	100 µl	100 µl	100 µl
<b>Incubate for exactly 15 min at room temperature in the dark</b>					
Stop Solution	100 µl	100 µl	100 µl	100 µl	100 µl
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

## NovaTec Immundiagnostica GmbH

### Technologie & Waldpark

Waldstr. 23 A6  
 D-63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760 Fax: +49 (0) 6074-487629

Email: [info@NovaTec-ID.com](mailto:info@NovaTec-ID.com)

Internet: [www.NovaTec-ID.com](http://www.NovaTec-ID.com)

[STRO0690engl,dt - 19122014-AL]