SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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
Tab.4: LightCycler® 480II real-time PCR Profil

Program	Cycles / Analysis Mode	Temperature targets			
		Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp rate [°C/s]
Initial Denat.	1 / none	95	none	00:00:30	4.4
Cycling	5 / Quantification	95	none	00:00:15	4.4
		60	single	00:00:30	2.2
TM-Analyse	1 / Melting Curves	95	none	00:00:01	4.4
		50	none	00:00:30	2.2
		70	continuous		0.14 (Acquisitions per °C = 1)
Cooling	1 / none	50	none	00:00:01	2.2

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The screenshot displays the LightCycler 480 Software interface with the following components:

- Header:** LightCycler® 480 Software release 1.5.0 SP1. Instrument: 25586 / Not Connected. Database: Research Database (Research). User: System Admin.
- Navigation:** Experiment, Subset Editor, Sample Editor, Analysis, Report, Sum.
- Setup:** Detection Format: R-Biopharm. Block Size: 96. Plate ID: [blank]. Reaction Volume: 20.
- Programs Table:**

Program Name	Cycles	Analysis Mode
Initial Denaturierung	1	None
Cycling	5	Quantification
TM	1	Melting Curves
Cooling	1	None
- Initial Denaturierung Temperature Targets Table:**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step Size (°C)	Step Delay (cycles)
95	None	00:00:30	4.4	0	0	0	0
- Overview Graph:** A line graph showing Temperature (°C) on the y-axis (30 to 100) versus Estimated Time (hh:mm:ss) on the x-axis (0:00 to 9:25). The profile shows a series of heating and cooling cycles, with a final ramp up to 70°C.
- Buttons:** Apply Template, End Program, + 10 Cycles, Start Run.
- Warning Message:** 10/21/2011 2:51:46 PM An error occurred when trying to broadcast to TIDFormeProperty: Invalid state change from 999 TPersistentClean as a result of invoking Dirty! Object is locked by user System Admin!
- Taskbar:** Windows taskbar showing the start button, system tray, and open applications (LightCycler® 480, Rotor-Gene Q Series, Dokument1 - Microsoft...).

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The screenshot shows the LightCycler 480 Software interface. The 'Subset Editor' button is highlighted with a red circle. The interface includes a 'Run Protocol' section with 'Detection Format' set to 'R-Biopharm'. Below this is a 'Programs' table with the following data:

Program Name	Cycles	Analysis Mode
Initial Denaturierung	1	None
Cycling	5	Quantification
TH	1	Melting Curves
Cooling	1	None

Below the programs table is an 'Initial Denaturierung Temperature Targets' table with the following data:

Target (°C)	Acquisition Mode	Hold (h:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step Size (°C)	Step Delay (cycles)
95	None	00:00:30	4.4	0	0	0	0

At the bottom of the interface is a graph titled 'Overview' showing Temperature (°C) on the y-axis (ranging from 30 to 100) and Estimated Time (h:mm:ss) on the x-axis (ranging from 0:00 to 9:25). The graph shows a series of temperature steps: a ramp up to 95°C, a hold, a ramp down to 60°C, a ramp up to 95°C, a hold, a ramp down to 60°C, and so on, ending with a final ramp up to 95°C.

A warning message is displayed at the bottom of the software window: "Warning 10/21/2011 2:51:46 PM An error occurred when trying to broadcast to TIDFormsProperty: Invalid state change from 999 TPersistentClean as a result of invoking Dirty! Object is locked by user System Admin!"

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