SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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



REF PG4915

1. Intended use

For research use only. Not intended for diagnostic procedures.

2. Test principle

The RIDA®GENE Monkeypox Virus RUO test is a multiplex real-time PCR for the detection of monkeypox virus (MPXV) DNA.

After DNA isolation, the specific gene fragment (if present) for the monkeypox virus (*A29L*) is amplified.

The target sequence to be amplified is detected using a hydrolysis probe that is labeled with a quencher at one end and a fluorescent reporter dye (fluorophore) at the other. The probes hybridize to the amplicon in the presence of a target sequence. During extension, **Taq-Polymerase** separates the reporter from the quencher. The reporter emits a fluorescent signal that is detected by the optical unit of a real-time PCR instrument. The fluorescent signal increases with the quantity of formed amplicons. The RIDA®GENE Monkeypox Virus RUO multiplex real-time PCR test contains an **Internal Control DNA** (ICD) to be able to control sample preparation and/or potential PCR inhibition.

3. Reagents provided

Tab. 1: Reagents provided (The reagents provided in the kit are sufficient for 100 determinations.)

Kit code	Reagent	Amount		Lid color
1	Reaction Mix	2 ×	1050 µL	yellow
2	Taq-Polymerase	1 ×	80 µL	red
D	Internal Control DNA	2 ×	1700 µL	orange
N	No Template Control	1 ×	450 µL	white
P	Positive Control	1 ×	200 µL	blue

4. Storage instructions

- Please follow the handling guidelines in Table 2 and store the kit directly after use according to the information specified.
- All reagents must be stored away from light at -16 °C to -28 °C and, if unopened, can be used until the expiration date printed on the label. The quality guarantee is no longer valid after the expiration date.
- All reagents should be carefully thawed prior to use (e.g., in a refrigerator at 2 °C - 8 °C).
- Repeated freezing/thawing of up to 5 times does not have an impact on the test property (if necessary, create aliquots after the first thaw and refreeze reagents immediately).
- Cool all reagents appropriately during PCR preparation (2 °C - 8 °C).

Tab. 2: Storage conditions and information

	Storage temperature	Maximum storage time
unopened	-16 °C to -28 °C	Can be used until the printed expiration date
opened	-16 °C to -28 °C	5 freeze-thaw cycles

5. Reagents required but not provided

5.1 Laboratory equipment

The following equipment is needed for performing the RIDA®GENE Monkeypox Virus RUO test:

Equipment
Extraction platform: <ul style="list-style-type: none">• RTP® Pathogen Kit (Invitex Molecular)• Maxwell® RSC (Promega)• MagNA Pure 96 instrument (Roche)
Real-time PCR instrument: <ul style="list-style-type: none">• RIDA®CYCLER (R-Biopharm AG)• LightCycler® 480 II (Roche)• ABI 7500 Fast Dx (Applied Biosystems)• CFX96™ Dx (Bio-Rad)• Rotor-Gene Q (QIAGEN) <p>Note: When using Rotor-Gene Q (QIAGEN), use only 0.1-mL reaction vials.</p>
RIDA®GENE Color Compensation Kit IV (PG0004) (R-Biopharm)*
Real-time PCR consumables (plates (low profile, white wells, clear frame), reaction vials, foils)
Centrifuge with rotor for plates
Vortexer
Pipettes (0.5 - 20 µL, 20 - 200 µL, 100 - 1000 µL)
Pipette tips with filters
Powder-free disposable gloves

* Necessary for applications using LightCycler® 480 II (Roche)

The RIDA®GENE Monkeypox Virus RUO kit can be used in conjunction with compatible workflows. Alternative nucleic acid extraction procedures and real-time PCR instruments must be verified/validated by the user. Please contact R-Biopharm AG at pcr@r-biopharm.de to check the compatibility.

6. Warnings and precautions for the users

For research use only.

This test must be carried out only by qualified laboratory personnel. The guidelines for working in medical laboratories must be observed.

Always adhere strictly to the operating manual when carrying out this test.

Do not pipette samples or reagents using your mouth. Avoid contact with broken skin and mucous membranes.

Wear personal protective equipment (appropriate gloves, lab coat, safety glasses) when handling reagents and samples, and wash hands after completing the test.

Do not smoke, eat, or drink in areas where samples are handled.

Separate rooms, special clothing, and instruments for extraction, PCR preparation, and PCR must be used to prevent cross-contamination and false-positive results.

Avoid contaminating the samples and components of the kit with microbes and nucleases (DNase/RNase).

Clinical samples must be viewed as potentially infectious and must be disposed of appropriately, like all reagents and materials that come into contact with potentially infectious samples.

Do not exchange or mix the components (Reaction Mix, Taq-Polymerase, Internal Control DNA, Positive Control, No Template Control) of one lot with the components of another lot.

Do not use the test kit after the expiration date. Users are responsible for the proper disposal of all reagents and materials after use. For disposal, please adhere to national regulations.

7. Collection and storage of samples

7.1 DNA preparation from human swab samples

A commercially available nucleic acid extraction kit (e.g., RTP[®] Pathogen Kit, Invitek Molecular) or nucleic acid extraction system (e.g., Maxwell[®] RSC, Promega) is recommended for DNA preparation from human swab samples. The manufacturer's instructions must be observed.

The RIDA[®]GENE Monkeypox Virus RUO test contains an **Internal Control DNA** that indicates potential PCR inhibition, checks the integrity of the reagents, and confirms successful nucleic acid extraction. The **Internal Control DNA** can be used either solely as an inhibition control or as a process control (extraction and inhibition control).

When the **Internal Control DNA** is used only as an inhibition control for amplification, 1 µL of the **Internal Control DNA** must be added to the Master Mix (see Tab. 4).

When the **Internal Control DNA** is used as an extraction control for sample preparation **and** as an inhibition control for amplification, then 20 µL of the **Internal Control DNA** must be used during extraction.

The **Internal Control DNA** should be added to the sample lysis buffer mix and should **not** be added directly to the sample material. It is recommended to pipette 1 µL of the **Internal Control DNA** per reaction into the PCR mix of the negative control and of the positive control.

8. Test procedure

8.1 Preparation of the Master Mix

The total number of the reactions needed for PCR (samples and control reactions) must be calculated. One positive control and one negative control must be included in each test run.

Adding an additional 10 % volume to the Master Mix is recommended in order to compensate for any pipetting loss (see Tab. 3). Prior to use, thaw the **Reaction Mix**, the **Taq-Polymerase**, the **Positive Control**, and the **No Template Control**, vortex (except for Taq-Polymerase), and briefly centrifuge. Always cool reagents appropriately during work steps (2 °C - 8 °C).

Tab. 3: Example of the calculation and production of the Master Mix for 10 reactions

Kit code	Components of the Master Mix	Quantity per reaction	10 reactions (plus 10 %)
1	Reaction Mix	19.3 µL	212.3 µL
2	Taq-Polymerase	0.7 µL	7.7 µL
	Total	20 µL	220 µL

Mix the Master Mix and then centrifuge for short time.

Tab. 4: Example of the calculation and production of the Master Mix for 10 reactions (ICD only as inhibition control)

Kit code	Components of the Master Mix	Quantity per reaction	10 reactions (plus 10 %)
1	Reaction Mix	19.3 µL	212.3 µL
2	Taq-Polymerase	0.7 µL	7.7 µL
D	Internal Control DNA	1.0 µL	11 µL
	Total	21.0 µL	231.0 µL

Mix the Master Mix and then centrifuge for short time.

8.2 Preparation of the PCR mix

Pipette 20 µL of the Master Mix into each reaction vial (vials/plates).

Negative control: Pipette 5 µL of the **No Template Control** into the pre-pipetted Master Mix.

Note: When the **Internal Control DNA** is used as an extraction control for sample preparation and as an inhibition control, it is recommended to add 1 µL of the **Internal Control DNA** to each PCR mix of the negative control.

Samples: Add 5 µL eluate to the pre-pipetted Master Mix.

Positive control: Add 5 µL of the **Positive Control** to the pre-pipetted Master Mix.

Note: When the **Internal Control DNA** is used as an extraction control for sample preparation and as an inhibition control, it is recommended to add 1 µL of the **Internal Control DNA** to each PCR mix of the positive control.

Seal the reaction vials or plates, briefly centrifuge at slow speed, and transfer into the real-time PCR instrument. Start PCR according to PCR instrument setup (see Tab. 5, Tab. 6, Tab. 7).

8.3 Device settings

8.3.1 Universal real-time PCR profile

To harmonize the RIDA®GENE assays, the RIDA®GENE Monkeypox Virus RUO assay was verified in the universal profile. This makes it possible to combine DNA and RNA assays with each other.

Tab. 5: Universal real-time PCR profile for LightCycler® series and RIDA®CYCLER

<u>Reverse transcription</u>	10 min, 58 °C
Initial denaturation	1 min, 95 °C
Cycles	45 cycles
<u>PCR</u> Denaturation	10 sec, 95 °C
Annealing/extension	15 sec, 60 °C
Temperature transition rate/ ramp rate	Maximum

Note: **Annealing and extension take place in the same step.**

Tab. 6: Universal real-time PCR profile for ABI 7500 Fast Dx, Rotor-Gene Q, and CFX96™ Dx

<u>Reverse transcription</u>	10 min, 58 °C
Initial denaturation	1 min, 95 °C
Cycles	45 cycles
<u>PCR</u> Denaturation	15 sec, 95 °C
Annealing/extension	30 sec, 60 °C
Temperature transition rate/ ramp rate	Maximum

Note: **Annealing and extension take place in the same step.**

Note: **The universal real-time PCR profile can also be used for RNA tests if RIDA®GENE DNA and RIDA®GENE RNA real-time PCR tests are combined in one run.**

8.4 Detection channel setting

Tab. 7: Selection of appropriate detection channels

Real-time PCR instrument	Detection	Detection channel	Note
R-Biopharm AG RIDA®CYCLER	MPXV	Green	-
	ICD	Yellow	
Roche LightCycler® 480 II	MPXV	465/510	RIDA®GENE Color Compensation Kit IV (PG0004) is required.
	ICD	533/580	
Applied Biosystems ABI 7500 Fast Dx	MPXV	FAM	-
	ICD	VIC	
Bio-Rad CFX96™ Dx	MPXV	FAM	-
	ICD	VIC	
QIAGEN Rotor-Gene Q	MPXV	Green	The gain settings must be set to 5 (factory default) for all channels.
	ICD	Yellow	

9. Quality control

Samples are evaluated using the analysis software of the respective real-time PCR instrument according to the manufacturer's instructions. Negative and positive controls must show the correct results (see Tab. 8).

The **Positive Control** is present in a concentration of 10^3 copies/ μ L. It is used in a total quantity of 5×10^3 copies in every PCR run.

Tab. 8: A valid PCR run must meet the following conditions:

Sample	Result	ICD Ct	Target gene Ct
Positive control	Positive	N/A *1	See Certificate of Analysis
Negative control	Negative	Ct > 20	0

*1 A Ct value for the ICD is not needed to obtain a positive result of the positive control.

If the positive control is not within the specified Ct range but the negative control is valid, all reactions need to be reanalyzed, including the controls.

If the negative control is not negative, but the positive control is valid, all reactions need to be reanalyzed, including the controls.

If the specified values are not met, check the following items before repeating the test:

- Expiration date of the reagents used
- Functionality of the equipment being used
- Correct test procedure

10. Result interpretation

The sample results are interpreted according to Table 9.

Tab. 9: Result interpretation*

Detection of		
MPXV	ICD	Result
+	+/-	Monkeypox virus detectable* ¹
-	+	Target gene not detectable
-	-	Invalid

* + = positive
- = negative

*¹ see Limitations of the method, No. 2

A sample is positive if the sample DNA and the **Internal Control DNA** show an amplification signal in the detection system.

A sample is also positive if the sample DNA shows an amplification signal, but no amplification signal can be seen for the **Internal Control DNA** in the detection system. Detecting the **Internal Control DNA** is not necessary in this case because high amplicon concentrations can result in a weak or absent signal of the **Internal Control DNA**.

A sample is negative if the sample DNA does not show an amplification signal, but an amplification signal is visible for the **Internal Control DNA** in the detection system. Inhibition of the PCR reaction can be ruled out by the detection of the **Internal Control DNA**.

A sample is invalid if the sample DNA and the **Internal Control DNA** do not show an amplification signal in the detection system. PCR inhibitors are present in the sample or an error occurred during the extraction process.

11. Limitations of the method

1. This test is used to detect monkeypox virus DNA.
2. The primer and probe system used can have cross-reactivities with other orthopoxviruses, such as variola virus, vaccinia virus, cowpox virus, camelpox virus, ectromelia virus, taterapox virus, Uasin Gishu disease viruses, raccoonpox virus, volepox virus, and skunkpox virus.
3. Improper specimen sampling, transport, storage, and handling or a pathogen load below the test's analytical sensitivity can lead to false-negative results.
4. The presence of PCR inhibitors can lead to non-evaluable results.
5. Mutations or polymorphisms in the primer or probe binding sites can interfere with the detection of new or unknown variants and can lead to false-negative results using the RIDA[®]GENE Monkeypox Virus RUO test.
6. As with all PCR-based tests, extremely low concentrations of the target sequences under the limit of detection (LoD) can be detected. The results obtained are not always reproducible.
7. A positive test result does not necessarily indicate the presence of viable organisms. A positive result indicates that the target gene (*A29L*) is present.
8. This assay should be performed in compliance with the regulation on good laboratory practice (GLP). Users must precisely follow the manufacturer's instructions when performing the test.

12. Performance characteristics

Not applicable

13. Version history






Version number	Section and designation
2022-05-25	Release version

14. Explanation of symbols

General symbols

	For research use only
	Follow instructions for use
	Batch number
	Use before
	Storage temperature
	Item number
	Number of tests
	Date of manufacture
	Manufacturer

Test-specific symbols

	Reaction Mix
	Taq Polymerase
	Extraction/inhibition control
	Negative control
	Positive control