SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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REF PG6835RUO

1. Intended use

For research use only. Not intended for diagnostic procedures.

The RIDA®GENE SARS-CoV-2 Lineage II RUO test is a multiplex real-time RT-PCR for the detection of SARS-CoV-2 mutations ((single) nucleotide polymorphisms).

2. Summary and explanation of the test

At the end of December 2019 in the Chinese metropolis of Wuhan, numerous cases of pneumonia of unknown cause occurred.¹ At the beginning of January 2020, Chinese authorities identified a novel coronavirus (SARS-CoV-2) as the cause of these illnesses.¹ The disease caused by SARS-CoV-2 was officially named COVID-19 (“coronavirus disease 2019”) and is transmitted from human to human.² Due to this pathogen’s novelty, the epidemic rapidly evolved into a pandemic.

There have already been 170,363,852 cases recorded worldwide (as of June 01, 2021).³ At the end of January 2020, the first cases were confirmed in Germany as well. Germany has had 3.681.126 cases so far (as of June 01, 2021).⁴ The WHO declared a Public Health Emergency of International Concern on 2020-01-31.¹

There has been an increase in COVID-19 cases in the UK since September, leading to increased epidemiological and virological investigations. Analysis of the viral genome sequences revealed that a large proportion of the cases belonged to a new phylogenetic cluster. The new variant is characterized by several spike protein mutations (deletion 69-70, deletion 144, N501Y, A570D, D614G, P681H, T716I, S982A, D1118H) as well as mutations in other genomic regions (ORF1ab, ORF8, N-Gene). There are initial indications (limited data) that this new SARS-CoV-2 variant B.1.1.7 (Alpha) is associated with increased case mortality. Furthermore, initial studies have shown that it is even easier transmissible from person to person than previously circulating variants. Infections with the above-mentioned variants were identified in Germany since December 2020. Some cases with the new variant have been reported in Denmark, the Netherlands and Belgium.^{5,6,7}

B.1.1.7 with E484K is a special form of the variant that has been detected several times in Great Britain, but is currently rarely found in Germany. This variant has an additional mutation in the S protein (E484K), which also occurs in the variants

B.1.351 (Beta) und P.1 (Gamma) and makes the virus less sensitive to neutralizing antibodies that have already been formed. Therefore, it is suspected that the currently available vaccines against this variant may have a lower efficacy.⁷

Another virus variant, B.1.351 (Beta), was first reported in South Africa in December 2020. Initial tests showed for this variant, too, a potentially higher transmissibility. Initial studies also suggest that protection by neutralizing antibodies against this variant could be reduced in people who were previously infected with the original variant or who received a vaccine based on the original variant.⁷

The SARS-CoV-2 variant P.1 (Gamma) circulated for the first time in the Brazilian state of Amazonas. Its changes are similar to the South African variant. An increased transmissibility and a possible reduction in the effectiveness of neutralizing antibodies in convalescent patients or vaccinated persons are discussed for this variant.⁷

Another virus variant B.1.617 (consisting of subvariants B.1.617.1 (Kappa), B.1.617.2 (Delta) and B.1.617.3) was first detected in the Indian state of Maharashtra and is spreading strongly there. Subvariant B.1.617.2 (Delta) is currently spreading strongly in the United Kingdom. B.1.617 is characterized by mutations associated with reduced efficacy of the immune response and mutations that may increase the transmissibility of the virus.⁷ The ecdc currently classifies variant B.1.617.2 (Delta) as a variant of concern. In contrast, variants B.1.617.1 & B.1.616.3 are currently classified as variants of interest.⁸

3. Test principle

RIDA®GENE SARS-CoV-2 Lineage II RUO is a multiplex real-time RT-PCR for the detection of SARS-CoV-2 mutations ((single) nucleotide polymorphisms).

Detection is done in a one-step real-time RT-PCR format: reverse transcription (RT) and subsequent PCR take place in one reaction vial. In the process, the isolated RNA is transcribed into cDNA with the help of a reverse transcriptase. Next, real-time PCR is used to amplify the specific gene fragments for the virus variants Delta B.1.617.2 (S gene T478K), Kappa B.1.617.1 & B.1.617.3 (S gene E484Q) and Beta B.1.351 (S gene E484K) (NCBI database). The amplified target sequences are detected using hydrolysis probes that are labeled with a quencher at one end and a fluorescent reporter dye (fluorophore) at the other. The probes hybridize to the amplicon in the presence of a target sequence. During extension, Taq-Polymerase separates the reporter from the quencher. The reporter emits a fluorescent signal that is detected by the optical unit of a real-time PCR device. The fluorescent signal increases with the quantity of formed amplicons.

4. Reagents provided

Table 1: Reagents provided (The reagents provided in the kit are sufficient for 100 determinations)

Kit code	Reagent	Amount		Lid color
1	Reaction Mix	2x	1050 µl	Yellow
2	Enzyme Mix	1x	80 µl	Red
N	No Template Control	1x	450 µl	White
P	Positive Control	1x	200 µl	Blue

5. Storage instructions

- All reagents must be stored away from light at -20 °C and, if unopened, can be used until the expiration date printed on the label. After the expiration date, the quality guarantee is no longer valid.
- All reagents should be carefully thawed prior to use (e.g., in a refrigerator at 2 °C to 8 °C).
- Repeated freezing/thawing of up to 5 times does not have an impact on the test property (if necessary, create aliquots after the first thaw and refreeze reagents immediately).
- Cool all reagents appropriately during PCR preparation (2 °C to 8 °C).

6. Reagents required but not provided

The RIDA®GENE SARS-CoV-2 Lineage II RUO multiplex real-time RT-PCR test can be used with the following extraction platforms and real-time PCR devices:

Table 2: Necessary equipment

Extraction platforms	
R-Biopharm	RIDA® Xtract
Promega	Maxwell® RSC
Roche	MagNA Pure 96
Real-time PCR devices	
R-Biopharm	RIDA®CYCLER
Roche	LightCycler® 480 II
Applied Biosystems	ABI 7500 Fast Dx
Bio-Rad	CFX96™ Dx
QIAGEN	Rotor-Gene Q

Note: When using Rotor-Gene Q (QIAGEN), use only 0.1 ml reaction vials.

Should you have to use other extraction procedures or real-time PCR devices, please contact R-Biopharm to check the compatibility at mdx@r-biopharm.de.

- RIDA®GENE Color Compensation Kit IV (PG0004) when using LightCycler® 480II
- Real-time PCR consumables (plates, reaction vials, films)
- Centrifuge with rotor for reaction vials or plates
- Vortexer
- Pipettes (0.5 to 20 µl, 20 to 200 µl, 100 to 1,000 µl)
- Pipette tips with filters
- Powder-free disposable gloves
- PCR water (nuclease-free)

7. Warnings and precautions for the users

For research use only.

- This test must be carried out only by trained laboratory personnel. The guidelines for working in medical laboratories must be followed.
- Always adhere strictly to the instructions for use when carrying out this test.
- Do not pipette samples or reagents using your mouth. Avoid contact with broken skin and mucous membranes.
- Wear personal protective equipment (appropriate gloves, lab coat, safety glasses) when handling reagents and samples, and wash hands after completing the test.
- Do not smoke, eat, or drink in areas where samples are handled.
- Ensure that the extraction, PCR preparation, and PCR are carried out in different rooms in order to avoid cross-contaminations.
- Clinical samples must be viewed as potentially infectious and must be disposed of appropriately, like all reagents and materials that come into contact with potentially infectious samples.
- Dispose of test kit once the expiration date has lapsed.
- Users are responsible for proper disposal of all reagents and materials after use. For disposal, please adhere to national regulations.

For further details, see the safety data sheets (SDS) on R-Biopharm's website.

8. Collection and storage of samples

8.1 RNA preparation from human respiratory samples

A commercially available nucleic acid extraction kit (e.g., RIDA® Xtract [R-Biopharm]) or nucleic acid extraction system (e.g., Maxwell® RSC [Promega]) is recommended for RNA preparation from human respiratory samples. The manufacturer's instructions must be observed.

9. Test procedure

9.1 Master Mix preparation

The total number of the reactions needed for PCR (samples and control reactions) must be calculated. One positive control and one negative control must be included in each test run.

Adding an additional 10% volume to the master mix is recommended in order to balance out the pipette loss (see Table 3, Table 4). Prior to use, thaw the **Reaction Mix**, the **Enzyme Mix**, the **Positive Control** and the **No Template Control**, mix thoroughly (except for the enzyme mix), and centrifuge for a short time. Always cool reagents appropriately during the work steps (2 °C to 8 °C).

Table 3: Example of the calculation and preparation of the master mix for 10 reactions

Kit code	Components of the master mix	Quantity per reaction	10 reactions (plus 10%)
1	Reaction Mix	19.3 µl	212.3 µl
2	Enzyme Mix	0.7 µl	7.7 µl
	Total	20 µl	220 µl

Mix the master mix and then centrifuge for short time.

9.2 Preparation of the PCR Mix

Pipette 20 µl of the master mix into each reaction vial (vial/plate).

Negative control: Pipette 5 µl of the **No Template Control** into the pre-pipetted master mix.

Samples: Add 5 µl eluate to the pre-pipetted master mix.

Positive control: Add 5 µl of the **Positive Control** to the pre-pipetted master mix.

Seal the reaction vials or plates, briefly centrifuge at slow speed, and transfer into the real-time PCR device. Start PCR according to PCR instrument set-up (see Table 5, Table 6, Table 7).

9.3 PCR instrument set-up

9.3.1 Universal real-time RT-PCR profile

Table 5: Universal real-time RT-PCR profile for LightCycler® series and RIDA®CYCLER

<u>Reverse transcription</u>	10 min, 58 °C
Initial denaturation	1 min, 95 °C
Cycles	45 cycles
<u>PCR</u> Denaturation	10 sec, 95 °C
Annealing/Extension	15 sec, 60 °C
Temperature transition rate/ ramp rate	Maximum

Note: **Annealing and extension take place in the same step.**

Table 6: Universal real-time RT-PCR profile for ABI 7500 Fast Dx, Rotor-Gene Q, and CFX96™ Dx

<u>Reverse transcription</u>	10 min, 58 °C
Initial denaturation	1 min, 95 °C
Cycles	45 cycles
<u>PCR</u> Denaturation	15 sec, 95 °C
Annealing/Extension	30 sec, 60 °C
Temperature transition rate/ ramp rate	Maximum

Note: **Annealing and extension take place in the same step.**

Note: **The universal real-time PCR profile can also be used for DNA tests if RIDA®GENE DNA and RIDA®GENE RNA real-time PCR tests are combined in one run.**

9.4 Detection channel setting

Table 7: Selection of appropriate detection channels

Real-time PCR device	Detection	Detection channel	Comment
R-Biopharm RIDA®CYCLER	SARS-CoV-2 S gene T478K	Green	-
	SARS-CoV-2 reference S gene	Yellow	
	SARS-CoV-2 S gene E484Q	Orange	
	SARS-CoV-2 S gene E484K	Red	
Roche LightCycler® 480 II	SARS-CoV-2 S gene T478K	465/510	RIDA®GENE Color Compensation Kit IV (PG0004) is required.
	SARS-CoV-2 reference S gene	533/580	
	SARS-CoV-2 S gene E484Q	533/610	
	SARS-CoV-2 S gene E484K	618/660	
ABI 7500 Fast Dx	SARS-CoV-2 S gene T478K	FAM	Stellen Sie den passiven Referenzfarbstoff ROX auf none
	SARS-CoV-2 reference S gene	VIC	
	SARS-CoV-2 S gene E484Q	ROX	
	SARS-CoV-2 S gene E484K	Cy5	
Bio-Rad CFX96™ Dx	SARS-CoV-2 S gene T478K	FAM	-
	SARS-CoV-2 reference S gene	VIC	
	SARS-CoV-2 S gene E484Q	ROX	
	SARS-CoV-2 S gene E484K	Cy5	
Qiagen Rotor- Gene Q	SARS-CoV-2 S gene T478K	Green	The gain settings must be set to 5 (factory default) for all channels.
	SARS-CoV-2 reference S gene	Yellow	
	SARS-CoV-2 S gene E484Q	Orange	
	SARS-CoV-2 S gene E484K	Red	

10. Quality control

Samples are evaluated using the analysis software of the respective real-time PCR device according to the manufacturer's instructions. Negative and positive controls must show the correct results (see Table 8).

The **Positive Control** is present in a concentration of 10^3 copies/ μ l. It is used in a total quantity of 5×10^3 copies in every PCR run.

Table 8: A valid PCR run must meet the following conditions:

Sample	Result	Target gene Ct
Positive control	Positive	See Quality Assurance Certificate
Negative control	Negative	Not detectable

The positive and negative controls are valid when they meet the conditions specified in the table. The Ct range for the positive control is specified on the Quality Assurance Certificate included with the product. If one of the two controls does not meet the conditions for a valid run, all the reactions need to be re-analyzed, including the controls.

If the specified values are not met, check the following items before repeating the test:

- Expiration date of the reagents used
- Functionality of the equipment being used
- Correct test procedure

11. Result interpretation

The results interpretation is done according to Table 9.

Table 9: Sample interpretation

Detection of				Result SARS-CoV-2 VOC & VOI ⁸
SARS-CoV-2 S gene T478K	SARS-CoV-2 reference S gene	SARS-CoV-2 S gene E484Q	SARS-CoV-2 S gene E484K	
negative	positive	negative	negative	SARS-CoV-2 detectable
positive	positive	negative	negative	SARS-CoV-2 B.1.617.2 (Delta) detectable
negative	positive	positive	negative	SARS-CoV-2 B.1.617.1 (Kappa) / B.1.617.3 detectable
negative	positive	negative	positive	SARS-CoV-2 B.1.1.7 + E484K / B.1.351 (Beta) / P.1 (Gamma) / B.1.525 (Eta) / P.3 (Theta) / B.1.620 / B.1.621
negative	negative	negative	negative	Target gene not detectable

VOC: Variants of concern; VOI: Variants of interest

A sample is positive if the sample RNA shows an amplification signal in the detection system.

A sample is negative if the sample RNA does not show an amplification signal.

12. Limitations of the method

1. This test is intended only for human respiratory samples.
2. Improper specimen sampling, transport, storage, and handling or a pathogen load below the test's analytical sensitivity can lead to false negative results.
3. The presence of PCR inhibitors can lead to non-evaluable results.
4. Mutations or polymorphisms in the primer or probe binding sites can interfere with the detection of new or unknown variants and can lead to false negative results using the RIDA[®]GENE SARS-CoV-2 Lineage II RUO test.
5. As with all PCR-based tests, extremely low concentrations of the target sequences under the limit of detection (LoD) can be detected. The results obtained are not always reproducible.

13. Performance characteristics


Not applicable

14. Version history


Version number	Section and designation
2021-06-01	Release version


15. Explanation of symbols


General symbols

	For research use only
	Consult instructions for use
	Batch number
	Use before
	Store at
	Item number
	Number of tests
	Date of manufacture
	Manufacturer

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