



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



# Instructions For Use

## PTA-IFU

205 South 600 West Logan, Utah 84323, U.S.A. – Tel. (800) 729-8350 – Tel. (435) 755-9848 – Fax (435) 755-0015 – www.scytek.com Rev. 4, 7/25/2022

### PTAH Stain Kit

(Phosphotungstic Acid Hematoxylin)

#### Description and Principle

The PTAH Stain Kit is intended for use in the histological visualization of collagen, striated muscle, glial fibers and collagen without using Zenker's Fixative with Mercuric Chloride as a mordant. This kit may be used on formalin-fixed, paraffin-embedded sections.

Zinc chloride acts to mordant the tissue for trichrome staining and is a safe and effective replacement for traditionally used mercuric chloride treatment. Aldehydes are formed by the oxidizing action of ferric ammonium sulfate. Hematoxylin and phosphotungstic form a blue dye lake staining fibrin, muscle, glial fibers, and nuclei. Collagen is dominantly stained by phosphotungstic acid resulting in lighter orange, salmon, or red/brown

#### Expected Results

Fibrin, Striated Muscle, Glial Fibers:	Blue to Purple
Collagen:	Light Orange/Salmon to Brownish/Red
Nuclei:	Blue to Purple

#### Kit Contents

1. Zinc Chloride Solution (10%)
2. Ferric Ammonium Sulfate
3. PTAH Solution

#### Storage

18-25° C  
18-25° C  
18-25° C

#### Suggested Controls (not provided)

Striated Muscle

#### Uses/Limitations

For In-Vitro Diagnostic use only.

Do not use if reagents become cloudy or precipitate

Do not use past expiration date.

Use caution when handling reagents.

Non-Sterile

Intended for FFPE sections cut at 5-10µm.

This procedure has not been optimized for frozen sections.

Frozen sections may require protocol modification.

#### Storage


Store kit and all components at room temperature (18-25°C).

#### Safety and Precautions

Please see current Safety Data Sheets (SDS) for this product and components GHS classification, pictograms, and full hazard/precautionary statements.

#### Procedure (Standard):

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Pour Zinc Chloride Solution (10%) into plastic staining jar and set in 60°C. water bath for 10 minutes to equilibrate temperature.
3. Place slide in warmed Zinc Chloride Solution (10%) and incubate for 20 minutes at 60°C.
4. During step 3, pour Ferric Ammonium Sulfate Aqueous Solution into a second plastic staining jar and set in 60° C water bath for 10 minutes to equilibrate temperature.
5. Rinse slide in running tap water for 1 minute.



Human Skeletal Muscle and Collagen at 10X stained with the PTAH Stain Kit (Phosphotungstic Acid Hematoxylin)

6. Rinse in distilled water for 1 minute.
  7. Place slide in warmed Ferric Ammonium Sulfate Aqueous Solution and incubate for 5 minutes at 60°C.
  8. During step 7, pour Phosphotungstic Acid Hematoxylin Solution into a third plastic staining jar and set in 60° C. water bath for 10 minutes to equilibrate temperature.
  9. Rinse slide in running tap water for 2 minutes.
  10. Rinse in distilled water for 1 minute.
  11. Place slide in warmed Phosphotungstic Acid Hematoxylin Solution and incubate for 60 minutes at 60°C.
  12. Differentiate section in 95% Reagent Alcohol. Check section using microscope for proper differentiation.
- Note:** Graded alcohols will remove some stain.
13. Dehydrate in 3 changes of Absolute Alcohol.
  14. Clear in 3 changes of fresh Xylene or Xylene Substitute, and mount in synthetic resin.

#### Procedure (Dropper):

**Equipment Needed:** 500 Watt Microwave Oven

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Place slide to fresh distilled water for 1 minute.
3. Fill a plastic staining jar with DI water. Place in microwave and heat until hot but not boiling.
4. Carefully lay slide across the top of the staining jar and apply 5 drops of Zinc Chloride Solution (10%). Heat in microwave for 10 seconds. Leave jar with slide in the microwave during the incubation period to maintain temperature.

5. Incubate slide for 15 minutes. **Note:** Water in staining jar will maintain reagent temperature during staining procedure.
6. Rinse slide in running tap water for 1 minute.
7. Rinse in distilled water for 1 minute.  
**Note:** During rinse step, reheat water in staining jar to hot but not boiling.
8. Carefully lay slide across the top of the staining jar and apply 5 drops of Ferric Ammonium Sulfate Solution and heat in microwabe for 10 seconds. Leave jar with slide in the microwave during the incubation period to better maintain temperature.
9. Incubate slide for 2 minutes.
10. Rinse slide in running tap water for 2 minutes.
11. Rinse in distilled water for 1 minute.  
**Note:** During rinse step, reheat water in staining jar to hot but not boiling.
12. Carefully lay slide across the top of the staining jar and apply 5 drops of Phosphotungstic Acid Hematoxylin Solution and heat in microwave for 10 seconds. Leave jar with slide in the microwave during the incubation period to better maintain temperature.
13. Incubate slide for 15 minutes.
14. Shake off Phosphotungstic Acid Hematoxylin Solution and repeat steps 13 and 14.
15. Differentiate section in 95% Reagent Alcohol. Check section using microscope for proper differentiation.
16. Dehydrate in 3 changes of Absolute Alcohol.
17. Clear in 3 changes of fresh Xylene or Xylene Substitute and mount in synthetic resin.

### **References**

1. Shapiro, S.H., Sohn, L.C.; Rapid Microwave Phosphotungstic Acid-Hematoxylin Stain for Paraffin and Glycol Methacrylate Sections; The Journal of Histotechnology; Volume 17, Number 2, June 1994, pages 125-126.



ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
435-755-9848  
U.S.A.



Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague, The Netherlands