



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ImmPRESS® HRP Goat Anti-Rat IgG (Mouse Adsorbed) Polymer Kit

Peroxidase

Cat. No. MP-7444

Storage Store reagents in original bottles at 2-8 °C. Do not freeze.

Description Instructions for Immunohistochemical staining using rat primary antibodies on mouse tissue.

The ImmPRESS Polymer Kit is an enzymatic, non-biotin amplification system that produces crisp, highly sensitive, specific staining with low background.

The reagents in the ImmPRESS Polymer Kit are ready-to-use—no mixing or titering is necessary to obtain optimal staining. Dilution of the reagents or changes in suggested incubation time may affect performance. The reagents are supplied in convenient dropper bottles.

A number of different wash buffers can be used with the ImmPRESS reagents. One of the most common is 10 mM sodium phosphate, pH 7.5, 0.9% saline (PBS). 0.1% Tween 20 detergent may be added to the wash buffer and is especially recommended for use with automated stainers.


Kit Components

Product Name	Volume
Normal Goat Serum, 2.5%	15 ml
ImmPRESS Goat Anti-Rat IgG (Mouse Adsorbed) Polymer Reagent	15 ml

ImmPRESS Polymer Kit (15 ml) will stain approximately 75-150 sections.

Staining Procedure

- For paraffin sections, deparaffinize and hydrate tissue sections through xylenes or other clearing agents and graded alcohol series.
For frozen sections or cell preparations fix with acetone or an appropriate fixative for the antigen under study, if necessary.
Wash for 5 minutes in tap water.
- If antigen unmasking is required, perform this procedure using a Vector® Antigen Unmasking Solution, Citrate-based, pH 6.0 (H-3300) or Tris-based, pH 9.0 (H-3301).



Human Tonsil (FFPE): CD45 (rat) detected with ImmPRESS Goat Anti-Rat IgG (Mouse Adsorbed) and ImmPACT® AEC Substrate (red). Counterstained with Vector Hematoxylin QS (blue).

- If quenching of endogenous peroxidase activity is required, incubate the sections in BLOXALL® Blocking Solution (SP-6000) for 10 minutes.
- Wash in buffer for 5 minutes.
- Incubate sections for 20 minutes with Normal Goat Serum, 2.5%.
- Tip off excess serum from sections.
- Incubate with rat primary antibody diluted in appropriate antibody diluent solution, such as diluted normal goat serum or BSA.
- Wash in buffer for 5 minutes.
- Incubate for 30 minutes with ImmPRESS Polymer Reagent.
- Wash for 2 x 5 minutes in buffer.
- Incubate in peroxidase substrate solution (not included) until desired stain intensity develops.
- Rinse sections in tap water.
- Counterstain (optional), clear and mount.

Detailed product listings, specifications, protocols and additional information is available on our website: vectorlabs.com