



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Vector® VIP

Peroxidase Substrate

Cat. No. SK-4600

Storage Store reagents in original bottles at 2-8 °C. Avoid storing reagents or working solution in strong direct light.

Description This kit contains all of the reagents necessary to prepare the substrate working solution. Vector VIP Substrate produces a purple reaction product in the presence of peroxidase (HRP) enzyme. Vector VIP Substrate is also suitable for darkfield and electron microscopy (EM).

Components

Product Name	Volume
Vector VIP Reagent 1	7 ml
Vector VIP Reagent 2	5 ml
Vector VIP Reagent 3	5 ml
Vector VIP Reagent 4	7.2 ml

Preparation of Substrate Working Solution

To 5 ml of PBS (10 mM sodium phosphate, 0.9% saline, pH 7.5)

- Add 3 drops (≈ 90 µl†) of Vector VIP Reagent 1
- Add 3 drops (≈ 75 µl†) of Vector VIP Reagent 2
- Add 3 drops (≈ 75 µl†) of Vector VIP Reagent 3
- Add 3 drops (≈ 120 µl†) of Vector VIP Reagent 4

Mix well before use. Use immediately.

† Drop volumes differ due to solvent compositions.


IMPORTANT: Little is known about the toxicity and carcinogenicity of the substrate kit components. Appropriate care should be exercised when using this reagent including gloves, eye protection, lab coats, and good laboratory procedures. Dispose in accordance with local regulations.

Instructions for Use

After incubation with a peroxidase (HRP) detection system, rinse sections in buffer. Incubate with the substrate working solution for 2-15 minutes. Optimal development times should be determined by the investigator.

Wash slides for 5 minutes in water.

Counterstain, if desired, with Vector® Hematoxylin, Vector Hematoxylin QS, or Vector Methyl Green. Note: Vector VIP is not compatible with all methyl green formulations.



Tonsil (Double Label): CD3 (m) detected with ImmPRESS® Reagent (HRP) Anti-Mouse IgG and Vector VIP Substrate (purple), Ki67 (m) detected with ImmPRESS Reagent (HR) Anti-Mouse IgG and Vector SG Substrate (blue-gray). No Counterstain.

For permanent, non-aqueous mounting: Dehydrate, clear, and coverslip using a non-aqueous mounting media, such as VectaMount® Mounting Medium (H-5000) or VectaMount® Express Mounting Medium (H-5700). Do not mount with aqueous mounting media.

Notes

We recommend using glass-distilled water in the preparation of substrate buffer. Deionized water may contain inhibitors of the peroxidase reaction. Solutions containing sodium azide or other inhibitors of peroxidase activity should not be used in diluting the peroxidase substrate.

Detailed product listings, specifications, protocols and additional information are available on our website: vectorlabs.com