



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

WestVision™ Peroxidase (HRP) Polymer

for Western Blot Detection

Cat. No.	WB-1000 Anti-Rabbit IgG WB-2000 Anti-Mouse IgG
Storage	2–8 °C
Unit Size	0.8 ml
Description	WestVision Peroxidase Polymer antibody conjugates are intended to be used in western blot applications to detect primary antibodies made in mouse (WB-2000) or rabbit (WB-1000) on nitrocellulose or PVDF membranes.

Reagents Not Provided

- Wash Buffer: Phosphate buffered saline with Tween® 20; 10 mM Na₂HPO₄, pH 7.5, 150 mM NaCl, and 0.1% Tween 20 (PBST).
- Blocking Solutions: We recommend WestVision™ Block and Diluent (SP-7000).
- Primary Antibody: Dilute primary antibody in appropriate blocking solution or PBST according to manufacturer's recommendation.
- Peroxidase Substrate: Chromogenic or chemiluminescent.


WestVision Reagent Preparation

- WestVision Working Solution. Dilute the WestVision reagent in PBST or the appropriate blocking solution. The diluent must not contain sodium azide. For chromogenic substrates dilute the WestVision reagent to 1:500–1:2,500. For chemiluminescent substrates dilute the WestVision reagent to 1:5,000–1:200,000.

Detection Protocol

Recommended volumes are based on the development of 100 cm² membranes. Volumes may be proportionally adjusted for blots of different sizes.

- Remove the blot from the transfer apparatus and block the membrane in 10 ml blocking solution for 30–60 minutes at room temperature with gentle agitation.
- Incubate the membrane in 10 ml of primary antibody solution for 30–60 minutes at room temperature or at 4°C overnight with gentle agitation (or for a time established to be optimal for the concentration of primary antibody used).
- Wash the membrane 3 x 5 minutes each in 10 ml PBST with gentle agitation.
- Incubate the membrane for 45 minutes at room temperature with 10 ml WestVision working solution with gentle agitation.



- Wash the membrane 3 x 5 minutes each in 10 ml PBST with gentle agitation.
- Prepare substrate working solution according to the substrate kit instructions. Following are protocols for signal development using either a chromogenic or chemiluminescent substrate.

Chromogenic signal development with TMB Substrate Kit, Peroxidase (SK-4400)

- Equilibrate membrane for 2 minutes in PBS in a clean vessel.
- Incubate membrane in the substrate working solution at room temperature with gentle agitation for 5 minutes or until suitable staining develops. Briefly rinse the membrane in PBS and air-dry.

Chemiluminescent signal development using DuoLuX® Chemiluminescent Substrate (SK-6604):

- Equilibrate membrane for 2 minutes in PBS in a clean vessel.
- Remove excess buffer by holding the membrane vertically and touching the edge of the membrane to absorbent paper.
- Place membrane target-side-up on plastic wrap on a level surface.
- Pipette 5 ml of DuoLuX Substrate working solution onto the membrane surface.
- Incubate for 5 minutes under subdued light. (If high background is present, rinse the membrane in PBS for a few seconds and remove excess as in Step 8.)
- Place the membrane between two pieces of plastic wrap or a clear sheet protector. Acquire image with the gel imager or expose the membrane to x-ray film for the appropriate time.

Detailed product listings, specifications, protocols and additional information is available on our website: vectorlabs.com