



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ZymoTaq™ DNA Polymerase

Cat. Nos. E2001 (50 Reactions)
E2002 (200 Reactions)

Storage: -20 °C



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Product Information

Features:

- Hot start DNA polymerase
- Allows for quick and easy setup at room temperature
- Robust product formation
- Reduces the occurrence of non-specific product formation
- Ideal for the amplification of *bisulfite-treated* DNA for methylation detection
- Compatible with real-time and quantitative PCR and suitable for TA-cloning

Description:

The ZymoTaq™ DNA Polymerase product contains all the reagents needed to perform “hot start” PCR. The inclusion of a heat-activated, thermal-stable DNA polymerase reduces primer dimer and non-specific product formation that can occur when performing conventional PCR. This unique product is specifically designed for the amplification of *bisulfite-treated* DNA for methylation detection, real-time and quantitative PCR that are SYBR Green and probe based. The product generates specific amplicons with little or no by-product formation. Simple and easy to use: heat at 95 °C for 10 minutes to initiate polymerization.

ZymoTaq™ DNA polymerase is a heat-activated, “hot start” polymerase that has 3'-terminal transferase activity. The addition of “A” overhangs to amplified DNA makes it ideal for use in TA-cloning.

Product Contents:

| | E2001 (50 Rxns.) | E2002 (200 Rxns.) | Storage Temp. |
|-----------------------------------|---------------------|----------------------|------------------|
| ZymoTaq™ DNA Polymerase (5 U/μl) | 20 μl | 80 μl | -20 °C |
| 2X Reaction Buffer | 2 x 625 μl | 8 x 625 μl | -20 °C |
| dNTP Mix (25 mM each dNTP) | 25 μl | 100 μl | -20 °C |
| DNase/RNase-Free H ₂ O | 2 x 1 ml | 5 x 1 ml | -20 °C |

Storage:

Store at -20 °C for up to 12 months. Avoid repeated freeze/thawing of reagents. Prolonged storage at -80 °C.

Enzyme Concentration:

Reaction conditions at 1X (50 μl total volume) will contain 2 U of ZymoTaq™ DNA polymerase. ZymoTaq™ DNA Polymerase is at a concentration of 5 U/μl.

Unit Definition:

One unit (U) enzyme is defined as the amount of enzyme required for the incorporation of 10 nM dNTPs into an acid-insoluble form in 30 minutes at 72 °C.

Proposed Reaction Setup (50 μl):


| Reagent | Volume | Final conc. |
|----------------------------------|----------|-------------------|
| 2X Reaction Buffer | 25 μl | 1X |
| dNTP Mix | 0.5 μl | 0.25 mM each dNTP |
| Forward Primer (10 μM) | Variable | 0.3 to 1 μM |
| Reverse Primer (10 μM) | Variable | 0.3 to 1 μM |
| Template DNA | Variable | < 200 ng/50 μl |
| ZymoTaq™ DNA Polymerase (5 U/μl) | 0.4 μl | 2 U/50 μl |
| ddH ₂ O | to 50 μl | – |
| Total volume | 50 μl | |

Note: The final concentration of MgCl₂ in the reaction (above) is 1.75 mM. If required, scale reaction reagent volumes accordingly to optimize the MgCl₂, primer, and/or template concentrations.

Suggested Conditions For PCR:

| | | |
|----------------------|---------------------|--------------------------|
| Initial denaturation | 95 °C | 10 min. |
| Denaturation | 94 to 96 °C | 30 sec. |
| Annealing | Variable | 30 to 40 sec. |
| Extension | 72 °C | 30 to 60 sec. for ≤ 1kb* |
| | 30-40 Cycles | |
| Final extension | 72 °C | 7 min. |
| Hold | 4 °C | > 4 min. |

***Note:** Add 15 to 30 seconds to the extension time for each additional kb over 1 kb. Make adjustments to the temperature if necessary.



Efficient PCR amplification of bisulfite treated DNA for methylation detection. The figure shows 200 bp product formation from bisulfite-treated DNA using ZymoTaq™ PreMix versus the polymerase systems from Suppliers R, Q, S, and I, respectively. In each case, equal amounts of bisulfite-treated DNA (EZ DNA Methylation-Direct™ Kit from Zymo Research) were used for each PCR, and the products separated in a 2.0% (w/v) agarose/TAE/EtBr gel. Also evident from the image, is the low occurrence of both primer dimerization and non-specific product formation.

Related Products:

| | E2003 (50 Rxns.) | E2004 (200 Rxns.) | Conc. | Storage Temp. |
|-----------------------------------|---------------------|----------------------|-------|------------------|
| ZymoTaq™ PreMix | 2 x 625 μl | 8 x 625 μl | 2X | -20 °C |
| DNase/RNase-Free H ₂ O | 2 x 1 ml | 5 x 1 ml | - | -20 °C |

Continued on reverse side...

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

Also Available:

| Product Name | Size | Catalog number |
|-------------------------------------------------------------|-----------|----------------|
| EZ DNA Methylation™ Kit | 50 | D5001 |
| | 200 | D5002 |
| | 2 x 96 | D5003 |
| | 2 x 96 | D5004 |
| EZ DNA Methylation-Gold™ Kit | 50 | D5005 |
| | 200 | D5006 |
| | 2 x 96 | D5007 |
| | 2 x 96 | D5008 |
| EZ DNA Methylation-Direct™ Kit | 50 | D5020 |
| | 200 | D5021 |
| | 2 x 96 | D5022 |
| | 2 x 96 | D5023 |
| EZ DNA Methylation-Startup™ Kit | 1 Kit | D5024 |
| EZ Bisulfite DNA Clean-up Kit™ | 50 | D5025 |
| | 200 | D5026 |
| | 2 x 96 | D5027 |
| | 2 x 96 | D5028 |
| Universal Methylated DNA Standard | 1 set | D5010 |
| Universal Methylated Human DNA Standard | 1 set | D5011 |
| Universal Methylated Mouse DNA Standard | 1 set | D5012 |
| Human HCT116 DKO Methylation Standards | 1 set | D5014 |
| Human HCT116 DKO Non-methylated DNA Standard | 5 µg | D5014-1 |
| Human HCT116 DKO Methylated DNA Standard | 5 µg | D5014-2 |
| Bisulfite Converted Universal Methylated Human DNA Standard | 1 set | D5015 |
| <i>E. coli</i> Non-methylated Genomic DNA | 5 µg | D5016 |
| Methylated-DNA IP Kit | 10 | D5101 |
| ChIP DNA Clean & Concentrator™ | 50 | D5201 |
| | 50 | D5205 |
| Anti-5-Methylcytosine Monoclonal Antibody (clone 10G4) | 50 µg | A3001-50 |
| | 200 µg | A3001-200 |
| CpG Methylase (M.SssI) | 200 units | E2010 |
| | 400 units | E2011 |

Trademarks and Disclaimers:

™ Trademarks of Zymo Research Corporation.

This product is for research use only and should only be used by trained professionals. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

The Polymerase Chain Reaction (PCR) process is covered by U.S. Patent: #4,683,195; 4,683,202 assigned to Hoffmann-La Roche. Patent pending in other countries. No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser by the purchase of Zymo Research's products. Further information on purchasing licenses to practice the PCR process can be obtained from the director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

Version 1.0.3

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com