SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Plasminogen Activator Inhibitor 1 (PAI-1)

REF: RT-38

Determination of the 4G/5G insertion/deletion on PAI-1 in Real Time PCR

INTRODUCTION AND PURPOSE OF USE

It has been stated that different polymorphisms result in different levels of PAI-1 concentrations. In case of the PAI-1-675 4G/5G, increased levels of PAI-1 in plasma were detected specially for patients carrying 4G/4G genotypes. Some studies demonstrated association of PAI-1 4G4G genotype with myocardial infarction. The Plasminogen Activator Inhibitor 1 (PAI-1) Kit is a qualitative test that allows the allelic discrimination, by means of *Real Time PCR*, of 4G/5G polymorphism on PAI-1 gene. The Procedure allows the amplification of Wild-type alleles and mutated alleles of PAI-1. Allelic discrimination is performed making a scatter plot of mutated allele's fluorescence (Y axis) versus wild-type allele's fluorescence (X axis), discriminating in this way the three possible genotype: Homozygote Wild-Type, Homozygote Mutated and Heterozygote Mutated. The analysis of the results is made by an instrument of Real Time PCR, composed by a thermal cycler with a system of fluorescence detection.

CONTENT

The kit contains reagents enough to perform 48 amplification tests:

	Quantity	Description
R1	3 x 220 µl	Amplification mMix dNTPs, Tris-HCl, KCl, MgCl ₂ , Taq Polymerase, Nuclease-free water, ROX
R2	3 x 130 µl	PAI-1 probes Mix Upstream primer, downstream primer, Target probes (FAM for Wild Type and VIC for Mutant), Nuclease-free water
R3	3 x 35 µl	Positive Control Homozygote Wild-Type cloned DNA corresponding to PAI-1 Wild-Type gene.
R4	3 x 35 µl	Positive Control Homozygote Mutated cloned DNA corresponding to PAI-1 Mutated gene.
R5	1 x 30 µl	Negative Control

Instruction for use: **ST. RT38-ENG.2**

MATERIALS AND STRUMENTATION REQUIRED BUT NOT SUPPLIED

Disposable latex powder-free gloves or similar material;
Bench microcentrifuge (12,000 - 14,000 rpm);
Micropipettes and Sterile tips with aerosol filter;
Vortex;

Plastic materials (microplate and optical adesive cover);
EZ1 ADV XL DSP DNA Blood Card (ref. 9018702)

Reagents

The PAI-1 kit was developed and validated to be used with the following extraction method:

Manual Extraction

Ref. 51304/51306 – *QiAmp DNA Mini kit*. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 50/250 samples. (QIAGEN)

Automatic extraction

Ref. 62124
EZ1 DSP DNA Blood kit. The kit allows the automatic DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 48 samples. (QIAGEN)

Strumentation

Automatic Extraction

Ref. 9001492.
EZ1 Advanced XL. Robotic Workstation. (QIAGEN) for the automatic nucleic acid purification up to 14 samples simultaneously (QIAGEN).

Real Time PCR

The PAI-1 kit was developed and validated to be used with the following real time PCR instruments:

- 7500 Fast* from Lifetechnologies
- StepOne Plus* from Lifetechnologies
- Rotor-Gene Q* MDx from QIAGEN
- Versant kPCR AD* from *Siemens* or *Stratagene MX3005P/MX3000P*

SAMPLES AND STORAGE

The Plasminogen Activator Inhibitor 1 (PAI-1) system must be used with extracted DNA from the following biological samples: Whole Blood EDTA. Collected samples must be shipped and stored at +2 - +8°C and used within 3 days from the collected data.

Store the sample at -20°C if it is used after 3 days.

PRECAUTIONS USE

This kit is for *in vitro* diagnostic (IVD), for professional use only and not for *in vivo* use.

After reconstitution, the amplification master mix must be used in one time (16 reactions). Do not re-frost already used material. Repeat thawing and freezing of reagents (more than twice) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.

At all times follow Good Laboratory Practice (GLP) guidelines. Wear protective clothing such as laboratory coats and disposable gloves while assaying samples. Avoid any contact between hands and eyes or nose during specimens collection and testing. Handle and dispose all used materials into appropriate bio-hazard waste containers. It should be discarded according to local law. Keep separated the extraction and the reagents preparation. Never pipette solutions by mouth. Avoid the air bubbles during the master mix dispensing. Eliminate them before starting amplification. Wash hands carefully after handling samples and reagents. Do not mix reagents from different lots. It is not infectious and hazardous for the health (see Material Safety data Sheet – MSDS). Do not eat, drink or smoke in the area where specimens and kit reagents are handled. Read carefully the instructions notice before using this test. Do not use beyond the expiration date which appears on the package label. Do not use a test from a damaged protective wrapper.

LIMIT OF THE METHOD

Use only extracted DNA from **Whole blood collected in EDTA**. Do not use contaminated DNA with extracted mucoproteins or haemoglobin; the latter inhibit the amplification reaction of nucleic acids and may cause invalid results . The extreme sensitivity of gene amplification may cause false positives due to cross-contamination between samples and controls. Therefore, you should:

- physically separate all the products and reagents used for amplification reactions from those used for other reactions, as well as from post-amplification products;
- use tips with filters to prevent cross-contamination between samples;
- use disposable gloves and change them frequently;
- carefully open test tubes to prevent aerosol formation;
- close every test tube before opening another one.

The proper functioning of the amplification mix depends on the correct collection, correct transportation, correct storage and correct preparation of a biological sample. As with any diagnostic device, the results obtained with this product must be interpreted taking in consideration all the clinical data and other laboratory tests done on the patient. As with any diagnostic device, with this product there is a residual risk of obtaining invalid, false positives or false negatives results.

STORAGE AND STABILITY

Store the product Plasminogen Activator Inhibitor 1 (PAI-1) at -20°C.

The *PAI-1* kit is shipped on dry ice. The kit components should be frozen. If one or more components are not frozen upon receipt or if the tubes have been compromised during transport, contact Clonit srl for assistance.

An intact and well stored product has a stability of 12 months from the date of production. Do not use beyond the expiration date which appears on the package label. Repeat thawing and freezing of reagents (more than twice) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.

ANALYTICAL PROCEDURE

Human DNA Extraction

Manual Extraction

Ref. 51304/51306
Follow the instructions inside the kit *QiAmp DNA Mini Kit*. Elute the sample in 50 µl of buffer AE.

Automatic extraction

Ref. 62124 - *EZ1 DSP DNA Blood kit* on *EZ1 Advanced XL instrument*.
Follow the instructions inside the kit *EZ1 DSP DNA Blood kit*. Start from 200 µl of sample and elute it in 50 µl of buffer AE.

SOFTWARE SETTING:

Lifetechnologies 7500 fast/StepOne plus

Turn the instrument and the computer on and open the control software. Click on "**Advance Setup**": by default the software will show the page "**experiment properties**". Write in the "**experiment name**" the file name, choose the type of instrument (**ABI7500** or **ABI7500fast / Step One Plus**), the type of reaction (**Genotyping**),

the type of used reagent (**Taqman®Reagents**) and the reaction time of analysis (**Standard ≈ 2 hours to complete a run**).

Open the page named "**page setup**".

In the window "**Assign SNP assay to the selected wells**" open "**Create new SNP Assay**" and set:

SNP Assay Name: PAI-1

	Reporter	Quencher
Allele 1 Name: PAI-1 WT	FAM	NFQ
Allele 2 Name: PAI-1 MUT	VIC	NFQ

In the page "**plate setup**", move on the area "**Assign Sample to the selected Wells**": set the name of the analysing samples, of positive controls and negative controls.

Choose an area of the plate where positive controls will be placed: select in the blank "**Assign SNP assay to the selected well**" and assign the PAI-1 SNP Assay. After set these tasks:

- task Positive control Allele1/Allele1** for PAI-1 Wild Type homozygote target;
- task Positive control Allele2/Allele2** for PAI-1 mutated homozygote target;

Choose an area of the plate where negative control will be placed: select "**Assign SNP assay to the selected well**" the "task Negative control" for PAI-1 SNP. Select an area of the plate where samples will be placed: select the wells and set PAI-1 SNP. Link every well to a sample, through the window "**Assign samples to selected wells**".

For each sample, select in the blank "**Assign SNP to selected wells**" the "**task UnKnown (U)**" for the PAI-1 SNP. Set ROX as passive reference, using it as normalizer of detected fluorescence.

Open "**Run Method**" (sheet **Graphic View**) and set the thermal cycling as follows:

cycles	Pre PCR Read	denaturation	Annealing extension	Post PCR Read
1	60°C 1min			
1		95°C 10min		
35		95°C 15sec	58°C 1min	
1				60°C1min

In the window "**Reaction volume plate per well**" set a volume of 25 µl.

After having prepared the plate, and correctly inserted in the instrument, press the button "**Start Run**".

ROTOR-GENE Q MDx

The experiments can be **set using the Quick Start Wizard** or the **Advanced Wizard**, which appears when the software is started.

Select the wizard "**Advanced**". As a first step, select the model "**Two Step Reaction**" with a double click in the "**New Run**".

In the next window, select the type of rotor installed on the instrument from the list that appears. Check the "**Locking Ring Attached**", check the checkbox and then click "**Next**".

Enter the name of the operator and the reaction volume of 25 µl, and then click "**Next**".

In the next window click on "**edit profile**". Set the following thermal cycle:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	58° C 1 min

Select the annealing / extension from the thermal profile and click on "**Acquiring A to cycling**."

In the next window, select yellow from the available channels and add it to acquiring channel along with the green channel and click "**OK**". In the next window click on "**OK**" and then click "**Next**".

click on "Edit Gain" button and set the following values for each channel:

Reporter	Gain
Green	6
Yellow	4

To begin the course, click on the button "**Start Run**". You can save the model before you begin your run by clicking on "**Save Template**". After clicking on the button "**Start Run**" window appears "**Save As**". The stroke can be saved in the desired position by the user.

Once the run started, the window "**Edit Samples**" allows you to set the name of samples and controls in the positions in which they were loaded on the instrument.

Select the locations where they were positioned the Wild Type and Mutant controls designate them as **PAI-1 WT positive CTR** and **PAI-1 Mutant Positive CTR**. Clicking on the box next to "**Type**" correspondent, in the dropdown menu "**Samples**" you can select the type of sample being analysed. Select "**Positive Controls**".

Select the location where you placed the Negative Control and name it as **Negative Control**. Clicking on the box next to "**Type**" correspondent, in the dropdown menu "**Samples**" you can select the type of sample being analysed. Select "**Negative Controls**".

Select the location of each sample and enter the name or code of the patient. Clicking on the box next to "**Type**" correspondent, in the dropdown menu "**Samples**" you can select the type of sample being analysed. Select "**UnKnown**".

At the end of the operation click "**OK**" in the "**edit samples**" and wait until the end of the race for the analysis (see "**Interpretation of Results**").

Analysis (see "**Interpretation of Results**").

Versant kPCR AD or Stratagene MX3005P/MX3000PMX3000P

Turn on the instrument, the computer and start the control software. Turn on the light at least 20 minutes before starting a new experiment. You can click on the lamp image for turning on the light from toolbar or you can select "**Lamp On**" from menù "**Instruments**".

Verify the correct setting of fluorescent reporters gains: In the setting menu choose and then "**Filter set gain setting**"

Reporter	Gain
FAM	8
JOE/HEX	2
ROX	1

From the major screen of the software it will open a window "**New option – Select experiment/Project type**": select "**Allelic Discrimination/SNPs Real Time**". the software automatically open the window **plate setup**.

Choose a plate zone where you can put the positive control **PAI-1** Wild Type and select in the bar on the right from the menu.

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
Pos. Control FAM	FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window "**well information**" here you can set the name of the calibrator: **PAI-1** - WT

It's possible to set the name of the dye near the name of the analysed target:

FAM	HEX
PAI-1-WT	PAI-1-MUT

Choose a plate zone where you can put the positive control **PAI-1** mutated and select in the bar on the right from the menu.

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
Pos. Control HEX	FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window "**well information**" here you can set the name of the calibrator: **PAI-1**- MUT

It's possible to set the name of the dye near the name of the analysed target:

FAM	HEX
PAI-1-WT	PAI-1-MUT

Choose a plate zone where you can put the unknown samples and select in the bar on the right from the menu.

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
UnKnown	FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window "**well information**" here you can set the name of the sample. Every well "**Unknown**" should be singularly named clicking on single name it will open the window "**well information**" in this window it is possible to insert the sample name in analysis.

Choose a plate zone where you can put the negative control and select in the right bar on the right from the menu.

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
NTC	FAM/HEX/ROX	ROX	None

Select all the well corresponding to occupied places and clicking right key of the mouse selecting the voice "**well information**": corresponding to **FAM** fluorophore set in the specific place (**Assay**): **PAI-1 WT**. Corresponding to fluorophore **HEX** set in the specific place (**Assay**): **PAI-1 MUT**.

Open page "**Thermal Profile setup**" and set the right thermal cycle:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	58° C 1 min

After preparing the plate and correctly insert it in the instrument, push the button "**Start Run**". the software will ask the file name for the saving data.

PREPARATION OF THE REACTIONS:

Defrost a tube of **Amplification mMix**;
Defrost a tube of **PAI-1 probes Mix**;
Mix carefully by vortex **210µl of Amplification mMix** and **126µl of PAI-1 probes Mix** (the mix as produced is enough to prepare **16 reactions** of amplification: **2 positive controls, 1 negative control and 13 samples**).

Distribute, in the amplification plate, **20 µl of just reconstituted mix** in chosen positions and already setted on the instrument software.

Distribute, in the negative control position, **5µl** of solution taken by the **negative control** vial.

Distribute, in chosen position for each sample, **5µl of corresponding sample**.

Distribute, in chosen positions for the positive controls, **5µl of Positive Control Homozygote Wild-Type** and **Positive Control Homozygote Mutated**.

Seal up accurately the plate using an optical adhesive film and verify that there aren't air bubbles in the mix, to avoid the amplification interferences.

Transfer the plate in the instrument and push the button "**Start Run**".

QUALITATIVE ANALYSIS


Lifetechnologies 7500 Fast, StepOne Plus.

At the end of the amplification reaction, the software automatically shows the obtained results in the "**Allelic Discrimination Plot**". The Lifetechnologies 7500Fast/StepOne Plus instrument automatically perform the genotyping of unknown samples by comparing controls versus homozygous wild-type and homozygous mutated contained in the kit. A proper results' analysis needs a correct settings of the instrumentation. For this aim, set:

Baseline fluorescence level from cycle 6;

	PAI-1-WT - FAM Threshold	PAI-1-MUT - VIC Threshold
7500 Fast	0.02	0.02
StepOne Plus	0.02	0.02

It is recommended to verify correct placement of each individual sample on the scatter plot. For viewing the report containing all data obtained during the analysis, click the sheet "**view table well**".

The position of the negative control in other plot region could be a contamination of the reaction mix. Be sure that in the negative control not be increasing of specific fluorescence by examining targets (FAM and VIC).

ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose it is necessary a correct setting of the software

Set the level of background fluorescence (Baseline) from cycle 6; Set the following threshold:

	PAI-1-WT - FAM Threshold	PAI-1-MUT - VIC Threshold
7500 Fast	0.02	0.02
StepOne Plus	0.02	0.02

Export data to Excel and set the formula for each sample and control:

Allele2 Ct (PAI-1-mut) – Allele1 Ct (PAI-1-wt)

See paragraph "INTERPRETATION OF RESULTS"

ROTOR - GENE Q MDx

At the end of the amplification reaction, click on **Analysis** from RotorGene Software. In the **Analysis window** select **Allelic Discrimination** sheet, click on **Cycling A green-Cycling A yellow** using **CTRL** button and click **“show”**. The amplification plot will appear.

Set in the **“CT calculation – Threshold”** the following value:

	PAI-1-WT - FAM Threshold	PAI-1-MUT - VIC Threshold
RotorGene-Q	0.05	0.05

Click **genotypes** button and set:

	Reacting Channel	Reacting Channel
Wild type	cycling A (green)	
Heterozigous	cycling A (green)	cycling A (yellow)*
None		
Mutant		cycling A (yellow)*

In the **Analysis window** select **“Scatter graph analysis”** sheet, click on **Cycling A green-Cycling A yellow** and click **“show”**. The scatter plot will appear with the 2 clouds: Wild Type (bottom right) and Mutant (high left). Is it possible to print the final report.

ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

At the end of the PCR run open the **"Analysis"** window. Select the **“Quantification”** sheet and click on **"cycling A (green)"**.

Select from the menu **“Dynamic Tube”** and subsequently **“Slope correct”**.

Set the correct setting of the threshold in the space provided **“CT calculation – Threshold”**.

	PAI-1-WT – Green (FAM) Threshold
RotorGene-Q	0.05

Open the **"Analysis"** window. Select the **“Quantification”** sheet and click on **"cycling A (yellow)"**.

Select from the menu **“Dynamic Tube”** and subsequently **“Slope correct”**.

Set the correct setting of the threshold in the space provided **“CT calculation – Threshold”**.

	PAI-1-MUT – Yellow (VIC) Threshold
RotorGene-Q	0.02

Export data to Excel, save the file as "Excel Analysis Sheet" and enter the following formula for each sample and control:

Yellow Ct (PAI-1-MUT) – Green Ct (PAI-1-WT)

See paragraph "INTERPRETATION OF RESULTS"

Versant kPCR AD or Stratagene MX3005P/MX3000PMX3000P

Click on button **“Analysis”** in the toolbar. The software immediately open sheet **“Analysis Term Setting”**. Activate the keys FAM and HEX in below screen and select the samples and disposable controls.

On the window **“Analysis Term Setting”** open the sheet **“Results”**. Select in the right screen the voice **“Amplification Plot”** in the area **“Area to Analysis”**.

Set the values in the area **“Threshold Fluorescence”** the values:

	PAI-1-WT - FAM Threshold	PAI-1-MUT - VIC Threshold
Versant kPCR	0.02	0.08

Select in the right screen the voice **“Dual Color Scatter Plot”** in the area **“Area to Analysis”**. Select **“Fluorescence”** in the area **“Display value for”** and choose from menu **“Rlast/Rfirst”**.

Select in the area **“Allele Association”**:

A Allele: **PAI-1 Wild Type**

B Allele: **PAI-1 Mutated**

At the end it will be possible obtain the detailed account, clicking **“Text Report”** in the area **“Area to Analysis”**.

Only for Versant kPCR AD or Stratagene MX3005P/MX3000PMX3000P

If the software can't discriminate the wild type and mutated positive controls click **“Show all genotypes”**. It will be now possible to

modify the genotyping windows, select the samples placed near the positive wild type and identify them as wild type. Select the samples placed near the positive mutated and identify them as mutated. Select the samples placed in an intermediate position between wild type positive control and mutated positive control and identify them as Heterozygote for the 4G/5G polymorphism

ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

Click "Analysis" in the toolbar. Click the sheet "Results"; and choose the analysis "Amplification plot". Check the correct setting of the threshold in the window "Threshold fluorescence" and set the following values:

	PAI-1-WT - FAM Threshold	PAI-1-MUT - HEX Threshold
MX3000P – MX3005P	0.02	0.02

From the Text Report window you can export the results by clicking on the main menu: file, export

Export data to Excel and set the formula for each sample and control:

HEX Ct (PAI-1-MUT) – FAM Ct (PAI-1-WT)

See paragraph "INTERPRETATION OF RESULTS"

INTERPRETATION OF RESULTS

The use of positive and negative control in each amplification session allows to verify the correct functioning of the amplification mix and the absence of any contamination.

The instrument software is able to analyze the fluorescences that are emitted by the specific probe for **PAI-1** Wild-Type (FAM) and by the specific probe for the **PAI-1** mutated (VIC/HEX).

Scatter Plot Analysis

A regular functioning of the amplification mix can be verified analysing the correct position of positive controls and negative controls on the scatter plot.

- Positive control Homozygote Wild-Type: horizontal position on X axis (down on the right) Ct< 30
- positive Control homozygote Mutated: vertical position on Y axis (up on the left) Ct< 30
- negative Control: placed at the origin of cartesian plane (down on the left)

Genotyping tests (allelic discrimination) are endpoint experiments: fluorescence data are collected at the end of the reaction (Post PCR Read) and subtracted to initial read fluorescence (Pre PCR Read).

The software makes a scatter plot with obtained results: Y axis is the normalized fluorescence of Mutated Allele, while X axis shows the normalized fluorescence of Wild-Type Allele. The diagnosis obtained with the comparison between unknown samples and Homozygote Wild-Type and Homozygote Mutate, given by the kit. Selecting the controls, we obtain their disposition on the scatter plot, depending of their relation between fluorescence emitted by two probes FAM (Wild-Type) and VIC/HEX (Mutated).

ΔCt Analysis

By analyzing ΔCt you can identify the correct genotype of the sample being analyzed by performing subtraction of the assigned Mutated allele Ct to assigned wild-type allele Ct. The genotype is determined by following the table below:

Genotype	ΔCt (Ct Mutated - Ct Wild Type)
Wild Type	ΔCt > 2
Mutated	ΔCt < - 2
Heterozygous	-2 > ΔCt < 2

In the amplification of each positive control (wild-type and Mutated), the Ct values of the allele-specific probe are used to validate the assay.

Make sure that the fluorescence emitted by the amplification allele has a correct Ct <30.

If the result of the allele-specific amplification of each control has a Ct > 30 or undetermined the session can not be considered valid and must be repeated.

In the amplification of each sample, the Ct values of the allele-specific probe are used to validate the assay from the extraction process up to the stage of detection.

Make sure that the fluorescence emitted by the allele-specific amplification of the sample, identified after analysis, has not a Ct > 30.

If a sample has a Ct > 30 means that there are problems in the extraction phase or in the amplification and therefore could be assigned to a sample genotype wrong. Repeat sample.

Sample	Wild Type Allele (FAM)	Mutated Allele (VIC)	Assay and Genotype
Wild Type	Ct > 30	Not relevant	Not Valid
Mutated	Not relevant	Ct > 30	Not Valid
Heterozygous	Ct > 30	Ct > 30	Not Valid

PERFORMANCES

Clinical Sensitivity:

For the purposes of this evaluation it is considered as clinical sensitivity the skill of determining true heterozygote and homozygote mutated samples in the totality of screened samples.. The test is performed following the method advices. Positive samples are confirmed with another CE disposable method.

Clinical sensitivity is 100% for material extracted from EDTA blood.

Diagnostic Specificity:

For the purposes of this evaluation is considered as diagnostic specificity the skill of the method of determining real negative samples. The diagnostic specificity of the system is valued analysing human genomic samples tested and confirmed as negative with another disposable system.

Diagnostic Specificity is 100% for material extracted from EDTA blood.

Analytical specificity:

Test's specificity is guaranteed by the use of specific primers for determining Plasminogen Activator Inhibitor 1 gene and of probes intentionally designed on 4G/5G polymorphism..

The alignment of the choose regions for specific primers' hybridization with available sequences of present in database, demonstrated: their conservation and the complete specificity for the analyzed targets. Samples that are defined as positive for a determined genotype as much must be recognized by the amplification system.

INTERFERENCES:

Verify that in DNA extracted from the sample there aren't nucleoproteins and haemoglobin, in way of exclude possible inhibition of PCR reactions. The interference due to contaminants can be highlighted through the spectrophotometric analysis and obtained data report at 260 nm (maximum absorbtion of Nucleic Acids) and 280 nm (maximum absorbtion of Proteins). A pure DNA might have a rate of approximately 1.8.

QUALITY CONTROL

It is therefore recommended to insert as quality control of every extraction session, amplification and detection of a negative sample and of a positive sample which have already tested before or referential material with known concentration. In accordance with the Clonit srl ISO EN 13485 Certified quality Management System, each lot of Factor V of Leiden is tested against predetermined specification to ensure consistent product quality

BIBLIOGRAPHY

Soudabeh Hosseini et al., **Genetic risk factors in patients with deep venous thrombosis, a retrospective case control study on Iranian population**, Thrombosis Journal (2015) 13:35 DOI 10.1186/s12959-015-0064-y


Tajinder Kumar Parpugga et al. , **The Effect of PAI-1 4G/5G Polymorphism and Clinical Factors on Coronary Artery Occlusion in Myocardial Infarction**, Disease Markers Volume 2015, Article ID 260101, 8 pages http://dx.doi.org/10.1155/2015/260101

Rosa Divella et al., **Circulating Levels of PAI-1 and SERPINE1 4G/4G Polymorphism Are Predictive of Poor Prognosis in HCC Patients Undergoing TACE**, Translational Oncology, Volume 8 Number 4 August 2015 pp. 273–278

TECHNICAL ASSISTANCE

For any question and support please contact our Technical support:

e-mail: info@clonit.it
phone: +39 02 56814413

	In vitro diagnostic device
	Read the instruction's manual
	Range of temperature
	Use within (dd/mm/yyyy: year-month)
	Lot (xxxx)
	Code
	Manufacturer
	Contains sufficient for <n> tests

EDMA code: 16010490
CND: W01060104

The **Plasminogen Activator Inhibitor 1** kit is CE marked diagnostic kit according to the European *in vitro* diagnostic directive 98/79/CE.

CLONIT S.r.l.

Headquarter: Via Varese 20 – 20121 Milano
Production Site: Via B. Quaranta 57 - 20139 Milano
Tel. + 39. (0)2.56814413 fax. +39. (0)2.56814515
www.clonit.it - info@clonit.it

Plasminogen Activator Inhibitor 1 (PAI-1)

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