**SZABO  
SCANDIC**

Part of Europa Biosite

## Produktinformation

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[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](http://linkedin.com/company/szaboscandic)

## MTHFR A1298C

REF: RT-55

### Determination of the A1298C mutation on MTHFR in Real Time PCR

#### INTRODUCTION AND PURPOSE OF USE

Heredity thrombophilia is a genetically determined tendency to develop thromboembolic disease at a young age with recurrent character. Environmental risk factors (smoking, alcohol, stress, diet, ...) and genetic mutations, such as the A1298C mutation localized on the MTHFR gene, interact resulting in the expression of this disease.

The MTHFR A1298C Kit is a qualitative test that allows the allelic discrimination, by means of *Real Time PCR*, of A1298C mutation associated to Hereditary Thrombophilia.

The Procedure allows the amplification of Wild-type alleles and mutated alleles of MTHFR A1298C. Allelic discrimination is performed making a scatter plot of mutated allele's fluorescence (Y axis) versus wild-type allele's fluorescence (X axis), discriminating in this way the three possible genotype: Homozygote Wild-Type, Heterozygote Mutated and Heterozygote Mutated.

The analysis of the results is made by an instrument of Real Time PCR, composed by a thermal cycler with a system of fluorescence detection.

#### CONTENT

The kit contains reagents enough to perform 48 amplification tests:

Quantity	Description
R1 3 x 220 µl	Amplification mMix: dNTPs, Tris-HCl, KCl, MgCl <sub>2</sub> , Taq Polymerase, Nuclease-free water (Blue cap)
R2 3 x 130 µl	MTHFR A1298C probes Mix Upstream primer, downstream primer, Target probe (FAM for Wild Type and VIC for Mutant), Nuclease-free water (Green cap)
R3 3 x 35 µl	Positive Control Homozygote Wild-Type cloned DNA corresponding to MTHFR A1298C Wild-Type gene.
R4 3 x 35 µl	Positive Control Homozygote Mutated cloned DNA corresponding to MTHFR A1298C Mutated gene.
R5 1 x 30 µl	Negative Control (Yellow cap)

Instruction for use: ST. RT55-ENG.7

#### MATERIALS AND STRUMENTATION REQUIRED BUT NOT SUPPLIED

Disposable latex powder-free gloves or similar material;  
Bench microcentrifuge (12,000 - 14,000 rpm);  
Micropipettes and Sterile tips with aerosol filter;  
Vortex;  
Plastic materials (microplate and optical adhesive cover);  
Heat block (only for extraction)  
Dry block shacker for 1.5ml conical tubes  
Magnetic rack for 1.5ml conical tubes  
EZ1 ADV XL DSP DNA Blood Card (ref. 9018702)

#### Reagents

The MTHFR A1298C kit was developed and validated to be used with the following extraction method:

#### Manual Extraction

Ref. 51304/51306

QIAamp DNA mini kit. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 50/250 samples. (QIAGEN)

#### Automatic extraction

Ref. 62124

EZ1 DSP DNA Blood kit. The kit allows the automatic DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 48 samples. (QIAGEN)

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10629800 - VERSANT® Sample Preparation 1.2 Reagents kit box 1. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 96 samples.

10629801 - VERSANT® Sample Preparation 1.2 Reagents kit box 2. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 96 samples.

#### Strumentation

##### Automatic Extraction

Ref. 9001492.

EZ1 Advanced XL. Robotic Workstation. (QIAGEN) for the automatic nucleic acid purification up to 14 samples simultaneously (QIAGEN).

#### Real Time PCR

The MTHFR A1298C kit was developed and validated to be used with the following real time PCR instruments:

- 7500 Fast from Lifetechnologies
- StepOne Plus from Lifetechnologies
- Rotor-Gene Q MDx from QIAGEN
- Versant kPCR AD from Siemens or Stratagene MX3005P/MX3000P
- CFX96 Real Time PCR System from BioRad

#### SAMPLES AND STORAGE

The MTHFR A1298C system must be used with extracted DNA from the following biological samples: Whole Blood EDTA. Collected samples must be shipped and stored at +2 - +8°C and used within 3 days from the collected data.

Store the sample at -20°C if it is used after 3 days.

#### PRECAUTIONS USE

This kit is for *in vitro* diagnostic (IVD), for professional use only and not for in vivo use.

**After reconstitution, the amplification master mix must be used in one time (16 reactions). Do not re-frost already used material. Repeat thawing and freezing of reagents (more than twice) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.**

At all times follow Good Laboratory Practice (GLP) guidelines. Wear protective clothing such as laboratory coats and disposable gloves while assaying samples.

Avoid any contact between hands and eyes or nose during specimens collection and testing.

Handle and dispose all used materials into appropriate bio-hazard waste containers. It should be discarded according to local law. Keep separated the extraction and the reagents preparation.

Never pipette solutions by mouth.

Avoid the air bubbles during the master mix dispensing. Eliminate them before starting amplification.

Wash hands carefully after handling samples and reagents.

Do not mix reagents from different lots.

It is not infectious and hazardous for the health (see Material Safety data Sheet - MSDS).

Do not eat, drink or smoke in the area where specimens and kit reagents are handled.

Read carefully the instructions notice before using this test.

Do not use beyond the expiration date which appears on the package label.

Do not use a test from a damaged protective wrapper.

Follow the instructions inside the kit QIAamp DNA Mini Kit. Elute the sample in 50 µl of buffer AE.

#### Automatic extraction

Ref. 62124 - EZ1 DSP DNA Blood kit on EZ1 Advanced XL instrument.

Follow the instructions inside the kit EZ1 DSP DNA Blood kit. Start from 200 µl of sample and elute it in 50 µl of buffer AE.

#### Manual extraction (SIEMENS)

Ref. 10629800 - VERSANT® Sample Preparation 1.2 Reagents kit box 1.

Ref. 10629801 - VERSANT® Sample Preparation 1.2 Reagents kit box 2.

Follow the instructions supplied by Siemens and elute it in 70 µl of Elution buffer. Transfer 55 µl of eluted sample to an appropriately size tube.

Sample can be stored at -20°C.

#### SOFTWARE SETTING:

##### Lifetechnologies 7500 fast/StepOne plus

Turn the instrument and the computer on and open the control software. Click on "Advance Setup": by default the software will shows the page "experiment properties". Write in the "experiment name" the file name, choose the type of instrument (ABI7500 o ABI7500fast), the type of reaction (Genotyping), the type of used reagent (Taqman®Reagents) and the reaction time of analysis (Standard ~ 2 hours to complete a run).

Open the page named "page setup".

In the window "Assign SNP assay to the selected wells" open "Create new SNP Assay" and set:

SNP Assay Name: MTHFR A1298C:

	Reporter	Quencher
Allele 1 Name: MTHFR A1298C WT	FAM	NFQ
Allele 2 Name: MTHFR A1298C MUT	VIC	NFQ

In the page "plate setup", move on the area "Assign Sample to the selected Wells": set the name of the analyzing samples, of positive controls and negative controls.

Choose an area of the plate where positive controls will be placed: select in the blank "Assign SNP assay to the selected well" and assign the MTHFR A1298C SNP Assay. After set these tasks:

- "task Positive control Allele1/Allele1" for MTHFR A1298C Wild Type homozygote target;
- "task Positive control Allele2/Allele2" for MTHFR A1298C mutated homozygote target;

Choose an area of the plate where negative control will be placed: select "Assign SNP assay to the selected well" the "task Negative control" for MTHFR A1298C SNP.

Select an area of the plate where samples will be placed: select the wells and set MTHFR A1298C SNP. Link every well to a sample, through the window "Assign samples to selected wells".

Verify the correct setting of fluorescent reporters gains: In the setting menu choose and then "Filter set gain setting"

Reporter	Gain
FAM	8
JOE/HEX	2
ROX	1

From the major screen of the software it will open a window "New option - Select experiment/Project type": select "Allelic Discrimination/SNPs Real Time". The software automatically open the window plate setup.

Choose a plate zone where you can put the positive control MTHFR A1298C Wild Type and select in the bar on the right from the menu.

Well type:	Collect	Reference	Replicate
Pos. Control FAM	FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window "well information" here you can set the name of the calibrator: MTHFR-A1298C WT

It's possible to set the name of the dye near the name of the analyzed target:

FAM	HEX
MTHFR A1298C-WT	MTHFR A1298C-MUT

Choose a plate zone where you can put the positive control MTHFR A1298C mutated and select in the bar on the right from the menu.

Well type:	Collect	Reference	Replicate
Pos. Control HEX	FAM/HEX/ROX	ROX	None

Clicking on every well it will appear the dialogue window "well information" here you can set the name of the calibrator: MTHFR-A1298C MUT

It's possible to set the name of the dye near the name of the analyzed target:

1	95°C 10 min
35	95°C 15 sec

Select the annealing / extension from the thermal profile and click on "Acquiring A to cycling."

In the next window, select yellow from the available channels and add it to acquiring channel along with the green channel and click "OK". In the next window click on "OK" and then click "Next".

click on "Edit Gain" button and set the following values for each channel:

Reporter	Gain
Green	8
Yellow	2

To begin the course, click on the button "Start Run". You can save the model before you begin your run by clicking on "Save Template".

After clicking on the button "Start Run" window appears "Save As". The stroke can be saved in the desired position by the user.

Once the run started, the window "Edit Samples" allows you to set the name of samples and controls in the positions in which they were loaded on the instrument.

Select the locations where they were positioned the Wild Type and Mutant controls designate them as **MTHFR A1298C WT positive CTR** and **MTHFR A1298C Mutant Positive CTR**. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analyzed. Select "Positive Controls".

Select the location where you placed the Negative Control and name it as **Negative Control**. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analyzed. Select "Negative Controls".

Select the location of each sample and enter the name or code of the patient. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analyzed. Select "Unknown".

At the end of the operation click "OK" in the "edit samples" and wait until the end of the race for the analysis (see "Interpretation of Results").

Analysis (see "Interpretation of Results").

#### Versant kPCR AD or Stratagene MX3005P/MX3000P</h

## QUALITATIVE ANALYSIS

### Lifetechnologies 7500 Fast, StepOne Plus.

#### Scatter Plot Analysis

At the end of the amplification reaction, the software automatically shows the obtained results in the "Allelic Discrimination Plot".


The Lifetechnologies 7500Fast/StepOne Plus instrument automatically perform the genotyping of unknown samples by comparing controls versus homozygous wild-type and homozygous mutated contained in the kit.

A proper results' analysis needs a correct settings of the instrumentation. For this aim, set:

Baseline fluorescence level from cycle 6;

	MTHFR-A1298C WT - FAM Threshold	MTHFR-A1298C MUT - VIC Threshold
7500 Fast	0.35	0.3
StepOne Plus	0.35	0.3

It is recommended to verify correct placement of each individual sample on the scatter plot. For viewing the report containing all data obtained during the analysis, click the sheet "view table well".

The position of the negative control in other plot region could be a contamination of the reaction mix.  
Be sure that in the negative control not be increasing of specific fluorescence by examining targets (FAM and VIC).

#### ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose it is necessary a correct setting of the software

Set the level of background fluorescence (Baseline) from cycle 6; Set the following threshold:

	MTHFR-A1298C WT - FAM Threshold	MTHFR-A1298C MUT - VIC Threshold
7500 Fast	0.25	0.25
StepOne Plus	0.25	0.25

Export data to Excel and set the formula for each sample and control:

### Allele2 Ct (MTHFR A1298C -MUT) – Allele1 Ct (MTHFR A1298C -WT)

See paragraph "INTERPRETATION OF RESULTS"

## ROTOR - GENE Q MDx

#### Scatter Plot Analysis

At the end of the amplification reaction, click on **Analysis** from RotorGene Software. In the **Analysis window** select **Allelic Discrimination** sheet, click on **Cycling A green-Cycling A yellow** using **CTRL** button and click "**show**". The amplification plot will appear.

Set in the "CT calculation – Threshold" the following value:

	MTHFR A1298C -WT - FAM Threshold	MTHFR A1298C -MUT - VIC Threshold
RotorGene-Q	0.2	0.2

Click **genotypes** button and set:

	Reacting Channel	Reacting Channel
Wild type	cycling A (green)	cycling A (green)
Heterozygous	cycling A (green)	cycling A (yellow)*
None		
Mutant		cycling A (yellow)*

In the **Analysis window** select "Scatter graph analisys" sheet, click on **Cycling A green-Cycling A yellow** and click "**show**". The scatter plot will appear with the 2 clouds: Wild Type (bottom right) and Mutant (high left). Is it possible to print the final report.

#### ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

### CFX96 real time PCR system

#### Scatter Plot Analysis

At the end of the PCR, select the "Allelic Discrimination" sheet.

On the bottom of the screen, set "Selected fluorophores": X = FAM and Y = VIC. Choose RFU from "Display Mode" and "Normalize data". Set the following Threshold:

At the end of the PCR run open the "Analysis" window. Select the "Quantification" sheet and click on "cycling A (green)".

Select from the menu "Dynamic Tube" and subsequently "Slope correct".

Set the correct setting of the threshold in the space provided "CT calculation – Threshold".

	MTHFR A1298C -WT - FAM Threshold	MTHFR A1298C -MUT - VIC Threshold
RotorGene-Q	0.05	

	MTHFR A1298C -WT - FAM Threshold	MTHFR A1298C -MUT - VIC Threshold
RotorGene-Q	0.05	

Open the "Analysis" window. Select the "Quantification" sheet and click on "cycling A (yellow)".

Select from the menu "Dynamic Tube" and subsequently "Slope correct".

Set the correct setting of the threshold in the space provided "CT calculation – Threshold".

	MTHFR A1298C -WT - FAM Threshold	MTHFR A1298C -MUT - VIC Threshold
RotorGene-Q	0.05	

Export data to Excel, save the file as "Excel Analysis Sheet" and enter the following formula for each sample and control:

### Yellow Ct (MTHFR A1298C -MUT) – Green Ct (MTHFR A1298C -WT)

See paragraph "INTERPRETATION OF RESULTS"

### Versant kPCR AD or Stratagene MX3005P/MX3000P

#### Scatter Plot Analysis

Click on button "Analysis" in the toolbar. The software immediatetly open sheet "Analysis Term Setting". Activate the keys FAM and HEX in below screen and select the samples and disposable controls.

On the window "Analysis Term Setting" open the sheet "Results". Select in the right screen the voice "Amplification Plot" in the area "Area to Analysis".

Set the values in the area "Threshold Fluorescence" the values:

	MTHFR A1298C -WT - FAM Threshold	MTHFR A1298C -MUT - VIC Threshold
Kpcr Versant	0.55	0.55

Select in the right screen the voice "Dual Color Scatter Plot" in the area "Area to Analysis".

Select "Fluorescence" in the area "Display value for" and choose from menu "Rlast/Rfirst".

Select in the area "Allele Association":

A Allele: MTHFR A1298C Wild Type

B Allele: MTHFR A1298C Mutated

At the end it will be possible obtain the detailed account, clicking "Text Report" in the area "Area to Analysis".

### Only for Versant kPCR AD or Stratagene MX3005P/MX3000P

If the software can't discriminate the wild type and mutated positive controls click "Show all genotypes". It will be now possible to modify the genotyping windows, select the samples placed near the positive wild type and identify them as wild type. Select the samples placed near the positive mutated and identify them as mutated. Select the samples placed in an intermediate position between wild type positive control and mutated positive control and identify them as Heterozygote for the A1298C Mutation.

#### ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

	MTHFR A1298C -WT - FAM Threshold	MTHFR A1298C -MUT - VIC Threshold
MX300P – MX3005P	0.1	0.1

From the Text Report window you can export the results by clicking on the main menu: file, export

Export data to Excel and set the formula for each sample and control:

### HEX Ct (MTHFR A1298C -MUT) – FAM Ct (MTHFR A1298C -WT)

See paragraph "INTERPRETATION OF RESULTS"

### CFX96 real time PCR system

#### Scatter Plot Analysis

At the end of the PCR, select the "Allelic Discrimination" sheet.

On the bottom of the screen, set "Selected fluorophores": X = FAM and Y = VIC. Choose RFU from "Display Mode" and "Normalize data". Set the following Threshold:

	MTHFR A1298C -WT - FAM Vertical Threshold	MTHFR A1298C -MUT - VIC Horizontal Threshold
CFX 96	0.2	0.3

You can export the report pushing the paper block figure on the top of the screen.

#### ΔCt Analysis

Further analysis can be performed with the ΔCt study of the results. For this purpose you need a different setting of analysis and a correct setting of the software:

At the end of the PCR, select the sheet "Quantitation". Analyzing a fluorophore at a time, open in the menu "Setting" the voice Baseline Threshold and set the following values of the Threshold in the "User Defined"

	MTHFR A1298C -WT - FAM Threshold	MTHFR A1298C -MUT - VIC Threshold
CFX 96	50	50

You can export the report by clicking on the figure notepad at the top of the screen.

Export data to Excel, save the file as "Excel Analysis Sheet" and enter the following formula for each sample and control:

### VIC Ct (MTHFR A1298C -MUT) – FAM Ct (MTHFR A1298C -WT)

See paragraph "INTERPRETATION OF RESULTS"

## INTERPRETATION OF RESULTS

The use of positive and negative control in each amplification session allows to verify the correct functioning of the amplification and the absence of any contamination.

The instrument software is able to analyze the fluorescences that are emitted by the specific probe for MTHFR A1298C Wild-Type (FAM) and by the specific probe for the MTHFR A1298C mutated (VIC/HEX).

#### Scatter Plot Analysis

A regular functioning of the amplification mix can be verified analyzing the correct position of positive controls and negative controls on the scatter plot.

- Positive control Homozygote Wild-Type: horizontal position on X axis (down on the right) Ct < 30
- positive Control homozigote Mutated: vertical position on Y axis (up on the left) Ct < 30
- negative Control: placed at the origin of cartesian plane (down on the left)

Genotyping tests (allelic discrimination) are endpoint experiments: fluorescence data are collected at the end of the reaction (Post PCR Read) and subtracted to initial read fluorescence (Pre PCR Read).

The software makes a scatter plot with obtained results: Y axis is the normalized fluorescence of Mutated Allele, while X axis shows the normalized fluorescence of Wild-Type Allele. The diagnosis obtained with the comparison between unknown samples and Homozygote Wild-Type and Homozygote Mutate, given by the kit. Selecting the controls, we obtain their disposition on the scatter plot, depending of their relation between fluorescence emitted by two probes FAM (Wild-Type) and VIC/HEX (Mutated).

#### ΔCt Analysis

By analyzing ΔCt you can identify the correct genotype of the sample being analyzed by performing subtraction of the assigned Mutated allele Ct to assigned wild-type allele Ct. The genotype is determined by following the table below:

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