# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation

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[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

CE

# **Fluorescence Immunoassay Analyser**

## **Product Instruction**

Model: RTR-FS100

---

**Xiamen United Medical Instruments Co., Ltd.**

## **Information of the Product Instruction:**

Version No.: 01

Formulated/Revised on: August 20, 2020

## **Product Information:**

Product Name: Fluorescence Immunoassay Analyser

Model: RTR-FS100

Reference number: 110400000100

Registration Certificate No. for Medical Device: MXZZ 20182220129

No. of Product Technology Requirements: MXZZ 20182220129

Manufacturing Enterprise License No.: MSYJXSCX No. 20120324

Xiamen United Medical Instruments Co., Ltd.

5#, No.2008, Wengjiao West Rd, Haicang District, 361026 Xiamen, Fujian

Province, PEOPLE'S REPUBLIC OF CHINA

Date of Manufacture: See product label.

Service Life: 5 years

## **Contact Us:**

Zip Code: 361022

E-mail: [umic@xmumic.com](mailto:umic@xmumic.com)

TEL: 0592-3132888

FAX: 0592-3132888

## Statement

Copyright of the Product Instruction belongs to Xiamen United Medical Instruments Co., Ltd.

This Product Instruction is provided for the users to operate and understand the product or to learn component repair. Anyone is not allowed to use or publicize the Product Instruction in any way without authorization, and any improper conduct is regarded as illegal. Xiamen United Medical Instruments Co., Ltd. assumes no liabilities for the resulting illegal events or problems involving third-party interest.

All information involved herein is confirmed. Xiamen United Medical Instruments Co., Ltd. assumes no responsibilities for the accidental injuries or life-threatening events directly or indirectly induced by misuse or misoperation of the product.

Protection offered by the product may be damaged, once the product is not used as specified by the manufacturer.

# Contents

<b>1. Product Overview</b> .....	<b>- 6 -</b>
<b>2. Scope of Application</b> .....	<b>- 6 -</b>
<b>3. Product Name and Model</b> .....	<b>- 6 -</b>
<b>4. Product Composition</b> .....	<b>- 6 -</b>
<b>5. Work Principle</b> .....	<b>- 6 -</b>
<b>6. Matters Needing Attention</b> .....	<b>- 6 -</b>
<b>7. Graphic Illustrations to Product Appearance and Parts</b> .....	<b>- 7 -</b>
<b>8. Product Type and Technical Parameters</b> .....	<b>- 8 -</b>
8.1 Product type.....	- 8 -
8.2 Normal work conditions:.....	- 8 -
8.3. Technical parameters.....	- 8 -
8.4. Performance Indicators.....	- 8 -
8.5. Testing methods.....	- 9 -
8.6. Calculation of examination result.....	- 9 -
8.7. Literature references.....	- 10 -
<b>9. Product Installation</b> .....	<b>- 10 -</b>
9.1 Inspection.....	- 10 -
9.2 Installation.....	- 10 -
<b>10. Operating Instructions and Graphics</b> .....	<b>- 11 -</b>
10.1 User login.....	- 11 -
10.2 Test Function.....	- 12 -
10.2.1 Immediate Test.....	- 13 -
10.2.2 Timing Test.....	- 15 -
10.2.3 In/Out.....	- 16 -
10.2.4 Upload.....	- 16 -
10.2.5 Print.....	- 16 -
10.3 Record Search.....	- 16 -
10.4 Setting.....	- 19 -
10.5 System.....	- 19 -
10.6 Item import.....	- 20 -
<b>11. Electromagnetic Compatibility</b> .....	<b>- 21 -</b>
<b>12. Explanation to Label-included Marks</b> .....	<b>- 22 -</b>
<b>13. Maintenance</b> .....	<b>- 23 -</b>
<b>14. Transportation, Storage and Scrap</b> .....	<b>- 23 -</b>
<b>15. Troubleshooting and Statement</b> .....	<b>- 23 -</b>
<b>16. Contraindication: None</b> .....	<b>- 24 -</b>
<b>17. After-sales Services</b> .....	<b>- 24 -</b>
17.1 Warranty Coverage.....	- 24 -
17.2 After-sales Service Hotline and Address.....	- 25 -
<b>Appendix A Packing List</b> .....	<b>- 26 -</b>

## **1. Product Overview**

Fluorescence Immunoassay Analyzer (hereinafter referred to as the Fluorometer) is a microprocessor-based device for fluorescence immunoassay and analysis. It can be used for fluorescence immunoassay, characterized in easy operation, fast test speed (completed in 8s) and good repeatability (CV less than 2%).

## **2. Scope of Application**

Fluorometer is applicable for the immunofluorescence quantitative determination of human body samples, used together with the fluorescence-immunochromatography-based dry reagents manufactured by Beijing Wantai Biological Pharmacy Enterprise Co., Ltd.

## **3. Product Name and Model**

Product Name: Fluorescence Immunoassay Analyser

Model: RTR-FS100

Software Name & Version: FIA Software, V1.0

## **4. Product Composition**

Fluorometer is composed of optical unit, mechanical unit, control unit, RF card module, thermal printer, display screen, power adapter and built-in software (Version No.: V1.0).

## **5. Work Principle**



The Fluorometer, based on the immunofluorescence chromatography and the principle of photoelectric sensing and collection, analyzes fluorescence reaction intensity to test the reaction concentration of strips. When a test strip is sent into the test bin, LED lamp will give off spectrum, exciting the strip to fluoresce, collected by the sensor system, and test results will be given after integral logic analysis of the central control system.

## **6. Matters Needing Attention**

- 1) Do not operate under excessive temperature or humidity, avoiding high-intensity magnetic field, strong vibration, serious impact and highly corrosive gases. Please leave enough space around the Fluorometer to ensure ventilation.
- 2) Please protect the Fluorometer from liquids, dust or objects.
- 3) Do not wipe the screen and buttons with organic solvent.
- 4) Do not use the Fluorometer at the places, where inflammable materials are stored, to prevent explosion.

- 5) Do not disassemble the Fluorometer to replace any part inside without the written authorization from Xiamen United Medical Instruments Co., Ltd. or sales representative.
- 6) The Fluorometer should use the exclusive power adapter provided by Xiamen United Medical Instruments Co., Ltd. only.
- 7) The Fluorometer should be well grounded in the process of using, and well grounded socket should be chosen for power connection.
- 8) Enough space should be left around the Fluorometer, greater than 10cm. As display screen and printer are set at the top, enough space should be left for convenient operation.

### 7. Graphic Illustrations to Product Appearance and Parts

- ① Power Switch                      ② Power Outlet                      ③ RS232 Port
- ④ USB Port                              ⑤ Network Port

## 8. Product Type and Technical Parameters

### 8.1 Product type

Machine Type: Mini Desktop

Protection against electric shock: Transient Overvoltage Class II

Pollution Degree: Grade 2

### 8.2 Normal work conditions:

Ambient temperature: 5°C-40°C;

Relative humidity: ≤80% (no condensate water);

Maximum altitude: 2,000m;

Power supply: AC100-240V 50/60Hz;

Placed: on dry, clean, flat horizontal surface, free from direct sunlight, mechanical vibration and corrosive gas, in well-ventilated room.

### 8.3. Technical parameters

Size: 300mm\*280mm\*160mm

Weight: ≤4.5kg

Adapter power supply: AC100V-240V 50/60Hz

Main unit power supply: 24V

Main unit power: 35VA

Repeatability: CV (coefficient of variation) ≤2% for the repeated measurements of standard cards.

### 8.4. Performance Indicators

Performance Indicators of Fluorometer are shown in the table.

Inspection Item	Performance Indicators
Reading of blank strip (negative)	<0.02
Reading repeatability (coefficient of variation, CV)	CV≤2%
LOD (Limit of Detection)	T value ≤2000
Measurement time	≤15s
Abnormal result prompt	the Fluorometer will give the prompt as "Item Error, please retest!"
Stability	Difference between 0h, 4h and 8h is ≤15%
Accuracy	relative difference ≤15%
Linear dependence	(   r   ) ≥0.99.

## 8.5. Testing methods

Fluorometer is tested with the following methods.

Inspection Item	Testing method
Reading of blank strip (negative)	Test blank Strip for 10 times, take the mean value as test result.
Reading repeatability (coefficient of variation, CV)	Test CV Strips for 10 times, take the mean value as test result, then calculate the CV of each
LOD (Limit of Detection)	Test LOD Strip for 10 times, and take the mean value of the T line.
Measurement time	Test CV Strip for 10 times, and count the total test time from clicking "Read" to displaying the results.
Abnormal result prompt	Test Abnormal Strip
Stability	Test CV Strips for 3 times at 0h, 4h and 8h, take the mean value as test result, then calculate the relative differences among different times.
Accuracy	Test CV Strips for 3 times, take the mean value as test result, calculate the relative deviation between measured value and theoretical values.
Linear dependence	Test linear Strips (CRP item) for 3 times, take the mean results, then calculate the correlation coefficient (r) between measured results and reference background.

## 8.6. Calculation of examination result

The software will start the result settlement after reading the fluorescence value of each reaction strip (such as T line, C line).

### 8.6.1 Qualitative items

- ① Cutoff value setting. Set the cutoff value of reagent according to T/C.
- ② Then import the lot information and cutoff value settings into the fluorometer.
- ③ The fluorometer will calculate the result of T/C according to the T and C value.
- ④ The sample is judged to be positive when  $T/C \geq \text{cutoff value}$ , and is judged to be negative

when  $T/C < \text{cutoff value}$ .

### 8.6.2 Quantitative items

- ① Formulate the curve equation based on the T/C value of reference product.
- ② Then import the lot information and curve equation into the fluorometer.
- ③ The fluorometer will calculate the result of T/C according to the T and C value.
- ④ The curve equation will calculate the protein content according to the result of T/C.

## **8.7. Literature references**

GB 4793.1-2007: Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use - Part 1: General Requirements (IEC 61010-1:2001 IDT)

GB 4793.9: Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use - Part 9: Particular Requirements for Automatic and Semi-automatic Laboratory Equipment for Analysis and Other Purposes (IEC 61010-2-081:2009, IDT)

YY 0648-2008: Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use - Part 2-101: Particular Requirements for in vitro Diagnostic (IVD) Medical Equipment (IEC 61010-2-101:2002, IDT)

GB/T 18268.1: Electrical Equipment for Measurement, Control and Laboratory Use - EMC Requirements - Part 1: General Requirements (IEC 61326-1:2005, IDT)

GB/T 18268.26: Electrical Equipment for Measurement, Control and Laboratory Use - EMC Requirements - Part 26: Particular Requirements - In vitro Diagnostic (IVD) Medical Equipment (IEC 61326-2-6:2005, IDT)

## **9. Product Installation**

### **9.1 Inspection**

Before unpacking, please check the packing box for damage and deformation caused during transportation.

Please unpack with care and check according to the packing list; then check whether the machine is intact and whether the symbols, marks and labels are clear and intact. Please contact Xiamen United Medical Instruments Co., Ltd. or local agent in case of any defects, referring to Appendix A for the packing list.

NOTE: please keep the packing box and materials properly for return transportation in the future.

### **9.2 Installation**

1) Power supply: voltage AC100V-240V, frequency 50Hz-60HZ, 2A+ three-pin socket. Please plug the power line of power grid into the power port of adapter, and then plug the DC line of adapter into the power port of Fluorometer.

2) Put the Fluorometer on a flat and stable platform, leaving around 20cm space at the back for the operation of power switch.

3) Add printing paper, referring to Appendix B.

## 10. Operating Instructions and Graphics

Operating requirements:

1. Operator should have some knowledge of biomedical testing and get familiar with the operation of Fluorometer after training. Please wear latex gloves while processing the samples and reagents to protect the samples or reagents from contacting skin.

2. In case of anything abnormal to the Fluorometer while using, please stop immediately and then restart to confirm. If it still doesn't work, please contact the after-sales maintenance personnel for confirmation.


3. When the Fluorometer is running, the door on front cover is open. Please pay attention not to put anything at the front cover. Do not put your hand or anything else into the door to avoid injury.

4. After testing, please discard the wastes, such as used consumables and samples, into the designated container for medical wastes.

### 10.1 User login

Login interface is shown as follows:

Operating process is presented below:

Step 1: Double click the input field to get soft keyboard. Both Chinese and English are available and you can press “EN” to switch between.

Step 2: Enter your account and password. account: admin, password: admin. After typing, press “Enter” to hide the soft keyboard and click “Login” to enter the main interface, as shown below:

## 10.2 Test Function

Click “Test Function” and the following interface will pop up to enter the number, name, sex, age and sample type.

No.	Name	Sex	Age	Type	Lot.	Item	Result	Unit	Ref
1		-		-					

## 10.2.1 Immediate Test


### 10.2.1.1 For a used Batch

Step 1:


Import the item information (refer to the section on item import for the method of importing item information)

Step 2:

Insert reagent strip into the Fluorometer and click “Read”. If no item information is found, the following box will pop up:

In case of abnormal reagent strip, the following prompt box will pop up:


### 10.2.1.2 For A new Batch.

Step 1:

Import the item information (refer to the section on item import for the method of importing item information)

Step 2:

Click “Lot.” to open pull-down menu and select item interface.

Choose the item name and item batch, then click “Confirm” button to pick this item batch.

Step 3:


Insert reagent strip into the Fluorometer and click “Read”.

### 10.2.2 Timing Test

Operating process is presented below:

Insert reagent strip into the Fluorometer and click “Timing Test”, a countdown box will pop up, as shown below. Countdown time can be set according to the item information.

### 10.2.3 In/Out

Click “In/Out” button to get in and out of the test bin.

### 10.2.4 Upload

Upload the test results to LIS system (if LIS system connected).

### 10.2.5 Print

Print the test results.

## 10.3 Record Search


Step 1:

Click the “Return” button to main interface. Click “Record” to enter the record search interface, and records can be searched according to “No.”, “Item”, “Lot.”, “Name” and “Sex”.


Step 2:

Enter name and click “Search”. If no record is found, the following prompt box will pop up:


Step 3:

Select the required data.

Double click the search results (as shown in the red box above), the result box will pop up and click the “Print” button to print the information.

Step 4:

Select the required data. Export the results to USB flash disk.

Click the “Export Results” button. If no USB flash disk is available, a prompt box will pop up as follows:

After plugging the USB flash disk, click the “Export Results” button and it shows:

#### 10.4 Setting

Please contact the after-sales maintenance personnel if need.

#### 10.5 System

Please contact the after-sales maintenance personnel if need.


## 10.6 Item import

The IC card for item import is provided with the matching reagent. Item name and batch number can be imported into the Fluorometer with this IC card. After logging in the test interface, place IC card at the right side of the Fluorometer, as shown below, and information import can be finished by swiping the card.

After item import completed, the newly imported item can be chosen with the “Lot.” input module on test interface.

After item information being imported, the following information will be given:

## 11. Electromagnetic Compatibility

### CAUTION

- This device meets the emission and immunity requirements as stipulated in GB/T 18268.26, referring to the table below.
- User is obliged to ensure electromagnetic compatibility, for the normal operation of the device.
- We suggest users to evaluate the electromagnetic environment before using the device.

### CAUTION

- Design and test of the device follow the Class A equipment as stipulated in GB 4824. It may arouse radio interference in domestic situations. Please take some protective measures.
- Do not use the device by the side of severe radiation source (such as unshielded RF source), or it may interfere with normal operation of the device.
- The device is intended to be used in industrial and medical sites. We suggest users to evaluate electromagnetic compatibility before its application in other situations.

**Table 1:**

Electromagnetic Emission	
Emission Test	Conformance
GB 4824 RF Emission	Group 1
GB 4824 RF Emission	Class A
GB 17625.1 Harmonic Emission	Not Applicable
GB 17625.2 Voltage Fluctuation/Scintillation Emission	Not Applicable


**Table 2:**

Electromagnetic Immunity			
Immunity Test	Basic Standard	Test Value	Conformance
Electrostatic Discharge (ESD)	GB/T 17626.2	Contact Discharge: $\pm 2\text{kV}$ , $\pm 4\text{kV}$ Air Discharge: $\pm 2\text{kV}$ , $\pm 4\text{kV}$ , $\pm 8\text{kV}$	B
RF Electromagnetic	GB/T 17626.3	3V/m, 80MHz~2.0GHz, 80%AM	A

Field			
Pulse Train	GB/T 17626.4	Power Line: $\pm 1\text{kV}(5/50\text{ns}, 5\text{kHz})$	B
Surge	GB/T 17626.5	Line-to-Ground: $\pm 2\text{kV}$ Line-to-Line: $\pm 1\text{kV}$	B
RF Transmission	GB/T 17626.6	Power Line: $3\text{V/m}$ , $150\text{kHz}\sim 80\text{MHz}$ , $80\%\text{AM}$	A
Power Frequency Magnetic Field (PFMF)	GB/T 17626.8	$3\text{A/m}$ , $50/60\text{Hz}$	A
Voltage Sags & Interruption	GB/T 17626.11	1 cycle 0%; 5/6 cycle 40%; 25/30 cycle 70%; 250/300 cycle 5%.	B C C C
<p>Performance Discrimination:</p> <p>A. While testing, it presents normal performance within the limits.</p> <p>B. While testing, it suffers temporary function or performance drop or loss, but could recover automatically.</p> <p>B. While testing, it suffers temporary function or performance drop or loss, but requires for operator intervention or system reset.</p>			

## 12. Explanation to Label-included Marks

Table 1 Explanation to Label-included Marks

Symbol	Notes
	Direct Current
SN	Serial Number
	CAUTION (referring to attachments)
	In Vitro Diagnostic Medical Device
	Date of Manufacture
(ON)	(Power) ON
○ (OFF)	(Power) OFF
	Biorisk
	Keep Upright

	<p>Handle With Care</p>
	<p>Caution Against Wet</p>
	<p>Storage Conditions</p>
	<p>Stacking Limit</p>

### 13. Maintenance

Outer surface of the device can be cleaned with the cleaning cloth dampened with 70% alcohol.

Please use dehaired gauze to clean the edges and corners of the device.

Please cover the device with dust cap when it is not used.

WARNING: After unpacking the device, Xiamen United Medical Instruments Co., Ltd. assumes no responsibilities if any part is damaged while cleaning.

### 14. Transportation, Storage and Scrap

#### Transportation

Transportation of the Fluorometer should follow the order contract signed. During transportation, please place the product as the marks on packing box and avoid collision, violent vibration or direct rain and snow.

#### Storage

The packed system should be stored in well-ventilated room free from corrosive gas, under -20°C~+55°C and humidity not greater than 93%.

#### Disposal

We suggest users to follow the rules for Class B electronic instruments, concerned with disposal of the Fluorometer.

### 15. Troubleshooting and Statement

Fault	Handling Methods
After powering on, LCD screen displays nothing and printer gives no green light.	Check the power for good connection and check the power port for looseness.
	Check the power switch for damage.
	Check the adapter for damage.
Program crashes and doesn't response.	Power off and then power on again.
Test bin doesn't work or is clamped while IN/OUT.	
Printer doesn't print or eat paper, or print nothing.	Press the "SEL" key, the green light is on, and the printer is online.
	Check the printing paper for wrong installation direction.
	Pull out the printing paper for a little.
Data backup failed	Check the USB flash disk for damage.

We hereby declare that normal operation of this product can be ensured only when you totally follow the Product Instruction published by manufacturer. We assume no liabilities for the indirect or corresponding damages in addition.

## **16. Contraindication: None**

## **17. After-sales Services**

### **17.1 Warranty Coverage**

Service life of the product is five years, warranty period of one year included. While ensuring to abide by the rules for transportation, storage and application, once the product cannot operate normally for quality problems, the manufacturer will repair or replace the product or parts free from charge within 12 months since delivery.

The following conditions are not included in warranty coverage.

1. The product is damaged for being used not according to the Product Instruction or the general rules for laboratory instruments;
2. The product is damaged due to disassembly without permission;

**Xiamen United Medical Instruments Co., Ltd.**  
Warranty Card

Product Name: \_\_\_\_\_  
Model: \_\_\_\_\_  
Product SN: \_\_\_\_\_  
Manufacture Date: \_\_\_\_\_

Purchaser: \_\_\_\_\_  
Purchase Date: \_\_\_\_\_  
Address: \_\_\_\_\_  
Quality Assurance Date: \_\_\_\_\_  
Warranty Card No.: \_\_\_\_\_

Figure of Warranty Card

## 17.2 After-sales Service Hotline and Address

Company Name: Xiamen United Medical Instruments Co., Ltd.

Residence: 5#, No.2008, Wengjiao West Rd, Haicang District, 361026 Xiamen, Fujian Province, PEOPLE'S REPUBLIC OF CHINA

Production Address: 5#, No.2008, Wengjiao West Rd, Haicang District, 361026 Xiamen, Fujian Province, PEOPLE'S REPUBLIC OF CHINA

Zip Code: 361022

Customer Service Hotline: 400-068-0279

E-mail: [umic@xmumic.com](mailto:umic@xmumic.com)

Website: <http://www.xmumic.com>

## Appendix A Packing List

### A.1 Packing List

Table A.1 Packing List

SN	Item	Quantity
1	Fluorescence Immunoassay Analyser	1 SET
2	Power Adapter/Power Line	1 SET
3	Product Instruction	1 COPY
4	Printing Paper	1 ROLL
5	Quality Certificate	1 COPY
6	Warranty Card	1 COPY

CAUTION: In case of any part missing or damage while unpacking the product, please contact Xiamen United Medical Instruments Co., Ltd. or local sales representative in time.

## Appendix B Instructions for Printer

### B.1 Control panel:

LED indicator:

After powering ON, if the indicator flickers twice, with interval of 1s,

it starts normally and the indicator may:

Flicker once: normal operation after self-check.

Flicker twice: printer not found.

Flicker 3 times: no printing paper.

Flicker 5 times: heater band of printer is overheated.

Flicker 10 times: chip of Chinese fonts not found.

Key pressed: feeding paper.

Open the top cover of printer

### B.2 Add printing paper:

1. Open the top cover of printer;
2. Put the paper roll in printer with reverse side up;
3. Close the top cover and leave a little printing paper outside.

NOTE: Before adding the printing paper, please tear off the gummed paper on the roll to protect print head.

CE

# **Fluoreszenz-Immunoassay-Analysator**

## **Produktanleitung**

Modell: RTR-FS100

---

**Xiamen United Medical Instruments Co., Ltd.**

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Xiamen United Medical Instruments Co., Ltd.

5#, No.2008, Wengjiao West Rd, Haicang District, 361026 Xiamen,

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Name des Unternehmens: Xiamen United Medical Instruments Co., Ltd.