



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



Datasheet

RealBlue Peroxidase Substrate Kit

Catalog Number: KA0565

Regulatory Status: For research use only (RUO)

Product Description: RealBlue Peroxidase Substrate Kit is 50-100 times more sensitive than DAB and is non-carcinogenic. Strongly recommended for the detection of least-abundant markers. RealBlue forms a blue chromogenic product, providing excellent contrast with DAB and other substrates for multiple labeling experiments.

Applications: IHC

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Instruction: Store the kit at room temperature.

References:

1. Origin and development of septoclasts in endochondral ossification of mice. Bando Y, Sakashita H, Taira F, Miyake G, Ogasawara Y, Sakiyama K, Owada Y, Amano O. Histochem Cell Biol. 2018 Feb 20. [Epub ahead of print]