



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### Live Cell Endoplasmic Reticulum Staining Kit (Blue Fluorescence)

**Catalog Number:** KA6033

**Regulation Status:** For research use only (RUO)

**Product Description:** Live Cell Endoplasmic Reticulum Staining Kit (Blue Fluorescence) uses the ER Tracer Blue as an ER marker. ER Tracer Blue stain is a cell-permeant fluorescent dye that is highly selective for ER. This stain consists of a blue fluorescent dye and ER binder that selectively bind to ER in most of cell types.

**Applications:** Func

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Instruction:** Store the kit at -20°C and avoid from light.