



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

Rapid Colorimetric APC-Maleimide Quantitation Kit

Catalog Number: KA7236

Regulation Status: For research use only (RUO)

Product Description: Rapid Colorimetric APC-Maleimide Quantitation Kit facilitates the quantification of maleimide groups on maleimide-activated APC proteins, which has an absorption peak at approximately 555 nm.

Absorbance (nm): 900 to 250

Regulation Status : For research use only (RUO)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Storage Instruction: Store in the dark at -20°C. Avoid
prolonged exposure to light.
Avoid repeated freezing and thawing.