



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

Inosine Assay Kit

Catalog Number: KA7270

Regulation Status: For research use only (RUO)

Product Description: Inosine Assay Kit uses a single Working Reagent that combines the enzyme reactions and color reaction in one step. The change in fluorescence intensity at $\lambda_{ex/em} = 530/585$ nm is directly proportional to inosine concentration in the sample.

Calibration Range: 1 to 25 μ M

Limit of Detection: 1 μ M

Suitable Sample : Biological samples

Sample Volume : 20 μ L

Detection Method : FL530/585nm

Regulation Status : For research use only (RUO)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Storage Instruction: Store at -20°C.
Avoid repeated freezing and thawing.